

## Contents

1. Description
  - 1.1 Principle of the MACS<sup>®</sup> Separation
  - 1.2 Background information
  - 1.3 Applications
  - 1.4 Reagent and instrument requirements
2. Protocol
  - 2.1 Sample preparation
  - 2.2 Magnetic labeling of non-B cells
  - 2.3 Magnetic separation: Depletion of non-B cells
  - 2.4 Magnetic labeling of IgG<sup>+</sup> memory B cells
  - 2.5 Magnetic separation: Positive selection of IgG<sup>+</sup> memory B cells
3. Example of a separation using the IgG<sup>+</sup> Memory B Cell Isolation Kit

## Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

## 1. Description

This product is for research use only.

<b>Components</b>	<p><b>2 mL Memory B Cell Biotin-Antibody Cocktail, human:</b></p> <p>Cocktail of biotin-conjugated monoclonal antibodies against CD2, CD14, CD16, CD36, CD43, and CD235a (Glycophorin A).</p> <p><b>2×2 mL Anti-Biotin MicroBeads:</b></p> <p>MicroBeads conjugated to monoclonal anti-biotin antibodies (isotype: mouse IgG1).</p> <p><b>1 mL Anti-IgG MicroBeads, human:</b></p> <p>MicroBeads conjugated to monoclonal anti-human IgG antibodies (isotype: mouse IgG1).</p>
<b>Capacity</b>	For 2×10 <sup>9</sup> total cells.
<b>Product format</b>	All components are supplied in buffer containing stabilizer and 0.05% sodium azide.
<b>Storage</b>	Store protected from light at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label.

### 1.1 Principle of the MACS<sup>®</sup> Separation

The isolation of IgG<sup>+</sup> memory B cells is performed in a two-step procedure. First, the non-B cells are indirectly magnetically labeled with a cocktail of biotin-conjugated antibodies, as primary labeling reagent, and Anti-Biotin MicroBeads, as secondary labeling reagent. The labeled cells are subsequently depleted by separation over a MACS<sup>®</sup> Column, which is placed in the magnetic field of a MACS Separator.

In the second step, the IgG<sup>+</sup> memory B cells are directly labeled with Anti-IgG MicroBeads and isolated by positive selection from the pre-enriched B cell fraction by separation over a MACS Column, which is placed in the magnetic field of a MACS Separator. After removing the column from the magnetic field, the magnetically retained IgG<sup>+</sup> memory B cells can be eluted as the positively selected cell fraction.

#### Human PBMCs: Depletion of non-B cells

1. Indirect magnetic labeling of non-B cells with Memory B Cell Biotin-Antibody Cocktail and Anti-Biotin MicroBeads.
2. Magnetic separation using an LD Column or an autoMACS Column (program "Depl025").

#### Pre-enriched B cells (flow-through fraction): Positive selection of IgG<sup>+</sup> memory B cells

1. Direct magnetic labeling of IgG<sup>+</sup> memory B cells with Anti-IgG MicroBeads.
2. Magnetic separation using two MS Columns or an autoMACS Column (program "Posseld2").

#### IgG<sup>+</sup> memory B cells

### 1.2 Background information

The IgG<sup>+</sup> Memory B Cell Isolation Kit has been developed for the isolation of IgG<sup>+</sup> memory B cells from human peripheral blood mononuclear cells (PBMCs). Memory B cells, defined as quiescent antigen-experienced B cells, are generated in response to T cell-dependent and T cell-independent antigens. They are able to react quickly to a recurrent antigenic challenge thereby providing serological immune protection.

IgG<sup>+</sup> memory B cells are isolated by depletion of unwanted non-B cells and subsequent positive selection with Anti-IgG MicroBeads. Unwanted cells, for example, T cells, NK cells, monocytes, dendritic cells, granulocytes, platelets, and erythroid cells are depleted using a cocktail of biotinylated antibodies against CD2, CD14, CD16, CD36, CD43, and CD235a (Glycophorin A), and Anti-Biotin MicroBeads. For evaluation of MACS Separations, staining with Anti-IgG-APC and CD19-PE antibodies is recommended after separation.

### 1.3 Applications

- Isolation of IgG<sup>+</sup> memory B cells from human PBMCs for further phenotypical or functional characterization.

### 1.4 Reagent and instrument requirements

- Buffer:** Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS<sup>®</sup> BSA Stock Solution (# 130-091-376) 1:20 with autoMACS<sup>®</sup> Rinsing Solution (# 130-091-222). Keep buffer cold (2–8 °C). Degas buffer before use, as air bubbles could block the column.  
**▲ Note:** EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). Buffers or media containing Ca<sup>2+</sup> or Mg<sup>2+</sup> are not recommended for use.
- MACS Columns and MACS Separators:** Depletion of non-B cells is performed on an LD Column. The subsequent positive selection of IgG<sup>+</sup> memory B cells is performed on two MS Columns. Depletion and positive selection can also be performed by using the autoMACS Pro Separator.

Column	Max. number of labeled cells	Max. number of total cells	Separator
<b>Depletion</b>			
LD	10 <sup>8</sup>	5×10 <sup>8</sup>	MidiMACS, QuadroMACS, VarioMACS, SuperMACS II
<b>Positive selection</b>			
MS	10 <sup>7</sup>	2×10 <sup>8</sup>	MiniMACS, OctoMACS, VarioMACS, SuperMACS II
<b>Positive selection or depletion</b>			
autoMACS	2×10 <sup>8</sup>	4×10 <sup>9</sup>	autoMACS Pro

**▲ Note:** Column adapters are required to insert certain columns into the VarioMACS<sup>™</sup> or SuperMACS<sup>™</sup> II Separators. For details refer to the respective MACS Separator data sheet.

- (Optional) Fluorochrome-conjugated antibodies for flow cytometric analysis, e.g., CD19-PE and Anti-IgG-APC. For more information about antibodies refer to [www.miltenyibiotec.com/antibodies](http://www.miltenyibiotec.com/antibodies).
- (Optional) Propidium Iodide Solution (# 130-093-233) or 7-AAD for flow cytometric exclusion of dead cells.
- (Optional) Dead Cell Removal Kit (# 130-090-101) for the depletion of dead cells.
- (Optional) Pre-Separation Filters (30 µm) (# 130-041-407) to remove cell clumps.

## 2. Protocol

### 2.1 Sample preparation

When working with anticoagulated peripheral blood or buffy coat, peripheral blood mononuclear cells (PBMCs) should be isolated by density gradient centrifugation, for example, using Ficoll-Paque<sup>™</sup>.

**▲ Note:** To remove platelets after density gradient separation, resuspend cell pellet in buffer and centrifuge at 200×g for 10–15 minutes at 20 °C. Carefully aspirate supernatant. Repeat washing step.

When working with tissues or lysed blood, prepare a single-cell suspension using standard methods.

For details refer to the protocols section at [www.miltenyibiotec.com/protocols](http://www.miltenyibiotec.com/protocols).

**▲ Dead cells may bind non-specifically to MACS MicroBeads.** To remove dead cells, we recommend using density gradient centrifugation or the Dead Cell Removal Kit (# 130-090-101).



### 2.2 Magnetic labeling of non-B cells

**▲ Work fast, keep cells cold, and use pre-cooled solutions.** This will prevent capping of antibodies on the cell surface and non-specific cell labeling.

**▲ Volumes for magnetic labeling given below are for up to 10<sup>8</sup> total cells.** When working with fewer than 10<sup>8</sup> cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10<sup>8</sup> total cells, use twice the volume of all indicated reagent volumes and total volumes).

**▲ For optimal performance it is important to obtain a single-cell suspension before magnetic labeling.** Pass cells through 30 µm nylon mesh (Pre-Separation Filters (30 µm), # 130-041-407) to remove cell clumps which may clog the column. Moisten filter with buffer before use.

**▲ The recommended incubation temperature is 2–8 °C.** Higher temperatures and/or longer incubation times may lead to non-specific cell labeling. Working on ice may require increased incubation times.

- Determine cell number.
- Centrifuge cell suspension at 300×g for 10 minutes. Aspirate supernatant completely.
- Resuspend cell pellet in 400 µL of cold buffer per 10<sup>8</sup> total cells.
- Add 100 µL of Memory B Cell Biotin-Antibody Cocktail per 10<sup>8</sup> total cells.
- Mix well and incubate for 10 minutes in the refrigerator (2–8 °C).
- Add 300 µL of cold buffer and 200 µL of Anti-Biotin MicroBeads per 10<sup>8</sup> total cells.
- Mix well and incubate for additional 15 minutes in the refrigerator (2–8 °C).
- Wash cells by adding 10 mL of buffer per 10<sup>8</sup> cells and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
- Resuspend up to 10<sup>8</sup> cells in 1 mL of cold buffer.  
**▲ Note:** For higher cell numbers, scale up buffer volume accordingly.
- Proceed to magnetic separation (2.3).



### 2.3 Magnetic separation: Depletion of non-B cells

▲ Choose an appropriate MACS Column and MACS Separator according to the number of total cells and the number of IgG<sup>+</sup> memory B cells. For details refer to table in section 1.4.

▲ Always wait until the column reservoir is empty before proceeding to the next step.

#### Depletion with LD Columns

1. Place LD Column in the magnetic field of a suitable MACS Separator. For details refer to LD Column data sheet.
2. Prepare column by rinsing with 2 mL of buffer.
3. Apply cell suspension onto the column.
4. Collect unlabeled cells that pass through and wash column with 2×1 mL of buffer. Collect total effluent; this is the unlabeled pre-enriched memory B cell fraction. Perform washing steps by adding buffer two times. Only add new buffer when the column reservoir is empty.
5. Proceed to 2.4 for the labeling of IgG<sup>+</sup> memory B cells.

#### Depletion with the autoMACS<sup>®</sup> Pro Separator

▲ Refer to the respective user manual for instructions on how to use the autoMACS<sup>®</sup> Pro Separator.

▲ Buffers used for operating the autoMACS Pro Separator should have a temperature of ≥10 °C.

1. Prepare and prime the instrument.
2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample and collection tubes into the Chill Rack.
3. For a standard separation choose the following program:

#### Depletion: Depl025

Collect negative fraction in row B of the tube rack.

4. Proceed to 2.4 for the labeling of IgG<sup>+</sup> memory B cells.



### 2.4 Magnetic labeling of IgG<sup>+</sup> memory B cells

▲ Volumes for magnetic labeling given below are for an initial starting cell number of up to 10<sup>8</sup> total cells. For higher initial cell numbers, scale up all volumes accordingly.

1. Centrifuge cell suspension at 300×g for 10 minutes. Aspirate supernatant completely.
2. Resuspend cell pellet in 50 µL of buffer.
3. Add 50 µL of Anti-IgG MicroBeads.
4. Mix well and incubate for 15 minutes in the dark in the refrigerator (2–8 °C).
5. Wash cells by adding 10 mL of buffer and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.

6. Resuspend up to 10<sup>8</sup> cells in 500 µL of buffer.  
▲ **Note:** For higher cell numbers, scale up buffer volume accordingly.
7. Proceed to magnetic separation (2.5).



### 2.5 Magnetic separation: Positive selection of IgG<sup>+</sup> memory B cells

#### Positive selection with MS Columns

▲ To achieve highest purities, perform two consecutive column runs.

▲ Always wait until the column reservoir is empty before proceeding to the next step.

1. Place MS Column in the magnetic field of a suitable MACS Separator. For details refer to MS Column data sheet.
2. Prepare column by rinsing with 500 µL of buffer.
3. Apply cell suspension onto the column. Collect flow-through containing unlabeled cells.
4. Wash column with 3×500 µL of buffer. Collect unlabeled cells that pass through and combine with the effluent from step 3.  
▲ **Note:** Perform washing steps by adding buffer aliquots as soon as the column reservoir is empty.
5. Remove column from the separator and place it on a suitable collection tube.  
▲ **Note:** To perform a second column run, you may elute the cells directly from the first onto the second, equilibrated column instead of a collection tube.
6. Pipette 1 mL of buffer onto the column. Immediately flush out the magnetically labeled cells (IgG<sup>+</sup> memory B cells) by firmly pushing the plunger into the column.
7. To increase purity of IgG<sup>+</sup> memory B cells, the eluted fraction must be enriched over a second MS Column. Repeat the magnetic separation procedure as described in steps 1 to 6 by using a new column.

#### Positive selection with the autoMACS<sup>®</sup> Pro Separator


1. Prepare and prime the instrument.
2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample and collection tubes into the Chill Rack.
3. For a standard separation choose the following program:

#### Positive selection: Posseld2

Collect positive fraction in row C of the tube rack. This is the enriched IgG<sup>+</sup> memory B cell fraction.

### 3. Example of a separation using the IgG<sup>+</sup> Memory B Cell Isolation Kit

IgG<sup>+</sup> memory B cells were isolated from human PBMCs using the IgG<sup>+</sup> Memory B Cell Isolation Kit, an LD Column and a MidiMACS™ Separator, two MS Columns and a MiniMACS™ Separator. Cells are fluorescently stained with Anti-IgG-APC and CD19-PE. Cell debris and dead cells are excluded from the analysis based on scatter signals and propidium iodide fluorescence.



Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

### Legal notices

#### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

#### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

#### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

#### Trademarks

autoMACS, MACS, MidiMACS, the Miltenyi Biotec logo, MiniMACS, OctoMACS, QuadroMACS, SuperMACS, and VarioMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Ficoll-Paque is a trademark of GE Healthcare companies.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.