
Miltenyi Biotec

MultiMACS™ M96 Separator

User manual

Copyright © 2022 Miltenyi Biotec and/or its affiliates. All rights reserved.

No part of this user manual may be reproduced, stored in a retrieval system, transmitted, published, or distributed in any form or by any means, electronically, mechanically, by photocopying, microfilming, recording, or otherwise, without the prior written consent of Miltenyi Biotec; however, notwithstanding the foregoing, the owners of the MultiMACS™ M96 Separator may make copies solely for purposes of training personnel in the use and servicing of the unit within their business or organization. Maximal care has been taken by Miltenyi Biotec in the preparation of this user manual. However, Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this document. The information in this document is provided "as is" without warranty of any kind, and is subject to change without notice. Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. μ MACS, MACS, the Miltenyi Biotec logo, and MultiMACS are either registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

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	86	Keeping flow-through and wash fractions: 96 NEG/POS
	92	96 POS+BEAD
	96	Creating a new user-defined process
	96	Selection of parameter types
	96	A. Selecting Fraction of Interest
	97	B. Selecting Elution Mode
	98	C. Selecting Plate Heights
	100	D. Advanced parameter types
	102	Saving the created process
	103	Running a user-defined process
	104	Viewing process parameters
	105	Editing a process
	107	A. Editing Fraction of Interest
	108	B. Editing Elution Mode
	109	C. Editing Plate Height
	112	D. Advanced parameter types
	114	Saving the edited process
	115	Testing Plate Height or Distance for Plate Removal
	116	Replacing a user-defined process
	117	Deleting a user-defined process
	119	Service and preferences
	119	Firmware version
	120	Preparing the MultiMACS M96 Separator for transport
	122	Rotating MultiMACS 96 Magnet and Column Holder
Chapter 7	127	Maintenance
Chapter 8	129	Troubleshooting
	129	Cross-contamination by drops at the column tips
	130	Collision detection
	131	Buffer temperature
	132	Malfunctions and possible remedies
	132	Error code s
	133	Fuse exchange
Chapter 9	135	Technical support
Chapter 10	137	Legal notes
	137	Limited warranty
	139	Trademarks

Please read all information contained in this user manual before use. Failure to read and follow these guidelines could lead to improper or incorrect use, handling or care of your instrument and could cause hazards to users, unpredictable results, device malfunction or damage, premature wear and reduced life time of the instrument, and may void your warranty.


Keep this user manual in a safe place, accessible for anyone using the MultiMACS M96 Separator.

This chapter describes the safety instructions and site requirements for your MultiMACS M96 Separator. The following warnings and cautions are provided to help you prevent injury to yourself or damage to the device.

Symbols and hazard levels

Setup of safety notices

Example

The safety notices inform the user about potential risks if warnings and precautions outlined below are not followed. The icon on the left side specifies the risk. The hazard level at the top classifies the hazard, as mentioned below. The level, type, and source of the hazard as well as potential consequences, prohibitions, and measures are pointed as follows.

Symbols and hazard levels

The following chart is an illustrated glossary depicting the symbols that are used in this user manual and on the MultiMACS M96 Separator.

Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury.

Indicates a hazardous situation, which, if not avoided, could result in death or serious injury.

Safety warning: The documentation must be consulted in all cases where this safety symbol is used on the instrument.

Hazard of crushing and shearing.
Risk of crushing and shearing of bodily parts due to mechanical hazards.

Strong magnetic field.
The magnetic field can interfere with magnetizable objects and electronic devices or damage magnetic information carriers. Risk of severe personal injury to persons carrying pacemakers or electronic medical implants.

Risk of contamination if biohazardous material is used. Indicates the risk of death, severe injury to the instrument operator, or equipment damage due to potentially dangerous biological material.

Indicates the risk of death or severe injury to the instrument operator due to hazardous voltage.

On (power).

Off (power).

Documentation must be consulted before proceeding with installation and operation of the system.

WEEE (Waste of Electrical and Electronic Equipment).

Fuse.

Persons carrying pacemakers or electronic medical implants must maintain distance.

Manufacturer.

Order number.

Type number.

Serial number.

NRTL certification mark: Product meets consensus-based standards of safety, required by the Occupational Safety/Health Administration (OSHA), determined by the Nationally Recognized Testing Laboratories (NRTL) TÜV SÜD.

European Conformity marking.

United Kingdom Conformity Assessed marking.

WARNING! Avoid serious injury. Read and understand operating instructions before installation or use. If manual is missing, obtain replacement from Miltenyi Biotec.

Warnings and precautions

The MultiMACS M96 Separator employs state-of-the-art technology. It is a computer-controlled device specifically developed for molecular biology applications in microtiter plate format. The MultiMACS M96 Separator is designed to operate safely after installation and when used by trained personnel according to general safety practices and the instructions set forth in this user manual. The guidelines in this section explain the potential risks associated with the operation of the instrument and provide important safety information in order to minimize these risks. By carefully following the instructions, you can protect yourself and the equipment from potential hazards and create a safe work environment. If this instrument is used in a manner not specified by the manufacturer, protection may be impaired.

Important: Please read and follow all operating instructions in this user manual and pay attention to all warnings displayed on the instrument. Retain this user manual and any other safety and operating instructions provided with the instrument in a place accessible to all users for future reference.

Important: The MultiMACS M96 Separator is intended for indoor use only. Do not use the instrument in areas classified as hazardous locations such as oxygen-laden environments.


Contact your local authority governing electrical power supply, building constructions, maintenance, or safety for more information regarding the installation of the equipment.

If you have a serious concern regarding the safe use of your instrument, contact your authorized Miltenyi Biotec service provider or call Miltenyi Biotec Customer Service.

General precautions

To reduce potential risks associated with operating the MultiMACS M96 Separator, please observe the following general precautions. Failure to observe these precautions could result in fire, bodily harm, and/or damage to the instrument.


Hazard of electric shock and spread of fire

Warning: Electrical devices pose the risk of an electric shock. To reduce the risk of an electric shock, do not open any cover of the MultiMACS M96 Separator nor any other accessory hardware supplied by Miltenyi Biotec. All covers of the device and accessory hardware are to be removed by authorized personnel only. Special care must be taken while handling fluids. Clean up spillages immediately. Do not allow fluids to enter the interior of the device. Unplug the power cable before manually cleaning the device.

A potential risk exists if an opened, dropped, or damaged MultiMACS M96 Separator is used, if liquids are spilled into the instrument, if an object has entered the instrument through the ventilation slots, or if an object has been dropped into the instrument. If flames or smoke appear immediately switch off the MultiMACS M96 Separator, unplug the instrument from the electrical outlet, and contact an authorized Miltenyi Biotec service provider or the Miltenyi Biotec Customer Service. Use of a damaged instrument or an instrument with a damaged power cable is expressly prohibited.


Strong magnetic field

Warning: The MultiMACS M96 Separator is equipped with an extremely powerful magnet in microtiter plate-format, the MultiMACS 96 Magnet. Keep any magnetic information carriers (such as credit cards, magnetic tapes, and storage media), any electronic equipment (such as hearing aids, pacemakers, measuring and control instruments, computers, and watches), and magnetizable tools and objects at a distance of at least 80 cm from the magnet cover. These items may be affected or damaged by the magnetic field.

The prism-shaped MultiMACS Display contains less powerful magnets; therefore, maintain a distance of at least 30 cm from the above listed magnetically sensitive material. The force of attraction between two magnets or between a magnet and magnetizable material increases strongly the closer they get. Thus, keep any magnetic material – including any other MultiMACS Magnet and Display – at a safe distance.

Hazard of crushing and shearing

Warning: Keep fingers etc. away from all moving parts of the MultiMACS M96 Separator to avoid crushing and shearing injuries or damage to the instrument.

If an object has been accidentally placed under the down-moving magnet, the physical resistance of the object causes the instrument to tilt backward, similar to a hydraulic jack. Do not place objects on the MultiMACS M96 Separator. Leave a space of at least 20 cm distance above and behind the MultiMACS M96 Separator to ensure that the instrument is able to tilt backwards, thus avoiding potential crushing and shearing injuries and minimizing instrument damage. Do not circumvent any safety device or measurement. Keep the proper respective Protective Screen attached to the MultiMACS M96 Separator when operating it manually.

Caution: For usage in a liquid-handling robot, removal of the Protective Screen may be necessary. When using the MultiMACS M96 Separator integrated into a compatible automated liquid-handling system, assure that the combined system, i.e., MultiMACS M96 Separator, MultiMACS 96 Magnet, and liquid-handling robot, meets all requirements of safety and function before starting automated processes. Thereby, particularly take the protection of the operator into account.

Position of safety notices on the instrument

Please notice the positions of the safety notices on the MultiMACS M96 Separator, and keep them in a completely readable state.

Secure installation

This section describes the requirements your site must meet for safe installation and operation of your MultiMACS M96 Separator. Read the instructions in this section and ensure that your site is properly prepared before you connect the instrument to its power source.

When planning your site layout and equipment locations, keep in mind the precautions described in this section to help avoid instrument failures and reduce the possibility of environmentally caused shutdowns.

Important: At all times, local working area safety instructions, laboratory policies, and standards regarding laboratory health and safety and prevention of accidents must be adhered to.

Mounting accessories

Do not place the MultiMACS M96 Separator on an unstable table, cart, stand, tripod, or bracket. As a consequence, the instrument might fall down. This may cause serious bodily harm and/or serious damage to the instrument. Use only on a table, cart, stand, tripod, or bracket sold with the instrument or recommended by Miltenyi Biotec. Do not place the MultiMACS M96 Separator within a built-in apparatus or a confined space such as a shelf rack unless the apparatus has been specifically designed to accommodate the instrument, proper ventilation is provided, and the mounting instructions for the instrument have been followed.

Air circulation

The instrument should not be placed next to radiators, heat registers, stoves, or other pieces of equipment (including amplifiers) that produce heat. Allow sufficient air circulation around the MultiMACS M96 Separator—at least 15 cm on all sides—during operation to ensure adequate cooling of the instrument. Prevent direct exposure of the instrument to sunlight. Slots and openings of the instrument are provided for ventilation and should never be blocked or covered, as these ensure reliable operation of the MultiMACS M96 Separator and protect the device from overheating. Never push a foreign object through an opening into the instrument.

Water and moisture

Do not use the instrument in a wet or damp location. Avoid high humidity or condensation and protect the instrument against water splashes.

Grounded (earthed) product

The instrument is equipped with a three-wire electrical grounding-type plug that has a third pin for grounding. This plug only fits into a grounded power outlet. This is a safety feature. Do not try to insert the plug into a non-grounded power outlet. If you cannot insert the plug into the outlet, contact your local electrician to replace the outlet.

Power sources

The instrument should only be operated from a power source indicated on the product's electrical ratings label. If you have questions about the type of power source to use, contact your authorized Miltenyi Biotec service provider or local power company. Do not use extension cords or power strips. Do not overload an electrical outlet. The overall system load must not exceed 80% of the branch circuit rating.

Accessibility

Make sure that the power switch as well as the connector for the power cable are easily accessible and located as close to the operator of the instrument as possible. If it is necessary to disconnect the power supply, unplug the cable from the power outlet.

Peripheral devices

Only connect peripheral devices that comply to EN/IEC 60950 to the RS232 connector labeled with “com”. In addition, only connect original MultiMACS Equipment to the connectors labeled “module 1” and “module 2”, and the MultiMACS Display to the connector labeled with “display”.

Secure operation, maintenance, transport, and disposal

Observe the following instructions to ensure secure operation, maintenance, transport, and disposal of your MultiMACS M96 Separator.

Important: At all times, local working area safety instructions, laboratory policies, and standards regarding laboratory health and safety and prevention of accidents must be adhered to.

Safe operation

If the instrument is not working properly and instructions or messages on the display screen advise to contact technical service, secure operation is no longer possible. Immediately switch off the MultiMACS M96 Separator, unplug the instrument from the power outlet, and contact an authorized Miltenyi Biotec service provider or Miltenyi Biotec Technical Support.

Servicing

Important: Unless otherwise specifically noted in this user manual or other Miltenyi Biotec documentation, do not service the MultiMACS M96 Separator yourself. Servicing and repair must be performed by qualified service personnel. Improper or incorrect servicing or repair of your MultiMACS M96 Separator can cause hazards to users, lead to unpredictable results, device malfunction or damage, premature wear and reduced life time of the instrument, and may void your warranty.

Inquire with your local Miltenyi Biotec representative about Miltenyi Biotec’s extensive instrument service and support arrangements, or refer to www.miltenyibiotec.com/support.

Important: When replacement or spare parts are required, make sure that the service provider uses only genuine Miltenyi Biotec parts or third-party parts specified and recommended by Miltenyi Biotec. Using unauthorized replacement or spare parts can cause malfunction of the device and impair results. Miltenyi Biotec does not honor any warranty or accept any responsibility for device failure or damages resulting from the use of inappropriate replacement or spare parts. After completing any service or repair work, have your authorized Miltenyi Biotec service provider perform all safety checks required by the repair procedure to ensure that the instrument is in proper operational condition.

Only use options and upgrades recommended by Miltenyi Biotec.

Cleaning

Unplug the MultiMACS M96 Separator from the power outlet before cleaning. Do not use liquid or aerosol cleaning agents; always use a damp cloth.

Hazardous material

If biohazardous material is or has been used, the operator shall choose and wear personal safety equipment in accordance with warnings and precautions for the used substances. Wear protective gloves, protective clothing, and safety glasses to prevent contact with skin and eyes. Defective or inadequate safety equipment might endanger the operator. The MultiMACS M96 Separator shall be operated in a safety hood if hazardous or unknown materials are processed. If hazardous material has been used or spilled, care must be taken to thoroughly decontaminate the system.

Columns, plates, tubes, and any other consumables that were in contact with biohazardous material shall be autoclaved prior to disposal. Liquid waste shall be autoclaved or decontaminated using a disinfectant that is appropriate for the specific pathogen, e.g., 10% bleach, isopropyl alcohol, or 70% ethanol.

Waste disposal must be in accordance with any local regulations.

Transport

The MultiMACS M96 Separator should be transported with care in packaging specified by Miltenyi Biotec. Internal damage can occur if it is subjected to excessive vibration or if it is dropped. If the instrument needs to be shipped back to the manufacturer for service, decontaminate the instrument from any hazardous material prior to shipment. If you have questions regarding proper decontamination or shipment, contact Miltenyi Biotec Technical Support for assistance.

Disposal

Waste of Electrical and Electronic Equipment (WEEE) customer information

Dispose of your end-of-life Miltenyi Biotec products in accordance with the applicable WEEE and hazardous waste disposal legislation, which may differ by country or region.

Electrical equipment may contain hazardous substances that may have a serious detrimental effect on the environment and/or human health. All equipment must be specifically collected and treated by designated waste facilities and by qualified WEEE compliance schemas. By ensuring that you dispose of your unwanted electrical and electronic equipment according to the applicable WEEE and hazardous waste disposal legislation, you are helping to preserve our natural resources and protect human health.

Miltenyi Biotec is committed to protecting the environment. Miltenyi Biotec offers product end-of-life return programs in many countries, and partners with licensed WEEE compliance schemes throughout the world. Miltenyi Biotec takes back your end-of-life Miltenyi Biotec equipment for recycling free of charge. The terms and availability of this offer vary by geography because of differences in regulatory requirements. Note that, depending on the type and use of your equipment, additional requirements may apply. Before shipping the instrument back to the manufacturer for disposal, decontaminate the instrument to remove any hazardous material.

For more information, or if you wish to dispose of your end-of-life Miltenyi Biotec equipment, contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support.

Electromagnetic compatibility

For US only. Changes or modifications of the equipment unless expressly approved by Miltenyi Biotec may void your authority to operate the equipment pursuant to 47 CFR §15.

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MultiMACS M96 Separator

The MultiMACS M96 Separator is a magnetic separator specifically developed for molecular biology applications, like nucleic acid purification or protein isolation, in microtiter plate format. Thus, Miltenyi Biotec's magnetic bead technology is easily applicable for automated high-throughput procedures.

The MultiMACS M96 Separator can be employed as a semiautomated bench-top instrument or integrated into fully automated liquid-handling platforms. A strong 96-well magnet is used to bind and purify magnetically labeled biomolecules with Multi-96 Columns or Multi-8 Columns. Then, target molecules are eluted.

The key components of the MultiMACS M96 Separator are shown in figure 1.

The MultiMACS 96 Magnet generates a high-gradient magnetic field to bind magnetically labeled molecules in the columns. The MultiMACS Column Holder does hold a MultiColumn frame with up to twelve Multi-8 Columns (up to 96 columns).

The MultiMACS 96 Magnet can be moved vertically (up and down) together with or independent from the MultiMACS Column Holder. Thereby, either lowering of the columns into the wells of the magnet or the lifting of the columns out of the magnet is performed. The vertical movements of the Magnet and the Column Holder are driven by a motor, which is placed inside of the MultiMACS M96 Separator. The motor is controlled either manually by utilizing the touchscreen or fully automated using the RS232-Interface. Note that the orientation of the MultiMACS 96 Magnet and the MultiMACS Column Holder can be changed from landscape to portrait and vice versa.

During a separation or purification process with the MultiMACS M96 Separator, liquids such as samples, wash and elution buffers are dispensed onto the columns either manually, e.g., by an 8-channel pipette or automated by the liquid-dispensing device of a robot. For manual processes, the plates are placed onto the Tip-Touch Plate, which can be moved horizontally back and forth to effectively deposit the last drop into the wells (tip-touch). In automated processes, the plates are moved by a gripper to accomplish this. The Protective Screen is used to prevent injuries to fingers that can

otherwise become trapped between MultiMACS 96 Magnet and MultiMACS Column Holder when the MultiMACS M96 Separator is operated manually.

The flow-through is collected in 96-well plates. The fraction containing the biomolecules of interest is collected in 96-well plates. For this purpose, the tips of the columns are inserted about 1 mm deep into the wells of the plates. Before removal of a plate, or separation of magnet and columns, any drop on the column tips has to be taken off. This avoids contamination of the magnet, cross-contamination between samples, and allows reproducible small elution volumes.


Figure 1: The MultiMACS M96 Separator.


Figure 2: Electrical connectors on the rear side of the MultiMACS M96 Separator.

Applications

The MultiMACS M96 Separator can be used with MACS MicroBeads for the parallel isolation of a wide variety of different biomolecules including:

- mRNA from cells, tissue, whole blood, and total RNA
- Biotinylated biomolecules
- Biomolecules interacting with biotinylated probes
- Proteins using specific antibodies (immunopurification)

Principle of MACS Technology

By labeling with target-specific MicroBeads, molecules such as DNA, RNA, cDNA, or proteins are magnetically purified. MACS MicroBeads are conjugated to Oligo(dT) or antibodies specific either to proteins such as Protein A, Protein G, Streptavidin, epitope tag, or to single-stranded oligonucleotides; therefore, allowing fast, high-purity isolation of biomolecules.

The small size of the magnetizable MicroBeads, about 50 nm in diameter, facilitates their binding to target molecules. First, the colloidal MicroBead suspension is added to the sample, then the MicroBead-labeled sample is applied to the Multi-8/96 Column. Placed in the strong magnetic field of the MultiMACS 96 Magnet, the labeled molecules are retained in the magnetized column matrix, and the non-target molecules are efficiently rinsed out of the column. As the column flow relies on gravity, tedious centrifugation steps as well as removal of wash buffers after bead separation are eliminated.

For elution of the target molecules an appropriate elution buffer, which disrupts the binding between the target molecule and its immobilized ligand, is applied onto the column. Thus, the target molecule is collected in the eluate, and the MicroBeads remain in the column. A second option is the release of the target molecule–MicroBead complex by removing the column from the magnetic field. In this case, the complete complex can be eluted outside of the MultiMACS M96 Separator using a buffer of choice.

The elution volume can be as low as 30 μ L and allows a highly concentrated target fraction to be collected – convenient for further, small-volume downstream applications.

Technical data

Technical data	Specification
Model	MultiMACS M96 Separator
Order number	130-091-937
Dimensions (w × d × h)	230 mm × 435 mm × 230 mm
Weight	11 kg
Input voltage	100–240 V ~, 50/60 Hz
Power consumption	200 W
Electric current	3.5 A max.
Fuses	2 x T3.15AH/250 V, 5 × 20 mm
RS232-Interface (labeled “com”)	Pin 1, 4, 6, 7, 8, 9 NC Pin 2 RXD Pin 3 TXD Pin 5 GND
DC-Output	Module 1 24 V DC / 8 A Module 2 24 V DC / 0.315 A
Working temperature	15–30 °C
Humidity	0–85%
Altitude	max. 2,000 m
Emission sound pressure level at workstation	<70 dB(A)

Table 3.1: Technical data of the MultiMACS M96 Separator.

The MultiMACS M96 Separator is labeled as a protection class I device and must be plugged into a grounded power outlet, refer to **General precautions on page 8**.

The MultiMACS M96 Separator is not specified for use in the cold room. The MultiMACS M96 Separator can process up to 96 samples in parallel. The MultiMACS M96 Separator has been investigated in accordance with the standards UL 61010-1 and CAN/CSA C22.2 No. 61010.1.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna,
- Increase the distance between the equipment and the receiver,
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected,
- Consult the dealer or an experienced radio/TV technician for help.

For other safety considerations, refer to the product label or visit **www.miltenyibiotec.com**.

Design and specifications are subject to change without notice.

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach
Germany

This declaration relates exclusively to the machinery in the state in which it was placed on the market, and excludes components which are added and/or operations carried out subsequently by the final user.

The declaration of conformity refers to the machinery identified as follows:

Description:	Laboratory equipment
Model:	MultiMACS M96 Separator

The machinery complies with all essential requirements of the following directives.

2006/42/EC Machinery
2011/65/EU Restriction of the use of certain hazardous substances in electrical & electronic equipment
2014/30/EU Electromagnetic compatibility

The machinery is in conformity with the following harmonized standards.

EN 61010-1:2010
EN 61326-1:2013

Person authorized to compile the relevant technical documentation:

Dr. Bernd Schröder
Global Head Regulatory Affairs
Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach
Germany

UK Declaration of Conformity

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach
Germany

This declaration relates exclusively to the machinery in the state in which it was placed on the market, and excludes components which are added and/or operations carried out subsequently by the final user.

The declaration of conformity refers to the machinery identified as follows:

Description:	Laboratory equipment
Model:	MultiMACS M96 Separator

The machinery complies with all essential requirements of the following legislations:

Supply of Machinery (Safety) Regulations 2008
The Restriction of the Use of Certain Hazardous Substances in
Electrical and Electronic Equipment Regulations 2012
Electromagnetic Compatibility Regulations 2016

The machinery is in conformity with the following UK designated standards:

EN 61010-1:2010
EN 61326-1:2013

Person authorized to compile the relevant technical documentation:

Dr. Bernd Schröder
Global Head Regulatory Affairs
Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach
Germany

Caution: Carefully read **Important safety information on page 5** before installation and strictly follow the guidelines listed in this chapter!

Safety functions

Protective Screen


The Protective Screen is used to prevent injuries to fingers that can otherwise become trapped between MultiMACS 96 Magnet and MultiMACS Column Holder. Do not remove the Protective Screen if the MultiMACS M96 Separator is operated manually in a semiautomated mode. Remove the Protective Screen only if used in liquid-handling platforms as it can interfere with gripping plates.

In case the Protective Screen is tilted while moving upwards, a light barrier detects the tilted position and the movement is stopped. An error message will appear on the screen.

Caution: Before starting automated processes with the MultiMACS M96 Separator integrated into a compatible automated liquid-handling system, assure that the combined system, i.e., MultiMACS M96 Separator and liquid-handling robot, meets all requirements of safety and function before starting automated processes. Thereby, particularly take the protection of the operator into account.


Emergency break

The MultiMACS M96 Separator contains a motor-driven magnet. If the magnet is moving, the following message appears on the display screen.

Caution: The MultiMACS 96 Magnet is motor-driven. Keep away from all moving parts! Do not reach for plates, columns, or any item between the Tip-Touch Plate and the magnet! Do not reach between the magnet and the Column Holder, while the magnet is moving!

To stop the movement and halt the process in case of danger, accident, or malfunction, tap anywhere on the Touch Display. Then the following screen appears.

Tap **OK** to continue the move. Choose **MOVE BACK** to move the magnet to the last position. Select **ESC** to abort the process.

Unpacking

Unpack the following equipment supplied with the MultiMACS M96 Separator:

- MultiMACS M96 Separator including MultiMACS Display, Column Holder, and MultiMACS 96 Magnet
- Tip-Touch Plate
- User manual
- Declaration of Conformity
- Null modem cable
- Country-specific power cable

Note: Ensure that air can circulate freely behind the rear side of the MultiMACS M96 Separator, and that the instrument is not exposed to heat source, e.g., a radiator or direct sunlight. Failure to do so may result in reduced performance.

Semiautomated bench-top use

1. Place the MultiMACS M96 Separator on a clean, stable, and planar laboratory bench. Tilt the MultiMACS M96 Separator by lifting the device slightly and attach the Tip-Touch Plate into the holes beneath the instrument as shown in figure 3.


Figure 3: Tip-Touch Plate has to be attached in the corresponding slots.

2. Connect the MultiMACS Touch Display via its cable to the RS232 connector on the rear side of the MultiMACS M96 Separator labeled with “display” (refer to figure 2).
3. Connect the power cable to the outlet on the rear side of the MultiMACS M96 Separator and to an appropriate electrical socket.

Automated use in liquid-handling platforms

The MultiMACS M96 Separator can be integrated into liquid-handling platforms for automated use. For further information, contact Miltenyi Biotec Technical Support.

Reagents, columns, and plates

MultiMACS Isolation Kits are specifically designed for the large-scale parallel isolation of molecules using MACS MicroBeads. MultiMACS Isolation Kits consist of a reagent box and a column box for usage with the MultiMACS M96 Separator.

MultiMACS Products	Isolations	Order no.
MultiMACS mRNA Isolation Kit	12 × 8 4 × 96	130-092-520 130-092-519
MultiMACS Streptavidin Kit	12 × 8 4 × 96	130-092-948 130-092-949
MultiMACS Protein A Kit	24 × 8 4 × 96	130-092-944 130-092-945
MultiMACS Protein G Kit	24 × 8 4 × 96	130-092-946 130-092-947
MultiMACS ViralVirus HIV Isolation Kit	12 × 8 4 × 96	130-092-806 130-092-807

µMACS Isolation Kits for manual processing can also be used on the MultiMACS M96 Separator when a column box is purchased separately.

µMACS Products	Isolations	Order no.
µMACS Epitope Tag Protein Isolation Kits		
Anti-HA MicroBeads	40	130-091-122
Anti-c-myc MicroBeads	40	130-091-123
Anti-GFP MicroBeads	40	130-091-125
Anti-GST MicroBeads	40	130-091-370
Anti-His MicroBeads	40	130-091-124

Column box product	Isolations	Order no.
Multi-8 Column, molecular	12 × 8	130-092-444
Multi-96 Column, molecular	4 × 96	130-092-445

Multi-8 Column, molecular (12×8) contains one waste plate (Deep Well Block, 2.5 mL, with sealing foil), one MultiColumn Frame, twelve sterile, single-packed Multi-8 Columns, and one Microtiter Plate with sealing foil. Multi-96 Column, molecular (4×96) contains four Deep Well Blocks, 2.5 mL, in a sterile package with Multi-96 Columns for collection of fractions, and four Microtiter Plates with sealing foil.

Additional Deep Well Blocks (2.5 mL) can be purchased separately.

Product	Quantity	Order no.
Deep Well Block, with sealing foil (2.5 mL)	6	130-092-549

Use microtiter plates with appropriate height dimensions and footprints (width × length). The height is 14 mm for the Microtiter Plate, and for the Deep Well Block, 2.5 mL, it is 44 mm. The Tip-Touch Plate of the MultiMACS M96 Separator is designed for the footprints of microtiter plates that comply to ANSI/SBS standards (length 127.76 mm ±0.5 mm, width 85.48 mm ±0.5 mm).

Attention: Set the plate height correctly if using a different plate height than that specified in the process, refer to **Editing a process on page 105**.

When using MultiMACS mRNA Isolation Kits and filtration of lysate is necessary, gravity flow filter columns can be purchased separately, refer to list below.


MultiMACS Filter product	Isolations	Order no.
Multi-8 Filter	12 × 8	130-092-546
Multi-96 Filter	4 × 96	130-092-547
Multi-8 Filter + Frame	12 × 8	130-092-548

A complete list of MACSmolecular products and protocols is available at www.miltenyibiotec.com.

- 8-channel pipette and tips (electronic 8-channel pipette volume: 15–1250 μL , Matrix Technologies Corporation, is recommended for transferring liquids without foaming)
- Appropriate reagent reservoirs
- If the target is eluted from MicroBeads with hot elution buffer:
 - Heatable reagent reservoir, e.g., heating block or a thermocycler in a corresponding format, e.g., 96-well.
 - Electronic 8-channel pipette programmable in 1 μL steps, with a maximum pipetting volume <300 μL , is recommended

Semiautomated bench-top use with the Touch Display

In this chapter; all functions of the Touch Display are explained which are necessary to carry out a manual separation process. The following scheme presents an overview of all process levels.

Buttons and symbols

^ Scroll up

The scroll function is only visible if there are more items to choose, which are not displayed. The button is used to scroll up the list.

∨ Scroll down


The scroll function is only visible if there are more items to choose, which are not displayed. The button is used to scroll down the list.

|← Back

The button is used to go to the last screen without execution of the indicated step. The magnet will not move during this step.

MOVE BACK

The button is used to go to the last screen without execution of the step indicated on the display.

Caution: The magnet will move during this step.

OK

The button is used to perform the indicated step and to go to the next screen.

ESC (escape)

The button is used to leave the current process or program without execution of the indicated step.

Caution: If ESC is touched during a process, the magnet will move and instructions will be given to remove any plates and columns.

YES

The button is used to confirm a change.

NO

The button is used to abort a change.

DEL (delete)

The button is used to delete a customized process.

Pre-defined processes to isolate mRNA, Protein A/G, and for streptavidin-based applications: 96 POS

Running 96 POS

Caution: Read the chapter **Safety functions** carefully before running a process. Then follow the instructions in chapter **Semiautomated bench-top use**. Always be sure to use MultiMACS 96 Magnet, MultiMACS Column Holder, and the plates in the same orientation as described in **Rotating MultiMACS 96 Magnet and Column Holder on page 122**.


Note: If using a MultiMACS isolation kit, select a pre-defined process according to the kit protocol and follow its instructions. Visit www.miltenyibiotec.com for available firmware-updates of pre-defined processes.

The scroll function (↕) is only visible if there are more items to choose, which are not displayed.

The MultiMACS M96 Separator contains a list of pre-defined separation programs to choose from, visit www.miltenyibiotec.com for updates. 96 POS is the standard process used for applications where the target is eluted from the beads during the elution step, for example MultiMACS mRNA Isolation Kit, MultiMACS Protein A/G kits, and several MultiMACS Streptavidin Kit applications.

If a process with different process parameters must be run, either create a new process, as described in **Creating a new user-defined process on page 95**. Or, if only a few parameters have to be changed, edit the parameters, as described in chapter **Editing a process on page 105**. For information on user-defined processes, refer to **Running a user-defined process on page 103**.


1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds for the Process Selection Screen.


The last process performed on the MultiMACS Instrument is displayed on the upper left (default: 96 POS). The one but last process is listed below.

If necessary scroll through the list of available processes by touching the symbol  or  until 96 POS is displayed.


2. Tap 96 POS to go to the Process Management Screen.


It is possible to check the process parameters: Refer to **Viewing process parameters on page 104**.

3. Tap  to run the process and the magnet will move to the start position. Follow the instructions given on the Touch Display.

4. Insert a MultiColumn Frame with up to twelve Multi-8 Columns or a pre-packed Multi-96 Column into the MultiMACS Column Holder. Tap OK and the magnet will move upwards and unlock the MultiMACS Column Holder. The following screen emerges.

Note: The instrument produces an acoustic signal when the framed fields are activated. Occasional clacking sounds belong to normal function of the MultiMACS M96 Separator.

5. Insert the waste plate, e.g., the Deep Well Block, 2.5 mL, included in the kit, onto the Tip-Touch Plate. If using a plate with a different height, adjust the height of the waste plate, as described in **Editing a process on page 105**.

Tap OK and the MultiMACS 96 Magnet will move downwards until column tips slightly immerse in the waste plate. The next screen appears.

RINSE, APPLY SAMPLE, WASH. IF REQUIRED: PRE-ELUTE	MOVE BACK	OK
	ESC	

6. Rinse columns with equilibration buffer, apply sample, and wash. A pre-elution might be required to exchange the void volume of the column, refer to the respective kit data sheet for detailed instructions. Tap OK for the following screen.

TIP-TOUCH COLUMNS IN PLATE	←	OK
	ESC	

7. Move the Tip-Touch Plate firmly back and forth once so that the inner walls of the wells of the Deep Well Block touch the tips of the Multi-8 Columns; thereby, any drop on the column tips that did not fall by gravity is removed.


Tap OK and the MultiMACS 96 Magnet will move upwards to allow removal of the waste plate.


INSERT ELUTION PLATE: MTP (14 MM)	MOVE BACK	OK
	ESC	


8. Remove the Deep Well Block. If some wells were not used, aspirate waste and store plate. Insert the elution plate (Microtiter Plate), included in the MultiMACS Kit. If using a plate with a different height, adjust the height of the elution plate (refer to **Editing a process on page 105**). Tap **OK** and the MultiMACS 96 Magnet will move downwards until column tips slightly immerse in the elution plate. The subsequent screen appears.


9. Immerse the pipette tips slightly (1 mm) into the reservoirs of the columns without touching the rims and apply elution buffer directly onto the Multi-8/96 Column matrix. Tap **OK** for the next screen.

10. Wait approximately one minute until buffer has passed through. Move the Tip-Touch Plate firmly back and forth once so that the inner walls of the wells of the elution plate touch the tips of the Multi-8 Columns; thereby, any drop on the column tips that did not fall by gravity is removed.

Tap **OK** and the MultiMACS 96 Magnet moves upwards to allow removal of the elution plate.

11. Remove elution plate and use samples for downstream application. Alternatively, seal the plate with adhesive foil and store it according to protocol. Tap **OK** and the MultiMACS 96 Magnet moves apart from the Column Holder.

REMOVE COLUMNS, TOUCH OK TO END PROCESS	MOVE BACK	OK
	ESC	

12. Remove MultiColumn Frame with Multi-8/96 Columns. If less than twelve Multi-8 Columns were utilized, remove the used Multi-8 Columns and store the MultiColumn Frame. Tap **OK** to finish the process.


Keeping flow-through and wash fractions: 96 NEG/POS

Caution: Read the chapter **Safety functions** carefully before running a process. Then follow the instructions in chapter **Semiautomated bench-top use**. Always be sure to use MultiMACS 96 Magnet, MultiMACS Column Holder, and the plates in the same orientation as described in **Rotating MultiMACS 96 Magnet and Column Holder on page 122**.

The MultiMACS M96 Separator contains a list of pre-defined separation programs to choose from (visit www.miltenyibiotec.com for updates). 96 NEG/POS is used when flow-through and wash fractions should be captured separately. It is not the standard process for MultiMACS mRNA Isolation Kit, MultiMACS Protein A/G Kits, etc., where the target is eluted from the beads during the elution step. Deep Well Blocks are included in the MultiMACS Kits, which can be used as waste plate or flow-through plate. For purchase of additional Deep Well Blocks, refer to **Reagents, columns, and plates on page 75**.

If a process with different process parameters has to be run, either create a new process, as described in **Creating a new user-defined process on page 95**, or, if only a few parameters have to be changed, edit the parameters as described in **Editing a process on page 105**. For further details on user-defined processes, refer to **Running a user-defined process on page 103**.


1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds for the Process Selection Screen.


Note: The scroll function is only visible if there are more items to choose, which are not displayed.

The last process performed on the MultiMACS M96 Separator is displayed on the upper left (default: 96 POS). The one but last process is listed below.
If necessary scroll through the list of available process names by


- touching the symbol  or  until 96 NEG/POS is displayed.
2. Tap 96 NEG/POS to go to the Process Management Screen.


If necessary, check the process parameters, as described in **Viewing process parameters on page 104**.


3. Tap  to run the process and the magnet will move to the start position. Follow the instructions given on the Touch Display.

4. Insert a MultiColumn Frame with up to twelve Multi-8 Columns or a pre-packed Multi-96 Column into the MultiMACS Column Holder. Tap **OK** and the magnet will move upwards and unlock the MultiMACS Column Holder. The following screen appears.

5. Insert the waste plate (one Deep Well Block, 2.5 ml, included in the kits) onto the Tip-Touch Plate. If using a plate with a different height, adjust the height of the waste plate, as described in **Editing a process on page 105**.


Note: If using a MultiMACS Isolation Kit, select a pre-defined process according to the respective kit data sheet and follow its instructions.

Visit
www.miltenyibiotec.com
for available firmware updates of pre-defined processes.

Tap **OK** and the MultiMACS 96 Magnet will move downwards until column tips slightly immerse in the waste plate.


6. Rinse columns. Tap **OK** for the next screen.


Move the Tip-Touch Plate firmly back and forth once so that the inner walls of the wells of the waste plate touch the tips of the Multi-8 Columns; thereby, removing any drop on the column tips that did not fall off by gravity.


7. Tap **OK** and the MultiMACS 96 Magnet will move upwards to allow insertion of the flow-through plate.

Insert flow-through plate. If using a plate with a different height adjust the height of the flow-through plate, as described in **Editing a process on page 105**.

8. Tap **OK** and the MultiMACS 96 Magnet will move downwards until column tips slightly immerse in the waste plate.

9. Apply sample and let buffer run through the column. Tap **OK** for the next screen.

TIP-TOUCH COLUMNS IN PLATE	←	OK
	ESC	

10. Move the Tip-Touch Plate firmly back and forth once so that the inner walls of the wells of the flow-through plate touch the tips of the Multi-8 Columns, and any drop on the column tips that did not fall off by gravity will be removed.

Tap **OK** and the MultiMACS 96 Magnet will move upwards to allow removal of the flow-through plate.

INSERT FLOWTHRU PL. 2 OF 2: DWB (44 MM)	MOVE BACK	OK
	ESC	

11. Insert flow-through plate. If using a plate with a different height, adjust the height of the flow-through plate, as described in **Editing a process on page 105**.

Tap **OK** and the MultiMACS 96 Magnet will move downwards until column tips slightly immerse in the waste plate.

WASH COLUMNS. IF REQUIRED: PRE-ELUTE	MOVE BACK	OK
	ESC	

12. Wash columns and let buffer run through. A pre-elution might be required to exchange the void volume of the column. Refer to the respective kit data sheet for detailed instructions. Tap **OK** for the following screen.

TIP-TOUCH COLUMNS IN PLATE	←	OK
	ESC	

13. Move the Tip-Touch Plate firmly back and forth once so that the inner walls of the wells of the flow-through plate touch the tips of the Multi-8 Columns, and any drop on the column tips that did not fall off by gravity will be removed.

Tap **OK** and the MultiMACS 96 Magnet will move upwards to allow removal of the flow-through plate.

INSERT ELUTION PLATE:	MOVE BACK	OK
HTP (14 MM)	ESC	

14. Remove Deep Well Block. If some wells were not used pour off waste liquid and store plate. Insert an elution plate (Microtiter Plate, included in MultiMACS Kits). If using a plate with a different height, adjust the height of the elution plate, as described in **Editing a process on page 105**.


Tap **OK** and the MultiMACS 96 Magnet will move downwards until column tips slightly immerse in the elution plate.

APPLY ELUTION BUFFER	MOVE BACK	OK
	ESC	


15. Immerse the pipette tips slightly (1 mm) into the reservoirs of the columns without touching the rims and apply elution buffer directly onto the Multi-8/96 Column matrix. Tap **OK** for the next screen.

TIP-TOUCH COLUMNS IN PLATE	←	OK
	ESC	

16. Wait approximately one minute until buffer has passed through. Move the Tip-Touch Plate firmly back and forth once so that the inner walls of the wells of the elution plate touch the tips of the Multi-8 Columns, and any drop on the column tips that did not fall off by gravity will be removed.


Tap **OK** and the MultiMACS 96 Magnet moves upwards to allow removal of the elution plate (Microtiter Plate).

17. Remove elution plate and use samples for downstream application. Alternatively, seal the plate with adhesive foil (included in the kit) and store it according to the kit's protocol.

Tap **OK** and the MultiMACS 96 Magnet moves apart from the Column Holder.

18. Take away the MultiColumn Frame with Multi-8/96 Columns. If less than twelve Multi-8 Columns were used, detach the used Multi-8 Columns and store the MultiColumn Frame. Tap **OK** to finish the process.


96 POS + BEAD

Caution: Read the chapter **Safety functions** carefully before running a process. Then follow the instructions in chapter **Semiautomated bench-top use**. Always be sure to use MultiMACS 96 Magnet, MultiMACS Column Holder, and the plates in the same orientation as described in **Rotating MultiMACS 96 Magnet and Column Holder on page 122**.

The MultiMACS M96 Separator contains a list of pre-defined separation programs to choose from, visit www.miltenyibiotec.com for updates. 96 POS + BEAD is used for applications where the target fraction is the molecule-MicroBead complex. The target is eluted outside the magnet.

If a process with different process parameters has to be run, either create a new process, as described in **Creating a new user-defined process on page 95**, or, if only a few parameters have to be changed, edit the parameters as described in **Editing a process on page 105**. For further details on user-defined processes, refer to **Running a user-defined process on page 103**.


1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds for the Process Selection Screen.


The last process performed on the MultiMACS M96 Separator is displayed on the upper left (default: 96 POS). The one but last process is listed below.

If necessary scroll through the list of available process names by touching the symbol ^ or v until 96 POS+BEAD is displayed.


2. Tap 96 POS+BEAD to go to the Process Management Screen.

If necessary check the process parameters, as described in **Viewing process parameters on page 104**.


3. Tap ▶ to run the process and the magnet will move to the start position. Follow the instructions given on the Touch Display.

4. Insert a MultiColumn Frame with up to twelve Multi-8 Columns or a pre-packed Multi-96 Column into the MultiMACS Column Holder.

Tap OK to move the magnet upwards and for the next screen.

5. Insert the waste plate, e.g., a Deep Well Block, onto the Tip-Touch Plate. If using a plate with a different height, adjust the height of the waste plate, as described in **Editing a process on page 105**.

Tap OK and the MultiMACS 96 Magnet will move downwards until column tips slightly immerse in the waste plate. The next screen appears.

RINSE, APPLY SAMPLE, WASH. IF REQUIRED: PRE-ELUTE	MOVE BACK	OK
	ESC	

6. Rinse the columns with equilibration buffer, apply sample, and wash.

Tap OK for the following screen.

TIP-TOUCH COLUMNS IN PLATE	←	OK
	ESC	

7. Move the Tip-Touch Plate firmly back and forth once so that the inner walls of the wells of the Deep Well Block touch the tips of the Multi-8 Columns. This process removes any drop on the column tips that did not fall by gravity.

Tap OK to move the MultiMACS 96 Magnet upwards and to allow removal of the waste plate.

REMOVE WASTE PLATE	MOVE BACK	OK
	ESC	

8. Remove the Deep Well Block.

Tap OK to move the MultiMACS 96 Magnet away from the Column Holder. The following screen appears.

REMOVE COLUMNS FOR EXTERNAL ELUTION TOUCH OK TO END	MOVE BACK	OK
	ESC	

Tap OK to finish the MultiMACS process. Remove the columns and proceed with the target elution as described in the respective datasheet.

Caution: Read the chapter **Safety functions** carefully before running a process. Then follow the instructions in chapter **Semiautomated bench-top use**. Always be sure to use MultiMACS 96 Magnet, MultiMACS Column Holder, and the plates in the same orientation as described in **Rotating MultiMACS 96 Magnet and Column Holder on page 122**.

1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds for the Process Selection Screen.

WELCOME TO THE MULTIMACS SEPARATOR BY MILTENYI BIOTEC		
96 POS	NEW	^
96 NEG/POS	SET UP	v

2. Tap NEW and the following screen appears.

ENTER THE NAME FOR YOUR CUSTOMIZED PROCESS	←	OK
	ESC	

Wait for a few seconds or tap OK to go to the next screen where a new process name has to be entered.

ENTER NAME: A_			←	OK
ABC	DEF	GHI	...	DEL

3. Tap the ABC box once to enter an A, twice to enter a B, and three times for a C.
Tap the ... box — multiple times if necessary — to switch to the triplets of the letters J to Z, symbols, and numerals.

Tap DEL to delete letters and OK to confirm the entered name.
If the process name already exists, a message will appear on the Touch Display that another name has to be chosen.

Selection of parameter types

After entering a name, the following parameter types for the new process have to be selected:

- A. Selecting **Fraction of Interest**
- B. Selecting **Elution Mode** (depends on parameter A)
- C. Selecting **Plate Heights** (for waste, flow-through, and elution plate height; also depending on parameter A)

At this stage, save the process or, if needed, enter advanced parameter types.

- D1. **Distance for Plate Removal**
- D2. **Number of Flow-through Plates**

Note: The scroll function is only visible if there are more items to choose, which are not displayed.

A. Selecting Fraction of Interest

First you have to select the desired fraction after magnetic labeling of samples. The positive labeled fraction, the negative labeled fraction (flow-through), or both fractions can be chosen.

SEL. FRACTION OF INTEREST:	←	^
COLLECT POS. FRACTION	ESC	✓

- If only the labeled fraction is to be captured, tap COLLECT POS. FRACTION and proceed with section B. Selecting Elution Mode.
- To collect only the unlabeled flow-through, scroll through the list by touching the symbol ^ or ✓ until COLLECT NEG. FRACTION appears.

Note: An elution with or without MicroBeads is not distinguished in this step.

SEL. FRACTION OF INTEREST:	←	^
COLLECT NEG. FRACTION	ESC	✓

Tap **COLLECT NEG. FRACTION** and proceed with section C. Selecting Plate Heights.

- In order to collect the labeled fraction as well as the unlabeled flow-through, scroll through the list by touching the symbol $\wedge$ or $\vee$ until **COLLECT BOTH FRACTIONS** appears.

SEL. FRACTION OF INTEREST:	←	^
COLLECT BOTH FRACTION	ESC	∨

Tap **COLLECT BOTH FRACTIONS** and proceed with section B. Selecting Elution Mode.

B. Selecting Elution Mode

If **COLLECT POS. FRACTION** or **COLLECT BOTH FRACTIONS** were chosen for the fraction of interest, two options for the elution of target molecules are possible.

SELECT ELUTION MODE:	←	^
ELUTION IN- SIDE MAGNET	ESC	∨

- To retain the MicroBeads within the column while the target molecule is eluted, tap **INSIDE MAGNET** and proceed with section C. Selecting Plate Heights.
- To elute the target molecule-MicroBead complex by removing the column from the magnetic field, scroll through the list by touching the symbol $\wedge$ or $\vee$ until **ELUTION OUTSIDE MAGNET** appears.

SELECT ELUTION MODE:	←	^
ELUTION OUT- SIDE MAGNET	ESC	∨

Tap **ELUTION OUTSIDE MAGNET** and proceed with section C. Selecting Plate Heights.

C. Selecting Plate Heights

The Tip-Touch Plate of the MultiMACS M96 Separator is designed for microtiter plates that comply to the footprint as described in the ANSI/SBS standard, refer to **Reagents, columns, and plates on page 75**.

The plate heights for waste, flow-through, and eluate depend on the selected process parameters, e.g., **Fraction of Interest**, and have to be determined individually.

Fraction of interest	WASTE PLATE to collect	FLOWTHRU PLATE to collect	ELUTION PLATE to collect
COLLECT POS FRACTION	Equilibration buffer, flow-through, and wash buffer	–	Eluate**
COLLECT NEG FRACTION	Equilibration buffer	Flow-through + wash buffer*	–
COLLECT BOTH	Equilibration buffer	Flow-through + wash buffer*	Eluate**

* Flow-through and wash buffer can be separately collected in flow-through plates by increasing the number of flow-through loops.

**Parameter Elution Plate does not appear if Elution Mode is OUTSIDE MAGNET.

1. Depending on the chosen fraction of interest, the heights of the required plates can be adjusted. For example, when the positive fraction was chosen, screens will be presented to specify the height of the waste plate and elution plate. Below all three screens are shown.

SELECT WASTE PLATE:

←

^

DWB (44 MM)

ESC

✓

SELECT FLOWTHRU PL.:

←

^

DWB (44 MM)

ESC

✓

SELECT ELUTION PLATE:

←

^

MTP (14 MM)

ESC

✓

2. Tap the displayed plate type at the lower left of the Touch Display to select the default value. If a height different to the default value should be selected, scroll through the list by touching the symbol $\wedge$ or $\vee$ until the following screens appear.

SELECT WASTE PLATE:	$\leftarrow$	$\wedge$
ENTER HEIGHT	ESC	$\vee$

SELECT FLOWTHRU PL.:	$\leftarrow$	$\wedge$
ENTER HEIGHT	ESC	$\vee$

SELECT ELUTION PLATE:	$\leftarrow$	$\wedge$
ENTER HEIGHT	ESC	$\vee$

3. Tap ENTER HEIGHT for the next screen.

30.0 MM	+0.5 MM	+5 MM	$\leftarrow$	OK
	-0.5 MM	-5 MM	ESC	TEST

Tap +0.5 MM or +5 MM to increase the displayed plate height by this amount (in mm) and -0.5 MM or -5 MM to decrease it. The minimum plate height is 8.0 mm, and the maximum plate height is 80.0 mm.

For further details, refer to **Testing Plate Height or Distance for Plate Removal on page 115**.

Tap OK to select the displayed plate height.

D. Advanced parameter types

For special applications, it may be necessary to change advanced process parameters, these include **Distance for Plate Removal** and **Number of Flow-through Plates**.

DO YOU WANT TO CHANGE ADVANCED PARAMETERS?	←	YES
	ESC	NO

If you press YES, the parameter **Distance for Plate Removal** and the **Number of Flow-through Plates** have to be selected. Alternatively, tap NO and go to **Saving the edited process on page 114**.

D1. Selecting Distance for Plate Removal

The Distance for Plate Removal is the distance the magnet moves up to create space to remove a plate and to insert a new one.

DISTANCE FOR PLATE REMOVAL:	←	^
PLATE REMOVAL 25.0 MM	ESC	v

Press PLATE REMOVAL 25.0 MM to select the default value. To change this parameter, scroll through the list by touching the symbol ^ or v until ENTER DISTANCE appears.

DISTANCE FOR PLATE REMOVAL:	←	^
ENTER DISTANCE	ESC	v

Tap ENTER DISTANCE for the following screen.

30.0 MM	+0.5 MM	+5 MM	←	OK
	-0.5 MM	-5 MM	ESC	TEST

Note: The Tip-Touch Insertion Distance is calculated by the MultiMACS M96 Separator automatically. Therefore, enter the actual height of the preferred plate.

Press +0.5 MM or +5 MM to increase the displayed distance for plate removal by this amount (in mm) and -0.5 MM or -5 MM to decrease it. The minimum height is 1.0 mm, maximum height is 45.0 mm. For details, refer to **Testing Plate Height or Distance for Plate Removal on page 115**.

Tap OK to select the displayed plate removal distance.

D2. Selecting Number of Flow-through Plates

1. If COLLECT NEG FRACTION or COLLECT BOTH FRACTIONS was chosen as the fraction of interest, the flow-through and each wash buffer fractions can be collected separately by increasing the number of flow-through plates.

NUMBER OF FLOWTHRU PLATES: 1	+1	←	^
		ESC	v

2. The default value is one plate to keep flow-through and all wash fractions in a single plate. Optionally, tap +1 to increase the number of plates, which enables collection of wash fractions separately.

NUMBER OF FLOWTHRU PLATES: 2	+1	←	^
	-1	ESC	v

3. Optionally, tap -1 to decrease the number of plates. Confirm by touching OK.

Saving the created process

1. Tap **SAVE** to confirm the addition of the process to the list. If there are already ten user-defined processes stored in the list, an obsolete process has to be replaced, refer to **Replacing a user-defined process on page 116**. After saving, the following screen appears.

2. Tap **YES** to start the process or **NO** to go back to the Welcome Screen.

Note: The scroll function is only visible if there are more items to choose, which are not displayed.

Note: If using a MultiMACS Isolation Kit, select a pre-defined process according to the kit protocol and follow the protocol instructions.

Running a user-defined process

Caution: Read the chapter **Safety functions** carefully before running a process. Then follow the instructions in chapter **Semiautomated bench-top use**. Always be sure to use MultiMACS 96 Magnet, MultiMACS Column Holder, and the plates in the same orientation as described in **Rotating MultiMACS 96 Magnet and Column Holder on page 122**.

1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds for the Process Selection Screen.

WELCOME TO THE
MULTIMACS SEPARATOR
BY MILTENYI BIOTEC

Note: The scroll function is only visible if there are more items to choose, which are not displayed.

96 POS	NEW	^
96 NEG/POS	SET UP	v

The last process performed on the MultiMACS is displayed in the upper left (default: 96 POS).

2. If necessary, scroll through the list of available process names by touching the symbol ^ or v until the desired process is displayed. Tap the process name to select it, e.g., MY_OWN_PROG, and the Process Management Screen appears.

MY_OWN_PROG			←	▶
VIEW	EDIT	DEL	ESC	

Note: When viewing process parameters, the name of the program will be shortened from two lines to one line.

CAUTION


Tray in motion during operation.
Hazard of crushing fingers.
Operate only with attached protection screen.
Device must be able to be tilted backwards.

3. Tap ▶ to run the process and the MultiMACS 96 Magnet will move to the start position. Follow the instructions given on the Touch Display.

Viewing process parameters


It is possible to control the parameters before running an already saved process. Changing process parameters is not possible in the view mode, this can only be done in the edit function.

- 1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds for the Process Selection Screen.

The last process performed on the MultiMACS is displayed in the upper left (default: 96 POS). The one but last process is listed below.

- 2. If necessary, scroll through the list of available process names by touching the symbol ^ or v until the desired process is displayed. Tap the process name to select it, e.g., 96 POS.

- 3. Tap VIEW for the View Process Screen.

- 4. The chosen parameter is displayed in the lower left of the Touch Display beginning with the parameter type fraction of interest. Scroll through the list with the symbol ^ or v to see the other process parameters. For a description of process parameters, refer to **Creating a new user-defined process on page 95**.

5. Tap the arrow |← to go back to the Process Management Screen.

Editing a process

Changing of process parameters of an existing process should be done in the edit mode. If several process parameters have to be changed, refer to **Creating a new user-defined process on page 95**.

1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds for the Process Selection Screen.

EDIT FRACTION OF INTEREST:	←	^
COLLECT POS. FRACTION	ESC	∨

96 POS	NEW	^
96 NEG/POS	ESC	∨

The last process performed on the MultiMACS M96 Separator is displayed in the upper left (default: 96 POS).

2. If necessary, scroll through the list of available process names by touching the symbol ^ or ∨ until the desired process is displayed. Press the process name to select it and to go to the Process Management Screen.

MY_OWN_PROG			←	▶
VIEW	EDIT	DEL	ESC	

Note: If the process name already exists, a message will appear on the Touch Display that another name has to be chosen.

3. Tap EDIT to edit process parameters. If a user-defined process is edited, the following screen appears.

SAVE AS A NEW PROCESS	←
JUST EDIT THE PARAMETERS	ESC

4. Press **SAVE AS A NEW PROCESS** to change process parameters and to save the new process. Thereby, the old process name remains with unchanged process parameters. Alternatively, tap **JUST EDIT THE PARAMETERS** to preserve the process name while changing process parameters.

If saving as a new process the following screen appears.

ENTER THE NAME FOR YOUR CUSTOMIZED PROCESS	←	OK
	ESC	

5. Wait for a few seconds or tap **OK** to go to the next screen where a new process name has to be entered.

ENTER NAME: A_			←	OK
ABC	DEF	GHI	...	DEL

6. Tap the **ABC** box once to enter an A, twice to enter a B, and three times for a C.

Tap the **...** box — multiple times if necessary — to switch to the triplets of the letters from J to Z, symbols, and numerals. Tap **DEL** to delete letters and **OK** to confirm the entered name.

EDIT PROCESS MY_OWN_PROG:	←	^
COLLECT POS. FRACTION	ESC	✓

The current process parameter name is displayed on the upper left of the Touch Display and the parameter is on the lower left, beginning with the parameter type **Fraction of Interest**. Depending on the parameter of the chosen process, several process parameters can be edited.

A. Editing Fraction of Interest

1. To change the parameter type **Fraction of Interest**, tap the displayed parameter, e.g., COLLECT POS FRACTION, and the next screen appears.

EDIT FRACTION OF INTEREST:	←	^
COLLECT POS. FRACTION	ESC	v

2. The parameter type now appears in the upper left of the screen. To change the parameter, scroll through the list of options by touching the symbol ^ or v until the desired parameter is displayed.

EDIT FRACTION OF INTEREST:	←	^
COLLECT BOTH FRACTION	ESC	v

3. Tap the parameter to select it, e.g., COLLECT BOTH FRACTIONS.

If COLLECT POS FRACTION is changed to COLLECT BOTH FRACTIONS or COLLECT NEG FRACTION, additional parameter types have to be selected. For details, refer to **C. Selecting Plate Heights on page 98**.

SAVE PROCESS TO LIST	←	
MORE CHANGES	ESC	

4. If no further parameter should be edited, tap SAVE PROCESS TO LIST and go to **Saving the created process on page 102**. Alternatively, edit other parameters types by touching MORE CHANGES.

B. Editing Elution Mode

1. This parameter type only appears if COLLECT POS FRACTION or BOTH FRACTIONS were chosen as the fraction of interest.

EDIT PROCESS MY_OWN_PROG:	←	^
COLLECT POS. FRACTION	ESC	v

2. Scroll through the list by touching the symbol ^ or v until the process parameter of the elution mode, ELUTION INSIDE MAGNET or ELUTION OUTSIDE MAGNET, is displayed. The current process parameter is displayed on the lower left of the Touch Display.

EDIT PROCESS MY_OWN_PROG:	←	^
ELUTION IN- SIDE MAGNET	ESC	v

3. To change the parameter of the parameter, type **Elution Mode**, tap the displayed parameter, e.g., ELUTION INSIDE MAGNET.

EDIT ELUTION MODE:	←	^
ELUTION IN- SIDE MAGNET	ESC	v

4. The parameter type now appears in the upper left screen. To change the parameter, scroll through the list of options by touching the symbol ^ or v until the desired parameter is displayed.

EDIT ELUTION MODE:	←	^
ELUTION OUT- SIDE MAGNET	ESC	v

5. Tap the changed parameter to select it, e.g., ELUTION OUTSIDE MAGNET, and the next screen appears.

SAVE PROCESS TO LIST	←
MORE CHANGES	ESC

If no further parameter should be edited, tap SAVE PROCESS TO LIST and go to page 115, Saving the edited process. Alternatively, edit other parameters types by touching MORE CHANGES.

C. Editing Plate Height

The Tip-Touch Plate of the MultiMACS M96 Separator is designed for microtiter plates that comply to the footprint as described in the ANSI/SBS standard, refer to **Reagents, columns, and plates on page 75**.

Fraction of interest	WASTE PLATE to collect	FLOWTHRU PLATE to collect	ELUTION PLATE to collect
COLLECT POS FRACTION	Equilibration buffer, flow-through and wash buffer	–	Eluate**
COLLECT NEG FRACTION	Equilibration buffer	Flow-through + wash buffer*	–
COLLECT BOTH	Equilibration buffer	Flow-through + wash buffer*	Eluate**

* Flow-through and wash buffer can be separately collected in plates by increasing number of flow-through loops.

** Parameter ELUTION PLATE does not appear if Elution Mode is OUTSIDE MAGNET.

Depending on selected process parameters (refer to table above) several plate heights for collecting waste, flow-through, and/or elution can be edited. For example, when COLLECT POS FRACTION was selected the heights of the plates for waste and eluate can be specified.

EDIT PROCESS MY_OWN_PROG#	←	^
COLLECT POS. FRACTION	ESC	v

1. Scroll through the list by touching the symbol $\wedge$ or $\vee$ until the process parameter of waste, flow-through, or elution plate is displayed. The current process parameter is displayed on the lower left of the Touch Display. Below all three screens are listed.

EDIT PROCESS MY_OWN_PROG:	$\leftarrow$	$\wedge$
WASTE PLATE DWB (44 MM)	ESC	$\vee$

EDIT PROCESS MY_OWN_PROG:	$\leftarrow$	$\wedge$
FLOWTHRU PLATE DWB (44 MM)	ESC	$\vee$

EDIT PROCESS MY_OWN_PROG:	$\leftarrow$	$\wedge$
ELUTION PLATE MTP (14 MM)	ESC	$\vee$

2. Press the displayed plate height to change it. The corresponding screen appears. Again, all three screens are shown.

EDIT WASTE PLATE:	$\leftarrow$	$\wedge$
DWB (44 MM)	ESC	$\vee$

EDIT FLOW- THRU PLATE:	$\leftarrow$	$\wedge$
DWB (44 MM)	ESC	$\vee$

EDIT ELUTION PLATE:	$\leftarrow$	$\wedge$
MTP (14 MM)	ESC	$\vee$

3. Tap the displayed plate type at the lower left side of the Touch Display to select the default value. If a height different to the default value should be selected, scroll through the list by touching the symbol $\wedge$ or $\vee$ until the corresponding screen is presented.

EDIT WASTE PLATE:	$\leftarrow$	$\wedge$
ENTER HEIGHT	ESC	$\vee$

EDIT FLOW- THRU PLATE:	$\leftarrow$	$\wedge$
ENTER HEIGHT	ESC	$\vee$

EDIT ELUTION PLATE:	$\leftarrow$	$\wedge$
ENTER HEIGHT	ESC	$\vee$

4. Tap ENTER HEIGHT.

30.0 MM	+0.5 MM	+5 MM	$\leftarrow$	OK
	-0.5 MM	-5 MM	ESC	TEST

5. Tap +0.5 MM or +5 MM to increase the displayed plate height by this amount (mm) and -0.5 MM or -5 MM to decrease it. The minimum plate height is 8.0 mm. The maximum plate height is 80.0 mm. Tap OK to select the displayed height.

SAVE PROCESS TO LIST	$\leftarrow$
MORE CHANGES	ESC

If no further parameter should be edited, go to **Saving the created process on page 102**. Alternatively, edit further parameters by touching MORE CHANGES.

D. Advanced parameter types

D1. Editing Distance for Plate Removal

The **Distance for Plate Removal** is the user-defined distance (default 25 mm) the magnet moves up to create space for plate removal and insertion.

EDIT PROCESS MY_OWN_PROG#	←	^
COLLECT POS. FRACTION	ESC	v

1. Starting from the screen above, scroll through the list by touching the symbol ^ or v until the process parameter PLATE REMOVAL 25.0 MM is displayed (displayed values may vary). The current process parameter is shown on the lower left of the Touch Display.

EDIT PROCESS MY_OWN_PROG#	←	^
PLATE REMOVAL 25.0 MM	ESC	v

2. To change the parameter, tap PLATE REMOVAL 25.0 MM (displayed value may vary). The next screen appears.

EDIT DISTANCE PLATE REMOVAL#	←	^
PLATE REMOVAL 25.0 MM	ESC	v

3. The parameter type now appears in the upper left screen. Press PLATE REMOVAL 25.0 MM to select default value, or, to change the parameter, scroll through the list by touching the symbol ^ or v until ENTER DISTANCE is displayed.

EDIT DISTANCE PLATE REMOVAL#	←	^
ENTER DISTANCE	ESC	v

4. Tap ENTER DISTANCE for the next screen.

30.0 MM	+0.5 MM	+5 MM	←	OK
	-0.5 MM	-5 MM	ESC	TEST

Adjust the desired distance by touching the corresponding buttons, e.g., +0.5 MM and +5 MM to increase and -0.5 MM and +5 MM to decrease distance.

5. Tap OK for the next screen.

SAVE PROCESS TO LIST	←
MORE CHANGES	ESC

If no further parameter should be edited, go to SAVE PROCESS TO LIST, refer to **Saving the edited process on page 114**.

Alternatively, edit other parameters types by touching MORE CHANGES.

D2. Editing Number of Flow-through Plates

If COLLECT BOTH FRACTIONS or COLLECT NEG FRACTIONS was chosen as fraction of interest, the flow-through fraction and different wash buffer fractions can be collected separately by increasing the number of flow-through plates.

EDIT PROCESS MY_OWN_PROG:	←	^
COLLECT POS. FRACTION	ESC	v

1. Starting from the previous screen above, scroll through the list by touching the symbol ^ or v until the process parameter NUMBER OF FLOWTHRU PL (plates) is displayed. The current process parameter is shown on the lower left of the Touch Display.

EDIT PROCESS MY_OWN_PROG:	←	^
NUMBER FLOW- THRU PL.: 1	ESC	v

2. To change the number of flow-through plates, tap the displayed parameter. The next screen appears.

EDIT NUMBER OF FLOWTHRU PLATES: 1	+1	←	OK
		ESC	

3. The default value 1 keeps flow-through and all wash fractions in one plate. Tap +1 to increase the number of plates; thereby, collecting wash fractions separately. Confirm by touching OK and to switch to the next screen.

SAVE PROCESS TO LIST	←
MORE CHANGES	ESC

4. If no further parameter should be edited, go to SAVE PROCESS TO LIST described below. Alternatively, edit other parameters by touching MORE CHANGES.

Saving the edited process

1. Save the process with the changed parameters by touching SAVE PROCESS TO LIST, thus leaving the edit mode. In case ten user-defined processes have already been saved in the list, an obsolete process has to be replaced, refer to **Replacing a user-defined process on page 116**. The following screen appears.

DO YOU WANT TO START THIS PROCESS NOW?		YES
	ESC	NO

2. Tap YES to start the process immediately after saving; alternatively, tap NO or ESC to get to the Welcome Screen.

Testing Plate Height or Distance for Plate Removal

For your convenience, plate heights or distance for plate removal can be tested before saving.


Note: The Tip-Touch insertion distance is calculated by the MultiMACS M96 Separator automatically. Therefore, enter the actual height of the preferred plate.

1. After touching ENTER HEIGHT or ENTER DISTANCE, the following screen appears.

30.0 MM	+0.5 MM	+5 MM	←	OK
	-0.5 MM	-5 MM	ESC	TEST

2. Tap TEST and follow the instructions to check the selected values.

INSERT MULTI-8/96 COLUMNS	MOVE BACK	OK
	ESC	

3. Insert the MultiColumn Frame with up to twelve Multi-8 Columns. Tap OK and the MultiMACS 96 Magnet moves up to the next position.

INSERT PLATE: HEIGHT 30.0 MM	MOVE BACK	OK
	ESC	

4. Insert the plate type to be tested. Tap OK and the MultiMACS 96 Magnet moves down to the next position.

TEST PLATE HEIGHT 30.0 MM AT M8/96 COLUMN	MOVE BACK	OK
	ESC	

5. Check whether the height is acceptable by moving the Tip-Touch

Plate back and forth. The column tips should be inserted slightly in the well (1 mm). Tap OK and the MultiMACS 96 Magnet moves up.

REMOVE PLATE	MOVE BACK	OK
	ESC	

6. Remove the plate. Tap OK, the MultiMACS 96 Magnet moves down, and the first screen reappears.

30.0 MM	+0.5 MM	+5 MM	←	OK
	-0.5 MM	-5 MM	ESC	TEST

7. Tap OK if the height was acceptable, otherwise adapt plate height and retry the test.

Return to the previous sections, such as **C. Selecting Plate Heights on page 98** or **D1. Selecting Distance for Plate Removal on page 100**.

Replacing a user-defined process

In case ten user-defined processes were already saved and another user-defined process should be created, or a process should be edited and saved as a new process, an obsolete process has to be removed from the list of user-defined processes.

REPLACE PROCESS:	←	^
MY_OWN_PROG:	ESC	✓

1. If necessary, scroll through the list by touching the symbol $\wedge$ or $\vee$ until the desired process name to be replaced is displayed. Touch the process name, e.g., MY_OWN_PROG, to replace it. The following screen emerges.

DO YOU WANT TO REPLACE THIS PROCESS: MY_OWN_PROG:	$\leftarrow$	YES
	ESC	NO

2. Confirm the replacement by touching YES or go back by touching NO or $\leftarrow$. The new process has now been saved and with the following screen you can choose to start this process immediately.

DO YOU WANT TO START THIS PROCESS NOW?		YES
	ESC	NO

3. Tap YES to start the process. Alternatively, tap NO or ESC to go to the Welcome Screen.

Note: Only user-defined processes can be deleted.

Deleting a user-defined process

1. To delete a process from the list of available processes, switch on the instrument and then touch the Welcome Screen or wait for a few seconds until the Process Selection Screen appears.

WELCOME TO THE MULTIMACS SEPARATOR BY MILTENYI BIOTEC		
MULTI-8/96 POS	NEW	$\wedge$
MULTI-8/96 NEG/POS	SET UP	$\vee$

The last process performed on the MultiMACS M96 Separator is displayed in the upper left (default: 96 POS). The one but last process is listed in the lower left.

2. If necessary, scroll through the list of available process names by

touching the symbol  or  until the desired process is displayed. Touch the desired process name to select it. The Process Management Screen appears.

3. Tap DEL to proceed and confirm.


4. Confirm the deletion of the process by touching YES. Tap NO to stop the deletion.

Firmware version

Before contacting the Miltenyi Biotec Technical Support, determine the firmware version on the MultiMACS M96 Separator.


- 1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds until the Process Selection Screen appears.


Note: SERVICE TOOLS are for service purposes and are only to be used by an authorized Miltenyi Biotec service provider.

The scroll function is only visible if there are more items to choose, which are not displayed.


- 2. Tap SET UP for the Setup Screen.

- 3. Tap PREFERENCES to display the Preference Category Screen.

- 4. Scroll down until SHOW VERSION OF FIRMWARE appears.


5. Tap this button and the current version of the firmware is presented in the next screen.

Preparing the MultiMACS M96 Separator for transport

The MultiMACS M96 Separator has to be locked before shipping. Contact Miltenyi Biotec before sending the MultiMACS M96 Separator back.

1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds until the Process Selection Screen appears.

2. Tap SET UP for the Setup Screen.


3. Tap PREFERENCES to display the Preference Category Screen.



4. Scroll down to LOCK UNIT FOR TRANSPORT.

5. Tap LOCK UNIT FOR TRANSPORT for the MultiMACS 96 Magnet and MultiMACS Column Holder to move to the lock position.

6. As soon as the text above appears, the instrument can be switched off. The MultiMACS M96 Separator is now ready for shipping.

Note: SERVICE TOOLS are for service purposes and are only to be used by an authorized Miltenyi Biotec service provider.

The scroll function is only visible if there are more items to choose, which are not displayed.


Rotating MultiMACS 96 Magnet and Column Holder

The MultiMACS 96 Magnet and Column Holder are delivered in landscape orientation, but they can also be used in portrait orientation. In rare cases, if 96-well magnets needs to be exchanged, follow the instructions outlined below.


Warning: The MultiMACS M96 Separator is equipped with an extremely powerful magnet, the MultiMACS 96 Magnet. Keep this magnet at a distance of at least 80 cm from any magnetic information carriers such as credit cards, magnetic tapes, and storage media, and from any electronic equipment, such as hearing aids, pacemakers, measuring and control instruments, computers, and watches. These items may be affected or damaged by the extremely powerful magnetic field. The force of attraction between two magnets or between a magnet and magnetically attractable material increases strongly the closer they get. Thus, keep any magnetic material – including any other MultiMACS 96 Magnet and Display – at a safe distance.

Note: The scroll function is only visible if there are more items to choose, which are not displayed.


- 1. Switch on the instrument and touch the Welcome Screen or wait for a few seconds for the Process Selection Screen.

- 2. Tap SET UP for the Setup Screen.


3. Tap PREFERENCES to go to the Preference Category Screen.

4. Tap EXCHANGE MAGNET and the MultiMACS 96 Magnet will move to the middle exchange position while the Column Holder is locked in the top position.

5. Wait until the following text appears.

6. Switch-off and unplug the MultiMACS M96 Separator.

Note: In landscape orientation, the safety notices on the MultiMACS 96 Magnet have to face the front. The notched corner is at the lower left.

Rotation of the MultiMACS 96 Magnet

1. Unscrew all four knurled screws of the MultiMACS 96 Magnet. For rotation of magnet from **landscape to portrait** orientation, turn magnet counterclockwise by 90°. Note that the 'notched' corner, marked red in figure 4 a + b, is now at the lower right. For rotation of the magnet from **portrait to landscape** orientation, turn the magnet clockwise by 90° (refer to figure 4 a + b).
2. Align the MultiMACS 96 Magnet with the metal brackets. The holes in the metal brackets must be aligned with the screw threads in the MultiMACS 96 Magnet. The holes can only be aligned with the screw threads if the notched corner of the magnet is to the left (landscape orientation, refer to figure 5a) or lower right (portrait orientation, refer to figure 5b).

3. Now fasten all four knurled screws.

4. Connect the power cable to the outlet on the rear side of the MultiMACS M96 Separator and to an appropriate electrical socket.

Note: Magnet and Column Holder must be always aligned in the same orientation (portrait or landscape).


Figure 4a: Landscape orientation.


Figure 4b: Portrait orientation.

Rotation of the Column Holder

1. Unscrew all four knurled screws of the MultiMACS Column Holder. For rotation of Column Holder from landscape to portrait orientation, rotate Column Holder counter clockwise by 90° (the same way as the MultiMACS 96 Magnet).

2. For rotation of Column Holder from portrait to landscape orientation, rotate Column Holder clockwise by 90° in the same way as the MultiMACS 96 Magnet.

3. Align the MultiMACS Column Holder with the metal brackets as shown in figure 5 a + b. The holes in the metal brackets have to be aligned with the screw threads in the MultiMACS Column Holder.

4. Now fasten all four knurled screws.


Figure 5a: Landscape orientation.


Figure 5b: Portrait orientation.

Caution: Unplug the power cable before cleaning the MultiMACS M96 Separator. The MultiMACS M96 Separator and the MultiMACS 96 Magnet must not be submerged in water or cleaned in a dish washer, ultrasonic cleaner, autoclave, or similar appliances.

Clean up spills immediately. Use a wash solution containing detergent, e.g., 1% SDS solution, and/or 70% ethanol solution for disinfection and cleaning.

The MultiMACS M96 Separator should not move while performing Tip-Touch. If this happens accidentally, clean the instrument feet with alcohol.

Cross-contamination by drops at the column tips

Always make a Tip-Touch before changing plates or allowing separation of the MultiMACS 96 Magnet and the MultiMACS Column Holder. To take off drops at the column tips, move Tip-Touch Plate so that the Multi-8/96 Columns move slightly in the MultiMACS Column Holder. Control whether all drops at the column tips have been removed. If drops are still present, tap MOVE BACK, repeat the Tip-Touch. Tap OK to return to the screen.

The MultiMACS M96 Separator should not move while performing Tip-Touch. If this happens, clean the instrument feet, refer to **Maintenance on page <MS>127**.

Collision detection

CAUTION
Tray in motion during operation.
Hazard of crushing fingers.
Operate only with attached protection screen.
Device must be able to be tilted backwards.

If a tilting of the MultiMACS M96 Separator is registered while the magnet is moving, the movement stops and the following screen appears.

MULTIMACS SEP. IS TILING! TOUCH MOVE BACK OR RETRY	MOVE BACK ESC	RE- TRY
---	---------------------	------------

If a blocking object is visible, tap **MOVE BACK** to set back the magnet to its former position. Take out any objects causing the stop of the movement. Subsequently, tap **OK** to restart the move. In case there are no blocking objects visible, tap **RETRY**.

If a tilting of the MultiMACS Column Holder is registered while the magnet is moving, the movement stops and the following screen appears.

COLUMN HOLDER IS TILING! TOUCH MOVE BACK OR RETRY	MOVE BACK ESC	RE- TRY
--	---------------------	------------

If a blocking object is visible, touch **MOVE BACK** to set back the magnet to its former position. Take out any objects causing the stop of the movement. Subsequently, tap **OK** to restart the move. In case there are no blocking objects visible, touch **RETRY**.

The use of cold buffer on a tempered column is not recommended due to the possibility of air bubble formation in the column matrix.

Malfunctions and possible remedies

Symptom	Cause	Possible remedy
No light on the display when power is turned on.	Power cord or display cord is not plugged in. Fuse is blown or missing.	Make sure power source matches power rating stamped on rear label. Check fuse holder in rear, refer to Fuse exchange on page 133 .
No message appears on the display or display is frozen.	Microprocessor is temporarily out of normal sequence.	Turn power off, wait 5 seconds, and then turn on again.

Table 8.1: Malfunctions, causes, and suggestions for fixing.

Error codes

All errors are coded. Refer to the following table for explanations.

Displayed error message	Explanation and remedy
STEPPING LOSS ERROR -101 INITIALIZING SYSTEM	Motor out of sequence. Tap OK for intermediate reinitialization of the system.
WRONG MAG POS ERROR -102 INITIALIZING SYSTEM	Incorrect positioning of the magnet. Tap OK for intermediate reinitialization of the system.
WRONG COL POS ERROR -106 INITIALIZING SYSTEM	Column Holder is locked in wrong position. Tap OK for intermediate reinitialization of the system. Column holder will be repositioned.
POSITION FAILED ERROR -111 INITIALIZING SYSTEM	Incorrect positioning of the magnet. Make sure device is not exposed to direct sunlight. Tap OK for intermediate reinitialization of the system.
CANNOT PICKUP ERROR -117 INITIALIZING SYSTEM	Locking of Column Holder did not work. Tap OK for intermediate reinitialization of the system.
WRONG SIGNAL ERROR -118 CALL TECHNICAL SERVICE	Electronic hardware problem. Call Miltenyi Biotec Technical Support.
LOCK FAILURE ERROR -123 CALL TECHNICAL SERVICE	Electronic or mechanical hardware problem. Call Miltenyi Biotec Technical Support.
MAGN. SENSOR ERROR -205	Electronic hardware problem. Call Miltenyi Biotec Technical Support.

Table 8.2: Explanation of error code and suggestion for fixing.

First disconnect both the power and the communication cable. The fuse holder is located on the rear side of MultiMACS M96 Separator at the power connector (refer to figure 6). Detach the fuse holder with a sharp tool, e.g., a screwdriver. Remove the defect fuse and insert a new fuse with the type and rating 5×20 mm, T2.5AH250V. Reinsert the fuse holder.

Caution: Replace fuses only with same type and rating to ensure protection against fire risk.


Figure 6: Fuse holder on the rear side of the MultiMACS M96 Separator.

For technical support, contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support at Miltenyi Biotec headquarters:

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-3803
technicalsupport@miltenyi.com

Visit **www.miltenyibiotec.com** for local Miltenyi Biotec Technical Support contact information.

Limited warranty

Except as stated in a specific warranty statement which may accompany this product or as otherwise agreed in writing by an authorized representative of Miltenyi Biotec, Miltenyi Biotec's warranty to you, the original purchaser and end user ("you" or "your"), with respect to the product accompanied by this limited warranty shall be subject to the following provisions and the general terms and conditions of sale of the company within the Miltenyi Biotec group which supplied the product in effect at the date of purchase. Those terms and conditions of sale may vary by country and region. Nothing in this document should be construed as constituting an additional warranty.

Miltenyi Biotec warrants that this product will operate or perform substantially in conformance with Miltenyi Biotec's published specifications and be free from material defects in material and workmanship, when subjected to normal, proper and intended usage by properly trained personnel, for the period of time set forth in the product documentation or package inserts accompanying the product (the "Warranty Period").

Miltenyi Biotec agrees, during the Warranty Period, to repair or replace, at Miltenyi Biotec's option, the defective product so as to cause the same to operate in substantial conformance with said published specifications; provided that you shall (a) promptly notify Miltenyi Biotec in writing upon the discovery of any nonconformity or defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) return the nonconforming or defective product to Miltenyi Biotec, freight prepaid, only after receipt of a Return Material Authorization ("RMA") from Miltenyi Biotec, which may include biohazard decontamination procedures and other product-specific handling instructions, if applicable.

Miltenyi Biotec shall have no obligation to make repairs, replacements or corrections to the product or any component thereof required, in whole or in part, as the result of (i) normal wear and tear, (ii) improper handling, installation, operation, storage, service, maintenance or repair, (iii) failure to follow the instructions, cautions, warnings, and notes set forth in the product documentation provided with the product or provided by Miltenyi Biotec from time to time, (iv) abnormal use, misuse, neglect, abuse, mishandling, misapplication, modification or alteration of the

product, (v) use of the product in a manner for which it was not designed, (vi) causes external to the product such as, but not limited to, power failure or electrical power surges, (vii) use of the product in combination with equipment, accessories, consumables or software not supplied or approved by Miltenyi Biotec, or (viii) accident, disaster or acts of God. ANY INSTALLATION, MAINTENANCE, REPAIR, SERVICE OR ALTERATION TO OR OF, OR OTHER TAMPERING WITH, THE PRODUCT PERFORMED BY ANY PERSON OR ENTITY OTHER THAN MILTENYI BIOTEC AUTHORIZED PERSONNEL WITHOUT MILTENYI BIOTEC'S PRIOR WRITTEN APPROVAL, OR ANY USE OF REPLACEMENT PARTS NOT SUPPLIED BY MILTENYI BIOTEC, SHALL IMMEDIATELY VOID AND CANCEL ALL WARRANTIES WITH RESPECT TO THE AFFECTED PRODUCT.

Miltenyi Biotec's warranty does not cover products sold AS IS or WITH ALL FAULTS, or which had its serial number defaced, altered or removed, or any consumables, or parts identified as being supplied by a third party.

Miltenyi Biotec must be informed promptly if a claim is made under this limited warranty. If a material or manufacturing defect occurs within the Warranty Period, Miltenyi Biotec will take the appropriate steps, at Miltenyi Biotec's option, to make repairs, replacements or corrections to the product or any component thereof as may be required to restore the full usability of your product. If Miltenyi Biotec determines that a product for which you have requested warranty services is not covered by the warranty hereunder, you shall pay or reimburse Miltenyi Biotec for all costs of investigating and responding to such request at Miltenyi Biotec's then prevailing time and materials rates. If Miltenyi Biotec provides repair services or replacement parts that are not covered by this warranty, you shall pay Miltenyi Biotec therefor at Miltenyi Biotec's then prevailing time and materials rates.

THE OBLIGATIONS CREATED BY THIS WARRANTY STATEMENT TO REPAIR OR REPLACE A DEFECTIVE PRODUCT ARE EXCLUSIVE AND SHALL BE THE SOLE REMEDY OF BUYER IN THE EVENT OF A DEFECTIVE PRODUCT. EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY STATEMENT, MILTENYI BIOTEC HEREBY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, ORAL OR WRITTEN, WITH RESPECT TO THE PRODUCT, INCLUDING WITHOUT LIMITATION ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. MILTENYI BIOTEC DOES NOT WARRANT THAT THE PRODUCT IS ERROR-FREE OR WILL ACCOMPLISH ANY PARTICULAR RESULT.

UNDER NO CIRCUMSTANCES WILL MILTENYI BIOTEC BE LIABLE FOR ANY LOSS OF USE, INTERRUPTION OF BUSINESS OR ANY INDIRECT, SPECIAL, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES OF ANY KIND (INCLUDING LOST PROFITS) REGARDLESS OF THE FORM OF ACTION WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE), STRICT PRODUCT LIABILITY OR OTHERWISE, EVEN IF MILTENYI BIOTEC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. IN NO EVENT SHALL THE TOTAL LIABILITY OF MILTENYI BIOTEC HEREUNDER EXCEED THE GREATER OF 100.00 EUROS OR THE AMOUNT YOU ACTUALLY PAID FOR THE PRODUCT GIVING RISE TO SUCH LIABILITY, REGARDLESS OF THE CAUSE OF ACTION, IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE. NOT ALL JURISDICTIONS ALLOW SUCH LIMITATIONS OF DAMAGES SO THE FOREGOING LIMITATIONS MAY NOT APPLY TO YOU.

This warranty statement gives you specific legal rights and you may have other rights, which may vary from jurisdiction to jurisdiction.

Trademarks

µMACS, MACS, the Miltenyi Biotec logo, and MultiMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

MultiMACS™ M96 Separator

User manual

Original instructions

Issued: 2022-10

140-002-070.04

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

Phone +49-2204-8306-0
Fax +49-2204-85197
macsde@miltenyi.com

www.miltenyibiotec.com

Read the chapter Important safety information and all other information in this user manual before using the instrument

Before using the MultiMACS M96 Separator, read the chapter **Important safety information** and all other information contained in this user manual, including any safety and operating instructions.

Pay special attention to all warnings displayed on the instrument. Failure to read and follow these guidelines could lead to improper or incorrect usage and result in damage to the instrument. Improper usage could also cause severe personal injury, death, unpredictable results, instrument malfunction, and premature wear to components shortening the lifetime of the instrument. Such actions may void your warranty. Keep the user manual and any other safety and operating instructions provided with the instrument in a safe place accessible to all users for future reference.

If you have a serious concern regarding the safe use of your instrument, contact your authorized Miltenyi Biotec service provider or call Miltenyi Biotec Technical Support. Documentation needs to be consulted before proceeding with installation and operation of the system.