

Contents

1. Description
 - 1.1 Important safety information
 - 1.2 Background information
 - 1.3 Technical specifications
 - 1.4 Applications
 - 1.5 Reagent and instrument requirements
2. Preparation of M Columns

1. Description

This product is for research use only.

For use by professional laboratory personnel only.

Components 10 M Columns for the isolation of magnetically labeled molecules; sterile packed.

Storage Store columns dry at +10 to +35 °C and protected from light. Do not use after expiration date.

1.1 Important safety information

- ⚠ WARNING** Contamination or infection could result in death or serious injury depending on the material used.
- All biological material must be considered potentially infectious
- Regulations for the treatment and disposal of infectious materials must be observed.

1.2 Background information

The patented MACS® Column Technology is based on the use of MACS MicroBeads, MACS Columns and MACS Separators. M Columns have been developed for the gentle isolation of MicroBead-labeled molecules. As MACS MicroBeads are extremely small, superparamagnetic particles, a high-gradient magnetic field is required to retain the labeled molecules. M Columns contain an optimized matrix to generate this strong magnetic field when placed in a permanent magnet such as the MiniMACS™, OctoMACS™, or SuperMACS™ II Separator. Washing and elution steps are performed by simply rinsing the column with an appropriate buffer as described in the individual µMACS™ MicroBeads data sheet or as tested experimentally.

1.3 Technical specifications

- For capacity of the M Column, refer to the individual µMACS MicroBeads data sheet.
- Columns are "flow stop" and do not run dry.
- Void volume: 80 µL.
- Recommended filling volume: 1 mL.
- Typical flow rate for PBS: 300 µL/minute.
- The columns are for single use only.

1.4 Applications

M Columns are used for molecular biology and protein biochemistry applications such as isolation of mRNA, immunoprecipitated protein, or epitope-tagged protein in combination with µMACS MicroBeads and a MiniMACS, OctoMACS, or SuperMACS II Separator.

▲ M Columns are for molecule isolation only. Do not use M Columns for cell separation.

▲ Do not use M Columns in combination with magnetic particles other than µMACS MicroBeads. Magnetic forces in the column are very high and may damage biological material if other beads are used.

▲ M Columns are not suitable for particles larger than 30 µm. To remove clumps and to prevent aggregates in sample, resuspend material carefully and precipitate clumps by centrifugation before applying the sample on the column.

▲ Samples or buffers with high viscosity might cause reduced column flow or column clogging.

1.5 Reagent and instrument requirements

- Equilibration buffer for column preparation: Buffer supplemented with 1% detergent, e.g. SDS or Ecosurf™ EH-9. Degas buffer before use, as air bubbles could block the column

- Separation buffer: Any buffer suitable for magnetic separation (for details, see individual µMACS MicroBeads data sheet). If the separation buffer contains 1% detergent, it can also be used for column preparation. Degas buffer before use, as air bubbles could block the column

▲ **Note:** Use degassed buffer only! Degas buffer by applying vacuum or sonification for ten minutes, preferentially with buffer at room temperature. Excessive gas in running buffer will form bubbles in the matrix during isolation. This is particularly important, when the applied buffer has a different temperature as the M Column, e.g. when using cold buffer on a column at room temperature. Air bubble formation in the M Column may lead to clogging of the column and decrease the quality of isolation.


- µMACS MicroBeads for magnetic labeling of target molecules
- MiniMACS Separator (# 130-042-102), OctoMACS Separator (# 130-042-109), or SuperMACS II Separator (# 130-044-104)
- MACS MultiStand (# 130-042-303) in combination with MiniMACS Separator or OctoMACS Separator

- Adapter for MS, LS, and LD Columns for use with SuperMACS II Separator
- MACS Acrylic Tube Rack (# 130-041-406) or OctoMACS Acrylic Tube Rack (# 130-090-448)

2. Preparation of M Columns

1. Insert M Column with the column wings to the front into MACS Separator according to A) or B).

A) Use with MiniMACS or OctoMACS Separator



Attach MiniMACS or OctoMACS Separator to the MultiStand and place M Column in the separator. Place a collection tube under the M Column.

▲ **Note:** Check that the ejection blocks in the gap of the magnet are attached before placing the MACS Column into the magnetic field of the MiniMACS or OctoMACS Separator.

B) Use with SuperMACS II Separator

For use of M Columns with the SuperMACS II Separator please refer to the respective data sheet.

2. Apply 250 μ L of degassed Equilibration Buffer on top of the column and let the solution run through.
3. If separation buffer is different from equilibration Buffer, apply 100 μ L of degassed separation buffer and let the solution run through. The M Column is now ready for magnetic separation. Perform the magnetic separation as indicated on the individual μ MACS MicroBeads data sheets.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

MACS, the Miltenyi Biotec logo, MiniMACS, OctoMACS, SuperMACS, and μ MACS are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2024 Miltenyi Biotec and/or its affiliates. All rights reserved.