

Contents

1. Description
 - 1.1 Principle of the Neurosphere Dissociation Kit (P)
 - 1.2 Background information
 - 1.3 Applications
 - 1.4 Reagent and instrument requirements
2. Protocol
 - 2.1 Reagent and instrument preparation
 - 2.2 Neurospheres dissociation protocol
3. Appendix

1. Description

This product is for research use only.

Components	Neurosphere Dissociation Kit (P) 6 vials, containing: 2.5 mL of Enzyme P 2×50 mL of Buffer X (sterile) 1.5 mL of Buffer Y (sterile) 1 vial (5–15 mg*) of Enzyme A 1 mL of Buffer A
Size	For 100 digestions of 1 mL.
Storage	Upon arrival immediately store Enzyme P in aliquots at –20 °C. Store all other components at 2–8 °C upon arrival. Reconstitute Enzyme A before the date indicated on the vial label. For information about reconstitution and storage after reconstitution of the lyophilized components refer to chapter 2.1. *The vials contain a consistent amount of units. Variability in ranges of mg within the vial is possible – according to variability of activity per mg

1.1 Principle of the Neurosphere Dissociation Kit (P)

The Neurosphere Dissociation Kit (P) provides a simple and fast procedure to enzymatically dissociate mouse, rat, or human neurospheres into single cells.

After harvesting, the neurospheres are digested enzymatically and further mechanically dissociated into a single-cell suspension by trituration.

1.2 Background information

The Neurosphere Dissociation Kit (P) has been designed for a gentle but rapid and efficient generation of single-cell suspensions from cultivated neurospheres. Neural stem cells isolated from whole mouse brain proliferate under special culture conditions forming clonally derived spheres, called neurospheres. Dissociated cells can be quantified and frequently subcultured to generate more neurospheres, as well as analyzed *in vitro* for phenotype distributions or other functional studies.

1.3 Applications

- Dissociation of neurospheres to single-cell suspensions for subsequent cultivation and expansion of neurospheres.
- Phenotyping, functional, and genetic studies of neural stem cells by flow cytometry or fluorescence microscopy.

1.4 Reagent and instrument requirements

- Hanks' Balanced Salt Solution (HBSS) with Ca²⁺ and Mg²⁺ (Sigma-Aldrich # 55037C), in the following referred to as HBSS (w)
- 2 mL, 15 mL, and 50 mL tubes
- Pre-Separation Filters (30 µm) (# 130-041-407)
- (Sterile) glass Pasteur pipettes
- MACSmix™ Tube Rotator (# 130-090-753) in combination with an incubation oven at 37 °C

2. Protocol

2.1 Reagent and instrument preparation

▲ Volumes given below are for the dissociation of 40 mL of neurosphere suspension. When working with less than 40 mL, use the same volumes as indicated. When working with more than 40 mL, scale up all reagent volumes and total volumes accordingly.

1. Resuspend the lyophilized powder in the vial labeled Enzyme A with 1 mL of Buffer A. Do **not** vortex. This solution should then be aliquoted and stored at -20°C for later use.
2. Fire-polish three glass Pasteur pipettes so that decreasing tip diameters are achieved. For details refer to 3. Appendix.

2.2 Neurospheres dissociation protocol

▲ All steps should be performed under sterile conditions.

▲ The MACSmix Tube Rotator is used with slow, continuous rotation.

1. Pre-heat 960 μL of Buffer X for up to 40 mL of neurosphere suspension at 37°C for 10–15 minutes before use.
2. Detach neurospheres by tapping the flask against the benchtop and pass a stream of media across the attached neurospheres. Harvest suspended neurospheres and transfer the suspension to an appropriate-sized sterile tube (50 mL).
3. Centrifuge at $300\times g$ for 10 minutes at room temperature and aspirate the supernatant carefully.
4. Resuspend the cell pellet in 960 μL of pre-heated Buffer X.
5. Add 25 μL of Enzyme P, 10 μL of Buffer Y, and 5 μL of Enzyme A to the neurosphere suspension and mix gently.
6. Transfer the neurosphere suspension to a 2 mL tube.
7. Incubate the closed tube for 10 minutes at 37°C under slow, continuous rotation using the MACSmix™ Tube Rotator.
8. Dissociate neurospheres mechanically using the wide-tipped, fire-polished Pasteur pipette by pipetting up and down 10 times slowly. Avoid forming air bubbles.

▲ **Note:** To remove cells stuck to the tube wall, rinse the tube wall once or twice with the sample solution.
9. Incubate at 37°C for 5 minutes using the MACSmix Tube Rotator.
10. Dissociate neurospheres mechanically using the other two fire-polished Pasteur pipettes of decreasing diameter. Pipette slowly up and down 10 times with each pipette, being careful to avoid the formation of air bubbles.
11. Apply the single-cell suspension to a Pre-Separation Filter (30 μm) placed on a 15 mL tube.
12. Apply 5 mL of HBSS (w) through the Pre-Separation Filter (30 μm).
13. Discard filter and centrifuge cell suspension at $300\times g$ for 10 minutes at room temperature. Aspirate the supernatant completely.
14. Cells should be processed immediately for further applications or resuspend in fresh media for subsequent cultivation after determination of the cell number.




Figure 1: Protocol overview

3. Appendix: Tips & hints

▲ For up-to-date information regarding antigen compatibilities of the Neurosphere Dissociation Kit (P) for subsequent MACS® Cell Separations, please refer to www.miltenyibiotec.com/130-095-943.

Production of appropriate Pasteur pipettes

For this protocol three Pasteur pipettes with openings of decreasing diameter are needed. The opening of the first pipette should be rounded without significant decrease in the size of the opening. The smallest opening should not be smaller than 0.5 mm so that the cells are not forced through with too much pressure. To produce openings that get progressively smaller, rotate the pipettes quickly in the flame to fire polish them for a few seconds. Production is easier if you apply the rubber sucker. Too much time may fuse the hole. The edges should be rounded.

Yield of viable cells is too low (dissociation is insufficient)

Make sure that the neurospheres are agitated sufficiently during the entire time of incubation and do not stick to the bottom of the tube. Flick the tube after adding the enzyme mix if it is necessary. During the working steps at 37 °C the MACSmix™ Tube Rotator is convenient for this purpose. Cell clumps will not get through the Pre-Separation Filter. Therefore, keep pipetting until the single-cell suspension looks homogeneous. Then apply suspension to the Pre-Separation Filter.

Single-cell suspension contains only dead cells

Make sure the openings of the Pasteur pipettes are not too small. Pipette more slowly and do not vortex the cells. Avoid forming bubbles. Follow the protocol non-stop.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

MACS, MACSmix, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.