Miltenyi Biotec

# gentleMACS™ Dissociator

User manual

Copyright © 2022 Miltenyi Biotec. All rights reserved.

No part of this user manual may be reproduced, stored in a retrieval system, transmitted, published, or distributed in any form or by any means, electronically, mechanically, by photocopying, microfilming, recording, or otherwise, without the prior written consent of Miltenyi Biotec; however, notwithstanding the foregoing, the owners of the gentleMACS Dissociator may make copies solely for purposes of training personnel in the use and servicing of the unit within their business or organization.

Maximal care has been taken by Miltenyi Biotec in the preparation of this user manual. However, Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein, or for incidental or consequential damages in connection with the furnishing, performance, or use of this document. The information in this document is provided “as is” without warranty of any kind, and is subject to change without notice.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. gentleMACS, MACS, and the Miltenyi Biotec logo are either registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

# gentleMACS™ Dissociator

## User manual

Original instruction

Issued: 2022-10

140-002-197.09

Miltenyi Biotec B.V. & Co. KG  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany

Phone +49 2204 8306-0  
Fax +49 2204 85197  
macsde@miltenyi.com

[www.miltenyibiotec.com](http://www.miltenyibiotec.com)

# Contents

<b>Chapter 1</b>	<b>7</b>	<b>Important information</b>
	<b>11</b>	<b>Wichtige Informationen</b>
	<b>17</b>	<b>Informations importantes</b>
	<b>23</b>	<b>Informazioni importanti</b>
	<b>29</b>	<b>Información importante</b>
<b>Chapter 2</b>	<b>34</b>	<b>Introduction</b>
	34	Intended use
<b>Chapter 3</b>	<b>35</b>	<b>gentleMACS Products and Spare Parts</b>
<b>Chapter 4</b>	<b>37</b>	<b>Storage and handling</b>
<b>Chapter 5</b>	<b>38</b>	<b>Installation</b>
<b>Chapter 6</b>	<b>39</b>	<b>Operating the gentleMACS Dissociator</b>
	<b>39</b>	<b>General information</b>
	<b>41</b>	<b>Description of the instrument</b>
	41	Layout of the keyboard
	41	Panel buttons
	42	gentleMACS Program Card
	42	Function of the LED
	42	Display
	42	Power on
	<b>43</b>	<b>Function of the menus</b>
	43	Select and run a program
	44	Copy a program
	45	Delete a program
	<b>45</b>	<b>Emergency shutdown procedure</b>
<b>Chapter 7</b>	<b>46</b>	<b>Performance and limitations of use</b>
<b>Chapter 8</b>	<b>47</b>	<b>Maintenance</b>
	<b>47</b>	<b>Cleaning</b>
	48	Cleaning procedure
	48	Disassembly for cleaning and reassembly
	<b>49</b>	<b>Safety</b>
	<b>50</b>	<b>Internal quality control</b>
	<b>50</b>	<b>Disposal information</b>
<b>Chapter 9</b>	<b>51</b>	<b>Troubleshooting</b>
<b>Chapter 10</b>	<b>52</b>	<b>Technical specifications</b>
<b>Chapter 11</b>	<b>53</b>	<b>Limited warranty</b>
<b>Chapter 12</b>	<b>55</b>	<b>Technical support</b>
<b>Chapter 13</b>	<b>56</b>	<b>Declaration of conformity</b>

## **Please read before use!**

Please read all information contained in this chapter before use. The use of an instrument for a task other than that for which it is intended, as well as improper or insufficient maintenance can greatly reduce the life time of the instrument and will invalidate the instrument's warranty. Incorrect handling and care as well as misuse can lead to premature wear or can cause hazards to users. The following chart is an illustrated glossary depicting the symbols that are used in this user manual and on the gentleMACS Dissociator.

## **Symbols**

### **Warning**

Triangles indicate general warnings. Warnings include the risk of damage to the equipment, severe personal injury, or loss of life. Consult the user manual for further instructions and proceed with caution.

### **Warning**

This symbol indicates the risk of loss of life or severe injury to the instrument operator due to high voltage. Consult the user manual for further instructions and proceed with caution.

### **Warning**

This symbol indicates the risk of loss of life, severe injury to the instrument operator, or equipment damage due to potentially dangerous chemical or biological material. Consult the user manual for further instructions and proceed with caution.

### **Caution**

This symbol indicates the risk of equipment damage if instructions in the user manual are not followed precisely. It is therefore strongly recommended that the instrument operator consult the corresponding chapters in the respective documentation.

### Protective conductor terminal

Symbol is attached on the inside of the instrument. Warning for service personnel only.

**On (supply)**

**Off (supply)**


**Fuse**

**WEEE** (Waste of Electrical and Electronic Equipment)

All operating instructions must be read before installing and operating the instrument.

**WARNING!** Avoid serious injury. Read and understand operating instructions before installing or use. If manual is missing, obtain replacement from Miltenyi Biotec.

European conformity marking

United Kingdom Conformity Assessed marking

NRTL certification mark: Product meets consensus-based standards of safety, required by the Occupational Safety/Health Administration (OSHA), determined by the Nationally Recognized Testing Laboratories (NRTL) TÜV SÜD

## Warnings and precautions

This chapter contains information for the safe operation of the gentleMACS Dissociator. Detailed instructions regarding safety are provided throughout this manual, marked with the symbols described above.

The operator is responsible for the correct use of the gentleMACS Dissociator following the instructions given in this manual and in full accordance with all local, state, and federal laws that may apply. The impeccable and safe operation of the instrument presumes that

any operator understands the user manual and pays attention to all safety information.

**Keep this user manual in a safe place, accessible for anyone using the gentleMACS Dissociator.**

**At all times, all safety instructions in this user manual, local working area safety instructions, laboratory policies, and standards regarding laboratory safety and health, and prevention of accidents must be adhered to.**

In the case of safety-relevant defects of the equipment, Miltenyi Biotec technical service must be contacted.

To ensure safe operation of the gentleMACS Dissociator all covers or parts that have been removed from the gentleMACS Dissociator must be re-installed prior to setting the instrument into operation. Do not allow fluids or objects to enter the interior of the instrument.

Incorrect handling of the gentleMACS Dissociator bears risks in the form of mechanical, electrical, chemical, and biological hazards.

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## Electrical hazards

Please be particularly aware of electrical hazards. Standard electrical safety precautions must be applied at all times. Only qualified personnel are permitted to perform electrical servicing. It is strictly forbidden to operate the instrument with disassembled or damaged parts.

The instrument is equipped with a three-wire electrical grounding-type plug. This is a safety feature. Do not insert the plug into a non-grounded power outlet. If you cannot insert the plug into the outlet or if you are not sure whether the outlet is grounded, contact your local electrician.

The instrument should only be operated from a power source indicated on the product's electrical ratings label. If you have questions about the type of power source to use, contact your authorized Miltenyi Biotec service provider or local power company.

Make sure that the mains switch as well as the connector for the power cable are easily accessible and located as close to the operator of the instrument as possible. If it is necessary to disconnect the power supply, unplug the cable from the power outlet.

After moving the gentleMACS Dissociator from a cold environment, for example, a cold room at 4 °C to room temperature, it is necessary to wait a minimum of 4 hours before setting the instrument into operation. Condensing liquid droplets can severely damage the instrument.

## Chemical and biological hazards

The proper handling of the gentleMACS Tubes minimizes exposure to the samples. However, potential hazards to the operator still exist. The operator must choose and wear personal safety equipment in accordance with warnings and precautions for the used substances. Defective or inadequate safety equipment might endanger the operator.

The gentleMACS Dissociator must be operated within a safety enclosure when hazardous or unknown materials are processed, to protect persons and the laboratory environment from contamination.

All liquid and solid waste must be considered biohazardous and, therefore, must be handled taking universal laboratory precautions. All clinical samples must be considered potentially infectious. Toxic or corrosive chemicals, poisonous, acidic, or radioactive substances may be present.

It is strictly prohibited for anyone other than the operator to handle contaminated accessories or contaminated parts of the gentleMACS Dissociator.

## Servicing

Unless otherwise mentioned in the Miltenyi Biotec documentation, do not service the equipment yourself. Service needed on components inside the instrument must be carried out only by personnel trained and authorized by Miltenyi Biotec.

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1

# 2 Introduction

The gentleMACS™ Dissociator is a benchtop instrument for the automated dissociation of tissues to obtain single-cell suspensions. Single-cell suspensions are a crucial prerequisite to achieve reliable results in cell separation, cell analysis, and cell culture. Moreover, the gentleMACS Dissociator can be used for the homogenization of cells and tissues to extract biomolecules for molecular biology or protein biochemistry applications.

Two types of gentleMACS Tubes have been specially designed: C Tubes for the gentle preparation of single-cell suspensions from tissues and M Tubes for the homogenization of tissues and cells. Both types of tubes are designed for safe and sterile sample handling.

A number of gentleMACS Programs have been specifically developed to offer standardized procedures for the gentle dissociation or homogenization of various tissues. A single sample or two samples in parallel can be processed.

A comprehensive gentleMACS Program overview and specific protocols for numerous applications are available for download at [www.gentlemacs.com](http://www.gentlemacs.com).

## Intended use

The gentleMACS™ Dissociator is used as a stand-alone, benchtop, research use only instrument for the generation of single cell suspensions from solid tissue or homogenization of tissues or cells to isolate subcellular material. A number of pre-defined programs offer a variety of processes for tissue dissociation and tissue homogenization. A single sample or two samples in parallel can be processed in gentleMACS Tubes.

The gentleMACS Dissociator is for research-use only in laboratory environment. It is to be used by laboratory professional users only.



## gentleMACS Products

Product	Description	Order no.
<b>gentleMACS Dissociator</b>	gentleMACS Dissociator 1 Power cord 1 gentleMACS Dissociator user manual	130-093-235
<b>C Tubes</b>	For the dissociation of tissues to obtain single-cell suspensions. 25 tubes, single-packed	130-093-237
	100 tubes, 4×25 tubes per bag	130-096-334
<b>M Tubes</b>	For the homogenization of tissues or cells to isolate subcellular material. 25 tubes, single-packed	130-093-236
	100 tubes, 4×25 tubes per bag	130-096-335
<b>M Tubes with Strainer</b>	50 tubes per bag (strainer: 600 µm mesh)	130-094-392

For more information, visit [www.gentlemacs.com](http://www.gentlemacs.com).

## gentleMACS Spare Parts

Product	Description	Order no.
<b>gentleMACS Sleeves</b>	2 adapters for attachment of gentleMACS Tubes	130-093-464
<b>gentleMACS Trough</b>	Safety container for spilled liquid	130-093-465
<b>gentleMACS Safety Shield</b>	Safety shield for protection	130-093-466
<b>gentleMACS Program Card</b>	Memory card containing gentleMACS Programs	130-093-462

Only use accessories recommended by Miltenyi Biotec. Failure to use recommended spare parts and disposables may compromise safe operation and may result in inaccurate results and reduced life time of the instrument. Miltenyi Biotec does not honor any warranty or accept any responsibility for damages resulting from the use of inappropriate accessories.

The gentleMACS Dissociator must be stored and operated under dry conditions (non-condensing humidity). Storage temperature must be between  $-15\text{ }^{\circ}\text{C}$  and  $+70\text{ }^{\circ}\text{C}$ . The instrument should be operated at a temperature between  $+4\text{ }^{\circ}\text{C}$  and  $+35\text{ }^{\circ}\text{C}$ .

# 5 Installation

Please unpack the instrument carefully and check for completeness and for any damage that may have occurred during transportation. If there is any damage, contact your local Miltenyi Biotec representative or Miltenyi Biotec technical service. Please keep the original packaging material.

The gentleMACS Dissociator is designed exclusively for indoor use.

## Components included in delivery

- 1 gentleMACS Dissociator
- 1 Power cord
- 1 gentleMACS Dissociator user manual

## Installation

The gentleMACS Dissociator requires a flat surface with the dimensions 250 mm × 400 mm and a power connection. To operate the instrument, the keyboard area and the top of the instrument must be accessible. The gentleMACS Dissociator is a ready-to-use instrument.

# Operating the gentleMACS Dissociator 6

The instrument must be used only with properly installed gentleMACS Tubes whose lids are tightly closed. Do not operate the instrument with gentleMACS Tube lids only.

## General information

The gentleMACS Dissociator is used as a stand-alone, benchtop instrument. A number of pre-defined programs offer a variety of processes for tissue dissociation and tissue homogenization. The specially designed gentleMACS Tubes are required for processing the samples. The gentleMACS Dissociator must be used exclusively with these Tubes.

C Tubes are designed for the dissociation of tissues to obtain single-cell suspensions. M Tubes are designed for homogenization of tissues or cells to isolate subcellular material. A single sample or two samples in parallel can be processed.


For details on the sample preparation, see the respective gentleMACS Protocols, which are available at [www.miltenyibiotec.com/protocols](http://www.miltenyibiotec.com/protocols).

Below, please find brief instructions for the use of the instrument.  
For detailed information, see the next sections.

1. Make sure that the gentleMACS Dissociator is switched on.  
The mains switch is located at the back panel of the instrument.
2. Pre-defined programs are provided by the internal gentleMACS Dissociator memory.  
Alternatively: Insert the gentleMACS Program Card to run a program from the Program Card or to copy a program from the Program Card to the internal memory of the instrument.
3. Transfer the samples into the Tubes. Make sure that the Tubes are securely closed.
4. Make sure to properly install the Tubes vertically and upside down onto the Sleeves. Proper installation will ensure that the Tubes are held in position against rotational and axial forces.
5. All functions of the gentleMACS Dissociator must be controlled via the display; see the following sections for details.
6. Select one of the programs offered by the instrument memory or by the gentleMACS Program Card.
7. Start the process by pressing the start button.
8. The display indicates the program name and the remaining time (seconds) for the current program to finish.
9. The end of the process is displayed for 5 seconds. Subsequently, the display indicates the program that has just been used.  
By pressing the start button, this program can be applied immediately to process the next sample.
10. Remove the Tubes containing the sample by pulling it vertically from the Sleeve.

### Layout of the keyboard

### Panel buttons

-  Selection of one of the three modes that allow the user to select a program for sample processing, or copy, or delete programs
-   Selection of a program
-  Confirmation of a selection or function
-   Run or stop a process

## gentleMACS Program Card

Insert the gentleMACS Program Card in the direction indicated by the arrow. The card will then automatically be selected to access programs.

If the gentleMACS Program Card is inserted, only the programs from the Program Card are available. The programs from the Program Card can be copied to the internal gentleMACS Dissociator memory. To select programs from the internal memory, remove the gentleMACS Program Card.

## Function of the LED

The LED indicates the instrument status for convenient monitoring.

<b>Green:</b>	System ready
<b>Green flashing:</b>	Process running
<b>Red:</b>	System error
<b>Red flashing:</b>	Error during process

## Display

Using the menu button, the operator can switch between **select**, **copy**, and **delete** modes.

<Program name> <Tube name>	Select mode
Copy <Program name>	Copy mode
Delete <Program name>	Delete mode


## Power on

gentleMACS V1.xx Miltenyi Biotec	Display during initialization.
<Program name> <Tube name>	After initialization, the green LED is on and the gentleMACS Dissociator is ready for operation.

## Select and run a program

To select a program from the internal gentleMACS Dissociator memory, press the menu button until the select menu is displayed. If the gentleMACS Program Card is inserted, only the programs from this card are available.


To select programs from the gentleMACS Dissociator memory, remove the gentleMACS Program Card. To copy programs from the gentleMACS Program Card to the gentleMACS Dissociator memory, see the section Copy a program.

# 6

## Copy a program

To copy a program from the gentleMACS Program Card to the gentleMACS Dissociator memory, insert the gentleMACS Program Card. Press the menu button until the copy menu is displayed. If the gentleMACS Program Card is not inserted, a program from the gentleMACS Dissociator memory will be selected.

Select the desired program from the gentleMACS Program Card using the arrow buttons. If a program is copied within the internal memory, "Copy memory" is displayed instead of "Copy card".

Press **ok**.

Remove the gentleMACS Program Card.

Select the destination on the gentleMACS Dissociator memory using the arrow buttons. If there is already a program stored at the selected destination, the name of this program is displayed. If the destination is empty, the line remains blank.

Press **ok** to store the program or press **menu** to cancel and select a different program.


If there is already a program stored at the selected destination, confirm "overwrite?" by pressing **ok** or cancel by pressing **menu**.

The program is being copied.

If an error occurs, confirm the error message by pressing **stop**.

## Delete a program

To delete a program from the gentleMACS Dissociator memory, press the menu button until the delete menu is displayed.

## Emergency shutdown procedure

1. Press the **stop** button.
2. Press the mains switch at the back panel of the instrument.
3. Unplug the instrument.

# 7 Performance and limitations of use

- A maximum of 99 programs can be stored on the internal gentleMACS Dissociator memory.
- Editing programs on the gentleMACS Dissociator is not possible.
- If an error occurs, the gentleMACS Dissociator displays a message that needs to be confirmed by the operator. Press the **stop** button in order to continue operation.
- Both motors run at the same time, theoretically at identical speed and acceleration. However, the friction and resistance of the sample will decrease the effective speed. The gentleMACS Dissociator will display an error message if a pre-set tolerance has been exceeded during the process.
- If one of the motors is blocked due to high resistance of the sample in the tube, the instrument will try four times to get started again. The message "Error: Motor 1 / 2 blocked" appears and the red LED flashes. The operator must confirm the incident by pressing the **stop** button.
- If a motor blockage can be resolved, the gentleMACS Dissociator continues the procedure from where it was interrupted. The message "Warning: x Retries" appears for 5 seconds during the incident and disappears automatically.
- If a motor blockage cannot be resolved, the power of the motor is turned off.

## Cleaning

In case of spillage or contamination, clean the instrument appropriately. For decontamination, use 70% ethanol or a 1% hypochlorite solution.


- Do not allow fluids or objects to enter the interior of the gentleMACS Dissociator.
- Make sure to clean the instrument exclusively with cleaning agents recommended by Miltenyi Biotec. Do not use acetone or any aggressive agents.
- Do not leave Tubes attached to the instrument while it is not in use.
- When processing infectious, radioactive, poisonous, or any other hazardous liquids, always abide by the necessary safety precautions.
- Before using organic solvents or aggressive chemicals in the dissociation or homogenization process, the compatibility with the gentleMACS Dissociator parts and Tubes must be tested.
- Large amounts of accidentally spilled liquid are automatically collected in the gentleMACS Trough. Use a pipette to aspirate the liquid before removing the Trough for cleaning and sterilizing.
- If servicing of the gentleMACS Dissociator is required, send the instrument back to a service provider authorized by Miltenyi Biotec. Make sure that all parts are decontaminated and sterilized. Use the original packaging material.
- Always use your personal safety equipment while removing or disassembling parts of the gentleMACS Dissociator for cleaning or sterilization.

## Cleaning procedure

For details on disassembly and reassembly, see the picture below.

1. Remove the Safety Shield.
2. Remove the Sleeves from the instrument for cleaning and/or sterilization.
3. Remove the Trough for cleaning and/or sterilization.
4. The remaining surface of the gentleMACS Dissociator should be cleaned with a cloth soaked with a detergent-containing aqueous solution. Carefully wipe all parts.
5. Repeat the procedure with another cloth soaked with distilled water.
6. If the homogenized substances are not water-soluble, it is recommended to wipe all parts with ethanol or an alternative organic solvent. Contact Miltenyi Biotec technical support for further advice.
7. The removed parts should be washed in a laboratory dishwasher. It is recommended to use a cleaning agent containing detergent and alkaline cleaning solution.
8. If necessary, exchange the Sleeves and the Trough.
9. Reassemble the parts and make sure that the removed parts snap into position to ensure the correct positioning.

## Disassembly for cleaning and reassembly

- Safety Shield, Sleeves, and Trough are removed by pulling precisely vertically with a slight force.
- The Trough can be detached easily by slightly pressing the rippled sides.
- Do not use any tools to detach the described parts.
- Reassemble the parts in reverse order.
- Correct assembly of all parts can be sensed by a slight snap as the part slips into its final position.

The gentleMACS Dissociator and the removable parts must be cleaned as described above. The lifetime of the Sleeves might be limited due to the mechanical strain during the homogenization procedure. Periodically, check the parts for wear.




Decontamination is within the responsibility of the operator of the instrument. Contaminated instruments or removable parts must not be given to any person.

After detachment and reassembly of any removable part, it is necessary to check that the part still functions properly.

## Applied materials

- 1 Housing middle part: PA 66/6 (with 30% glass fibers)
- 2 Housing front part: PA 66/6 (with 30% glass fibers)
- 3 Housing rear part: PA 66/6 (with 30% glass fibers)
- 4 Keyboard: Polyester film
- 5 gentleMACS Safety Shield: PP
- 6 gentleMACS Trough: PA 66/6 (with 30% glass fibers)
- 7 gentleMACS Sleeve: PA 66/6 (with 30% glass fibers)

## Internal quality control

During the lifetime of the instrument no calibration is necessary. To check the quality of the instrument, the gentleMACS Dissociator must be sent to Miltenyi Biotec.

## Disposal information

The gentleMACS Dissociator shall be disposed of according to the biosafety level of its use. Potentially contaminated parts must be sterilized before being disposed of, according to local laboratory regulations.

Decontamination is within the responsibility of the operator of the instrument. Contaminated instruments or contaminated parts of the instrument must not be sent or given to any person.

### **Waste of Electrical and Electronic Equipment (WEEE) customer information**

Please dispose of your end-of-life Miltenyi Biotec products in accordance with the applicable WEEE and hazardous waste disposal legislation, which may differ by country or region.

Electrical and electronic equipment may contain hazardous substances which may have a serious detrimental effect on the environment and/or human health. That is why all equipment must be specifically collected and treated by designated waste facility centres and by qualified WEEE compliance schemes. By ensuring that you dispose of your unwanted electrical and electronic equipment according to the applicable WEEE and hazardous waste disposal legislation, you are helping to preserve our natural resources and protect human health.

Miltenyi Biotec is committed to protecting the environment. Miltenyi Biotec offers product end-of-life return programs in many countries, and partners with licensed WEEE compliance schemes throughout the world. Miltenyi Biotec lets you recycle your end-of-life Miltenyi Biotec equipment free of charge. The terms and availability of this offer vary by geography because of differences in regulatory requirements. Please note that, depending on the type and use of your equipment, additional requirements may apply.

For more information, or if you wish to dispose of your end-of-life Miltenyi Biotec equipment, please contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support.

Error	Source of error	User action
<b>Liquid droplets in Trough</b>	Tube is not closed properly.	Close Tube tightly beyond the first resistance.
	Content of Tube is jammed between cap and Tube.	The top area of cap and Tube must be perfectly clean for proper functionality of the seal.
	Tube seal is leaky.	Replace Tube. Clean the instrument as described in chapter 8.
<b>Tube vibrates during run</b>	Sleeve is not installed correctly.	Ensure correct installation of the Sleeve. Make sure that the Sleeve snaps into position.
	Inappropriate sample material.	Refer to the respective datasheet for technical specifications of the gentleMACS Tubes.
<b>Tube is not tightly held in position during run</b>	Tube does not fit properly into the snap-fit of the Trough.	Ensure correct installation of the Tube.
	Tube is not closed properly.	Close Tube tightly beyond the first resistance.
<b>gentleMACS Dissociator is without power</b>	Power cord is not connected.	Check power cord.
	Instrument is not switched on.	Check mains switch at the rear panel of the instrument.
	Fuse is blown.	Check fuses at the rear panel of the instrument. Replace fuses only with those of the same type and rating.
<b>gentleMACS Dissociator does not read the Program Card</b>	Dust deposit on Program Card, or Program Card is damaged.	Check Program Card.

If any of the problems persist, please contact Miltenyi Biotec technical support.

# 10 Technical specifications

Technical specifications		
<b>Size</b>	Length × Width × Height	347 mm × 215 mm × 195 mm
<b>Weight</b>		4.7 kg
<b>Main power</b>	Power consumption	max. 320 VA
	Supply voltage	100–240 V~, 50/60 Hz
	Fuse	2× T4AH250V
<b>Motor</b>	Rotational speed	200–4,000 rpm
	Rotation	clockwise or counter clockwise
	Rotational acceleration	100 r/s <sup>2</sup>
	Torque	minimum 150 mNm. Overload restriction over the range of 200–4,000 rpm
<b>Conditions of operation</b>	Working temperature:	+4 °C to +35 °C
	Humidity	15–80 %, relative, non-condensing
	Altitude	maximum 2000 m
	Storage temperature	–15 °C to +70 °C
<b>USB interface</b>	For service purposes only	Pin 1: USB +5 V
		Pin 2: USB –
		Pin 3: USB +
		Pin 4: USB GND

The gentleMACS Dissociator is labeled as a protection class I device and must be plugged into a grounded power outlet.

Supply voltage fluctuations up to +/- 10% of the nominal voltage.  
 Transient overvoltages present on the mains supply: category II.  
 The instrument is suitable for rated pollution degree 2.  
 Emission sound pressure level at workstation < 70 dB(A).  
 Design and specifications are subject to change without notice.

The gentleMACS Dissociator is designed to meet the requirements of the following standard:

- FCC 47 CFR Ch.1 part 15, class B (electromagnetic compatibility)

For other safety considerations, refer to the product label, or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com).

Except as stated in a specific warranty statement which may accompany this product or as otherwise agreed in writing by an authorized representative of Miltenyi Biotec, Miltenyi Biotec's warranty to you, the original purchaser and end user ("you" or "your"), with respect to the product accompanied by this limited warranty shall be subject to the following provisions and the general terms and conditions of sale of the company within the Miltenyi Biotec group which supplied the product in effect at the date of purchase. Those terms and conditions of sale may vary by country and region. Nothing in this document should be construed as constituting an additional warranty.

Miltenyi Biotec warrants that this product will operate or perform substantially in conformance with Miltenyi Biotec's published specifications and be free from material defects in material and workmanship, when subjected to normal, proper and intended usage by properly trained personnel, for the period of time set forth in the product documentation or package inserts accompanying the product (the "Warranty Period").

Miltenyi Biotec agrees, during the Warranty Period, to repair or replace, at Miltenyi Biotec's option, the defective product so as to cause the same to operate in substantial conformance with said published specifications; provided that you shall (a) promptly notify Miltenyi Biotec in writing upon the discovery of any nonconformity or defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) return the nonconforming or defective product to Miltenyi Biotec, freight prepaid, only after receipt of a Return Material Authorization ("RMA") from Miltenyi Biotec, which may include biohazard decontamination procedures and other product-specific handling instructions, if applicable.

Miltenyi Biotec shall have no obligation to make repairs, replacements or corrections to the product or any component thereof required, in whole or in part, as the result of (i) normal wear and tear, (ii) improper handling, installation, operation, storage, service, maintenance or repair, (iii) failure to follow the instructions, cautions, warnings, and notes set forth in the product documentation provided with the product or provided by Miltenyi Biotec from time to time, (iv) abnormal use, misuse, neglect, abuse, mishandling, misapplication, modification or alteration of the product, (v) use of the product in a manner for which it was not designed, (vi) causes external to the product such as, but not limited to, power failure or electrical power surges, (vii) use of the product in combination with equipment, accessories, consumables or software not supplied or approved

by Miltenyi Biotec, or (viii) accident, disaster or acts of God. ANY INSTALLATION, MAINTENANCE, REPAIR, SERVICE OR ALTERATION TO OR OF, OR OTHER TAMPERING WITH, THE PRODUCT PERFORMED BY ANY PERSON OR ENTITY OTHER THAN MILTENYI BIOTEC AUTHORIZED PERSONNEL WITHOUT MILTENYI BIOTEC'S PRIOR WRITTEN APPROVAL, OR ANY USE OF REPLACEMENT PARTS NOT SUPPLIED BY MILTENYI BIOTEC, SHALL IMMEDIATELY VOID AND CANCEL ALL WARRANTIES WITH RESPECT TO THE AFFECTED PRODUCT.

Miltenyi Biotec's warranty does not cover products sold AS IS or WITH ALL FAULTS, or which had its serial number defaced, altered or removed, or any consumables, or parts identified as being supplied by a third party. Miltenyi Biotec must be informed promptly if a claim is made under this limited warranty. If a material or manufacturing defect occurs within the Warranty Period, Miltenyi Biotec will take the appropriate steps, at Miltenyi Biotec's option, to make repairs, replacements or corrections to the product or any component thereof as may be required to restore the full usability of your product. If Miltenyi Biotec determines that a product for which you have requested warranty services is not covered by the warranty hereunder, you shall pay or reimburse Miltenyi Biotec for all costs of investigating and responding to such request at Miltenyi Biotec's then prevailing time and materials rates. If Miltenyi Biotec provides repair services or replacement parts that are not covered by this warranty, you shall pay Miltenyi Biotec therefor at Miltenyi Biotec's then prevailing time and materials rates.

THE OBLIGATIONS CREATED BY THIS WARRANTY STATEMENT TO REPAIR OR REPLACE A DEFECTIVE PRODUCT ARE EXCLUSIVE AND SHALL BE THE SOLE REMEDY OF BUYER IN THE EVENT OF A DEFECTIVE PRODUCT. EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY STATEMENT, MILTENYI BIOTEC HEREBY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, ORAL OR WRITTEN, WITH RESPECT TO THE PRODUCT, INCLUDING WITHOUT LIMITATION ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. MILTENYI BIOTEC DOES NOT WARRANT THAT THE PRODUCT IS ERROR-FREE OR WILL ACCOMPLISH ANY PARTICULAR RESULT.

UNDER NO CIRCUMSTANCES WILL MILTENYI BIOTEC BE LIABLE FOR ANY LOSS OF USE, INTERRUPTION OF BUSINESS OR ANY INDIRECT, SPECIAL, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES OF ANY KIND (INCLUDING LOST PROFITS) REGARDLESS OF THE FORM OF ACTION WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE), STRICT PRODUCT LIABILITY OR OTHERWISE, EVEN IF MILTENYI BIOTEC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. IN NO EVENT SHALL THE TOTAL LIABILITY OF MILTENYI BIOTEC HEREUNDER EXCEED THE GREATER OF 100.00 EUROS OR THE AMOUNT YOU ACTUALLY PAID FOR THE PRODUCT GIVING RISE TO SUCH LIABILITY, REGARDLESS OF THE CAUSE OF ACTION, IN CONTRACT, TORT, STRICT LIABILITY OR OTHERWISE.

NOT ALL JURISDICTIONS ALLOW SUCH LIMITATIONS OF DAMAGES SO THE FOREGOING LIMITATIONS MAY NOT APPLY TO YOU. This warranty statement gives you specific legal rights and you may have other rights, which may vary from jurisdiction to jurisdiction.

For technical support, contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support at Miltenyi Biotec headquarters:

Miltenyi Biotec B.V. & Co. KG  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany  
Phone: +49 2204 8306 3803  
technicalsupport@miltenyi.com

Visit **[www.miltenyibiotec.com](http://www.miltenyibiotec.com)** for local Miltenyi Biotec Technical Support contact information.

# 13 EC Declaration of conformity

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Miltenyi Biotec B.V. & Co. KG  
Friedrich-Ebert-Str. 68  
51429 Bergisch Gladbach  
Germany

This declaration relates exclusively to the machinery in the state in which it was placed on the market, and excludes components which are added and/or operations carried out subsequently by the final user.

The declaration of conformity refers to the machinery identified as follows:

Description: Laboratory equipment  
Model: gentleMACS Dissociator

The machinery complies with all essential requirements of the following directives:

2006/42/EC Machinery  
2011/65/EU Restriction of the use of certain hazardous substances in electrical & electronic equipment  
2014/30/EU Electromagnetic compatibility

The machinery is in conformity with the following harmonized standards:

EN 61010-1:2010  
EN 61010-2-051:2015  
EN 61326-1:2013

Person authorized to compile the relevant technical documentation:

Dr. Bernd Schröder  
Global Head Regulatory Affairs  
Miltenyi Biotec B.V. & Co. KG  
Friedrich-Ebert-Str. 68  
51429 Bergisch Gladbach  
Germany

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Miltenyi Biotec B.V. & Co. KG  
Friedrich-Ebert-Str. 68  
51429 Bergisch Gladbach  
Germany

This declaration relates exclusively to the machinery in the state in which it was placed on the market, and excludes components which are added and/or operations carried out subsequently by the final user.

The declaration of conformity refers to the machinery identified as follows:

Description: Laboratory equipment  
Model: gentleMACS Dissociator

The machinery complies with all essential requirements of the following legislations:

Supply of Machinery (Safety) Regulations 2008  
The Restriction of the Use of Certain Hazardous Substances in  
Electrical and Electronic Equipment Regulations 2012  
Electromagnetic Compatibility Regulations 2016

The machinery is in conformity with the following UK designated standards:

EN 61010-1:2010  
EN 61010-2-051:2015  
EN 61326-1:2013

Person authorized to compile the relevant technical documentation:

Dr. Bernd Schröder  
Global Head Regulatory Affairs  
Miltenyi Biotec B.V. & Co. KG  
Friedrich-Ebert-Str. 68  
51429 Bergisch Gladbach  
Germany

## Miltenyi Biotec

### Germany/Austria

Miltenyi Biotec B.V. & Co. KG  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany  
Phone +49 2204 8306-0  
Fax +49 2204 85197  
macsde@miltenyi.com

### USA/Canada

Miltenyi Biotec, Inc.  
2303 Lindbergh Street  
Auburn, CA 95602, USA  
Phone 800 FOR MACS  
Phone +1 866 811 4466  
Fax +1 877 591 1060  
macsus@miltenyi.com

### Australia

Miltenyi Biotec  
Australia Pty. Ltd.  
Unit 11, 2 Eden Park Drive  
Macquarie Park, NSW 2113  
Australia  
Phone +61 2 8877 7400  
Fax +61 2 9889 5044  
macsau@miltenyi.com

### Benelux

Miltenyi Biotec B.V.  
Sandifortdreef 17  
2333 ZZ Leiden  
The Netherlands  
macsnl@miltenyi.com

### Customer service

#### The Netherlands

Phone 0800 4020120  
Fax 0800 4020100

#### Customer service Belgium

Phone 0800 94016  
Fax 0800 99626

#### Customer service Luxembourg

Phone 800 24971  
Fax 800 24984

### China

Miltenyi Biotec Technology &  
Trading (Shanghai) Co., Ltd.  
Room A401, 4/F  
No. 1077, Zhangheng Road  
Pudong New Area  
201203 Shanghai, P.R. China  
Phone +86 21 6235 1005-0  
Fax +86 21 6235 0953  
macscn@miltenyi.com.cn

### France

Miltenyi Biotec SAS  
10 rue Mercoeur  
75011 Paris, France  
Phone +33 1 56 98 16 16  
macsfr@miltenyi.com

### Hong Kong

Miltenyi Biotec Hong Kong Ltd.  
Unit 301, Lakeside 1  
No. 8 Science Park West Avenue  
Hong Kong Science Park  
Pak Shek Kok, New Territories  
Hong Kong  
Phone +852 3751 6698  
Fax +852 3619 5772  
macshk@miltenyi.com.hk

### Italy

Miltenyi Biotec S.r.l.  
Via Paolo Nanni Costa, 30  
40133 Bologna  
Italy  
Phone +39 051 6 460 411  
Fax +39 051 6 460 499  
macsit@miltenyi.com

### Japan

Miltenyi Biotec K.K.  
NEX-Eitai Building 5F  
16-10 Fuyuki, Koto-ku  
Tokyo 135-0041, Japan  
Phone +81 3 5646 8910  
Fax +81 3 5646 8911  
macsjp@miltenyi.com

### Nordics and Baltics

Miltenyi Biotec Norden AB  
Medicon Village  
Scheeleorget 1  
223 81 Lund  
Sweden  
macsse@miltenyi.com  
**Customer service Sweden**  
Phone 0200 111 800  
Fax +46 280 72 99

#### Customer service Denmark

Phone 80 20 30 10  
Fax +46 46 280 72 99

#### Customer service

**Norway, Finland, Iceland,  
and Baltic countries**  
Phone +46 46 280 72 80  
Fax +46 46 280 72 99

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[www.miltenyibiotec.com](http://www.miltenyibiotec.com)

Miltenyi Biotec provides products and services worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. gentleMACS, MACS, and the Miltenyi Biotec logo are either registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. Copyright © 2022 Miltenyi Biotec and/or its affiliates. All rights reserved.

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Miltenyi Biotec Swiss AG
Miltenyi Biotec Ltd.