

Contents

1. Description
 - 1.1 Principle
 - 1.2 Background information
 - 1.2 Applications
2. Protocol: Dissolving of lyophilizate
3. Example
4. References

1. Description

Components	20 µg StemMACS™ Cas9 Nuclease mRNA encoding Cas9 (CRISPR associated protein 9), a RNA-guided DNA endonuclease from <i>Streptococcus pyogenes</i> A20. The coding sequence has been modified with N- and C-terminal nuclear localization signals (NLS) for efficient transport of the Cas9 protein into the nucleus ¹ .
	1 mL Double-distilled Water , RNase-free
Specifications	<i>In vitro</i> transcribed, polyadenylated and capped mRNA that has been modified with pseudouridine and 5-methyl-cytidine to reduce the innate antiviral response to single-stranded RNA.
Formulation	Lyophilized from a filtered (0.2 µm) solution.
Storage	Store the lyophilized product at -20 °C. The expiration date is indicated on the label. After reconstitution, the product can be stored at -70 °C for up to 3 month.
Quality control	mRNA size has been verified on an Agilent Bioanalyzer System. Cas9 protein expression and function after transfection was confirmed by eGFP gene knock out in mouse ES cells.

1.1 Principle

The transient expression of key developmental regulators, recombinases or markers via mRNA transfection is a powerful tool for modulating cell fate. StemMACS mRNAs are highly pure, *in vitro*-transcribed mRNAs that have been carefully optimized and validated to ensure high level expression after transfection.

1.2 Background information

Cas9 originates from the CRISPR (clustered regularly interspersed palindromic repeats) adaptive immunity system of *Streptococcus pyogenes*.

Cas9 nuclease can be used to introduce site-directed double strand breaks into the genomic DNA. These will subsequently be repaired through one of two repair pathways: non-homologous end joining (NHEJ) or homology directed repair (HDR). NHEJ repair is error prone and will frequently result in nucleotide deletions or insertions (indel). If such insertions/deletions result in a frame shift or introduction of a premature stop codon, they will disrupt the open reading frame of the target gene. Therefore, this approach can be used to generate a gene knock out.

Alternatively, the HDR pathway can be used to introduce specific changes into the gene of interest. By providing a DNA repair template that contains the desired modification flanked by regions that are homologous to the genomic target sequence, HDR enables genomic integration of the edited target sequence.

The position of the Cas9 cleavage site is determined by two factors: the target sequence recognized by the guide RNA and the presence of a protospacer adjacent motif (PAM, in the case of *Streptococcus pyogenes* Cas9 the nucleotide sequence NGG). The Cas9/guide RNA complex will recognize the conjunction of both sequence elements and generate a double strand break 3 base pairs upstream of the PAM².

1.3 Applications

- Site-specific gene knockout
- Gene editing by homologous recombination

2. Protocol: Reconstitution of lyophilizate

▲ RNA is susceptible to degradation by exogenous ribonucleases. Wear gloves, use RNase-free reagents, tubes, and pipette tips.


1. Dissolve StemMACS Cas9 Nuclease mRNA in 200 µL of Double-distilled Water. Vortex thoroughly. The final concentration will be 0.1 µg/µL.
2. Briefly centrifuge to collect the content at the bottom of the tube.
3. Prepare aliquots and store at -70 °C to -80 °C. Do not subject aliquots to more than two freeze-thaw cycles.

For satisfactory transfection results, use a protocol that is optimized for your specific cell type. StemMACS eGFP mRNA or StemMACS Nuclear eGFP mRNA allow easy evaluation of transfection efficiency and are recommended as positive controls.

3. Example

HM1 mouse embryonic stem cells carrying a GFP reporter gene were transfected with 750 ng StemMACS Cas9 Nuclease mRNA and 250 µg plasmid encoding the guide RNA using Lipofectamine® 2000. Cells were analyzed by flow cytometry 96 hours after transfection. 17.04% of all cells were GFP-negative indicating successful eGFP knockout.

Transfection efficiency was >90% as determined by independent electroporation of mCherry mRNA into control HM1 cells (not shown).



4. References

1. Cong, L. *et al.* (2013) Multiplex genome engineering using CRISPR/Cas systems. *Science* 339 (6121): 819–823.
2. Jinek, M. *et al.* (2012) A programmable dual-RNA-guided DNA endonuclease in adaptive bacterial immunity. *Science* 337 (6096): 816–821.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

The Miltenyi Biotec logo and StemMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Lipofectamine is a trademark of Life Technologies Corporation.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.