

## Contents

1. Description
  - 1.1 Important safety information
  - 1.2 Background information
  - 1.3 Technical specifications
  - 1.4 Applications
  - 1.5 Reagent and instrument requirements
2. Use of MS Columns
  - 2.1 Preparation of MS Columns
  - 2.2 Magnetic separation using MS Columns

## 1. Description

**This product is for research use only.**  
**For use by professional laboratory personnel only.**

**Components**      **MS Columns (# 130-042-201):**  
25 MS Columns and plungers, sterile packed.  
or  
**MS Columns plus tubes (# 130-122-727):**  
25 MS Columns and plungers (# 130-042-201),  
sterile packed, and 75× 5 mL Tubes for  
MS Columns (# 130-091-598), sterile packed as  
3× 25 tubes.

**Storage**            Store columns dry at +10 to +35 °C and  
protected from light. The expiration date is  
indicated on the box label. Do not use after this  
date.

### 1.1 Important safety information

**⚠ WARNING**      Contamination or infection could result in  
death or serious injury depending on the  
material used.

All biological material must be considered  
potentially infectious

- Regulations for the treatment and  
disposal of infectious materials must be  
observed.

### 1.2 Background information

The patented MACS® Column Technology is based on the use of MACS MicroBeads, MACS Columns, and MACS Separators. MS Columns have been developed for the gentle isolation of MicroBead-labeled cells. As MACS MicroBeads are extremely small, superparamagnetic particles, a high-gradient magnetic field is required to retain the labeled cells. MS Columns contain an optimized matrix to generate this strong magnetic field when placed in a permanent magnet such as the MiniMACS™ Separator, OctoMACS™ Separator, or SuperMACS™ II Separator.

### 1.3 Technical specifications

	Max. number of labeled cells	Max. number of total cells
Manual use	1×10 <sup>7</sup>	2×10 <sup>8</sup>

▲ **Note:** Column capacity may decrease when separating cells larger than lymphocytes. Please refer to the respective MACS Cell Separation Reagent data sheet for column capacity of other cells than lymphocytes.

- Recommended sample size for leukocytes: 10<sup>4</sup>–10<sup>7</sup> labeled cells in 10<sup>6</sup>–2×10<sup>8</sup> total cells. Sample concentration: up to 10<sup>8</sup> leukocytes/500 µL cell suspension.
- Typical enrichment rate: 50-fold to up to 1,000-fold, depending on the strength and specificity of the magnetic labeling. Up to 10,000fold enrichment can be achieved by separation over two sequential columns.
- Columns are “flow stop” and do not run dry.
- Void volume: 60 µL. Reservoir volume: 3.5 mL.
- Typical flow rate for phosphate-buffered saline (PBS) containing 0.5% bovine serum albumin (BSA): 0.35–0.55 mL/min.
- MS Columns are for single use only.

### 1.4 Applications

MS Columns have been developed for positive selection of human and animal cells, especially rare cells, out of a heterogeneous cell suspension in combination with a MACS Separator. MS Columns can also be used for depletion of cells which strongly express the magnetically labeled surface antigen. They can also be used to separate other biological material such as plant cells, bacteria, viruses, protozoa, cell organelles, etc.

▲ Do not use MS Columns in combination with magnetic particles other than MACS MicroBeads. Magnetic forces in the column are very high and may damage biological material if other beads are used.

▲ MS Columns are not suitable for particles larger than 30 µm. To remove clumps and to prevent aggregates in the sample, resuspend material carefully and pass through 30 µm nylon mesh (Pre-Separation Filters (30 µm), # 130-041-407) before separation.

▲ Samples or buffers with high viscosity might cause reduced column flow or column clogging.

### 1.5 Reagent and instrument requirements

- **Buffer:** Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (+2 to +8 °C). Degas buffer before use, as air bubbles could block the column.

▲ **Note:** The recommended buffer is PBS supplemented with EDTA and BSA. The suitability of other buffers has to be tested experimentally.

▲ **Note:** Use degassed buffer only! Degas buffer by applying vacuum, preferentially with buffer at room temperature. Excessive gas in running buffer will form bubbles in the matrix during separation. This may lead to clogging of the column and decrease the quality of separation.


- MACS MicroBeads for magnetic labeling of cells.
- MiniMACS Separator (# 130-042-102), OctoMACS Separator (# 130-042-109), or SuperMACS II Separator (# 130-044-104)
- MACS MultiStand (# 130-042-303) in combination with MiniMACS Separator or OctoMACS Separator
- Adapter for MS, LS, and LD Columns for use with SuperMACS II Separator
- MACS Acrylic Tube Rack (# 130-041-406) or OctoMACS Acrylic Tube Rack (# 130-090-448)
- (Optional) Pre-Separation Filters (30 µm) (# 130-041-407) to remove cell clumps.

## 2. Use of MS Columns

### 2.1 Preparation of MS Columns

1. Insert MS Column with the column wings to the front into a MACS Separator according to A) or B).

#### A) Use with MiniMACS or OctoMACS Separator



Attach MiniMACS or OctoMACS Separator to the MultiStand and place MS Column in the separator. Place a collection tube under the MS Column.

#### B) Use with SuperMACS II Separator

For use of MS Columns with the SuperMACS II Separator please refer to the respective data sheet.

2. Prepare MS Column by rinsing with buffer: apply 500 µL of degassed buffer on top of the column and let the buffer run through. MS Columns are “flow stop” and do not run dry.
3. Discard effluent and change collection tube. The MS Column is now ready for magnetic separation.

▲ **Note:** Use column immediately after filling to avoid formation of air bubbles caused by warming up. Do not store columns after filling.

▲ **Note:** The time for filling the column with buffer is dependent on the storage conditions, temperature, and humidity. Therefore, the time may vary from a few seconds to several minutes. This filling time has no influence on the quality of the separation.



## 2.2 Magnetic separation using MS Columns

▲ For details on magnetic labeling, refer to the MACS Cell Separation Reagent data sheets.

▲ Always wait until the column reservoir is empty before proceeding to the next step.

1. Resuspend up to  $10^8$  total cells in 500 µL of degassed buffer.
  - ▲ **Note:** For higher cell numbers, scale up buffer volume accordingly.
  - ▲ **Note:** When working with fresh anticoagulated blood or buffy coat, dilute before separation 1:2 with buffer. Alternatively, use the StraightFrom® Whole Blood MicroBeads or StraightFrom Buffy Coat MicroBead Kits in combination with Whole Blood Columns.
  - ▲ **Note:** To remove clumps, pass cells through Pre-Separation Filters (30 µm).
2. Apply cell suspension onto the prepared MS Column. Collect flow-through containing unlabeled cells.
3. Wash MS Column with  $3 \times 500$  µL degassed buffer. Collect unlabeled cells that pass through and combine with the flow-through from step 2.
  - ▲ **Note:** Perform washing steps by adding buffer aliquots as soon as the column reservoir is empty.
4. Remove MS Column from the separator and place it on a new collection tube.
5. Pipette 1 mL buffer onto the MS Column. Immediately flush out fraction with the magnetically labeled cells by firmly applying the plunger supplied with the column.
6. (Optional) To increase the purity of the magnetically labeled fraction, the eluted fraction can be enriched over a second MS Column. Repeat the magnetic separation procedure as described in steps 2 to 5 by using a new column.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for local Miltenyi Biotec Technical Support contact information.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

autoMACS, MACS, the Miltenyi Biotec logo, MiniMACS, OctoMACS, StraightFrom, and SuperMACS are registered trademarks or trademarks Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2024 Miltenyi Biotec and/or its affiliates. All rights reserved.