Miltenyi Biotec

CliniMACS® CD8 Reagent

MD

CE 0123

REF 200-070-264

English | en

Contents

QTY

One vial contains 7.5 mL CliniMACS® CD8 Reagent in a sterile solution.

Composition

Each vial contains 7.5 mL of a colloidal solution of iron-dextran beads conjugated to murine monoclonal CD8 antibody in PBS/EDTA buffer stabilized with Poloxamer 188 (0.03% w/v).

Intended purpose

The CliniMACS CD8 Reagent is intended for *in vitro* magnetic labeling of CD8 positive human cells to enable the enrichment or depletion of human CD8 positive cells with a CliniMACS System only.

Limitation

Miltenyi Biotec as the manufacturer of the CliniMACS System does not give any recommendations regarding the use of separated cells for therapeutic purposes and does not make any claims regarding a clinical benefit.

Side-effects

When the CD8 enriched or depleted cells are infused or injected into patients, they may receive traces of murine antibody and iron-dextran.

Iron-dextran beads and/or murine antibodies may cause allergic or anaphylactic reactions in patients. Intensive care equipment and medication should be available.

Precautions

The instructions for use including those described in the corresponding CliniMACS User Manuals must be followed.

For the manufacturing and use of target cells in humans the national legislation and regulations – e.g. for the EU the Directive 2004/23/EC (human tissues and cells) or the Directive 2002/98/EC (human blood and blood components) – must be followed. Thus, any clinical application of the target cells is exclusively within the responsibility of the user of a CliniMACS System.

The separation of CD8 positive cells using the CliniMACS CD8 Reagent must be performed by trained operators only.

Operators using the CliniMACS CD8 Reagent should have experience in the separation of cells from bone marrow or peripheral blood.

All materials which have come into contact with blood or blood products must be treated as infectious material. Regulations for the treatment of infectious material must be observed.

Regulatory and legal notes

In Canada, any clinical application of the output product must be performed in accordance with applicable Canadian legislation and regulations that pertain to cellular therapies (e.g. for advanced cellular therapies, the applicable sections of the Food and Drugs Act and the Food and Drug Regulations). Refer to the CliniMACS Prodigy T Cell

Transduction Brochure for information on the performance and cell separation procedures of the CliniMACS CD8 Reagent as part of the CliniMACS Prodigy T Cell Transduction System.

Warnings

The CliniMACS Reagents and Biotin Conjugates are intended for *in vitro* use only and are not designated for therapeutic use or direct infusion into patients.

The CliniMACS CD8 Reagent is not recommended for use with patients known or suspected to have sensitivity against mouse immunoglobulins or iron-dextran.

Patients may develop human anti-mouse antibodies (HAMA).

See also "Side-effects".

Do not use after the use-by date printed on the vial label.

Do not use if package is damaged. Use reagent only if vial is undamaged and sealed.

Storage

The CliniMACS CD8 Reagent is shipped refrigerated and must be stored at +2 °C to +8 °C (+36 °F to +46 °F) immediately after receipt. Do not freeze the reagent.

Further information

Endotoxin content ≤ 2 EU/mL as determined by kinetic Limulus Amebocyte Lysate (LAL) assay (Pharmacopoeia Europaea (Ph. Eur.))

STERILE A

Sterile. Manufactured aseptically, sterile filtered, filled aseptically.

The depicted symbol indicates the single sterile barrier of the product.

Any serious incident that has occurred in relation to this product should be reported to Miltenyi Biotec B.V. & Co. KG – using the contact information provided – and the competent authority of the member state in which the user of this product is established.

Performance

The performance of the product depends on the chosen application. For information on the different applications, refer to the corresponding CliniMACS User Manual.

Handling instructions

The use of the product is described in the corresponding user manual of the chosen application.

Refer to the CliniMACS User Manual of the corresponding application for the validated combinations of instruments, separation programs, reagents, tubing sets, and buffer.

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