

## Index

1. Description
  - 1.1 Principle and product applications
  - 1.2 Recommended reagent dilution
  - 1.3 Reagent requirements
2. General protocol for MACS Comp Reagent
3. Examples of compensation with MACS Comp Reagent

## Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

## 1. Description


**This product is for research use only.**

<b>Product format</b>	1 mL MACS Comp Reagent, mouse: Anti-leukocyte CD45 monoclonal antibody (rat IgG2b) conjugated to anti-fluorochrome monoclonal antibody (mouse IgG1).  Antibody conjugates are supplied in a solution containing stabilizer and 0.05% sodium azide.
<b>Product size</b>	100 tests.
<b>Storage</b>	Store protected from light at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label.

### 1.1 Principle and product applications

The MACS Comp Reagent is a dual antibody system that allows for compensation of PE- or APC-conjugated antibodies and several tandem fluorochrome conjugates. For a list, please see the table in section 1.2.

In contrast to bead-based compensation systems, the cell-based MACS Comp Reagent system takes into account the autofluorescence of cells during compensation and avoids overcompensation. Equal quantities of stained and unstained cells are mixed and compensation is performed using a lymphocyte gate.



The MACS Comp Reagent employs a dual antibody system that selectively targets CD45-expressing cells (all leukocytes) as well as the fluorochrome moiety of the antibody to be compensated for. Mouse cells are incubated with the MACS Comp Reagent for labeling via CD45 after which the fluorochrome-conjugated antibody in question is added. Thus, all CD45-positive cells within the sample are fluorescently labeled with the target antibody—irrespective of the antibody's specificity. Unlabeled cells are spiked into the sample to serve as an internal negative control. Cells are analyzed by flow cytometry and compensation is performed using the fluorescence signal of the target antibody.

## Product applications

- Efficient compensation for flow cytometric analysis of rare cells.
- Efficient compensation of tandem-fluorochrome conjugates.
- Efficient compensation for flow cytometric analysis of weakly stained cells.

### 1.2 Recommended reagent dilution

For labeling of mouse cells.

	MACS Comp Reagent – PE <sup>a</sup>	MACS Comp Reagent – APC <sup>a</sup>
<b>Reagent dilution for flow cytometry</b>	1:11	1:11
<b>Compatible fluorochrome conjugates</b>	PE PE-Cy5 PE-Cy5.5 PE-Cy7 PE-Alexa Fluor® 700 PE-Alexa Fluor 750 PE-Alexa Fluor 647	APC APC-Cy7 APC-Alexa Fluor 750

a) Suitable for use with fresh or fixed cells.

▲ **Note:** Though the MACS Comp Reagent has been tested to function with a broad range of tandem fluorochrome conjugates, Miltenyi Biotec cannot guarantee that the product will work with the above listed tandem fluorochrome conjugates from all suppliers.

### 1.3 Reagent requirements

- **Buffer:** Prepare a solution containing phosphate-buffered saline (PBS) pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (4–8 °C). Degas buffer before use, as air bubbles could block the column.
- ▲ **Note:** EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). BSA can be replaced by other proteins such as mouse serum albumin, mouse serum, or fetal calf serum. Buffers or media containing Ca<sup>2+</sup> or Mg<sup>2+</sup> are not recommended for use.
- A minimum of 10<sup>6</sup> freshly isolated mouse cells, e.g. splenocytes.

## 2. General protocol for MACS Comp Reagent

▲ Work fast, keep the cells cold, use pre-cooled solutions which will prevent capping of antibodies on the cell surface and a non-specific cell labeling.


▲ Working on ice requires increased incubation times. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling.

1. Resuspend  $5 \times 10^5$  mouse splenocytes in 100  $\mu$ L of buffer.
2. Add 10  $\mu$ L of MACS Comp Reagent.
3. Mix well and refrigerate for 10 minutes in the dark (4–8 °C).
4. Wash cells by adding 1–2 mL of buffer and centrifuge at  $300 \times g$  for 10 minutes. Aspirate supernatant completely.
5. Resuspend cells in 100  $\mu$ L of buffer.
6. Add the antibody conjugate to be compensated for at the recommended titer.
7. Mix well and refrigerate for 10 minutes in the dark (4–8 °C).
8. Wash cells by adding 1–2 mL of buffer and centrifuge at  $300 \times g$  for 10 minutes. Aspirate supernatant completely.
9. Resuspend cells in 500  $\mu$ L of buffer.
10. **Immediately** before measurement, add a further  $5 \times 10^5$  unstained splenocytes in 500  $\mu$ L of buffer to serve as an internal negative control.
11. Proceed to the compensation procedure. For automatic compensation, follow the instructions provided by the manufacturer of your flow cytometer. Manual compensation can be performed as outlined in the following steps.
12. Set a lymphocyte gate. In general, only cells with identical autofluorescent characteristics should be viewed.
13. Make sure that the voltage adjustment in each channel is optimized for unstained cells. Create a dot plot displaying the fluorescence channel of the used fluorochrome. Set the opposing axis to the appropriate fluorescence channel to eliminate the fluorescence overlap (e.g. FITC fluorescence channel for PE-conjugated antibodies). Create separate regions for the analysis of positive and negative cells.
14. Open the statistics window to display the median fluorescence intensity of both populations. Adjust the compensation values in the channel to be corrected so that the median fluorescence intensities of the stained cell population equals that of the unstained cells.
15. Repeat steps 12–14 for compensation of additional antibody conjugates/fluorescence channels.


## 3. Examples of compensation with MACS Comp Reagent

Mouse splenocytes were labeled with MACS Comp Reagent followed by the addition of the fluorochrome-conjugated antibody to be compensated for. Unlabeled cells were added directly before measurement. Data depicting the optimal compensation of CD25-PE (# 130-091-013) (a) and CD19-APC (# 130-092-039) (b) are shown.

- (a) Optimized compensation of CD25-PE using the MACS Comp Reagent-PE and mouse splenocytes.



- (b) Optimized compensation of CD19-APC using the MACS Comp Reagent-APC and mouse splenocytes.



Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for local Miltenyi Biotec Technical Support contact information.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

autoMACS, MACS, and the Miltenyi logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only.

Cy is a registered trademark of GE Healthcare UK Limited.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.