

Contents

1. Description
 - 1.1 Background information
 - 1.2 Applications
2. Protocol
 - 2.1 Preparation of complete medium
 - 2.2 Reconstitution of StemMACS mRNA Reprogramming Cocktail, human
 - 2.3 Reprogramming of fibroblasts

1. Description

This product is for research use only.

Components	1× StemMACS mRNA Reprogramming Box, human containing: 2×25 µg StemMACS mRNA Reprogramming Cocktail, human 2×1 mL Double-distilled water 1×130 µL StemMACS mRNA Transfection Reagent 1×25 mL StemMACS mRNA Transfection Buffer 1× StemMACS Repro-Brew XF, human containing: 1×500 mL StemMACS Repro-Brew Basal Medium XF, human 1×10 mL StemMACS Repro-Brew Supplement XF (50×), human
Specifications	StemMACS Repro-Brew XF, human: pH: 7.2–7.6 Osmolality: 300–330 mOsmol/kg
Quality control	Functionality-tested for reprogramming of human fibroblasts into iPSCs.

Storage

Upon arrival, store lyophilized StemMACS mRNA Reprogramming Cocktail, human at –20 °C. After reconstitution, the cocktail can be stored at –70 °C to –80 °C for up to three months.

Store StemMACS Repro-Brew Basal Medium XF, human protected from light at 2–8 °C. Do not freeze.

Store StemMACS Repro-Brew Supplement XF (50×), human at –20 °C. Avoid repeated freeze-thaw cycles.

Store StemMACS mRNA Transfection Reagent and StemMACS mRNA Transfection Buffer at 2–8 °C. Do not freeze.

The expiration dates are indicated on the labels.

1.1 Background information

The StemMACS iPSC mRNA Reprogramming Kit, human has been designed for the rapid and efficient reprogramming of human fibroblasts into induced pluripotent stem cells (iPSCs). With the StemMACS iPSC mRNA Reprogramming Kit, human iPSCs can be generated within 14 days under feeder-free conditions without the need for a conditioned medium or the use of B18R protein. In contrast to traditional viral- or DNA-based reprogramming methods, mRNA reprogramming eliminates the risk for genomic integration and the safety concerns connected with using viral-based vectors. StemMACS mRNA Reprogramming Cocktail, human contains transcripts of the following genes: *POU5F1* (*OCT3/4*), *SOX2*, *KLF4*, *MYC*, *NANOG*, *LIN28A*, and *SOCS1*.

1.2 Applications

- mRNA-based reprogramming of human fibroblasts into iPSCs

2. Protocol

2.1 Preparation of complete medium

▲ Before StemMACS Repro-Brew XF, human can be used in cell culture, the two components need to be mixed according to the following protocol to obtain the complete medium.

1. Thaw StemMACS Repro-Brew Supplement XF (50×), human at room temperature prior to use. If a precipitate is visible at the bottom of the vial, resuspend thoroughly or swirl the vial in a 37 °C waterbath for 1–2 minutes until the supplement is clear.
2. Prepare appropriate aliquots and store them at –20 °C. Avoid repeated freeze-thaw cycles.
3. To obtain the complete medium, dilute StemMACS Repro-Brew Supplement XF (50×), human 1:50 with StemMACS Repro-Brew Basal Medium XF, human, e.g., add 2 mL StemMACS Repro-Brew Supplement XF (50×), human to 98 mL StemMACS Repro-Brew Basal Medium XF, human. The complete medium can be used for up to 14 days when stored at 2–8 °C.

2.2 Reconstitution of StemMACS mRNA Reprogramming Cocktail, human

▲ RNA is susceptible to degradation by exogenous ribonucleases. Wear gloves and use RNase-free reagents, tubes, and pipette tips.

1. Dissolve 25 µg StemMACS mRNA Reprogramming Cocktail, human in 250 µL of double-distilled water. Vortex thoroughly. The final concentration is 0.1 µg/µL.
2. Briefly centrifuge to collect the content at the bottom of the tube.
3. Prepare aliquots and store at –70 °C to –80 °C. Avoid repeated freeze-thaw cycles. Thawed aliquots can be stored at 2–8 °C for up to five days.

2.3 Reprogramming of fibroblasts

▲ The detailed protocol is available at www.miltenyibiotec.com/reprogramming-fibroblasts.

The reprogramming protocol includes plating of fibroblasts in complete medium and allowing them to adapt. mRNA transfections are performed twice a day for five consecutive days. After several days of medium changes iPSC colonies can be isolated for further expansion.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Limited use license

This product was developed by Miltenyi under the license of intellectual property licensed by iPS Academia Japan, Inc. Under this license Miltenyi is entitled to sell this product for non-commercial research use only (whether the buyer is an academic or a for-profit entity) under a non-transferable, limited-use license. Non-commercial research use only excludes the following: clinical use (i.e. administration of this product or any material made by using this product to humans or animals for therapeutic, diagnostic, or prophylactic purposes in any form), administration or implantation of any material made using this product into an animal for purposes other than basic pre-clinical research applications, any use for commercial purposes (i.e. any use of this product or any material made using this product for the manufacture of product, e.g. culture medium and equipment, for providing services or information, for screening pharmaceutical products or for sales or transfer), transfer of this product or any material made using this product to a for-profit entity. Purchasers wishing to use the product or any material made using this product for purposes other than non-commercial research use should contact Miltenyi's Licensing Office below.

Miltenyi's Licensing Office: licensing@miltenyibiotec.de

Trademarks

The Miltenyi Biotec logo and StemMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2023 Miltenyi Biotec and/or its affiliates. All rights reserved.