

## Contents

1. Description
  - 1.1 Background information
  - 1.2 Technical specifications
  - 1.3 Reagent and instrument requirements
2. Use of the MACSwell Imaging Frames
  - 2.1 Sample preparation
  - 2.2 Loading the MACSwell Imaging Frame into the MACSima™ Imaging System

## 1. Description


This product is for research use only.

**Components** Sterile and single-packed MACSwell Imaging Frames, 10 pieces per box.

Product	Order no.
MACSwell One	130-124-673
MACSwell One Small	130-126-794
MACSwell Two	130-124-675
MACSwell Four	130-124-676

**Storage** Store MACSwell Imaging Frames dry and protected from light at room temperature (19–25 °C). The expiration date is indicated on the box label. Do not use after this date.

### Parts of the MACSwell Imaging Frames



**Figure 1:** Parts of the MACSwell Imaging Frames.

## 1.1 Background information

The MACSwell Imaging Frames have been designed for the use with the MACSima Imaging System. The MACSima Imaging System is a fully automated instrument capable of staining hundreds of markers on one sample using the MICS (MACSima Imaging Cyclic Staining) technology. The frames allow for the use of microscope slides in the MACSima Imaging System. They form a liquid-tight contact with the microscope slide containing the sample, which enables analysis on the MACSima Imaging System.

## 1.2 Technical specifications

- For single use only.
- For professional laboratory use only.
- Frame material: cyclo olefin polymer (COP).
- Slider material: Polycarbonate.
- Frames are not compatible with acetone.
- Recommended working temperature between 4 °C and 40 °C.
- Size: 95 mm × 31 mm × 14.5 mm (length × width × height).

	MACSwell One	MACSwell One Small or MACSwell Two	MACSwell Four
Weight [g]	11.3	11.3	11.3
Number of wells	1	1 or 2	4
Well area [mm <sup>2</sup> ]	931	437	190
Maximum filling volume [mL]	12.5	5.9	2.7
Working volume [μL] (staining)	1000	500	250
Maximum sample size [mm]	47.25×17	21.25×17	8×17

**Table 1:** Data of MACSwell Imaging Frames.

▲ **Note:** MACSwell One Small has one centered well, MACSwell Two has two wells.


## 1.3 Reagent and instrument requirements

- MACSima Imaging System (# 130-121-164)
- MACSima Running Buffer (#130-121-565)
- Standard microscope slide (75–76 mm × 25–26 mm × 1 mm) or one variant of MACSwell Micro Slide, e.g., MACSwell Micro Slide (17 μm) (# 130-124-678)

## 2. Use of the MACSwell Imaging Frames

- Please refer to the MACSima Imaging System user manual for detailed information on using the instrument.
- For cell culture experiments prior to microscopy experiment, all steps should be performed under sterile conditions.
- Wet microscope slides can be used. Ensure that the sealing of the frame is free from any contamination before assembling the sample carrier.
- If the sealing is in contact with a written part of the slide, the related well is not usable.

### 2.1 Sample preparation



**Figure 2:** Use of the MACSwell Imaging Frames.

#### 2.1.1 When working with fixed tissue samples

- Every box is supplied with a black positioning template for the MACSwell Imaging Frame model ordered. Arrange the template underneath the microscope slide in use. The hole in the template indicates the optimal position of the sample on the microscope slide.

1. Prepare a microscope slide with a fixed tissue sample. For

positioning the sample(s) on the microscope slide refer to figure 3. Maximum sample sizes are listed in table 1 in chapter 1.2.

▲ **Note:** When working with the MACSwell Two, up to two samples can be processed in parallel, with the MACSwell Four up to four samples in parallel.

2. Remove frame including the lid from the sliding frame. Turn the frame upside down. The sealing is now on top.
 

▲ **Note:** The lid is only used for subsequent cell culture experiments.
3. Place the microscope slide with a fixed sample on the sealing of the frame with the sample orientated to the sealing (figure 2, A.1).
4. Ensure that the sample is not placed underneath the sealing (figure 2, A.2).
 

▲ **Note:** To adjust the position move the microscope slide carefully. Avoid fingerprints on the microscope slide.
5. Place the sliding frame on the frame. Start inserting the sliding hooks in the wider spot of the gaps, away from the MACSima Logo (figure 2, B.1).
 

▲ **Note:** Make sure the sliding hooks of the sliding frame are correctly inserted in the gaps of the frame.
6. Adjust the sliding frame in the direction of the MACSima Logo. The end position is defined by the stopping edge on the sliding frame (figure 2, B.2). Stop when the sliding frame can't be shifted anymore without applying too much pressure to avoid breakage.
7. Turn the MACSwell Imaging Frame back up (figure 2, C).
8. Add MACSima Running Buffer to each well for a total volume of:
  - 1900  $\mu$ L for MACSwell One,
  - 950  $\mu$ L for MACSwell One Small,
  - 950  $\mu$ L for MACSwell Two, or
  - 475  $\mu$ L for MACSwell Four Imaging Frames.
9. Proceed with loading the MACSwell Imaging Frames into the MACSima Imaging System (figure 2, D).

#### 2.1.2 When working with cell suspension samples

- MACSwell Micro Slides are compatible with all variants of the MACSwell Imaging Frames. If using MACSwell Micro Slides, refer to the respective data sheet.

1. Remove the frame including the lid from the sliding frame.
 

▲ **Note:** The lid is only used for subsequent cell culture experiments.
2. Turn the frame upside down. The sealing is now on top. Place the empty MACSwell Micro Slide onto the sliding frame. The microcavities, where the sample will later be applied (figure 2, A.1), should be orientated towards the sealing. Ensure that the microcavities will not be placed under the sealing (figure 2, A.2).
 

▲ **Note:** To adjust the position move the MACSwell Micro Slide carefully. Avoid fingerprints on the MACSwell Micro Slide.
3. Place the sliding frame on the frame. Start inserting the sliding hooks in the wider spot of the gaps, away from the MACSima Logo (figure 2, B.1).
 

▲ **Note:** Make sure the sliding hooks of the sliding frame are correctly inserted in the gaps of the frame.

- Adjust the sliding frame in the direction of the MACSima Logo. The end position is defined by the stopping edge on the sliding frame (figure 2, B.2).
- Turn the MACSwell Imaging Frame back up.
- Add the previously fixed cell suspension to the cavity of the imaging frame and wait the appropriate time for the cells to settle in the cavities.  
▲ **Note:** An initial equilibration step might be needed.
- After seeding the cells, add MACSima Running Buffer to each well for a total volume of:  
1900 µL for MACSwell One,  
950 µL for MACSwell One Small,  
950 µL for MACSwell Two, or  
475 µL for MACSwell Four Imaging Frames.
- Proceed with loading the MACSwell Imaging Frames into the MACSima Imaging System (figure 2, D).

### 2.1.3 When working with adherent cell samples

- First mount the frame on a microscope slide (figure 2, A–C).
- (Optional) Coat the slide.
- Seed cells and follow the culture protocol.
- Fix cells according to the protocol. Wash sample thoroughly to remove the fixating agent.  
▲ **Note:** Do not use acetone fixation!
- Add MACSima Running Buffer to each well for a total volume of:  
1900 µL for MACSwell One,  
950 µL for MACSwell One Small,  
950 µL for MACSwell Two, or  
475 µL for MACSwell Four Imaging Frames.
- Proceed with chapter 2.2 (figure 2, D).

## 2.2 Loading the MACSwell Imaging Frame into the MACSima Imaging System

- For more information please refer to the MACSima Imaging System manual.

- Start the MACSima Imaging System and software.
- Follow the instructions of the instrument. When the software asks to scan the barcode: Scan the 2D barcode on the MACSwell Imaging Frame with the barcode scanner of the MACSima Imaging System.  
▲ **Note:** Barcodes of MACSwell Imaging frames can be read even if still in primary packaging.
- When the MACSima Imaging System door is open, insert the MACSwell Imaging Frame into the MACSima Imaging System so that the MACSima Logo is facing upwards on the far side of the cavities.  
▲ **Note:** Make sure to always remove the lid before loading the MACSwell Imaging Frame into the MACSima Imaging System.
- Follow the instructions of the instrument.
- After the MICS experiment has finished, the MACSwell Imaging Frame can be removed from the instrument.

Follow the instructions of the instrument. Disassemble the MACSwell Imaging Frame by removing the sliding frame and subsequently the microscope slide. Discard the MACSwell Imaging Frame.




Figure 3: Sample area of the MACSwell Imaging Frames.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience

which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

## **Licenses**

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

## **Trademarks**

MACSima, MACSwell, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.