

## Contents

1. Description
  - 1.1 Principle of the MACS LENTiCheck Kit
  - 1.2 Background information
  - 1.3 Reagent and instrument requirements
2. Protocol
  - 2.1 Reagent preparation
  - 2.2 Protocol procedure
3. Data acquisition and analysis
4. Recommended acceptance criteria
5. Recommended criteria for a positive result
6. Troubleshooting

## 1. Description

This product is for research use only.

<b>Components</b>	<p><b>MACS LENTiCheck Kit Primers and Probe Mix</b> (64 µM per primer, 40 µM of the probe): The probe is labeled with FAM dye.</p> <p><b>MACS Check Kit Positive Control (pLTG2713):</b> Plasmid (gag'-VSVg-RRE-PTBP2) (10<sup>8</sup> copies/µL)</p>
<b>Capacity</b>	For 96 reactions (1–2 test samples).
<b>Product format</b>	<p>MACS LENTiCheck Primers and Probe Mix is supplied in a buffered solution of 10 mM TRIS-Cl and 1 mM EDTA, at pH 8.0, using nuclease-free water.</p> <p>MACS LENTiCheck Kit Positive Control (pLTG2713) is supplied in a buffered solution of 5 µM TRIS-Cl and 0.5 µM EDTA, at pH 8.0, using nuclease-free water.</p>
<b>Storage</b>	Store protected from light at –20 °C. The expiration date is indicated on the vial labels.

### 1.1 Principle of the MACS LENTiCheck Kit

The replication competent lentivirus (RCL) assay is intended to detect the lentiviral envelope gene sequence (vesicular stomatitis virus G glycoprotein (VSVg) gene) of potential RCLs that might arise, e.g., during unintended recombination of plasmids used in vector production or due to recombination of the viral vector with viral sequences in cells after transduction.

The MACS LENTiCheck Kit provides specific primers to amplify an amplicon of the lentiviral envelope gene sequence VSVg. The envelope gene (VSV-G) is thought to be the most suitable target for detecting RCLs as incorporation of an envelope gene sequence is required to generate a replication competent virus.

A VSVg-specific probe comprising a fluorescent dye (FAM) and quencher moiety is additionally used. During the PCR

amplification, the probe binds to the VSVg gene and gets hydrolyzed by the 5'-3'-exonuclease activity of the PCR polymerase. Thereby, the fluorescent dye is separated from the quencher leading to the emission of a fluorescence signal. The threshold cycle (Ct) value determined during the quantitative PCR is related to the amount of the respective gene in the sample.

For reliable quantification, the plasmid standard MACS Check Kit Positive Control (pLTG2713) is used to determine a standard curve for the amplification of the VSVg gene. Additional positive and negative controls are used to ensure proper detection of the VSVg gene. No signals should be gained for the viral vector or for the transduced cells under investigation.

Please note that the primers included in the MACS LENTiCheck Kit, human are only compatible with the lentiviral vectors provided by Lentigen and Miltenyi Biotec.

### 1.2 Background information

During the production of lentiviral vectors as well as in cells transduced with lentiviral vectors, there is the risk of an interplay between the viral vector and viral properties gained from other potentially present viruses or virus genes. Such an interplay might give rise to RCLs, also named replication competent retroviruses (RCR). The MACS LENTiCheck Kit has been developed to detect the gene encoding for the pseudotyping envelope protein VSV-G. The VSV-G gene is used to express the envelope protein during the production of lentiviral vectors to enable the packaging of the vector. However, VSV-G must not be packed into the viral vector. Therefore, the viral vector as well as the transduced cells are intended to be free of VSV-G gene. The VSV-G gene can be determined using the MACS LENTiCheck Kit. The absence of the VSV-G gene in *ex vivo* transduced cells according to the PCR assay can be used for lot release testing, in particular, when time constraints are present.

### 1.3 Reagent and instrument requirements

- Genomic DNA purified from 1×10<sup>6</sup> cells
- Nuclease-free water
- Taqman™ Fast Advanced Master Mix (Thermo Fisher Scientific Cat. No. 4444556) or similar
- Multicolor real-time PCR detection system (thermal cycler) for the detection of FAM (excitation maximum 495 nm, emission maximum 520 nm)
- Depending on the thermal cycler, a reference standard dye, such as ROX Reference Dye (Thermo Fisher Scientific Cat. No. 12223-012) or similar
- 1.5 mL microcentrifuge tubes
- 96-well PCR plates for fast PCR protocols as recommended for the multicolor real-time PCR detection system by the respective supplier
- Adhesive sealing films
- 2× PCR workstations (1 for DNA-free work and 1 for handling template DNA)

- 2× 1000 µL pipettes (DNA-free and for handling template DNA) with disposable filter tips
- 2× 100 µL pipettes (DNA-free and for handling template DNA) with disposable filter tips
- 2× 20 µL pipettes (DNA-free and for handling template DNA) with disposable filter tips
- Vortex mixer
- Microcentrifuge for 1.5 mL tubes
- Microcentrifuge tube racks

## 2. Protocol

### 2.1 Reagent preparation

▲ Mix all solutions briefly with a vortex mixer before use and perform a short spin before opening the tube.

▲ All components of the MACS LENTiCheck Kit (Primers and Probe Mix and Positive Control (pLTG2713)) are ready-to-use. Bring all reagents to room temperature (19–25 °C) before use.

### 2.2 Protocol procedure

▲ When mixing DNA solutions with 2× master mix: Always work with solutions of lower DNA concentration first, followed by solutions of successively higher DNA concentrations.

1. Prepare a blank solution to serve as no-template control (NTC) in the DNA-free PCR workstation: Add 50 µL nuclease-free water to a new 1.5 mL microcentrifuge tube.
2. Prepare standards for a standard curve by serial dilution of the MACS Check Kit Positive Control (pLTG2713) standard stock solution in the DNA-free PCR workstation according to table 1. “ST” is the abbreviation of “standard”. Mix well.

Vial name	Target standard concentration	Component 1		Component 2	
		Name	Volume (µL)	Name	Volume (µL)
ST 0	1×10 <sup>7</sup> copies/µL	Stock (1×10 <sup>8</sup> copies/µL)	10	Nuclease-free water	90
ST 1	1×10 <sup>6</sup> copies/µL	ST 0 (1×10 <sup>7</sup> copies/µL)	10	Nuclease-free water	90
ST 2	1×10 <sup>5</sup> copies/µL	ST 1 (1×10 <sup>6</sup> copies/µL)	10	Nuclease-free water	90
ST 3	1×10 <sup>4</sup> copies/µL	ST 2 (1×10 <sup>5</sup> copies/µL)	10	Nuclease-free water	90
ST 4	1×10 <sup>3</sup> copies/µL	ST 3 (1×10 <sup>4</sup> copies/µL)	10	Nuclease-free water	90
ST 5	1×10 <sup>2</sup> copies/µL	ST 4 (1×10 <sup>3</sup> copies/µL)	10	Nuclease-free water	90
ST 6	1×10 <sup>1</sup> copies/µL	ST 5 (1×10 <sup>2</sup> copies/µL)	10	Nuclease-free water	90
ST 7	1 copy/µL	ST 6 (1×10 <sup>1</sup> copies/µL)	20	Nuclease-free water	180
ST 8	0.5 copies/µL	ST 7 (1 copy/µL)	50	Nuclease-free water	50

**Table 1:** Dilution scheme for standard curve.

3. (Optional) Prepare additional quality control samples by serially diluting the standard 1 (ST 1) of the MACS Check Positive Control (pLTG2713) in the template DNA PCR workstation according to table 2. “QC” is the abbreviation of “quality control”. Mix well.

Vial name	Target standard concentration	Component 1		Component 2	
		Name	Volume (µL)	Name	Volume (µL)
QC 1	5×10 <sup>5</sup> copies/µL	ST 0 (1×10 <sup>7</sup> copies/µL)	10	Nuclease-free water	190
QC 2	5×10 <sup>4</sup> copies/µL	QC 1 (5×10 <sup>5</sup> copies/µL)	10	Nuclease-free water	90
QC 3	5×10 <sup>3</sup> copies/µL	QC 2 (5×10 <sup>4</sup> copies/µL)	10	Nuclease-free water	90
QC 4	5×10 <sup>2</sup> copies/µL	QC 3 (5×10 <sup>3</sup> copies/µL)	10	Nuclease-free water	90
QC 5	5×10 <sup>1</sup> copies/µL	QC 4 (5×10 <sup>2</sup> copies/µL)	10	Nuclease-free water	90
QC 6	5 copies/µL	QC 5 (5×10 <sup>1</sup> copies/µL)	10	Nuclease-free water	90

**Table 2:** Dilution scheme for quality control samples.

4. Discard QC 1, QC 2, and QC 5. QC 3, QC 4, and QC 6 will serve as quality controls.
5. Prepare genomic DNA from 1×10<sup>6</sup> cells. Elute the genomic DNA in 300 µL nuclease-free water. This starting material will be called “test sample” within the text.

▲ **Note:** Two replicates of the test samples should be analyzed with five replicate qPCR measurements for each replicate.

To be able to control the assay performance, it is recommended to prepare additional control samples:

- A negative control sample (non-transduced cells),
- a DNA extraction negative control sample, and
- a non-transduced cell sample with the MACS Check Positive Control (pLTG2713) spiked to the cells before genomic DNA extraction (DNA extraction positive control sample). Spiking 300 copies per 1×10<sup>6</sup> cells is equivalent to 50 copies per µg genomic DNA, assuming 6 µg of genomic DNA per 10<sup>6</sup> cells.

Sample description	Cell number/ pellet	Spike	Volume of PBS for pellet resuspension (µL)
Test sample (transduced T cells)	1×10 <sup>6</sup>	No.	200
Negative control sample (non-transduced cells)	1×10 <sup>6</sup>	No.	200
DNA extraction negative control sample	0 (PBS only)	No.	200
DNA extraction positive control sample	1×10 <sup>6</sup>	Yes. Add 30 µL of ST 6 (10 copies/ µL MACS Check Positive Control (pLTG2713))	170

**Table 3:** Overview of the test sample and the corresponding controls.

6. Prepare and aliquot a 2× master mix solution for the duplex assay in the DNA-free PCR workstation according to table 4, e.g., for one sample use 1188 µL Taqman Fast Advanced Master Mix and 12 µL MACS LENTiCheck Primers and Probe Mix. Mix well.

Number of test samples incl. controls	Taqman Fast Advanced Master Mix (μL)	MACS LENTiCheck Primers and Probe Mix (μL)
1	1188	12
2	1584	16

**Table 4:** Preparation of 2× master mix solution.

7. (Optional) In case a reference standard dye, e.g. ROX, is required for the multicolor real-time PCR detection system, the respective amounts of standard dyes should be added to the 2× master mix. Typically, the reference standard dye is provided as 100× or 50× concentrated stock solution (table 5).

Number of test samples incl. controls	Taqman Fast Advanced Master Mix (μL)		MACS LENTiCheck Primers and Probe Mix (μL)	Reference standard dye stock solution (μL)	
	Concentration of reference standard dye stock solution			Concentration of reference standard dye stock solution	
	100x	50x		100x	50x
1	1176	1164	12	12	24
2	1568	1552	16	16	32

**Table 5:** Preparation of 2× master mix solution including a reference standard dye with 50× or 100× concentrated stock solution.

8. In the DNA-free PCR workstations add 50 μL 2× master mix to the 1.5 mL microcentrifuge tubes containing the no-template control (NTC).
9. In the DNA-free PCR workstations add 50 μL 2× master mix to new 1.5 mL microcentrifuge tubes for the standard curve (8 tubes) and for the three quality control samples (3 tubes). A loading scheme is shown in table 6.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ST 8 (0.5 copies/ μL)			ST 7 (1 copy/μL)			ST 6 (1×10 <sup>1</sup> copies/μL)			ST 5 (1×10 <sup>2</sup> copies/μL)		
B	ST 4 (1×10 <sup>3</sup> copies/μL)			ST 3 (1×10 <sup>4</sup> copies/μL)			ST 2 (1×10 <sup>5</sup> copies/μL)			ST 1 (1×10 <sup>6</sup> copies/μL)		
C	DNA extraction negative control sample 1	Negative control sample 1 (non-transduced cells)		Replicate 1 of test sample 1 (transduced T cells)		Replicate 2 of test sample 1 (transduced T cells)		DNA extraction positive control sample 1		DNA extraction negative control sample 2		DNA extraction positive control sample 2
D												
E												
F												
G												
H	NTC			QC 6 (5 copies/ μL)			QC 4 (5×10 <sup>2</sup> copies/μL)			QC 3 (5×10 <sup>3</sup> copies/μL)		

**Table 6:** Example of a loading scheme of a 96-well PCR plate.

10. In the DNA-free PCR workstations add 75 μL 2× master mix to new 1.5 mL microcentrifuge tubes for each of the two replicates of each test sample, the respective negative control sample (non-transduced cells), the DNA extraction negative

control sample, and the DNA extraction positive control sample (5 tubes). A loading scheme is shown in table 6.

11. In the template DNA PCR workstation add 50 μL of each standard (table 1) to the 50 μL 2× master mix prepared in step 9. Mix well.
12. (Optional) In the template DNA PCR workstation, add 50 μL of the three quality control samples QC 3, QC4, and QC 6 (table 2) to the 50 μL 2× master mix prepared in step 9. Mix well.
13. In the template DNA PCR workstation add 75 μL of each of the two replicates of each of the test sample, the respective negative control sample (non-transduced cells), the DNA extraction negative control sample, and the DNA extraction positive control sample (table 3) to the respective 75 μL of 2× master mix prepared in step 10. Mix well.
14. Briefly before use, vortex all prepared solutions thoroughly.
15. In the template DNA PCR workstation, add 25 μL of each prepared solution into the 96-well PCR plate according to table 6. Each replicate of each test sample and the respective controls will be applied and analyzed in quintuple.
16. Seal the plate by using an adhesive sealing film and spin down before loading into the instrument. Run plate in a multicolor realtime PCR detection system under following conditions (table 7).


Temperature (°C)	Time (seconds)	Data collection	Number of cycles
50	120	No	1
95	20	No	1
95	5	No	40
60	30	Yes (plate read)	

**Table 7:** Thermal cycling conditions.

### 3. Data acquisition and analysis

Perform the primary data analysis according to the recommendations for the used multicolor real-time PCR detection system. Typically, analysis settings, e.g., for baseline and threshold values, need to be defined. The primary data analysis will result in threshold cycle or crossing point (Ct) values for each dye of each sample.

Use the data of the plasmid standard (table 1) to generate a standard curve for VSV-G (FAM). Plot the Ct values against the log(10) of the number of copies used for the standard.



**Figure 1:** Example of a standard curve depicting the measured Ct value versus the log(10) of the copies per μL of the MACS Check Kit Positive Control (pLTG2713) standard dilution (refer to table 2).

The mean Ct value for each triplicate sample is used to determine the log(10) of the number of copies according to the respective standard curve. The determined value to the power of 10 gives the copy number of VSV-G in the respective sample.

The linear equation of the standard is:

$$Ct = \text{slope} \times \log(\text{copies}/\mu\text{L}) + \text{Y-intercept}$$

Calculation example from figure 1:

$$Ct = (-3.333) \times \log(\text{copies}/\mu\text{L}) + 38.331$$

Accordingly, the copy number of VSV-G in the respective sample can be calculated by transforming the equation:

$$\text{Copies}/\mu\text{L} = 10^{(Ct - \text{Y-intercept}) \div \text{slope}}$$

Calculation example from figure 1:

$$\text{Copies}/\mu\text{L} = 10^{(Ct - 38.331) \div (-3.333)}$$

#### 4. Recommended acceptance criteria

Determine if the assay met system suitability acceptance criteria:

- The slope of the standard curve must be between -3.1 and -3.6.
- The R2 for the standard curve must be greater than or equal to 0.990.
- At least two wells containing the NTC must return an average Ct value of "N/A".
- The DNA extraction negative control sample and the negative control sample (non-transduced cells) must have no signal in five wells.
- The quality control sample QC 3 must be within 0.3 logs of the known quantity.
- The quality control samples QC 4 and 6 must be within 0.5 logs of their known quantities.

Valid Options:

- It is valid to omit up to three wells from the standard curve that are greater than 0.5 Ct from the rest of the replicates.

All acceptance criteria must be met for the assay to be considered valid. If the assay fails to meet the acceptance criteria, the assay is invalid and must be repeated.

#### 5. Recommended criteria for a positive result

▲ Results < 40 Ct are positive.

DNA extraction positive control sample:

- Sample will be deemed positive if at least three out of five wells are positive.
- If less than three wells shows amplification, it is recommended to repeat the assay.

Test sample:

- Sample will be deemed positive if at least three out of five wells are positive in duplicate samples.
- If less than three wells are positive for a replicate, that replicate will be deemed negative.
- If both replicates are negative, then the sample will be deemed negative.
- If one replicate is negative and one replicate is positive, it is recommended to repeat the assay.
  - If in the repeat assay both replicates are negative, then the sample will be reported as negative.
  - If in the repeat assay one or both replicates are positive, then the sample will be reported as positive.

#### 6. Troubleshooting

##### Signal in negative control sample

Signals in negative control sample typically indicate contamination of the used material with MACS Check Positive Control (pLTG2713) plasmid or PCR products from previous runs. It is recommended to set up the PCR reaction in a different area separated from PCR reactions or analysis. Dedicated equipment (pipettes, tips, etc.) should be used for preparing the PCR mixes. To reduce the risk of cross-contaminations, filter tips should be used. Avoid aerosols that might be generated, e.g., by vortexing, rigorous pipetting, or during centrifugation. The area and the equipment for preparing the PCR reactions should be regularly cleaned.

##### Deviation of replicate samples

Uniform mixtures are important to gain reproducible results. During set up of the PCR reactions, the components should be mixed well, e.g., by vortexing for 5 – 10 seconds. Pipetting errors can give rise to increased variation. It is recommended to regularly calibrate the pipettes for PCR set up. Uneven temperature or data collection on the multicolor real-time PCR detection system can also cause inconsistent results. The used instrument should be checked regularly.

##### No signal detected

MACS LENTiCheck Primers and Probe Mix might have not been added when preparing the 2× master mix solution or they were exposed to excessive light or temperature for a prolong period. Store kit components protected from light, minimize exposure to room temperature, and mix well before master mix preparation.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for local Miltenyi Biotec Technical Support contact information.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

MACS and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Taqman is a trademark of Roche Molecular Systems, Inc.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.