

Contents

1. Description
 - 1.1 Important safety information
 - 1.2 Background information
 - 1.3 Technical specifications
 - 1.4 Applications
 - 1.5 Reagent and instrument requirements
2. Use of Large Cell Columns
 - 2.1 Preparation of Large Cell Columns
 - 2.2 Magnetic separation using Large Cell Columns

1. Description

This product is for research use only.

For use by professional laboratory personnel only.

Components	25× Large Cell Column, 25 Large Cell Columns and plungers, sterile packed 25× flow resistor (23G)
Storage	Store columns dry at +10 to +35 °C and protected from light. The expiration date is indicated on the box label. Do not use after this date.

1.1 Important safety information



Contamination or infection could result in death or serious injury depending on the material used.

All biological material must be considered potentially infectious

- Regulations for the treatment and disposal of infectious materials must be observed.

1.2 Background information

The patented MACS® Column Technology is based on the use of MACS MicroBeads, MACS Columns, and MACS Separators. Large Cell Columns have been developed for the gentle isolation of MicroBead-labeled cells. As MACS MicroBeads are extremely small, superparamagnetic particles, a high-gradient magnetic field is required to retain the labeled cells. Large Cell Columns contain an optimized matrix to generate this strong magnetic field when placed in a permanent magnet such as the MiniMACS™ Separator or OctoMACS™ Separator.

1.3 Technical specifications

	Max. number of labeled cells	Max. number of total cells
Manual use	1×10 ⁷	2×10 ⁸

▲ **Note:** Column capacity may decrease when separating cells larger than lymphocytes. Please refer to the respective MACS Cell Separation Reagent data sheet for column capacity of other cells than lymphocytes.

- Recommended sample size for leukocytes: 10⁴–10⁷ magnetically labeled cells in 10⁶–2×10⁸ total cells. Sample concentration: up to 10⁸ leukocytes/500 µL cell suspension.
- Typical enrichment rate: 50-fold to up to 1,000-fold, depending on the strength and specificity of the magnetic labeling.
- Void volume: 80 µL. Reservoir volume: 3.5 mL.
- Typical flow rate without flow resistor for phosphate-buffered saline (PBS) containing 0.5 % bovine serum albumin (BSA): 1.2–2.0 mL/minute.
- Large Cell Columns are for single use only.
- Never remove plastic cover from flow resistor.

1.4 Applications

Large Cell Columns have been developed for positive selection of human and animal cells, especially large cells, out of heterogeneous cell suspensions in combination with the MiniMACS Separator and OctoMACS Separator. They can also be used to separate other biological material such as plant cells and protozoa.

▲ Positive selection of large cells, e.g. megakaryocytes, labeled with MACS MicroBeads from up to 2×10⁸ total cells.

1.5 Reagent and instrument requirements

- Buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold +2 to +8 °C). Degas buffer before use, as air bubbles could block the column.

▲ **Note:** The recommended buffer is PBS supplemented with EDTA and BSA. The suitability of other buffers has to be tested experimentally.


▲ **Note:** Use degassed buffer only! Degas buffer by applying vacuum, preferentially with buffer at room temperature. Excessive gas in running buffer will form bubbles in the matrix during separation. This may lead to clogging of the column and decrease the quality of separation.

- MACS MicroBeads for magnetic labeling of cells
- MiniMACS Separator (# 130-042-102) or OctoMACS Separator (# 130-042-109) in combination with the MACS MultiStand (# 130-042-303)
- (Optional) Pre-Separation Filters (30 µm) (# 130-041-407) to remove cell clumps.

2. Use of Large Cell Columns

2.1 Preparation of Large Cell Columns

1. Attach the MACS Separator to the MACS MultiStand and place Large Cell Separation Column with the column wings to the front in the separator. Place a collection tube under the column.



2. Attach the flow resistor (23G needle) to the Large Cell Separation Column.
3. Apply 500 μ L of degassed buffer on top of the column and let the buffer run through. Then discard effluent and change collection tube.



2.2 Magnetic separation using Large Cell Columns

▲ For details on magnetic labeling, refer to the MACS Cell Separation Reagent data sheets.

▲ Always wait until the column reservoir is empty before proceeding to the next step.

1. Resuspend up to 10^8 total cells in 500 μ L of degassed buffer.
 - ▲ **Note:** Use column immediately after filling to avoid formation of air bubbles caused by warming up of the buffer in the column.
 - ▲ **Note:** For higher cell numbers, scale up buffer volume accordingly.
 - ▲ **Note:** To remove clumps, pass cells through Pre-Separation Filters (30 μ m).
2. Apply cell suspension onto the prepared Large Cell Column. Collect flow-through containing unlabeled cells.
3. Wash Large Cell Column with 3×500 μ L degassed buffer. Collect unlabeled cells that pass through and combine with the flow-through from step 2.
 - ▲ **Note:** Perform washing steps by adding buffer aliquots as soon as the column reservoir is empty.
4. Remove Large Cell Column from the separator. Remove flow resistor and place it on a new collection tube.
5. Pipette 1 mL buffer onto the Large Cell Column. Immediately flush out fraction with the magnetically labeled cells by firmly applying the plunger supplied with the column.
6. (Optional) To increase the purity of the magnetically labeled fraction, the eluted fraction can be enriched over a second Large Cell Column. Repeat the magnetic separation procedure as described in steps 2 to 5 by using a new column.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, the Miltenyi Biotec logo, and MiniMACS are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2024 Miltenyi Biotec and/or its affiliates. All rights reserved.