

Contents

1. Description
 - 1.1 Background
 - 1.2 Applications
2. How to dilute

1. Description

Components	6×75 mL MACS® BSA Stock Solution
Capacity	For preparation of 6×1.5 L of MACS Separation Buffer.
Product format	MACS BSA Stock Solution contains phosphate-buffered saline (PBS) and 10% bovine serum albumin (BSA). This product is preservative-free.
Storage	Store protected from light at 2–8 °C. The expiration date is indicated on the bottle label. After opening the MACS BSA Stock Solution should be stored refrigerated and used within 3 days.

1.1 Background

The MACS BSA Stock Solution is used for the preparation of a preservative-free MACS Separation Buffer. When diluted in autoMACS™ Rinsing Solution (# 130-091-222) an optimal buffer for all MACS Separations is obtained.

1.2 Applications

The MACS Separation Buffer is recommended for any labeling and washing step during the MACS Separation procedure:

- Magnetic labeling of cells with MACS MicroBeads or labeling with other antibody conjugates.
- Cell washing steps, for example, during the separation procedure.

▲ **Note:** This MACS Separation Buffer can be used as an alternative, preservative-free Running Buffer solution for MACS Separations on the autoMACS Pro Separator or the autoMACS Separator. For more information on autoMACS Buffers please refer to the autoMACS Running Buffer (# 130-091-221) and the autoMACS Rinsing Solution (# 130-091-222) data sheets.

2. How to dilute

For the preparation of a preservative-free MACS Separation Buffer empty one bottle of MACS BSA Stock Solution into one bottle of autoMACS Rinsing Solution (# 130-091-222, 6×1.45 L). Mix well by inverting the autoMACS Rinsing Solution bottle several times. The MACS Separation Buffer contains a final BSA concentration of 0.5%.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Unless otherwise specifically indicated, all Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.