

Contents

1. Description
 - 1.1 Principle of the B Cell Expansion Kit
 - 1.2 Applications
 - 1.3 Reagent and instrument requirements
2. Protocol
 - 2.1 Reconstitution of lyophilized reagents
 - 2.2 B cell isolation
 - 2.3 B cell expansion
3. Example of B cell expansion using the B Cell Expansion Kit
4. Appendix: Culture plate sizes for B cell expansion

1. Description

This product is for research use only.

Components **1000 µg Human CD40-Ligand Multimer Kit:**
Recombinant human CD40-Ligand, premium grade (2×500 µg) and Cross-Linking Antibody (2×1 mL).

5 µg Human IL-4, premium grade
500 mL StemMACS™ HSC Expansion Medium XF, human

or

200 µg Human CD40-Ligand Multimer Kit:
Recombinant human CD40-Ligand, premium grade (2×100 µg) and Cross-Linking Antibody (2×0.2 mL).

5 µg Human IL-4, premium grade
100 mL StemMACS™ HSC Expansion Medium XF, human

Biological activity

CD40-Ligand:
The ED₅₀ of multimerized Human CD40-Ligand is ≤250 ng/mL corresponding to a specific activity of ≥4×10³ U/mg. The biological activity of the specific cytokine lot is indicated on the Certificate of Analysis (<https://www.miltenyibiotec.com/certificates>).

▲ **Note:** The ED₅₀ of multimerized Human CD40-Ligand is determined by proliferation assay using enriched CD19⁺ B cells in the presence of the cross-linking antibody and 50 IU/mL interleukin 4 according to Spriggs, M.K. *et al.*¹ Activity was also shown by activation assay using enriched CD19⁺ B cells in the presence of 50 IU/mL interleukin 4.

IL-4:

The ED₅₀ is ≤0.2 ng/mL corresponding to an activity of ≥5×10⁶ IU/mg. The biological activity of the specific cytokine lot is indicated on the Certificate of Analysis (<https://www.miltenyibiotec.com/certificates>).

▲ **Note:** The ED₅₀ is determined by proliferation assay using TF-1 cells according to Kitamura *et al.*¹ The proliferation assay was calibrated with the international standard for human IL-4 (NIBSC code 88/656) provided by the WHO/National Institute for Biological Standards and Control.

Product format All components are supplied in azide-free buffer containing stabilizer. Low endotoxin.

Storage The Cross-Linking Antibody should be stored at 2–8 °C. Lyophilized reagents and medium should be stored at –20 °C. The expiration date is indicated on the vial label. For information about reconstitution of the lyophilized reagents and storage after reconstitution refer to chapter 2.1.

1.1 Principle of the B Cell Expansion Kit

The B Cell Expansion Kits have been developed for the activation and expansion of human B cells. The CD40-Ligand Multimer mimics T cell-dependent activation of B cells. The B cell expansion is achieved by culturing and restimulation at day 7 and 10 of culture.

1.2 Applications

- Activation and expansion of human B cells.

1.3 Reagent and instrument requirements

- CD19 MicroBeads, human (# 130-050-301)
- AB serum
- Humidified incubator
- (Optional) Fluorochrome-conjugated antibodies for flow cytometric analysis, for example, CD80-PE, CD86-FITC, and CD19-APC. For more information about fluorochrome-conjugated antibodies refer to www.miltenyibiotec.com/antibodies.
- (Optional) Propidium Iodide Solution (# 130-093-233) or 7-AAD Staining Solution (# 130-111-568) for the flow cytometric exclusion of dead cells.

2. Protocol

▲ All steps in the protocol have to be performed under sterile conditions.

2.1 Reconstitution of lyophilized reagents

It is recommended to reconstitute lyophilized Human CD40-Ligand with deionized sterile filtered water to a final concentration of 0.5 mg/mL in a volume of 200 μ L (for 100 μ g size) or 1000 μ L (for 500 μ g size) and store at -20°C .

To increase the biological activity, we highly recommend multimerization of Human CD40-Ligand and Cross-Linking Antibody. Therefore, equal volumes of both components (e.g. 40 μ L of the reconstituted CD40-Ligand and 40 μ L of the Cross-Linking Antibody) should be incubated together for 30 minutes at room temperature ($19-25^{\circ}\text{C}$) prior to media application. Reconstitute Human IL-4 with deionized sterile-filtered water to a final concentration of 2.5×10^5 IU/mL, e.g. if the lot has a biological activity of 8×10^6 IU/mg, 160 μ L water should be used for reconstitution of 5 μ g. Mix by pipetting up and down until resuspended. To avoid repeated freeze-thaw cycles, reconstituted reagents should be aliquoted and stored at -20°C or below until use.

2.2 B cell isolation

Isolate B cells using CD19 MicroBeads, human (# 130-050-301). For details refer to the respective data sheet.

2.3 B cell expansion

This protocol is optimized for the activation and expansion of B cells, using the Human CD40-Ligand Multimer Kit. B cells can either be expanded in low density (0.15×10^6 cells/mL) or in high density (1×10^6 cells/mL) cell culture (refer to table in section 4). The low density protocol will increase B cell expansion by up to 10-fold.

▲ Volumes given below are for 10 mL of final expansion medium. When working with larger volumes please scale up accordingly.


▲ The final expansion medium should be prepared directly before the stimulation. Long-term storage is not recommended.

1. Prepare 10 mL StemMACS™ HSC Expansion Media XF with 2 μ L reconstituted Human IL-4, and 5% AB serum.
2. Pre-incubation of the components of the CD40-Ligand Multimer Kit: Mix 40 μ L reconstituted Human CD40-Ligand with 40 μ L Cross-Linking Antibody. Incubate for 30 minutes at room temperature. Subsequently, add the multimerized Human CD40-Ligand (80 μ L) to 10 mL prepared StemMACS HSC Expansion Media XF (final expansion medium).
3. Resuspend B cells at a density of 10^6 cells (high density) or 0.15×10^6 cells (low density) in 1 mL of the final expansion medium and culture in suitable cell culture plate. For details refer to chapter 4. Appendix: Culture plate sizes for B cell expansion.
4. Incubate at 37°C , 5% CO_2 .
▲ **Note:** In high density culture exchange half of the medium on day 5 with fresh medium supplemented with 5% AB serum and 50 units IL-4/mL (0.2 μ L IL-4 per mL media).
5. On day 7 and 10 harvest B cells and repeat step 1–4.
6. On day 14 harvest B cells and proceed to downstream applications.

3. Examples of B cell expansion using the B Cell Expansion Kit

For the expansion of B cells using the B Cell Expansion Kits 80 μ L multimerized Human CD40-Ligand are used per 10 mL prepared StemMACS HSC Expansion Media XF. Please note, that optimal concentration and ratio will depend on the application as well as the used assay.

The components of the CD40-Ligand Multimer Kit were pre-incubated for 30 minutes at room temperature (as described in 2.4.2). B cells were expanded with a starting density of 0.15×10^6 cells/mL (low density) according to the protocol, as described in section 2.4. On day 7 and 10 cells were harvested, counted, and restimulated. All cells were harvested, counted, and analyzed on day 14. With the low density protocol a B cell expansion of up to 10-fold was achieved.



4. Appendix: Culture plate sizes for B cell expansion

For B cell expansion the cells should be resuspended in culture medium at 1×10^6 cells/mL or 0.15×10^6 cells/mL.

The following table lists culture plate sizes suitable for different cell numbers. It also indicates the appropriate amount of medium to be added.

Total cell number (high density)	Total cell number (low density)	Final medium volume	Culture plate
0.5×10^6	0.075×10^6	0.5 mL	48 well
1.0×10^6	0.15×10^6	1 mL	24 well
2.0×10^6	0.3×10^6	2 mL	12 well
4.0×10^6	0.6×10^6	4 mL	6 well

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

MACS, the Miltenyi Biotec logo, StemMACS, and QuadroMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Ficoll-Paque is a trademark of GE Healthcare companies

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.