

Contents

1. Description

- 1.1 Background information
- 1.2 Applications
- 1.3 Reagent requirements

2. Protocol

- 2.1 Freezing of primary cells
- 2.2 Thawing of primary cells
- 2.3 Freezing of tissues
- 2.4 Thawing of tissues

1. Description

This product is for research use only.

Components 50 mL MACS Freezing Solution

Specifications pH: 7.4–7.8
Contains 10% DMSO as cryoprotectant.

Storage Store protected from light at 2–8 °C. The expiration date is indicated on the vial label.

1.1 Background information

The MACS Freezing Solution is a ready-to-use animal component-free media formulation designed for the xeno- and serum-free cryopreservation of primary cells and solid tissues, such as human and mouse tumor samples, peripheral blood mononuclear cells (PBMCs), and cells from dissociated solid tissues (e.g. cells from dissociated tumors). Cells and tissues cryopreserved in MACS Freezing Solution show high preservation of viability and good recovery after thawing.

1.2 Applications

- Cryopreservation of primary cells for downstream applications (e.g. cell analysis, cell culture)
- Cryopreservation of solid tissues for downstream applications (e.g. dissociation into single-cell suspensions)

1.3 Reagent requirements

- Cell culture medium
- Cryogenic vials
- Freezing container
- (Only for tissue pieces) MACS SmartStrainers (70 µm) (#130-098-462)
- (Optional) MACS Tissue Dissociation Kits (e.g. Tumor Dissociation Kit, human (# 130-095-929)) in combination with the gentleMACS[™] Octo Dissociator with Heaters (#130-096-427) for the preparation of single-cell suspensions from solid tissues

2. Protocol

▲ For the handling of the freezing container refer to the instructions of the provider.

2.1 Freezing of primary cells

1. Determine cell number.
2. Centrifuge cell suspension at 300×g for 10 minutes. Aspirate supernatant completely.
3. Resuspend cell pellet in MACS Freezing Solution to reach a concentration of 10⁷ cells per mL.
4. Quickly transfer the cell suspension into cryogenic vials.
5. Place the vials in a freezing container and immediately store at –80 °C.
6. After 24 hours transfer the cryogenic vials with the cells into a liquid nitrogen tank for long-term storage.

2.2 Thawing of primary cells

▲ Work quickly to avoid loss of cells.
▲ Pre-heat water bath at 37 °C before taking cells out of liquid nitrogen.

1. Take a vial with cells out of the liquid nitrogen tank.
2. Incubate the vial in a pre-heated water bath at 37 °C until only a little lump of ice is left.
3. Quickly transfer cell suspension into a 15 mL conical tube and dropwise add 5 mL of cell culture medium.
4. Wash cryogenic vial with 1 mL of cell culture medium and add to conical tube.
5. Centrifuge cells for 5 minutes at 300×g. Aspirate supernatant carefully.
6. Resuspend the cell pellet in appropriate buffer according to downstream application.

2.3 Freezing of tissues

1. Determine tissue weight.
2. Cut tissue into small pieces as recommended for downstream application (e.g. tissue dissociation).
▲ **Note:** If using the MACS Tissue Dissociation Kits in combination with the gentleMACS Octo Dissociator with Heaters for downstream tissue dissociation, cut the tissue as indicated in the respective data sheet.
3. Transfer up to 200 mg of tissue pieces per 1 mL of capacity into the cryogenic vials.
4. Add 1 mL of MACS Freezing Solution to up to 200 mg of tissue pieces. For larger tissue amounts scale up accordingly.
5. Close the cryogenic vials and resuspend tissue pieces by inverting the vial several times until pieces float in MACS Freezing Solution.

6. Place the vials in a freezing container and immediately store at -80°C .
7. After 24 hours transfer the cryogenic vials with the tissue pieces into a liquid nitrogen tank for long-term storage.

2.4 Thawing of tissues

▲ Work quickly to avoid loss of cells.

▲ Pre-heat water bath at 37°C before taking tissues out of liquid nitrogen.

1. Place a MACS SmartStrainer ($70\ \mu\text{m}$) on a 50 mL conical tube.
2. Rinse the MACS SmartStrainer ($70\ \mu\text{m}$) with 1 mL of cell culture medium.
3. Take a vial with tissue pieces out of the liquid nitrogen tank.
4. Incubate the vial in a pre-heated water bath at 37°C until only a little lump of ice is left.
5. Pour tissue pieces onto the MACS SmartStrainer ($70\ \mu\text{m}$).
6. Wash cryogenic vial with 1 mL of cell culture medium and pour it onto the MACS SmartStrainer ($70\ \mu\text{m}$).
7. Wash tissue pieces with additional 15 mL of cell culture medium.
8. Pieces are ready to be used for downstream application, e.g., tissue dissociation using the MACS Tissue Dissociation Kits in combination with the gentleMACS Octo Dissociator with Heaters.
9. (Optional) To recover cells remaining in the medium, centrifuge the medium for 5 minutes at $300\times g$.
10. (Optional) Aspirate supernatant and resuspend cell pellet in appropriate buffer according to downstream application.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, the Miltenyi Biotec logo, and gentleMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.