

# Complete workflows allow comprehensive tumor microenvironment analysis of cell subsets from limited tumor patient samples

Anne Richter, Ramona Siemer, Elvira Criado-Moronati, Anna Baranska, Philipp Gert, David Agorku, Olaf Hardt, Daniela Vorholt, Aparajita Singh, Ruth Kläver, Christian Dose, Bianca Heemskerk, and César Evaristo. Miltenyi Biotec GmbH, Bergisch Gladbach, Germany

## Introduction

therapeutic interventions on tumor-infiltrating leukocytes (TILs) and other cell populations within the tumor microenvironment.

Immunotherapy approaches that engage T cells to attack tumors However, the phenotypic and functional analysis of TILs is techhave proven clinical efficacy and tremendous potential in multiple nically challenging and labor intensive. The number of TILs can be cancers. However, there are several mechanisms of resistance to very low and small subpopulations might escape analysis as they current therapeutic strategies. Therefore, clinical benefit is expe- are lost in the background noise. Importantly, flow cytometryrienced by only a subset of patients, and often only partial responses based phenotyping of different cell populations requires dividing are achieved. Thus, additional research is necessary to improve the limited tumor material into multiple samples for analysis using immunotherapeutic strategies. In particular, it is important to ana-population-specific antibody panels, which dramatically reduces lyze steady-state anti-tumor immunity and monitor the effects of the number of cells available for analysis. Therefore, it is fundamental to use tools and workflows that can maximize the amount and quality of the data obtained from limited tumor patient samples.

## Results

#### Characterization of tumor microenvironment by immunohistochemistry, using the MACSima™ Imaging Platform

One tissue section (approx. 1 cm wide) was snap-frozen and in incomplete in its incomplete in its incomplete in its incomplete in the incomplete in its inco using the MACSima™ Imaging Platform. The underlying MICS Imaging Platform in a completely automated manner. (multiparameter imaging cell screen) technology enables staining of hundreds of markers on a single sample. It uses the principle of an iterative staining process comprising three main steps,

later fixed with acetone and analyzed by immunohistochemistry fluorescence signal, all of which are conducted by the MACSima

Figure 2 shows the expression of CD326 (EpCAM, in red) overlaid with the indicated marker (in white) within the same tumor section.

body fragments instead of antibodies to label specific cell surface REAlease CD4/CD8 (TIL) MicroBead Kit, human. Frequencies of markers. The antibody fragments are engineered to have a low T cells after gating on i) live cells, ii) FSC/SSC, iii) single cells was affinity for markers when present as monomers. However, when increased by more than 10-fold in the isolated T cell fraction the fragments are multimerized as a complex they bind markers compared to the unseparated cells (fig. 4B). (C) Within the isolated with high avidity. REAlease Technology can control the multimer/ population, critical tumor-specific subpopulations such as CD279 monomer state of the fragments and thus allows for a controlled release, where monomerized antibody fragments dissociate from the cell surface. Therefore, the technology enables users to obtain cells that are free from antibody fragments and magnetic labe. (fig. 4A).

REAlease® Technology relies on recombinantly engineered anti- T cells were isolated from the dissociated tumor sample, using the (PD1)hiCD366 (TIM-3)+CD39+CD8+ T cells and CD279 (PD1)hiCD137 (4-1BB)+CD39+CD4+ T cells could be identified easily (fig. 4C). Using enriched target cells, the number of events that needed to be acquired for proper flow cytometry analysis was 10-fold lower than with unseparated bulk tumor cells. This greatly reduced the time required for analysis.

Target cell with antigen

**REAlease Complex** 


Anti-Biotin MicroBead

♥ ◆ Non-target cells

## T cell isolation increases sensitivity of single-cell TCR sequencing

Cryopreserved dissociated tumor sample was thawed and filtered using a MACS SmartStrainer (70 µm). Cell viability in the thawed sample was below 30%. By eliminating dead cells using the (gray bars), and the corresponding number of cells was plotted on Dead Cell Removal Kit, the viability rate was increased to 83% the y-axis. The number of cells with the same TCRβ CDR3 clono-(fig. 5A). T cells were magnetically isolated using the REAlease CD4/ type in the isolated T cell population is displayed in black bars CD8 (TIL) MicroBead Kit, human to a purity of 80%. Single-cell TCR (fig. 5B). Conversely, the top 50 TCRβ CDR3 clonotypes identified sequencing was then performed to assess clonality of tumorinfiltrating T cells, either using unseparated bulk cells or isolated Γ cells. Target cell number for sequencing was 10,000/sample. bulk sample is displayed in gray bars (fig. 5C).

The top 50 TCRB CDR3 clonotypes identified in the bulk sample were ranked by order of abundancy from left to right on the x-axis in the isolated T cell population were ranked by order of abundancy (black bars). The number of cells with the same clonotype in the



## Materials and methods

### Isolation and characterization of TIL subsets from limited tumor samples

tissue was automated and optimized for epitope preservation performed on unseparated bulk cell samples. using the gentleMACS™ Octo Dissociator with Heaters and the Tumor Dissociation Kit, human. One aliquot of the dissociated tu- For single-cell sequencing on the 10x Genomics Platform, 17,000 were phenotyped by flow cytometry to assess expression of in table 1. activation and exhaustion markers. The remaining cell suspension was cryopreserved in vials containing  $1\times10^7$  cells.

At a later time point,  $1\times10^7$  cells were thawed, and dead cells were removed using the Dead Cell Removal Kit. T cells were isolated using the REAlease CD4/CD8 TIL MicroBead Kit and used for single-cell TCR sequencing using the 10x Genomics Platform. To


Resected ovarian carcinoma was stored and shipped overnight in T cell-depleted fraction for sequential B cell isolation using CD19 MACS® Tissue Storage Solution that was shown to maintain cell MicroBeads, human and StraightFrom Whole Blood and Bone viability and phenotype up to 48 h after collection. One small Marrow CD138 MicroBeads, human. Isolated B cells were subjected tissue section was analyzed by immunohistochemistry, using the to single-cell BCR sequencing using the 10x Genomics Platform. MACSima™ Imaging Platform. Dissociation of the remaining tumor For comparison, single-cell TCR and BCR sequencing was also

mor sample was stained with the 8-Color Immunophenotyping cells per sample were loaded in order to obtain 10,000 cells within Kit, human to determine the immune cell composition. T cells were droplets. Only T and B cells were sequenced, and so the number magnetically isolated from 1×10<sup>7</sup> total dissociated cells using the of analyzed cells differed greatly between the unseparated bulk REAlease® CD4/CD8 (TIL) MicroBead Kit, human. Isolated T cells cell samples and the samples of isolated cells. An example is shown

Sample	Number of analyzed cells
Bulk (T cells)	616
Isolated T cells	2,049
Bulk (B cells)	268
Isolated B cells	1,588
Table 1	

Single-cell sequencing


maximize the number of cells available for analysis, we used the An overview of the workflow is shown in figure 1.




#### TIL immunophenotyping by flow cytometry reveals low frequency of T cells within the tumor microenvironment


Dissociated tumor was stained with the 8-Color Immunophenotyping Kit, human and analyzed by flow cytometry using the (fig. 3A) and their frequency among total living cells (fig. 3B) are MACSQuant® Analyzer 16 to determine the immune cell composi-

tion. The gating strategy used to identify the indicated cell types






Effective isolation of tumor-infiltrating T cells using the REAlease® CD4/CD8 (TIL) MicroBead Kit, human enables improved analysis



#### B cell isolation increases sensitivity of single-cell BCR sequencing

To maximize the number of cells available for analysis, we used the T cell-depleted fraction (obtained during T cell isolation, fig. abundancy from left to right on the x-axis (gray bars), and the 5A) for sequential B cell isolation. B cells were isolated using CD19 MicroBeads, human and StraightFrom Whole Blood and Bone number of cells with the same BCR identified in the isolated Marrow CD138 MicroBeads, human, to obtain all possible B cells B cell population is displayed in black bars (fig. 6A). Conversely, and plasma cells. Single-cell BCR sequencing was then performed to assess clonality of tumor-infiltrating B cells, either using unseparated bulk cells or isolated B cells. Target cell number for cells with the same BCR in the bulk sample is displayed in gray sequencing was 10,000/sample. The top 25 heavy-chain BCR bars (fig. 6B).

clonotypes identified in the bulk sample were ranked by order of corresponding number of cells was plotted on the y-axis. The the top 25 BCR clonotypes identified in the isolated B cell population were ranked by order of abundancy (black bars). The number of



### Conclusion and outlook

- We developed novel tools that allow the comprehensive analysis of the tumor microenvironement.
- The MACSima Imaging Platform allows identification and phenotyping of multiple cell types while maintaining their spatial distribution in the tumor sample.
- T cell isolation improved phenotypic analysis by flow cytometry and, due to REAlease Technology, can be used to isolate further subpopulations of interest.
- Sequential isolation of the different cell types maximized the number of cells of each lineage available for analysis.
- T and B cell isolation greatly increased the sensitivity of single-cell immunoprofiling.
- These workflows greatly reduce experimental time and allow the performance of more complex experimental setups. We believe the use of these innovative tools and workflows can significantly increase the quality of the data obtained in immuno-oncology and immunotherapy research.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. gentleMACS, MACS, the MACS logo, MACSima, MACSQuant, REAfinity, REAlease, StraightFrom, Vio, Viobility, VioBlue, VioBright, and VioGreen are registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide. Copyright © 2019 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.

