

Contents

1. Description
 - 1.1 Important safety information
 - 1.2 Background information
 - 1.3 Technical specifications
 - 1.4 Applications
 - 1.5 Reagent and instrument requirements
2. Use of D Columns
 - 2.1 Preparation of D Columns
 - 2.2 Magnetic separation using D Columns

1. Description

This product is for research use only.
For use by professional laboratory personnel only.

Components	5× D Columns, sterile packed 5× flow resistor (20G) 5× flow resistor (21G) 5× flow resistor (22G) 5× 3-way stopcock
Storage	Store columns dry at +10 to +35 °C and protected from light. The expiration date is indicated on the box label. Do not use after this date.

1.1 Important safety information

- ⚠ WARNING** Contamination or infection could result in death or serious injury depending on the material used.
- All biological material must be considered potentially infectious
- Regulations for the treatment and disposal of infectious materials must be observed.

1.2 Background information

The patented MACS® Column Technology is based on the use of MACS MicroBeads, MACS Columns, and MACS Separators. The matrix of D Columns is composed of ferromagnetic fibers covered with a cell-friendly coating. When placed in the magnetic field of a MACS® Separator, the fibers amplify the magnetic field, thus inducing a high gradient within the column. This is crucial for isolation of cells which are only minimally labeled with MACS MicroBeads, leaving enough epitopes free for concurrent antibody staining. The space between the fibers is several times larger than primary and most cultured cells. This allows the cells to freely flow through the column. Magnetically labeled cells are held in suspension within the column and do not actually “bind” the column matrix. This suspension minimizes stress on the cells and allows for efficient sterile washing by avoiding cell aggregation.

1.3 Technical specifications

	Max. number of labeled cells	Max. number of total cells
Manual use	1×10 ⁹	1×10 ¹¹

- Recommended sample size for leukocytes: 10⁷–10⁹ magnetically labeled cells in 10⁹–10¹¹ total cells. Sample concentration: up to 10⁸ leukocytes/500 µL cell suspension.
- Void volume: 43 mL. Reservoir volume: 20 mL.
- Typical flow rates: refer to table on next page.
- D Columns are for single use only.
- Never remove plastic cover from flow resistor.

1.4 Applications

D Columns have been developed for depletion of large amounts of magnetically labeled human and animal cells out of a heterogeneous suspension using SuperMACS II Separator and MACS MicroBeads. They can be used to separate different biological material including plant cells, bacteria, viruses, protozoa, or cell organelles.

- ▲ Do not use D Columns in combination with magnetic particles other than MACS MicroBeads. Magnetic forces in the column are very high and may damage biological material if other beads are used.
- ▲ To remove clumps and to prevent aggregates in the sample, resuspend material carefully and pass through 30 µm nylon mesh (Pre-Separation Filters (30 µm), # 130-041-407) before separation.
- ▲ Samples or buffers with high viscosity might cause reduced column flow or column clogging.
- ▲ Do not store column after filling.
- ▲ If the buffer does not flow well, there may be an air bubble in the capillary of the flow resistor. Switch the stopcock to “rinse”, flush the needle with buffer from the side syringe, switch the stopcock back to “run” and continue the run. The needle may also become blocked by cell clumps, which means the needle has to be replaced.

1.5 Reagent and instrument requirements


- Buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (+2 to +8 °C). Degas buffer before use, as air bubbles could block the column.
 - ▲ **Note:** The recommended buffer is PBS supplemented with EDTA and BSA. The suitability of other buffers has to be tested experimentally.
 - ▲ **Note:** Use degassed buffer only! Degas buffer by applying vacuum, preferentially with buffer at room temperature. Excessive gas in running buffer will form bubbles in the matrix during separation. This may lead to clogging of the column and decrease the quality of separation.
- MACS MicroBeads for magnetic labeling of cells
- SuperMACS II Separator (# 130-044-104)

- Adapter for CS and D Columns for use with SuperMACS II Separator
- 70% ethanol in double-distilled water
- 50 mL syringe
- MACS Acrylic Tube Rack (# 130-041-406)
- (Optional) Pre-Separation Filters (30 μ m) (# 130-041-407) to remove cell clumps.


2. Use of D Columns

2.1 Preparation of D Columns

Positions of the 3-way stopcock



Assembly of the D Column



1. Remove yellow cap from separation column and attach 3-way stopcock to column at port A.
2. Fill a syringe with 70% ethanol and attach to port B of the 3-way stopcock.
3. Turn 3-way stopcock to position "fill".
4. Insert column into the mounted Column Adapter for CS and D Columns and move Column Adapter into the magnetic field of the SuperMACS II Separator by turning the handle (for details refer to respective data sheet).
5. Fill the column upright from the bottom with 70% ethanol from the syringe until the solution reaches the reservoir.
6. Turn the 3-way stopcock to position "run" and rinse column by filling from the top with buffer. Allow 70% ethanol to run into the column until it reaches the white filter device. Then, add fresh buffer. Rinse with 500 mL of buffer. Do not let the column "run dry" at any time during procedure.
7. Choose a flow resistor 22G, 21G, or 20G (refer to table below) and cut off the tip of the plastic sheath with the pliers supplied with the SuperMACS II Separator. Leave the plastic sheath in place for safety and attach the flow resistor to port C of the 3-way stopcock.

▲ **Note:** Do not remove plastic sheath from flow resistor.

▲ **Note:** The lower the flow rate, the higher is the depletion efficiency of magnetically labeled cells.

Flow resistor	Buffer flow rates in mL/min
22G	3.5
21G	4.0
20G	8.0

8. Fill the syringe with buffer and attach to port B of the 3-way stopcock. Leave the syringe attached during separation, except when refilling.
9. Turn the stopcock to position "rinse" and flush the air from the flow resistor. Turn the stopcock in position "closed". The column is now ready for separation.

2.2 Magnetic separation using D Columns

1. Apply magnetically labeled cell suspension in appropriate volume of buffer (up to 10^8 cells per 500 μ L) onto the column with the flow resistor attached. Turn the stopcock to position "run" and allow cell suspension to penetrate the matrix.
2. Wash the column with 200 mL of buffer. Collect effluent as negative fraction that is depleted of the magnetically labeled cells.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, the Miltenyi Biotec logo, and SuperMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2024 Miltenyi Biotec and/or its affiliates. All rights reserved.