

## Contents

1. Description
  - 1.1 Background information
  - 1.2 Applications
  - 1.3 Reagent and instrument requirements
2. Protocol

## Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

## 1. Description

This product is for research use only.

<b>Components</b>	<p><b>60 µL Oct3/4 Antibody, anti-human/mouse, APC, REAfinity</b> (clone: REA622, isotype: recombinant human IgG1).</p> <p><b>60 µL FoxA2 Antibody, anti-human, APC, REAfinity</b> (clone: REA506, isotype: recombinant human IgG1).</p> <p><b>60 µL OTX2 Antibody, anti-human, Vio® B515, REAfinity</b> (clone: REA1178, isotype: recombinant human IgG1).</p> <p><b>60 µL PAX-6 Antibody, anti-human, PE, REAfinity</b> (clone: REA507, isotype: recombinant human IgG1).</p> <p><b>60 µL TTF-1 Antibody, anti-human, Vio B515, REAfinity</b> (clone: REA1090, isotype: recombinant human IgG1).</p> <p><b>60 µL Sox1 Antibody, anti-human, PE, REAfinity</b> (clone: REA698, isotype: recombinant human IgG1).</p> <p><b>25 mL Fixation/Permeabilization Solution 1</b></p> <p><b>2×40 mL Fixation/Permeabilization Solution 2</b></p> <p><b>40 mL Permeabilization Buffer (10×)</b></p>
<b>Capacity</b>	30 tests and up to 1×10 <sup>6</sup> total cells/test.
<b>Product format</b>	Antibodies are supplied in buffer containing stabilizer and 0.05% sodium azide.
<b>Storage</b>	Store antibodies and Fixation/Permeabilization Solution 1 and 2 protected from light at 2–8 °C. Do not freeze. Store Permeabilization Buffer (10×) at room temperature. The expiration date is indicated on the vial label.

### 1.1 Background information

Pluripotent stem cell (PSC) differentiation is a core aspect of PSC research. In the last decade, protocols for the generation of PSC-derived midbrain dopaminergic (mDA) neurons have been established and optimized. Nevertheless, differences in the manufactured cell lots remain and cannot be avoided. Therefore, a quality control step that allows reliable *in vitro* evaluation of the differentiation efficiency is pivotal.

The PSC-mDA Neuron Analysis Kit, anti-human, REAfinity has been developed for *in vitro* phenotyping the identity and purity of mDA neurons. The kit allows the detection of specific regional markers that are expressed already early during differentiation. It can assess the identity of the cells, the cell number of the different sub-populations, and detect non-differentiated cells that might contaminate the culture. Cellular controls specific to other brain regions, such as ventral and dorsal forebrain cells, ventral hindbrain cells, as well as for remnant pluripotent stem cells are necessary for reliable gate definition and detection of contaminating cell populations.

The kit applies recombinant engineered REAfinity Antibodies. REAfinity Antibodies are recombinant antibodies that provide superior lot-to-lot consistency and purity compared to mouse or rat hybridoma-derived, monoclonal antibodies. They have been recombinantly engineered to produce highly specific antibodies that require no FcR blocking step. Additionally, they all have the same IgG1 isotype, requiring less isotype controls.

In summary, this flow cytometric assay results in a quick and reliable qualitative and quantitative analysis during the differentiation process.

### 1.2 Applications

- Flow cytometry-based quality control for *in vitro* phenotyping of mDA neurons derived from human PSCs.

### 1.3 Reagent and instrument requirements

- Buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS® BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (2–8 °C).
  - ▲ **Note:** EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). BSA can be replaced by other proteins such as human serum albumin, human serum, or fetal bovine serum (FBS). Buffers or media containing Ca<sup>2+</sup> or Mg<sup>2+</sup> are not recommended for use.
- Deionized or distilled water.
- Flow cytometer equipped with a red (640 nm) and a blue (488 nm) laser, e.g., MACSQuant Analyzer 10 (# 130-096-343) or MACSQuant Analyzer 16 (# 130-109-803).
- (Optional) Propidium Iodide Solution (# 130-093-233) for detection of dead cells during flow cytometric analysis.

- (Optional) MACS Comp Bead Kit, anti-REA (# 130-104-693) for optimal compensation of the fluorescence spillover from fluorochrome-conjugated antibodies.

## 2. Protocol

▲ For a detailed step-by-step protocol, as well as a comprehensive list of reagent and instrument requirements refer to the application protocol “Flow cytometry-based quality control assay for PSC-derived midbrain dopaminergic neurons” on [www.miltenyibiotec.com/PSC-mDA-neurons](http://www.miltenyibiotec.com/PSC-mDA-neurons). The application protocol provides information on cell culture, generation of cellular controls, and detailed analysis.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for local Miltenyi Biotec Technical Support contact information.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery (“Product Warranty”). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order (“Product Warranty”). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

MACS, the Miltenyi Biotec logo, REAfinity, and Vio are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2022 Miltenyi Biotec and/or its affiliates. All rights reserved.