




Introduction

bind to its target, and can be used to enumerate CAR T cells during manufacturing and immunomonitoring. Those detection reagents are used in our flow cytometric assays which were developed for 1) in-process control, QC release testing, and concomitant research during cell manufacturing, and 2) for determination of CAR T cell persistence and cell phenotyping during patient immunomonitoring.

For all flow assays mentioned above so-called Express Modes have been programmed, that allow an automated acquisition and analysis of stained samples on MACSQuant[®] Analyzers. These Express Modes feature predefined experiment settings, automate the sample measurement, and apply a fully automated gating strategy with computational gate adjustment for each individual data file. This allows for a high standardization by reducing operator variability, full reproducibility of data analysis, and integration into automated workflows.




1 Flow cytometry assays for analysis of CART cells during cell manufacturing and immunomonitoring

Table 1											
		Cellular composition	Transduction efficiency	Hematopoietic cells	Malignant B cells	T stem cell memory	Differentiation	Proliferative ability	Exhaustion	Activation	Staining control
Channel	Fluorochrome										
V1	ViobBlue®	CD45	CD45	–	CD45	CD197	CD45RA	CD27	CD223	CD154	–
V2	Viogreen™	CD4	CD4	–	Ig kappa	CD4	CD4	CD4	CD4	CD4	CD4
B1	FITC	CD3	CD3	CD45	CD19	CD3	CD3	CD3	CD3	CD3	CD3
B2	PE	CD16/CD56	CAR-DR	CD133/2	CD34	CD86	CAR-DR	CAR-DR	CAR-DR	CAR-DR	CAR-DR
B3		7-AAD	7-AAD	7-AAD	7-AAD	7-AAD	7-AAD	7-AAD	7-AAD	7-AAD	7-AAD
B4	PE-Vio® 770	CD19	–	–	CD20	CD62L	CD62L	CD279	CD279	CD25	–
R1	APC	CD14	CD14	CD34	CD5/CD10	CD45RO	CD45RO	CD127	CD366	CD137	–
R2	APC-Vio® 770	CD8	CD8	–	Ig lambda	CD8	CD8	CD8	CD8	CD8	CD8
Express Mode		Immune_Cell_Compotion_h_02	CAR_T_Cell_Transduction_h_02		B_Malignant_Cells_h_01	CAR_T_Cell_Tscm_h_01	CAR_T_Cell_Differentiation_h_01	CAR_T_Cell_Proliferative_Ability_h_01	CAR_T_Cell_Exhaustion_h_01	CAR_T_Cell_Activation_h_01	CAR_T_Cell_Staining_Control_h_01


exhaustion status during cell manufacturing and immunomonitoring. All stainings have been established and tested on leukapheresis and cultured CAR-T cells for cell manufacturing and on whole blood for immunomonitoring.

1 Enumeration of immune cell subsets in leukapheresis products




the starting material (fig. 2). Using the same panel on the final cellular product allowed assessment of cellular impurities and quantification of T cells.

2 Assessment of transduction efficiency and viability of CART cells




Cell_Transduction_h_02" enabled the determination of CART cell transduction efficiency and viability (fig. 3).

3 Analysis of CART cell persistence and differentiation during immunomonitoring



"Persistence_h_01" was used as a model for the analysis of CART cell persistence and differentiation after adoptive transfer during immunomonitoring (fig. 4).

4 Evaluation of donor-dependent variations and automated Express Mode gating



fully automated gating procedure of the Express Mode adapts to each data file via pre-defined algorithms, which yields a high level of standardization across samples, lacks operator variability, and ensures full reproducibility.

patient immunomonitoring. This will help with establishing complex individualized therapies and will enable us to understand in greater detail the phenotypic changes occurring throughout the life time of a CART cell.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. The ClinMACS® System components, including Reagents, Tubing Sets, Instruments, and PBS/EDTA Buffer, are designed, manufactured and tested under a quality system certified to ISO 13485.

In the EU, the ClinMACS System components are available as CE-marked medical devices for their respective intended use, unless otherwise stated. The ClinMACS Reagents and Biotinylated Conjugates are intended for *in vitro* use only and are not designated for therapeutic use or direct infusion into patients. The ClinMACS Reagents in combination with the ClinMACS System components are not intended for use in humans. The ClinMACS System components are not intended for use in humans for any of the separated cells for therapeutic purposes and does not make any claims regarding a clinical benefit. For the manufacturing and use of target cells in humans the national legislation and regulations – e.g. for the EU the Directive 2004/23/EC (human tissues and cells), or the Directive 2002/98/EC (human blood and blood components) – must be followed. Thus, any clinical application of the target cells is exclusively within the responsibility of the user of a ClinMACS System.

In the US, the ClinMACS CD34 Reagent System, including the ClinMACS Plus Instrument, ClinMACS CD34 Reagent, ClinMACS Tubing Sets TS and LS, and the ClinMACS PBS/EDTA Buffer, are FDA approved for all products of the ClinMACS Product Line are available for use only under an approved Investigational New Drug (IND) application or Investigational Device Exemption (IDE) for use in the US. ClinMACS MicroSort and the ClinMACS MicroSort Plus are not for human therapeutic or diagnostic use.

In the US, the ClinMACS Prodigy™ T-Cell Transduction Process is available for research use only.

In the US, ClinMACS Prodigy, MACS, the MACS logo, MACSQuant, ViO, ViObio, and ViOgreen are registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide. Copyright © 2019 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.