
Miltenyi Biotec

MACSQuant[®] Instrument

User manual

MACS
QuantVYB

MACS
Quant16

MACS10
Quant

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No part of this user manual may be reproduced, stored in a retrieval system, transmitted, published, or distributed in any form or by any means, electronically, mechanically, by photocopying, microfilming, recording, or otherwise, without the prior written consent of Miltenyi Biotec; however, notwithstanding the foregoing, the owners of the MACSQuant Instrument may make copies solely for purposes of training personnel in the use and servicing of the unit within their business or organization. Maximal care has been taken by Miltenyi Biotec in the preparation of this user manual. However, Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein, or for incidental or consequential damages in connection with the furnishing, performance, or use of this document. The information in this document is provided "as is" without warranty of any kind, and is subject to change without notice. Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. autoMACS, MACS, MACSima (MACSQuant/MACSima Storage Solution), MACSQuant, MACSQuantify, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide.

MACSQuant[®] Instrument

User manual

Original instructions

Issued: 2025-08

140-005-358-006.06

Miltenyi Biotec B.V. & Co. KG

Friedrich-Ebert-Straße 68

51429 Bergisch Gladbach

Germany

Phone +49 2204 8306-0

Fax +49 2204 85197

macsde@miltenyi.com

www.miltenyibiotec.com

Content

Important safety information	7
Wichtige Sicherheitsinformationen	17
Información importante de seguridad	29
Informations de sécurité importantes	41
Importanti istruzioni di sicurezza	53
1 Introduction	65
1.1 Intended use	65
1.2 Instrument description	65
1.2.1 Key features of the MACSQuant Instrument	65
2 Installation	69
2.1 Components included in the delivery	69
2.2 Consumables	70
2.3 Accessories	71
2.4 Service contracts	71
2.5 Assembly	71
2.5.1 Installing the fluid bottles	71
2.5.2 Installing hydrophobic air filters	73
2.5.3 Installing the Single tube rack	73
2.5.4 Connecting the instrument to the power supply	74
2.5.5 Installing the MACSQuant Column	74
2.5.6 Installing the needle guard	76
2.5.7 Installing the MACS MiniSampler Plus	76
2.5.8 Installing the Chill Rack and Universal Reagent Rack	79
2.5.9 Installing a webcam	81
2.5.10 Installation checklist	81
3 Operating the instrument	83
3.1 Switching on	83
3.2 Priming the instrument	83

3.3 Switching off	83
3.3.1 Manual shutdown	84
3.3.2 Automatic shutdown	84
4 Calibration	85
4.1 Calibrating the sample uptake unit	85
4.1.1 Calibrating the needle arm	85
4.2 Calibrating racks and plates	87
4.2.1 Calibrating the Single tube rack	87
4.2.2 Calibrating the Chill 5 Rack	87
4.2.3 Calibrating the Chill 96 Rack	88
4.2.4 Calibrating the MACSPlex Filter Plate	88
4.2.5 Calibrating the Universal Reagent Rack	89
4.2.6 Calibrating the custom 96 racks	89
4.3 Calibrating the rack detection	90
4.3.1 Calibrating the Single tube rack detection	90
4.3.2 Calibrating the rack detection	90
4.4 Initial calibration	91
5 Instrument monitoring	93
5.1 Hardware monitor	93
5.1.1 Fluidics	93
5.1.2 Sample uptake unit	94
5.1.3 Lasers and detectors	94
5.2 Instrument status LEDs	95
6 MACSQuant Live Support	97
7 Maintenance	99
7.1 Surface disinfection	99
7.2 Disinfection of the fluidics system	99
7.3 Exchanging buffer and solution bottles	99
7.4 Exchanging the waste bottle	100

7.5 Exchanging hydrophobic air filters	100
7.6 Exchanging the MACSQuant Column	100
7.7 Cleaning the instrument	101
7.7.1 Cleaning programs	101
7.7.2 Cleaning recommendations	102
7.7.3 Removing a potential blockage	102
7.7.4 Cleaning the washing station	104
7.7.5 Cleaning the uptake needle	104
7.8 Long-term storage of the instrument	104
8 Troubleshooting	105
8.1 Leakage	108
8.1.1 Replacing the check valve	108
8.1.2 Syringe leakage	109
8.1.3 Column leakage	109
8.2 Air in system	109
8.2.1 Loose connections	109
8.2.2 Releasing air from the sheath particle filter	109
9 Technical data and specifications	111
9.1 Technical data of the MACSQuant Instrument	111
9.2 Excess volume of the MACSQuant Analyzer 10, VYB, and Analyzer 16	115
9.3 FCC compliance statement	115
9.4 EC/EU Declaration of Conformity	116
9.5 UK Declaration of Conformity	117
9.6 Technical data and specifications of the MACS MiniSampler Plus	119
9.7 EC/EU Declaration of Conformity	120
9.8 UK Declaration of Conformity	121
10 Technical support	123
11 Legal notes	125
11.1 Limited warranty	125

11.2 Trademarks	126
12 Index	127

Read the user manual first

WARNING

Before using the MACSQuant Instrument (MACSQuant Analyzer 10, VYB, or Analyzer 16), read the chapter **Important safety information** and all other information contained in this user manual, including any safety and operating instructions. Pay special attention to all warnings displayed on the instrument. Failure to read and follow these guidelines could lead to improper or incorrect usage and result in damage to the instrument. Improper usage could also cause severe personal injury, death, unpredictable results, instrument malfunction, and premature wear to components shortening the lifetime of the instrument. Such actions may void your warranty. Keep the user manual and any other safety and operating instructions provided with the instrument in a safe place accessible to all users for future reference.

If you have a serious concern regarding the safe use of your instrument, contact your authorized Miltenyi Biotec service provider or call Miltenyi Biotec Technical Support.

General safety instructions

THIS CHAPTER DESCRIBES POTENTIALLY HAZARDOUS SITUATIONS ASSOCIATED WITH THIS INSTRUMENT AND PROVIDES IMPORTANT SAFETY INFORMATION TO MINIMIZE THE RISKS AND PROTECT YOURSELF AND OTHERS.

This chapter provides important information for your personal safety and the correct use of the instrument. Read and observe all instructions carefully before proceeding with the installation and use of the instrument. Observe general safety practices in addition to this user manual.

- Use this instrument only as indicated in this user manual to avoid personal injury and property damage.
- Observe local working area safety instructions and laboratory policies, as well as standards for health, safety, and prevention of accidents.
- Keep this user manual in a place that is always accessible to all users.
- In case of severe accidents, damages to the instrument, or if smoke or flames appear, cut the power supply immediately.
- To entirely disconnect the instrument from the power supply, unplug the power cable.
- Ensure that the power switch and the plug of the power cable of the instrument are easily accessible.

Hazard levels

Signal words are used to warn against hazardous situations and property damages. The following signal words are used in this user manual.

⚠ WARNING or **WARNING!** indicates a hazardous situation that, if not avoided, could result in death or serious injury.

⚠ CAUTION or **CAUTION!** indicates a hazardous situation that, if not avoided, could result in minor or moderate injury. It is also used to warn against unsafe practices.

NOTICE or **NOTICE** indicates information considered important but not hazard related (e.g. messages relating to property damage).

Symbols

SAFETY WARNING: THE DOCUMENTATION MUST BE CONSULTED IN ALL CASES WHERE THIS SAFETY SYMBOL IS USED ON THE INSTRUMENT

BIOLOGICAL HAZARD

MOVING PARTS

STRONG MAGNETIC FIELD

PERSONS WEARING A PACEMAKER OR OTHER ELECTRONIC IMPLANTS MUST MAINTAIN DISTANCE

ON (POWER ON)

OFF (POWER OFF)

READ THE USER MANUAL BEFORE USING THE INSTRUMENT

FUSE

WEEE (WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT)

ORDER NUMBER

SERIAL NUMBER

TYPE NUMBER

MANUFACTURER

EUROPEAN CONFORMITY MARKING

NRTL CERTIFICATION MARK: PRODUCT MEETS CONSENSUS-BASED STANDARDS OF SAFETY, REQUIRED BY THE OCCUPATIONAL SAFETY/HEALTH ADMINISTRATION (OSHA), DETERMINED BY THE NATIONALLY RECOGNIZED TESTING LABORATORIES (NRTL) TÜV SÜD

UNITED KINGDOM CONFORMITY ASSESSED MARKING

Safety labels

Notice the hazard points and safety symbols of the instrument.

- Keep safety labels and safety markings clean and legible.
- Inspect the safety labels and safety markings regularly and replace them if they are not legible or identifiable from a safe viewing distance.
- Contact Miltenyi Biotec for replacement labels.


Figure 1: Hazard areas and safety symbols on the front side of the instrument


Figure 2: Hazard areas and safety symbols on the rear side of the instrument


Figure 3: Hazard areas and safety symbols on the left side of the instrument

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Protective housing

The housing of the instrument reduces the risk of an electric shock and short circuits. A short circuit can be caused by disassembled parts or objects that have been dropped into the instrument. Emitted sparks, caused by a short circuit, might ignite combustible vapors or material. The protective housing prevents the spread of fire. Do not use the instrument in areas classified as hazardous, for example oxygen-laden environments.

- Do not remove or penetrate any cover of the housing except for the front access covers.
- Only Miltenyi Biotec service personnel are allowed to remove any other cover of the instrument.
- Do not use the instrument if it is opened or disassembled.
- Do not use the instrument if it has been dropped or is damaged.
- Do not use the instrument if an object entered the inside of the instrument.

Overheating and fire hazards

Protect the instrument from overheating.

- Leave enough space around the instrument (at least 15 cm on all sides).
- Do not place the instrument next to heat sources, for example heaters.
- Do not cover the ventilation slots and openings of the instrument.
- Ensure adequate air circulation during operation.

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- Protect the instrument against spilled liquids.
- Clean up spilled liquids immediately.

- Do not use the instrument if liquids have entered the inside of the instrument.
- Unplug the instrument before cleaning it.
- Use only small amounts of cleaning agents on a soft cloth to wipe the instrument. Do not spray or pour liquid cleaning agents onto or into the instrument.

Humid and dusty environments

Humid and dusty environments can cause short circuits.

- Do not use the instrument in a wet location or areas with high humidity or condensation.
- After moving the instrument from a cold environment to room temperature, wait for the instrument to dehumidify before using it.

Cables and power supplies

Using the instrument with other than the supplied cables is potentially hazardous. The instrument has a three-wire electrical grounding plug with a third pin for grounding. This is a safety feature.

- Only insert the plug into a grounded electrical outlet. Do not try to insert the plug into a non-grounded electrical outlet.
- Do not use the instrument if the power cable is damaged.
- Only use the included power cable.
- Do not use extension cables or power strips.
- Do not overload an electrical outlet.

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Biological material

All biological material must be considered potentially infectious.

- Do not open the front access covers while the instrument is in operation. Exceptions are maintenance procedures that are explicitly described in this manual.
- Keep away from the robotic needle arm while the instrument is in operation.
- Wear personal protective equipment (such as gloves, safety glasses, etc.) as indicated in the safety data sheet for the particular substance.
- Droplets might leak from the system. Operate the instrument in a biological safety cabinet suitable for the used specimen if hazardous or potentially infectious materials are processed.
- Decontaminate the instrument after processing hazardous or potentially infectious material.
- If hazardous or potentially infectious material has been spilled or leaked from the system, decontaminate the affected area.
- Do not continue to use contaminated accessories or parts of the instrument.
- Run the **Clean** program first before maintenance work on any part of the fluidics system.
- Do not load or unload samples or racks while the instrument is in motion.
- Leave enough space around the MACS MiniSampler Plus (at least 15 cm on the left and right side).

Waste

All liquid and solid waste must be considered hazardous.

- Prefill the empty waste bottle with an appropriate disinfectant according to the specification of the manufacturer.
- Immediately replace the waste bottle after unmounting and fasten the bottle closure to the new bottle.
- Always have an empty waste bottle available.
- To avoid spillages, install only one bottle at a time.
- Autoclave or alternatively decontaminate waste with an appropriate disinfectant.
- Follow the general laboratory safety regulations when handling liquid and solid waste.
- Observe local regulations regarding waste disposal.

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- Always inspect the fluidics system and check for leakages before using the instrument.
- Exchange hydrophobic air filters once a year to avoid clogging through dust deposits.
- Exchange hydrophobic air filters if they came into direct contact with any liquid to avoid clogging of the filters and to prevent contamination of liquids.

Chemical hazards

WARNING

Substances and reagents can be hazardous.

- All safety measures in section **Biological hazards** also apply to any hazardous substances and reagents that may be present in the sample.
- Operate the instrument in a fume hood if hazardous substances and reagents are processed.
- Use any substances and reagents only as stated in the respective safety data sheet.

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Risk of crushing or cutting.

Robotic needle arm

The robotic needle arm moves while the instrument is in operation.

- Do not operate the instrument without the needle guard.
- Keep away from the robotic needle arm while the instrument is in operation.
- Do not obstruct the movement of the robotic needle arm.

MACS MiniSampler Plus

The carriage of the MACS MiniSampler Plus moves.

- Keep away from the MACS MiniSampler Plus while the instrument is in operation.
- Always stop or abort a process before handling accessory hardware or consumables.
- Do not load or unload samples or racks while the instrument is in motion.

Fluidics system

Syringe pumps are part of the fluidics system. They move while the instrument is in operation.

- Keep away from the fluidics system while the instrument is in operation.
- Do not open the front access covers while the instrument is in operation.

Optical radiation hazards

CAUTION

Exposure to optical radiation could result in eye injury.

Bottle illumination

Powerful LEDs are used to illuminate the bottles. Radiation of disassembled units may lead to eye injury.

According to the international standard IEC 62471, this lamp system is in excess of the Exempt Risk Group and there are risks dependent on the use of the instrument. The exposure hazard value (EHV) measured at a distance of 20 cm from the bottle holders is 0.91, the hazard distance is 19 cm for Risk Group 1.

- Do not look directly at LED radiation.
- Do not disassemble, modify, or remove the installed LED radiation sources or their mounting brackets. LED radiation sources do not automatically stop emitting when disassembled.
- Do not remove the bottle holders unless integrating the instrument into a liquid handling system.

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- Do not look at the VCSEL port through optical instruments such as lenses, magnifying glasses, etc.

- Keep a distance of 100 mm.
- Avoid that the path of the light beam is at the same height as the human eye during operation.

Laser

The instrument has three Continuous-Wave Lasers. The lasers are classified as laser class 3B. The lasers are secured by a protective housing.

- Do not remove the protective housing.

	Output power	Pulse time	Wavelength
MACSQuant Analyzer 10			
Rack detection	3.3 mW	215 μ s	850 nm
Internal laser 1	40 mW	Continuous	405 nm
Internal laser 2	30 mW	Continuous	488 nm
Internal laser 3	20 mW	Continuous	640 nm
MACSQuant Analyzer 16			
Rack detection	3.3 mW	215 μ s	850 nm
Internal laser 1	63 mW	Continuous	405 nm
Internal laser 2	50 mW	Continuous	488 nm
Internal laser 3	70 mW	Continuous	640 nm
MACSQuant VYB			
Rack detection	3.3 mW	215 μ s	850 nm
Internal laser 1	40 mW	Continuous	405 nm
Internal laser 2	50 mW	Continuous	488 nm
Internal laser 3	100 mW	Continuous	561 nm

Table 2.1: Laser output power, pulse time, and wavelength

Ergonomic hazards

WARNING

Heavy instrument.

Risk of tearing or straining muscles.

- Two people are required to lift the instrument.
- Grip the instrument at the base of the orange bottle baskets located at both sides of the instrument.

Servicing and transportation

Servicing

Improper servicing or repair of the instrument or use of unauthorized parts can cause malfunction of or damage to the instrument. This can cause hazards to users.

- Unless otherwise specifically stated in this user manual or other documentation by Miltenyi Biotec, do not service the instrument yourself.
- Servicing and repair must be performed by Miltenyi Biotec certified and qualified service personnel.

- If the instrument needs servicing, decontaminate the instrument to remove any hazardous material.
- If you have questions regarding proper decontamination or shipment, contact Miltenyi Biotec Technical Support for assistance.
- Only use accessories and upgrades recommended by Miltenyi Biotec.
- Inquire with your local Miltenyi Biotec representative about Miltenyi Biotec's extensive instrument service and service contracts, or refer to www.miltenyibiotec.com/support.

Transportation

The instrument must be transported with care in packaging specified by Miltenyi Biotec. Internal damage could occur if the instrument is subjected to excessive vibration or if it is dropped. If the instrument needs to be shipped back to the manufacturer for service, contact Miltenyi Biotec for instructions and packaging materials.

Disposal

WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT (WEEE) CUSTOMER INFORMATION

Dispose of your end-of-life Miltenyi Biotec products in accordance with the applicable WEEE and hazardous waste disposal legislation, which may differ by country or region. Electrical equipment may contain hazardous substances that may have a serious detrimental effect on the environment and/or human health. All equipment must be specifically collected and treated by designated waste facilities and by qualified WEEE compliance schemes. By ensuring that you dispose of your unwanted electrical and electronic equipment according to the applicable WEEE and hazardous waste disposal legislation, you are helping to preserve our natural resources and protect human health. Miltenyi Biotec is committed to protecting the environment. Miltenyi Biotec offers product end-of-life return programs in many countries and partners with licensed WEEE compliance schemes throughout the world. Miltenyi Biotec takes back your end-of-life Miltenyi Biotec equipment for recycling free of charge. The terms and availability of this offer vary by geography because of differences in regulatory requirements. Note that, depending on the type and use of your equipment, additional requirements may apply. Before shipping the instrument back to the manufacturer for disposal, decontaminate the instrument to remove any hazardous material. For more information or if you wish to dispose of your end-of-life Miltenyi Biotec equipment, contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support.

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This user manual applies to following instruments running with MACSQuantify Software version 3.1 or higher.

- MACSQuant Analyzer 10 (# 130-096-343)
- MACSQuant VYB (# 130-096-116)
- MACSQuant Analyzer 16 (# 130-109-803)

In the following, the name MACSQuant Instrument refers to all three instruments. Unless further specified, the stated information applies to all instrument models.

MACSQuant Instruments are compact flow cytometers that combine multisample and multiparameter sample analysis with unrivaled ease of use. With processing rates of up to 15,000 events per second, they facilitate also absolute quantification of cell populations. MACS Antibodies and Cell Separation Reagents are designed for easy and convenient handling with MACSQuant Instruments. Miltenyi Biotec is leading provider of comprehensive solutions for flow cytometry.

1.1 Intended use

The MACSQuant Instrument is a digital benchtop flow cytometer equipped with three lasers for cellular analysis. Units are equipped with an enrichment unit for the pre-analysis enrichment of rare cells. It is intended for both phenotypic analysis with fluorescent antibodies and functional analysis of fluorescently labeled cells. The MACSQuant Instrument is intended for research use only. Only laboratory professional users are allowed to use the MACSQuant Instrument. The instrument is intended to be used only as stated in the provided instructions.

1.2 Instrument description

1.2.1 Key features of the MACSQuant Instrument

- An integrated computer controlled by a touchscreen runs the MACSQuantify Software, which is used for sample acquisition as well as for data analysis.
- A fully automated computer-controlled robotic needle arm delivers samples to the sample injection port, which is located in the washing station. The uptake needle as well as the sample injection port are automatically washed during and after cell processing.
- The Single tube rack is used for the analysis of individual samples.
- The optional MACS MiniSampler Plus allows fully automated multiple sample processing. It holds different racks that can be automatically recognized by a 2D barcode reader. The MACS MiniSampler Plus can also be combined with the Universal Reagent Rack.

- A color-coded and easy to access fluidics system supplies the MACSQuant Instrument with all buffers and solutions required during and after operation. The sheath particle filter prevents any debris from entering the fluidics system.
- Bottles for MACSQuant Running Buffer, MACSQuant Washing Solution, and MACSQuant/MACSimu Storage Solution are placed in bottle holders. Liquid waste is collected in a separate waste bottle. All fluid levels are constantly monitored and visualized by the illumination of the fluid bottles.
- The MACS Cell Enrichment Unit with the optional MACSQuant Column allows easy analysis of particularly challenging cell samples, such as rare cell populations in high volumes of buffer, or precious low volume samples.


Figure 1.1: Front view of the MACSQuant Analyzer 10, VYB, and Analyzer 16

- | | | | |
|----------|--------------------------------|-----------|------------------------|
| 1 | Touchscreen | 6 | Rack detector |
| 2 | Front cover | 7 | 2D barcode reader |
| 3 | Bottle holder | 8 | LED status light |
| 4 | Main switch | 9 | Hydrophobic air filter |
| 5 | Slot for MACS MiniSampler Plus | 10 | Robotic needle arm |


Figure 1.2: Front view of the MACSQuant Analyzer 10, VYB, and Analyzer 16 without front cover

- | | | | |
|----------|------------------------------|----------|---------------------------|
| 1 | Fluidics (tubings not shown) | 3 | MACSQuant Column |
| 2 | Syringe pumps | 4 | MACS Cell Enrichment Unit |


Figure 1.3: Rear view of the MACSQuant Analyzer 10, VYB, and Analyzer 16

- | | | | |
|----------|------------------------|-----------|---------------------------|
| 1 | CFast-card slot | 7 | Power CAN |
| 2 | COM port | 8 | Power cable port |
| 3 | USB ports (4×) | 9 | External CAN |
| 4 | LAN ports (2×) | 10 | RS232/BCR |
| 5 | Display ports (2×) | 11 | Bottle sensor port |
| 6 | Sheath particle filter | 12 | Bottle distribution block |

2

Installation

Unpacking and installation of the MACSQuant Instrument should be only done by qualified Miltenyi Biotec representatives.

⚠ WARNING

Electric shock, short circuit, overheating, fire, and explosion could result in death or serious injury.

- Do not use the instrument if it has been dropped or is damaged.

The instrument is designed to fit onto benchtops, into standard size laminar flows, into safety cabinets, or in connection to fully automated liquid handling systems. It should be installed on a stable, flat, and vibration-free surface. The operating environment should be dust-free, sufficiently ventilated, and free from sources of electromagnetic radiation. To ensure a flat surface in laminar flow hoods, the MACSQuant Instrument can be placed on an autoMACS Pro Laminar Hood Plate.

2.1 Components included in the delivery

Component
MACSQuant Instrument
Fluid sensor cable module
Power cable
6 × Fluidics tubing
4 × Bottle closure with sensors
Single tube rack
User manual
Empty bottles with lids
Keyboard
Hydrophobic air filters
MACSQuant Advanced Calibration Beads
MACSQuant Starting Buffer Kit

Table 2.1: Components included in the delivery

2.2 Consumables

Component	Description	Order no.
MACSQuant Columns	For use with the MACS Cell Enrichment Unit	130-094-458
Hydrophobic Air Filter	To vent fluid bottles	130-090-385
MACSQuant Advanced Calibration Beads	For instrument calibration	130-134-208
MACSQuant Running Buffer	6 × 1.5 L	130-092-747
MACSQuant Running Buffer (16x)	1 × 1.25 L	130-111-562
MACSQuant Running Buffer (16x)	6 × 1.25 L	130-111-747
MACSQuant Washing Solution	6 × 1.5 L	130-092-749
MACSQuant/MACSimas Storage Solution	6 × 1.5 L	130-092-748
MACSQuant Starting Buffer Kit	4 × 1.5 L MACSQuant Running Buffer 1 × 1.5 L MACSQuant Washing Solution 1 × 1.5 L MACSQuant/MACSimas Storage Solution	130-094-190
MACS MultiSpecies Comp Beads	One-bottle solution for compensation of the fluorescence spillover from fluorochrome-conjugated antibodies from multiple species	130-137-398
MACS All-in-One Comp Bead Kit	Single product for compensation of the fluorescence spillover from amine-reactive live/dead discriminating dyes and fluorochrome-conjugated antibodies from multiple species	130-138-867
MACS Comp Bead Kit, anti-REA	For compensation of the fluorescence spillover from fluorochrome conjugated REAfinity antibodies	130-104-693
MACS Comp Bead Kit, anti-mouse Igk	For compensation of the fluorescence spillover from mouse Igk fluorochrome-conjugated antibodies	130-097-900
MACS Comp Bead Kit, anti-rat Igk	For compensation of the fluorescence spillover from rat Igk fluorochrome-conjugated antibodies	130-107-755
MACS Comp Bead Kit, anti-human Igk	For compensation of the fluorescence spillover from human Igk fluorochrome-conjugated antibodies	130-104-187

Table 2.2: Consumables

2.3 Accessories

Complementary components can be ordered separately.

Component	Description	Order no.
MACS MiniSampler Plus	Use with the Universal Reagent Rack and one of the available Chill Racks for automated multisample processing	130-105-745
Chill 5, 15, 50 Rack Set	Set of cooling racks	130-097-038
Chill 5 Rack	Cooling rack for 24 × 5 mL tubes	130-092-951
Chill 15 Rack	Cooling rack for 15 × 15 mL and 5 × 5 mL tubes	130-092-952
Chill 50 Rack	Cooling rack for 6 × 50 mL, 3 × 15 mL, and 3 × 5 mL tubes	130-092-953
Chill 96 Rack	Cooling rack for 96-well plates	130-094-459
MACSQuant Buffer Supply Station	To connect larger fluid containers	130-101-841
Universal Reagent Rack	Rack for autolabeling reagents	130-115-722
MACS Reagent Rack 4	Rack for autolabeling reagents	130-094-574

Table 2.3: Accessories available for the MACSQuant Instrument

2.4 Service contracts

For information about available service contracts for the MACSQuant Instrument, visit www.miltenyibiotec.com.

2.5 Assembly

2.5.1 Installing the fluid bottles

WARNING

Biohazardous waste.

Contamination or infection may lead to severe personal injury or death, depending on the material used.

- Prefill the empty waste bottle with an appropriate disinfectant according to the specification of the manufacturer.

CAUTION

Powerful LEDs are used to illuminate the bottles.

Exposure to optical radiation may cause eye injury.

- Do not remove the bottle holders.

Operating the instrument requires MACSQuant Running Buffer, MACSQuant Washing Solution, and MACSQuant/MACSiMa Storage Solution. An additional bottle is required for waste collection. Always operate the instrument with ready-to-use MACSQuant buffers and solutions only. The bottles, closures, tubings, and fluid sensor cables are color-coded for easier handling.


Figure 2.1: Connect fluid bottles

- 1 Install only one fluid bottle at a time. Place the new bottle (3) into the orange bottle holder (4). Note the corresponding color coding.
- 2 Unscrew the lid of the fluid bottle and replace it with the appropriate bottle closure delivered with the instrument.
- 3 Remove the caps from the bottle distribution block and connect the tubing to the appropriate color-coded fluid port of the bottle distribution block on the rear side of the instrument. See **Figure 1.3**.
- 4 Note the color coding of the tubes and connect them to the respective bottle closure (2).
- 5 Connect the 4-end fluid sensor cable to the port labeled **Bottle Sensor** at the rear side of the instrument. See **Figure 1.3**. Fasten the screws of the sensor cable plug securely.
- 6 Note the color coding and connect each fluid sensor cable to the respective bottle closure (1).

NOTICE! Do not screw in the sensor cable.

- 7 Proceed with **Installing hydrophobic air filters on the facing page**.

	Running Buffer	Washing Solution	Storage Solution	Waste container
Icon				
Color of sensor cables and tubes	Blue	Green	Black	Red

Table 2.4: Bottle symbols and their respective color-coding

2.5.2 Installing hydrophobic air filters

The MACSQuant Instrument is equipped with hydrophobic air filters to allow air flow and to prevent contamination. Use only hydrophobic air filters from Miltenyi Biotec. See **Consumables on page 70**.

- 1 Remove the Chill Rack.
- 2 Open the front cover.
- 3 Connect hydrophobic air filters via Luer connectors. See **Table 2.5** for an overview about the locations. Avoid any contact with liquids because this could cause clogging of the filter.

Location	Amount
On fluidic bottles	4×
Behind the front cover	1×
On the rear side of the instrument	1×

Table 2.5: Locations of hydrophobic air filters on the MACSQuant Analyzer 10, VYB, and Analyzer 16

2.5.3 Installing the Single tube rack

The instrument is delivered with the Single tube rack. The Single tube rack is compatible with all standard 5 mL flow cytometry tubes and can also be used with 1.5 mL and 2 mL reaction tubes.

- 1 Gently insert the Single tube rack (1) into the corresponding slots (2) located at the front of the instrument. The rack should click into place.
- 2 To remove the Single tube rack, hold it between two fingers and pull gently.


Figure 2.2: Installation of the Single tube rack

2.5.4 Connecting the instrument to the power supply

WARNING

Electric shock, short circuit, overheating, fire, and explosion could result in death or serious injury.

- Do not use the instrument if the power cable is damaged.
- Only use the included power cable.

- 1 Connect the power cable to the power port at the rear side of the instrument. The main switch is located on the right-hand side of the instrument. See **Figure 1.1**. Ensure that the main switch is in position O (off) before connecting the power cable.
- 2 Attach the cable coming from the 2D barcode reader to the port labeled **RS232/BCR** at the rear side of the instrument and fasten securely.
- 3 Plug in the power cable to a grounded electrical outlet.
- 4 Connect the keyboard included in the delivery to a USB port.
- 5 Use the main switch on the right-hand side on the instrument to turn the instrument into standby mode.
- 6 Tap the screen to start the MACSQuantify Software. Wait until the MACSQuantify Software shows the sign-in window.

2.5.5 Installing the MACSQuant Column

WARNING

The instrument has a powerful magnet.

Magnetizable objects can suddenly move towards the magnet.

- Keep all magnetic storage devices, electronic equipment, and magnetizable objects at a distance of at least 30 cm from the instrument.

WARNING

The instrument has a powerful magnet.

Strong magnetic fields can influence the functioning of pacemakers or electronic medical implants.

If wearing pacemakers or electronic medical implants, keep a distance of at least 30 cm from the instrument.

The MACSQuant Instrument is delivered with a column substitute, which must be replaced by a MACSQuant Column before using the MACS Cell Enrichment Unit. The MACS Cell Enrichment Unit is integrated in the instrument. Magnetically and fluorochrome-labeled cells are pre-enriched in the MACSQuant Column prior to cell analysis allowing flow cytometric analysis of rare cell populations in a fraction of time. Up to 5×10^6 magnetically labeled cells can be enriched and analyzed in a single step. The entire process of labeling, pre-enrichment, and flow cytometric analysis is fully automated. The MACSQuant Column can be used for up to three months before it needs to be replaced. Keep the column substitute, it can be re-installed if the MACSQuant Column is no longer needed.

Follow the instructions below to replace the column substitute or to exchange an old column with a new column.

- 1 Open the front cover of the instrument to access the MACS Cell Enrichment Unit with the column substitute or MACSQuant Column.

- 2 Hold the top and bottom of the column (1) with both hands and pull gently but firmly to remove the column from the MACS Cell Enrichment Unit (2).


Figure 2.3: MACSQuant Column removal

- 3 Hold the column (4) in one hand and gently unscrew the top column connector (3) counterclockwise. Tilt the column downwards to empty any liquid into a paper towel.


Figure 2.4: MACSQuant Column overview

- 4 Unscrew the bottom column connector (5) to remove the column.
- 5 Screw the bottom column connector into the bottom end of the new column. Note the correct orientation of the column. The top end of the MACSQuant Column has a clearly visible spacer, while the metal beads reach further to the bottom end of the column. Insert the column only with the spacer in upward position.

- 6 Point the column towards the top of the MACS Cell Enrichment Unit (2) and screw in the top column connector.
- 7 Align the column so that the top column connector sits on the guide of the magnet cover. Push the column into the slot until you hear a click. Verify that the column is placed in the center of the magnet cover.
- 8 Close the front cover of the instrument.
- 9 Log in as MQ Administrator.
- 10 To prime the new column, go to the **Tools** tab in the side pane of the MACSQuantify Software.
- 11 Select **Prime Column**. Refer to the MACSQuantifySoftware user manual for further information about Pre-enrichment of cells.

2.5.6 Installing the needle guard

- 1 Check if the needle guard (1) is attached to the robotic needle arm (2).
- 2 If not, install as depicted in **Figure 2.5**.


Figure 2.5: Attachment of the needle guard

2.5.7 Installing the MACS MiniSampler Plus

If using the instrument with a MACS MiniSampler Plus, follow the instructions below. The MACS MiniSampler Plus can be installed to either the autoMACS Pro or the MACSQuant Instrument (MACSQuant Analyzer 10, VYB, Analyzer 16).

Use a MACS MiniSampler Plus only with one type of instrument. After a MACS MiniSampler Plus was used with a MACSQuant Instrument, it cannot be used with the autoMACS Pro anymore.

- 1 Switch off the instrument before installing the MACS MiniSampler Plus.
- 2 Remove the transparent protection foil from the lens of the rack detector. Note the position of the MACS MiniSampler Plus slot (1) located at the front of the instrument and the MACS MiniSampler Plus guiding (2).
- 3 Tilt the MACS MiniSampler Plus.


Figure 2.6: Install the MACS MiniSampler Plus

- 4 Slide the guiding into the receiving slot. Keep the MACS MiniSampler Plus in a tilted position until resistance is met.
- 5 Lower the MACS MiniSampler Plus to a horizontal position to lock it.


Figure 2.7: MACS MiniSampler Plus in the horizontal position

- 6 Open the front door of the instrument and fix the MACS MiniSampler Plus (4) with the screw (3).


Figure 2.8: Fix MACS MiniSampler Plus in installation position

- 7** Close the front cover of the instrument and ensure that the MACS MiniSampler Plus is properly set and secured.
- 8** Guide the MACS MiniSampler Plus cable underneath the instrument and connect it to the port labeled **External CAN** at the rear side of the instrument.

2.5.8 Installing the Chill Rack and Universal Reagent Rack

- 1 Install the Universal Reagent Rack (1) onto the MACS MiniSampler Plus (3) into the left recess by snapping the engagement hook (2) into the undercut.


Figure 2.9: Install the Universal Reagent Rack

- 2 Place a MACS Chill Rack (4) onto the MACS MiniSampler Plus (3) into the right-hand recess, ensuring that the rack bar code is facing the instrument. If automatic recognition of the tube rack fails, the instrument shows a window for manual selection of the tube rack. Ensure the rack is placed correctly into the recess before confirming the selection of the rack.


Figure 2.10: Install the MACS Chill Rack

- 3 Note the position of the lid guiding (5) at both sides of the MACS MiniSampler Plus and attach the protective cover (6). Keep the protective cover attached and closed while the instrument is in operation.


Figure 2.11: Attach the protective cover

2.5.9 Installing a webcam

Connect a webcam to a USB port at the rear side of the instrument. The webcam is recognized by the software automatically during a MACSQuant Live Support session.

Unplug the webcam during normal operation of the instrument. Use it only during a MACSQuant Live Support session.

2.5.10 Installation checklist

Ensure the correct installation of the instrument by the following checklist:

- Note the correct positioning of each fluid bottle and the waste bottle, recognizable by the color code and the symbols. This is crucial for successful analysis.
- Ensure that the correct tubing and fluid sensor cable is attached to the corresponding bottle, and that all tubing connections as well as the bottle closures are fastened. If necessary, tighten loose connections.
- Ensure that the fluid bottles are filled and installed correctly. For safe operation of the instrument, all fluid bottles must contain at least 150 mL of the respective solution. In order to prevent contamination of the fluidics system the use of non-sterile buffers and solutions is not recommended. Only use buffers and solutions supplied by Miltenyi Biotec for operation of the instrument. A reproducible and optimal performance of the instrument cannot be guaranteed when the instrument is operated using home-made buffers and/or solutions purchased from another manufacturer.
- Make sure that the waste bottle is empty and installed correctly.

WARNING! When working with biohazardous samples, it is recommended to fill the waste bottle with 100 mL of disinfectant before use. Empty bottles from the Running Buffer, Washing Solution, or Storage Solution may be used as waste bottles only with the appropriate safety label on it. For proper disposal, follow local regulation.

- Open the front cover and check that the column substitute or the MACSQuant Column is installed correctly. Ensure that the tubes are securely fastened to the column and that no part of the visible tubing is pinched or obstructed.
- Ensure that the needle guard is attached to the instrument.

WARNING! There is a hazard of crushing, shearing, or puncturing bodily parts. Operate only with attached needle guard.

When all points of this installation checklist are fulfilled, close all covers.

3

Operating the instrument

3.1 Switching on

- 1 Use the main switch on the right-hand side of the instrument to switch the instrument to standby mode. A red LED in the lower left corner of the monitor indicates that the instrument is in standby mode.
- 2 Tap the touchscreen to start the MACSQuantify Software. A sign-in window appears after a short initialization phase.
- 3 Optional: See **Priming the instrument below** to switch to acquisition mode.

For a detailed instruction how to use the MACSQuantify Software, refer to the MACSQuantify Software user manual.

3.2 Priming the instrument

WARNING

Defective or inadequate equipment can cause a biological hazard. Contamination or infection may lead to severe personal injury or death, depending on the material used.

- Always inspect the fluidics system and check for leakages before using the instrument.

- 1 To prime the instrument (power up the lasers, activate fluidics), click the **Power** button in the upper right corner of the software and select **Acquisition mode**.
- 2 Let the optical bench warm up after priming for at least 30 minutes.
- 3 Go to **View > Hardware > Lasers and detectors** to check the current laser temperature. See **Lasers and detectors on page 94**.

While the instrument is warming up, it is recommended to perform a **Flush** cycle followed by a **Clean** program. See **Cleaning programs on page 101**.

3.3 Switching off

Before switching off the instrument, start the shutdown procedure. During the shutdown procedure, the fluidics system is cleaned and filled with MACSQuant/MACSima Storage Solution.

The shutdown procedure can be started manually. Alternatively, set the instrument to automatic shutdown after a defined idle time. Refer to the MACSQuantify Software user manual for details.

3.3.1 Manual shutdown

- 1 Click the **Power** button at the upper right corner of the screen to shut down manually.
- 2 Select **Instrument off** to switch the instrument to standby mode, or select **Data analysis** to proceed with data analysis after the shutdown procedure. After a seven-minute washing procedure (default setting), the instrument is in standby or data analysis mode.
- 3 From data analysis mode select **Instrument off** to switch to standby mode.
- 4 From standby, use the main switch on the right-hand side of the instrument to switch off completely.

3.3.2 Automatic shutdown

Automatic shutdown is enabled by default after a defined idle time. A MQ Administrator can change the timer settings. Refer to the MACSQuantify Software user manual.

4

Calibration

Before using the MACSQuant Instrument for the first time, the hardware must be calibrated. Hardware calibration can only be performed by a MQ Administrator. To calibrate the hardware the user must be familiar with the MACSQuantify Software. Refer to the MACSQuantify Software user manual.

4.1 Calibrating the sample uptake unit

4.1.1 Calibrating the needle arm

 **CAUTION**

Moving robotic needle arm.

Risk of crushing or cutting.

- Do not obstruct the movement of the robotic needle arm.
- Keep away from the robotic needle arm while the instrument is in operation.

During calibration, the needle arm moves between the samples and the washing station with its sample injection port and washing port.



Figure 4.1: The washing station, top view (needle guard not shown)

- | | | | |
|---|-----------------------|---|--------------------|
| 1 | Sample injection port | 3 | Washing port |
| 2 | Waste port | 4 | Robotic needle arm |

- Go to the **Tools** tab in the side pane of the MACSQuantify Software and select the **Calibrate uptake unit** box. A red closed circle indicates that the respective uptake component is not calibrated.


- Click **Start calibration** in the section **Needle arm calibration relative to dropping station**. The robotic needle arm moves towards the washing station before being inserted into the sample injection port.

If the uptake needle pops out of the robotic needle arm holder, a system failure is caused and a prompt screen to reinitialize appears. Re-insert the uptake needle, adjust the robotic needle arm appropriately, and click **OK** to confirm.

- Manually adjust the robotic needle arm to place the uptake needle over the center of the sample injection port.
- Gently lower the robotic needle arm manually into the sample injection port until it makes first contact with the bottom of the orifice. The uptake needle must not pop out of the needle holder and always remain perpendicular to the horizontal plane. Do not bend the needle or insert it diagonally.
- In the field Z-offset, enter 30.
- Click **Save calibration**.

- The closed green circle indicates that the calibration is completed. Click **Test** to confirm successful calibration.

4.2 Calibrating racks and plates

4.2.1 Calibrating the Single tube rack

- Go to the **Tools** tab in the side panel of the MACSQuantifySoftware and select the **Calibrate uptake unit** box. Red closed circles indicate that the uptake components are not calibrated.
- Attach the Single tube rack with an empty sample tube to the instrument. See **Installing the Single tube rack on page 73**.
- Click **Start calibration** in the section **Single tube rack**. The robotic needle arm moves automatically forward.
- Adjust the uptake needle to the center of the tube (i.e. equidistant from the tube edges). Move the needle arm forward or backward.
- Lower the robotic needle arm carefully until it almost touches the tube bottom (i.e. only a fraction of a millimeter from the bottom of the tube). To check the uptake needle position, gently wiggle the tube to ensure that there is a small amount of movement.
- Click **Save calibration**. The closed green circle indicates that the calibration is completed. Click **Test** to confirm that the correct coordinates are saved.

4.2.2 Calibrating the Chill 5 Rack

For the correct usage of the MACS MiniSampler Plus, calibrate the Chill 5 Rack.

- Attach the MACS MiniSampler Plus to the MACSQuant Instrument.
- Go to the **Tools** tab in the side pane of the MACSQuantify Software and select the **Calibrate uptake unit** box. Red closed circles indicate that the uptake components are not calibrated.

- 3 Place the Chill 5 Rack onto the MACS MiniSampler Plus and place an empty 5 mL tube in position D6.
- 4 Click **Start calibration** in the section **Chill 5 rack**. The robotic needle arm automatically moves forward.
- 5 Adjust the uptake needle to the center of the tube in position D6 (i.e. equidistant from the tube edges). Move the needle arm forward or backward. Move the MACS MiniSampler Plus for movement to the right and left side.
- 6 Lower the robotic needle arm carefully until it almost touches the tube bottom (i.e. only a fraction of a millimeter from the bottom of the tube). To check the uptake needle position, gently wiggle the tube to ensure that there is a small amount of movement.
- 7 Click **Save calibration**. The closed green circle indicates that the calibration is completed.
- 8 Select test position in the rack window in the side pane.
- 9 Click **Test** to confirm successful calibration. The MACSQuant Instrument tests selected positions.

4.2.3 Calibrating the Chill 96 Rack

- 1 Go to the **Tools** tab in the side pane of the MACSQuantify Software and select the **Calibrate uptake unit** box. Red closed circles indicate that the uptake components are not calibrated.
- 2 Place a Chill 96 Rack onto the MACS MiniSampler Plus and load an empty 96-well plate.
- 3 Click **Start calibration** in the section **Chill 96 rack**. The robotic needle arm automatically moves forward.
- 4 Adjust the uptake needle to the center of the well H12 (i.e. equidistant from the well edges). Move the needle arm forward or backward. Move the MACS MiniSampler Plus for movement to the right and left side.
- 5 Lower the robotic needle arm carefully until it almost touches the well bottom (i.e. only a fraction of a millimeter from the bottom of the well). To check the uptake needle position, gently wiggle the plate to ensure that there is a small amount of movement.
- 6 Click **Save calibration**. The closed green circle indicates that the calibration is completed.
- 7 Select test position in the rack window in the side pane.
- 8 Click **Test** to confirm successful calibration. The MACSQuant Instrument tests selected positions.

4.2.4 Calibrating the MACSPlex Filter Plate

- 1 Go to the **Tools** tab in the side pane of the MACSQuantify Software and select the **Calibrate uptake unit** box. Red closed circles indicate that the uptake components are not calibrated.
- 2 Place a Chill 96 Rack onto the MACS MiniSampler Plus and load a MACSPlex Filter Plate.
- 3 Click **Start calibration** in the section **MACSPlex Filter Plate**. The robotic needle arm automatically moves forward.
- 4 Adjust the uptake needle to the center of the well H12. Move the needle arm forward or backward. Move the MACS MiniSampler Plus for movement to the right and left side.
- 5 Lower the robotic needle arm carefully until it almost touches the well bottom (i.e. only a fraction of a millimeter from the bottom of the well).
- 6 Go to the field **Z-offset** and enter a value of -7.
- 7 Press enter.

Be careful not to damage the membrane of well H12 when calibrating the MACSPlex Filter Plate.

- 8 Click **Save calibration**. The closed green circle indicates that the calibration is completed.
- 9 Select test position in the rack window in the side pane.
- 10 Click **Test** to confirm successful calibration. The MACSQuant Instrument tests selected positions.

4.2.5 Calibrating the Universal Reagent Rack

- 1 The MACSQuant Instrument can be used in combination with both the Universal Reagent Rack and the Reagent Rack 4. The instructions for calibration apply to both.
- 2 Go to the **Tools** tab in the side pane of the MACSQuantify Software and select the **Calibrate uptake unit** box. Red closed circles indicate that the uptake components are not calibrated.
- 3 Place the Universal Reagent Rack onto the MACS MiniSampler Plus.
- 4 Place an empty reagent vial on position 4 (front-most position) of the Universal Reagent Rack.
- 5 Remove the cap from the reagent vial.
- 6 Click **Start calibration** in the section **Reagent Rack 4**.
- 7 Select a reagent rack from the drop-down menu. A dialog box will prompt you to ensure that the Universal Reagent Rack is installed correctly and that the Single tube rack is removed.
- 8 Confirm the dialog box by clicking **OK**. The robotic needle arm automatically moves forward.
- 9 Adjust the uptake needle to the deepest point of the vial on position 4.
- 10 Lower the robotic needle arm carefully until it almost touches the vial bottom (i.e. only a fraction of a millimeter from the bottom of the vial). To check the uptake needle position, gently wiggle the vial to ensure that there is a small amount of movement.
- 11 Click **Save calibration**. The closed green circle indicates that the calibration is completed.
- 12 Click **Test** to confirm successful calibration. The MACSQuant Instrument tests all positions.

4.2.6 Calibrating the custom 96 racks

If 96-well plates of different format (e.g. deep well or round bottom) are used, they might have a different distance from the bottom of the rack to the bottom of the well. To pick the samples correctly, individual custom racks can be defined for each required rack format.

- 1 Go to the **Tools** tab in the side pane of the MACSQuantify Software and select the **Calibrate uptake unit** box. Red closed circles indicate that the uptake components are not calibrated.

If the rack type is not shown, use the next or previous window button in the **Calibrate uptake unit** window to switch to another page.

- 2 Click **Start calibration** in the section **Custom 96 Well Rack 1**. The robotic needle arm automatically moves forward.

If several types of 96 well plates are used, it is recommended to calibrate each type as a custom rack.

- 3 Proceed as described for standard 96 well plates. See **Calibrating the Chill 96 Rack on the previous page**.

4.3 Calibrating the rack detection

CAUTION


Direct exposure to laser beam could result in eye or even skin injury.

The instrument has four vertical cavity surface emitting lasers (VCSEL) for automated rack detection (Class 1M). The radiation is not visible.

- Keep a distance of 100 mm.

4.3.1 Calibrating the Single tube rack detection

- 1 Go to the **Tools** tab in the side panel of the MACSQuantify Software and select the **Calibrate rack detection** box. Red closed circles indicate that the components are not calibrated.
- 2 Click **Start calibration** in the section **Single tube rack detection**.


- 3 Follow the instructions on the screen.
- 4 The closed green circle indicates that the calibration is completed.

4.3.2 Calibrating the rack detection

This step is required for the automatic detection of Chill 5, 15, and 50 Racks by the MACSQuant Instrument if using the MACS MiniSampler Plus.

- 1 Attach the MACS MiniSampler Plus to the MACSQuant Instrument.
- 2 Go to the **Tools** tab in the side pane of the MACSQuantify Software and select the **Calibrate rack detection** box. Red closed circles indicate that the components are not calibrated.
- 3 Place a Chill 5 Rack onto the MACS MiniSampler Plus and click **Start calibration** in the section **Rack detection**.

5

Instrument monitoring

5.1 Hardware monitor

The hardware monitor can be accessed under the **View** menu. With the hardware monitor the real time status of the MACSQuant Instrument can be assessed, for example, to provide additional information in the case of error messages appearing on the screen. The serial number of the MACSQuant Instrument and photomultiplier tube (PMT) information can be found in the hardware monitor.

5.1.1 Fluidics

The **Fluidics** tab in the hardware monitor displays the pumps and valves as well as the status of the fluid bottles in live status.

The fluid bottles in the hardware monitor and on the instrument are illuminated in red when empty. The waste bottle is illuminated in red when full.

The waste (W), air (A), and fill pump symbols (F) are illuminated in green when active. Active valves 1-3 are indicated in green, defective valves in red. Valve status is indicated by o for open and c for closed.

The air pressure should be 300 mbar after priming when the instrument is in acquisition mode.

The reservoir indicates the level of MACSQuant Running Buffer in the reservoir. The level should be high or normal during acquisition mode.


Figure 5.1: Real-time hardware monitor of the fluidics components. Specifications may vary between instrument models.

Component	Further information
Fluid container	Displays liquid levels of bottles in real time (full/empty).
Pumps	Displays the status of the waste (W), air (A), and fill (F) pumps.
Separation unit	Displays the status of the MACSQuant Column and MACS Cell Enrichment Unit
Valves 1–3	Displays the position of the valves for the MACSQuant Column: C = closed, O = open, green = in use.
Valves 4–6	Displays the position of the general fluidics system valves.
Syringe drive	Displays the position of the plunger of the syringes.
Sensors	Displays the general system pressure and level of the sheath fluid reservoir. Air pressure refers to pressure applied when instrument is in acquisition mode, for example, 300 mbar.

Table 5.1: Panels of the **Fluidics** tab

5.1.2 Sample uptake unit

The **Sample uptake unit** tab indicates the status of the robotic needle arm and if the MACS MiniSampler Plus is connected. The relative position of the robotic needle arm is listed in the lower left box and if the Single tube rack is calibrated.


Figure 5.2: Real-time hardware monitor of the sample uptake unit

5.1.3 Lasers and detectors

The **Lasers and detectors** tab displays the status of the optical bench. The temperature, fan speed, PMT voltage, filter settings, and annotated path of each laser is shown. The operating temperature for the main laser should be between 30 °C and 50 °C. The main laser is the blue laser for MACSQuant Analyzer 10, and 16, and the yellow laser for the MACSQuant VYB. The side lasers are kept at a temperature of 25 °C. When in acquisition mode, colored lines from the lasers to the detectors indicate that lasers are switched on. The bench temperature is kept between 31 °C and 37 °C. The fan speed is regulated automatically depending on the room temperature and internal temperature of the MACSQuant Instrument. In case of errors, contact Miltenyi Biotec Technical Support.


Figure 5.3: Real-time hardware monitor of the MACSQuant Analyzer 16 optical bench. T indicates the current temperature, P the power of the main laser. Specifications may vary between instrument models.

5.2 Instrument status LEDs

The MACSQuant Instrument is equipped with LEDs, illuminating each bottle to indicate the status of the instrument in acquisition mode. If the fluid bottles are not illuminated, the MACSQuant Instrument is in data analysis mode and the lasers are off.

Bottle illumination	Description
Green	Instrument is in acquisition mode and ready to measure
Blue	Instrument processes a sample and is busy
Red	Critical bottle level or general instrument error. The blinking red light indicates which bottle is affected. A message in the Instrument status bar specifies the instrument error.

Table 5.2: Color coding of the fluid bottle illumination on the MACSQuant Instrument

6

MACSQuant Live Support

MACSQuant Live Support is a real-time diagnostic service provided by Miltenyi Biotec Technical Support. The Live support is available for all users without any costs 24 h a day from Monday through Friday. Highly trained MACSQuant specialists can be reached to assist with any queries. If required, the MACSQuant Instrument can be remotely controlled for diagnosis or file transfer. A webcam can be connected to a free USB port at the rear side of the instrument. See **Installing a webcam on page 81**. The MACSQuant Instrument must have network access to the internet to use the MACSQuant Live Support

For details on how to set up and use the MACSQuant Live Support, refer to the MACSQuantify Software user manual.

7

Maintenance

WARNING

Electric shock, short circuit, overheating, fire, and explosion could result in death or serious injury.

- Do not remove or penetrate any cover of the housing except for the front access covers.
- Only Miltenyi Biotec service personnel are allowed to remove any other cover of the instrument.

WARNING

Biohazardous material.

Contamination or infection could result in death or serious injury depending on the material used.

- Run the **Clean** program first before maintenance work on any part of the fluidics system.

Instrument parts wear over time and require regular maintenance and/or replacement. Contact Miltenyi Biotec Technical Support regarding spare parts. Miltenyi Biotec also offers service contracts for the MACSQuant Instruments. Visit www.miltenyibiotec.com for further information.

7.1 Surface disinfection

The uptake needle of the needle arm and the surface of the instrument must be decontaminated upon contact with biohazardous samples. Use 70% ethanol for disinfection of the instrument surface. Use 70% ethanol or isopropyl alcohol for cleaning of the needle arm.

WARNING! Liquids inside the instrument can cause short circuits. Unplug the instrument before cleaning it. Use only small amounts of cleaning agents on a soft cloth to wipe the instrument. Do not spray or pour liquid cleaning agents onto or into the instrument.

7.2 Disinfection of the fluidics system

The fluidics system must be decontaminated by running the standard shutdown sequence if the instrument was used to process biohazardous samples. Additionally, use the **Clean** program before starting the standard shutdown sequence. See **Cleaning programs on page 101**.

7.3 Exchanging buffer and solution bottles

- 1 Install one bottle at a time. Remove the bottle to be exchanged from the bottle holder.
- 2 Place a new, sealed bottle into the orange bottle holder. Note the corresponding color-coding.
- 3 Unscrew the lid of the new fluid bottle.

- 4 Remove the bottle closure from the old bottle and install it on the new bottle. Do not disconnect the color-coded tubing and sensor cable from the bottle closure. Take care not to invert the bottle closures as this could cause liquid to enter the hydrophobic air filter.

7.4 Exchanging the waste bottle

Before exchanging the waste bottle, wait for ongoing measurements or washing cycles to be finished to avoid spillage.

- 1 Place an empty, safety labeled (biohazard) bottle next to the orange bottle holder.
- 2 Unscrew the cap of the empty bottle.
- 3 Remove the bottle closure from the full waste bottle in the bottle holder, and install the bottle closure on the empty bottle next to the bottle holder.
- 4 Close the full waste bottle and remove it from the bottle holder.
- 5 Place the empty bottle into the bottle holder.

7.5 Exchanging hydrophobic air filters

WARNING

Defective or inadequate equipment can cause a biological hazard.

Contamination or infection could result in death or serious injury depending on the material used.

- Exchange hydrophobic air filters once a year to avoid clogging through dust deposits.
- Exchange hydrophobic air filters if they came into direct contact with any liquid to avoid clogging of the filters and to prevent contamination of liquids.

The MACSQuant Instrument is equipped with hydrophobic air filters to allow air flow and to prevent contamination. Use only hydrophobic air filters from Miltenyi Biotec. See **Consumables on page 70**.

- 1 Remove the Chill Rack.
- 2 Open the front cover.
- 3 Unscrew the hydrophobic air filters. For an overview about the location of hydrophobic air filters, see **Table 2.5**.
- 4 Replace with a new hydrophobic air filter.

7.6 Exchanging the MACSQuant Column

WARNING

The instrument has a powerful magnet.

Magnetizable objects can suddenly move towards the magnet.

- Keep all magnetic storage devices, electronic equipment, and magnetizable objects at a distance of at least 30 cm from the instrument.

⚠ WARNING

The instrument has a powerful magnet.

Strong magnetic fields can influence the functioning of pacemakers or electronic medical implants.

If wearing pacemakers or electronic medical implants, keep a distance of at least 30 cm from the instrument.

The MACSQuant Column should be exchanged every three months or replaced with the column substitute after three months.

- 1 Open the front cover of the instrument to access the MACS Cell Enrichment Unit with the column substitute or MACSQuant Column.
- 2 Hold the top and bottom of the column with both hands and pull gently but firmly to remove the column from the MACS Cell Enrichment Unit.
- 3 Hold the column in one hand and gently unscrew the top column connector counterclockwise. Tilt the column downwards to empty any liquid into a paper towel.
- 4 Unscrew the bottom column connector to remove the column.
- 5 Screw the bottom column connector into the bottom end of the new column. Note the correct orientation of the column. The top end of the MACSQuant Column has a clearly visible spacer while the metal beads reach further to the bottom end of the column. Insert the column only with the spacer in upward position.
- 6 Point the column towards the top of the MACS Cell Enrichment Unit and screw in the top column connector.
- 7 Align the column so that the top column connector sits on the guide of the magnet cover. Push the column into the slot until you hear a click. Verify that the column is placed in the center of the magnet cover.
- 8 Close the front access cover.
- 9 To prime the new column, go to the **Tools** tab in the side pane of the MACSQuantify Software.
- 10 Select **Prime Column**. Refer to the MACSQuantify Software user manual for further information about Pre-enrichment of cells.

7.7 Cleaning the instrument

7.7.1 Cleaning programs

Cleaning programs	Description
Rinse	Rinses the needle with MACSQuant Running Buffer. Duration 2 minutes.
Clean	Cleans the needle and flow cell with 0.25 mL 1% hypochlorite solution. Duration 12 minutes.
Flush	System backflush rinse if a blockage of the fluidics system occurs. Duration 16 minutes.

Table 7.1: Cleaning programs available for the MACSQuant Instrument

To start a cleaning program, select one of the following options:

- Click the **Rinse** button in the status bar to start a needle **Rinse**.
- Right-click **Rinse** in the status bar. Select **Clean** to start a **Clean** program. Follow the instructions on the screen.
- Right-click **Rinse** in the status bar. Select **Flush** to start a **Flush** program. Follow the instructions on the screen.

Figure 7.1: Rinse button in the status bar with the options **Clean** and **Flush**

7.7.2 Cleaning recommendations

Indication	Cleaning programs
After each user	Clean program
After sticky material, such as bone marrow or tumor samples	Clean program followed by a Flush program
During laser warm up in acquisition mode	Flush program followed by a Clean program
Debris or high noise	Clean program, hot bleach
Blockage of the fluidics system can also be used as weekly cleaning routine	Hot bleach, see Removing a potential blockage below

Table 7.2: Cleaning recommendations for the MACSQuant Instrument


See **Cleaning programs on the previous page** and **Removing a potential blockage below** for details how to start the cleaning programs.

7.7.3 Removing a potential blockage

Clogs or blockages in the sample tube or flow cell can cause greatly reduced event rates. Optimized sample preparation can prevent most clogs. Use a hot bleach cycle as described below to remove blockages. Contact Miltenyi Biotec Technical Support if the problem persists.

The following protocol can also be used as a weekly cleaning routine.

- 1 Go to the **Experiment** tab.
- 2 Select the Chill 5 Rack.
- 3 Highlight column 1.

- 4 Enter Hot bleach in the **Sample ID** field.
- 5 Highlight column 2.
- 6 Enter Running buffer in the **Sample ID** field.
- 7 Highlight all wells.
- 8 Set the flow rate to high, the uptake volume to 450 μL , and the sample volume to 500 μL . It is not necessary to open an analysis window.
- 9 Optional: Enter weekly clean in the **Project** field.

The screenshot shows the 'Experiment' settings window with the following configurations:

- Experiment:** Rack: Chill 5 rack; File: adm2023-05-03.0001; Project: weekly clean; Sample ID: Running buffer; Description: (empty)
- Flow rate:** Low, Med, High (selected)
- Pickup and measure:** Mix sample: Off; Mode: Standard; Uptake volume: 450 μl ; Sample volume: 500 μl

- 10 Optional: Check all settings via **View > Experiment**.

Sample ID	Description	Flow rate	Auto flow rate on	Mix sample	Mode	Uptake volume	Sample volume
A1	Hot bleach	High	No	Off	Standard	450	500
B1	Hot bleach	High	No	Off	Standard	450	500
C1	Hot bleach	High	No	Off	Standard	450	500
D1	Hot bleach	High	No	Off	Standard	450	500
A2	Running buffer	High	No	Off	Standard	450	500
B2	Running buffer	High	No	Off	Standard	450	500
C2	Running buffer	High	No	Off	Standard	450	500
D2	Running buffer	High	No	Off	Standard	450	500

- 11 Save this as an experiment file called **weekly cleaning** in the Public folder, so that the file can be opened by any user.
- 12 Open the **weekly cleaning** experiment file.
- 13 Prepare 2 mL hot bleach solution. Mix 1 mL of 1% hypochlorite solution with 1 mL boiled distilled water. Use a microwave oven or water kettle to boil the water.
- 14 Place four tubes containing 0.5 mL hot bleach solution into wells A1-D1 of the Chill 5 Rack.

Use a Chill 5 Rack stored at room temperature.

- 15 Place four tubes containing 0.5 mL MACSQuant Running Buffer into wells A2-D2.
- 16 Place the Chill 5 Rack onto the MACS MiniSampler Plus.
- 17 Click the **Start** button.

7.7.4 Cleaning the washing station

WARNING

Electric shock, short circuit, overheating, fire, and explosion could result in death or serious injury.

Liquids inside the instrument can cause short circuits.

- Unplug the instrument before cleaning it.
- Use only small amounts of cleaning agents on a soft cloth to wipe the instrument. Do not spray or pour liquid cleaning agents onto or into the instrument.

- 1 Switch off the instrument.
- 2 Remove the Single tube rack and other racks.
- 3 Open the cover of the washing station to the left side.

WARNING! Potentially contaminated liquid may spill out of the orifice of the washing station and the tubing.

- 4 Clean the washing station by wiping it with a tissue and an appropriate disinfectant, for example, 70% ethanol or isopropyl alcohol.
- 5 Rinse with distilled water.
- 6 Close the cover of the washing station.

7.7.5 Cleaning the uptake needle

Clean the uptake needle regularly to prevent contamination or clogging.

Be careful not to bend or damage the needle while cleaning.

- 1 Switch off the instrument.
- 2 Take the needle at the top and pull it out of the magnetic holder.
- 3 Clean the needle by wiping it with a tissue and an appropriate disinfectant, for example, 70% ethanol or isopropyl alcohol.
- 4 Wipe the needle holder with a tissue and an appropriate disinfectant, for example, 70% ethanol or isopropyl alcohol.
- 5 Clean the needle and the needle holder with a tissue with distilled water.
- 6 Reinstall the needle.
- 7 Calibrate the needle arm. See **Calibrating the needle arm on page 85**.

7.8 Long-term storage of the instrument

The instrument should be used every two weeks to prevent clogging. If the instrument is not in use, start the instrument every two weeks and switch it into acquisition mode. See **Priming the instrument on page 83**. During the priming procedure, the fluidics system is flushed. After the instrument has been primed, shut down the instrument. See **Manual shutdown on page 84**.

8

Troubleshooting

 **WARNING**

Biohazardous material.

Contamination or infection could result in death or serious injury depending on the material used.

- If hazardous or potentially infectious material has been spilled or leaked from the system, decontaminate the affected area.
- Do not continue to use contaminated accessories or parts of the instrument.

If you encounter an issue with the MACSQuant Instrument, see **Table 8.1** for general issues and to **Table 8.2** for issues with error messages.

Issue	Causes	Solution	Reference
Leakage	MACSQuant Column	Check MACSQuant Column.	Column leakage on page 109
	Loose syringe	Tighten the 0.5 and 5 mL syringe.	Loose connections on page 109
	Loose tubings	Check connections of buffer tubings.	
	Broken check valve	Replace check valve.	Replacing the check valve on page 108
Air in system	Loose syringe	Tighten the 0.5 and 5 mL syringe.	Loose connections on page 109
	Loose tubings	Check connections of buffer tubings.	
	Air in sheath particle filter	Release air from sheath particle filter.	Releasing air from the sheath particle filter on page 109
No events	Loose syringe	Tighten the 0.5 and 5 mL syringe.	Loose connections on page 109
	Wet air filter	Exchange the air filter.	Exchanging hydrophobic air filters on page 100
	Blockage	Clean the instrument.	Removing a potential blockage on page 102
	Incorrect needle arm calibration	Recalibrate the needle arm.	Calibrating the needle arm on page 85
Unusual pattern in dot plot	Wet air filter	Exchange the air filter.	Exchanging hydrophobic air filters on page 100
	Blockage	Clean the instrument.	Removing a potential blockage on page 102
Incomplete PMT calibration	Not enough events	Provide enough MACSQuant Advanced Calibration Beads, vortex beads thoroughly before use.	
		Repeat the PMT calibration.	
	Loose syringe	Tighten the 0.5 and 5 mL syringe.	Loose connections on page 109
Failed PMT calibration	High noise (CV) or low SI	Clean the instrument.	Cleaning recommendations on page 102
		Tighten the 0.5 and 5 mL syringe.	Loose connections on page 109
		Check the connections of the buffer tubings.	
	Release the air from the sheath particle filter.	Releasing air from the sheath particle filter on page 109	
High debris or high noise	Air in system	Tighten 0.5 and 5 mL syringe.	Loose connections on page 109
		Check connections of buffer tubings.	
		Release air from sheath particle filter.	Releasing air from the sheath particle filter on page 109
	Dirt in fluidics system	Clean the instrument.	Cleaning recommendations on page 102
Gap in event rate (HDR-T plot), uptake of air	Incorrect calibration	Recalibrate the plates.	Calibrating racks and plates on page 87
	Incorrect sample or uptake volume	Provide sufficient sample volume. Check the uptake volume. Consider excess volume.	Excess volume of the MACSQuant Analyzer 10, VYB, and Analyzer 16 on page 115

Table 8.1: Issues on the MACSQuant Instrument without error message

Error message	Causes	Solution	Reference
Plunger overload during acquisition DeviceError: 0x0309: Plunger overload	Clogging	Clean the instrument.	Cleaning recommendations on page 102
Error on MACS MiniSampler Plus Device Error: 0x0901: Position not allowed or parameter-range invalid	Obstruction	Ensure enough space around the MACS MiniSampler Plus during acquisition.	
	Wrong calibration	Restart the instrument.	
Error on uptake unit Y (movement back and forth) DeviceError: 0x0101: Position not allowed or parameter-range invalid	Obstruction	Ensure that the needle arm can move freely.	
	Wrong calibration	Calibrate the needle arm.	Calibrating the needle arm on page 85
Error on uptake unit Z (movement up and down) DeviceError: 0x0101: Position not allowed or parameter-range invalid	Obstruction	Ensure that the needle arm can move freely.	
	Dirty needle holder	Clean needle holder.	Cleaning the uptake needle on page 104
	Wrong calibration	Calibrate the needle arm.	Calibrating the needle arm on page 85
		Calibrate rack or plates.	Calibrating racks and plates on page 87
Laser overheating DeviceError: 0x005b : Laser rec. Timeout	Room conditions	Leave 15 cm space around the instrument. Do not place the instrument near radiators, heaters, or into direct sunlight. Keep the room temperature below 25°C.	
Error on valve DriverError: 0x04000000: Module init flag lost DeviceError: 0x0025: Unknown error 5	Valve blocked	Restart the instrument.	
	Wrong (homemade) buffer used	Use only MACSQuant Running buffer to avoid wearing of valves.	
Error on valve DeviceError: 0x001c: Rotate failed	Valve blocked	Restart the instrument.	
	Wrong (homemade) buffer used	Use only MACSQuant Running Buffer to avoid wearing of valve.	
System pressure control failed	Wet air filter	Exchange the air filter.	Exchanging hydrophobic air filters on page 100
	Loose syringe	Tighten the 0.5 and 5 mL syringe.	Loose connections on page 109
	Loose tubing	Check connections of buffer tubings.	
	Air in system	Tighten 0.5 and 5 mL syringe. Check connections of the buffer tubings.	Loose connections on page 109
		Release air from the sheath particle filter.	Releasing air from the sheath particle filter on page 109

Table 8.2: Error messages of the MACSQuant Instrument

8.1 Leakage

WARNING

Defective or inadequate equipment can cause a biological hazard.

Contamination or infection may lead to severe personal injury or death, depending on the material used.

- If hazardous or potentially infectious material has been spilled or leaked from the system, decontaminate the affected area.
- Do not continue to use contaminated accessories or parts of the instrument.

8.1.1 Replacing the check valve

Follow the instructions below to replace the check valve.

- 1 Switch off the instrument.
- 2 Remove the Single tube rack or remove Chill Racks from the MACS MiniSampler Plus.
- 3 Open the front cover.
- 4 Remove the MACSQuant Column from its holder. Do not disconnect the column from the connectors.
- 5 Slide the cover of the washing station to the left.
- 6 Slide the small door above the MACSQuant Column to the right. The purple part is the check valve.
- 7 Unscrew the check valve from the black Luer connectors.
- 8 Connect the new check valve to the Luer connectors. Fasten the connections carefully but tightly. Proceed with step 11.
- 9 Optional: It is possible to use the MACSQuant Instrument short term without a check valve. Proceed with step 10. Otherwise continue with step 11.
- 10 Connect the black Luer connectors without a check valve.

If the instrument is used without a check valve, make sure that the waste bottle level is lower than the MACSQuant Running Buffer bottle level when the instrument is switched off or in data analysis mode.

Note that the absolute cell count might be affected if the instrument is used without a check valve .

- 11 Start and prime the instrument.
 - 12 Check if the leakage has stopped during priming. If the instrument is used without a check valve and the leakage persists, the Luer connectors might be broken. Contact Miltenyi Biotec Technical Support.
 - 13 Close the door above the MACSQuant Column and the cover of the washing station.
 - 14 Put the MACSQuant Column back into its holder.
 - 15 Close the front cover.
 - 16 Attach the Single tube rack or place a Chill Rack on the MACS MiniSampler Plus.
- If the problem persists, contact Miltenyi Biotec Technical Support.

8.1.2 Syringe leakage

WARNING

Defective or inadequate equipment can cause a biological hazard.

Contamination or infection may lead to severe personal injury or death, depending on the material used.

- If hazardous or potentially infectious material has been spilled or leaked from the system, decontaminate the affected area.
- Do not continue to use contaminated accessories or parts of the instrument.

A loose syringe can lead to a leakage. See **Loose connections below**.

8.1.3 Column leakage

If a recently exchanged MACSQuant Column shows signs of leakage, do the following:

- 1 Check if the column is installed properly. The column should be inserted precisely into the column connector and fastened to the point of resistance. If the column is not installed correctly, loosen the column connector, insert the column precisely, and tighten the connector again. See **Installing the MACSQuant Column on page 74**.
- 2 Check if the column is broken or cracked.
- 3 Check if the leakage persists. If so, unscrew the column and check if the connectors of the columns are damaged. If necessary, exchange the leaking column with a new MACSQuant Column.
- 4 Check if the column connector is fastened properly.

If the problem persists, contact Miltenyi Biotec Technical Support.

8.2 Air in system

8.2.1 Loose connections

If connections are loose, air can enter the fluidics system. Follow the instructions below to check syringes and buffer tubing connections. Use your fingers to tighten the connections. Do not use any tools.

- 1 Open the front cover.
- 2 Tighten the glass syringes. Turn the metal part on top of the syringe counterclockwise.
- 3 Tighten the tubing connections on the buffer bottles clockwise.
- 4 Check the tubing connections coming from the bottles to the rear side of the instrument. Tighten them clockwise.

If the problem persists, contact Miltenyi Biotec Technical Support.

8.2.2 Releasing air from the sheath particle filter

Check the fluid level of the sheath particle filter regularly. If air enters the fluidics system or fluid evaporates, air can be trapped in the sheath particle filter. Follow the instructions below to release air from the sheath particle filter.

- 1 Prime the instrument. Release of air is only possible with a pressurized system.
- 2 Get a paper towel to collect liquid.

- 3** Carefully unscrew the white cap of the sheath particle filter until air has been fully released. Do not remove the cap entirely.
 - 4** Close the cap of the sheath particle filter when liquid comes out. Collect the liquid with a paper towel.
- If the problem persists, contact Miltenyi Biotec Technical Support.

9

Technical data and specifications

The MACSQuant Instrument is classified as protection class I device and must be plugged into a grounded electrical outlet.

The MACS MiniSampler Plus is labeled as a protection class III device and must be connected to the main instrument.

Note that for USA and Canada, the main power supply cord and plug of the instrument shall comply with the following specifications: UL listed and KAM cord, minimum type SJ, minimum 18 AWG, 3 conductors. Rated for a minimum temperature of 60 °C. Provided with grounding-type (NEMA 5-15P) attachment plug, rated 125 VAC, 10 A. Opposite end terminates in IEC 320-style connector, rated 125 VAC, 10A.

Supply voltage fluctuations up to +/- 10% of the nominal voltage. The instrument is suitable for pollution degree 2 environments.

For other safety considerations, refer to the product label or visit www.miltenyibiotec.com.

Design and specifications are subject to change without notice.

9.1 Technical data of the MACSQuant Instrument

MACSQuant Instrument	
Model	MACSQuant Analyzer 10 (# 130-096-343) MACSQuant Analyzer 16 (# 130-109-803) MACSQuant VYB (# 130-096-116)
Dimensions (width × depth)	669 mm × 400 mm (26.34 × 15.75 in)
Dimensions (width × depth) with MACS MiniSampler Plus	669 mm × 500 mm (26.34 × 19.69 in)
Footprint (width × depth)	385.5 mm × 284.5 mm (15.18 × 11.2 in)
Height (adjustable touchscreen)	394-553.5 mm (15.51–21.79 in)
Weight	50 kg (110 lbs)
Power requirements	100-240 V~, 50/60 Hz
Power consumption (standard use)	450 W
Fuses	2 × 8 AT, 250 V
Rack detection	850 nm, 3.3 mW, pulse time 215 μs
RS232 Interface	labeled "RS232/AUX" - not in use
AC Output	Labeled "Bottle Sensor"

MACSQuant Instrument

RS232 Interface + DC-Output	Labeled "RS232/BCR"
CAN Bus + DC-Output	Labeled "External CAN"
CAN Bus + DC-Output	Labeled "Power Can"
Monitor	15,6" LCD touchscreen
Working temperature	+20 °C to +25 °C
Working environment	Indoor use only
Humidity	0% to 85% relative humidity, non-condensing
Altitude	Max. 2,000 m
Emission sound pressure level at workstation	<61 dB(A)

Table 9.1: Technical data of the MACSQuant Instrument

PC

RS232 Interface	Labeled "COM 1"
USB ports	4× USB 3.0 ports 2× USB 3.0 ports (at touchscreen)
Display Port	Labeled D Port (2×)
Ethernet ports (GbE)	Labeled LAN (2×)
CFast-card slot	Inserted card contains Miltenyi Biotec's proprietary MBCore operating system, based on Linux®
RAM	2× 16 GB DDR4 (SO-DIMM)
Mass storage	1 TB SSD M.2
Mainboard battery	CR2032, button cell, lithium metal

Table 9.2: Technical data of the integrated PC

Optics

Model	Laser excitation
MACSQuant Analyzer 10	405 nm, 40 mW
	488 nm, 30 mW
	640 nm, 20 mW
MACSQuant Analyzer 16	405 nm, 63 mW
	488 nm, 50 mW
	640 nm, 70 mW
MACSQuant Analyzer VYB	405 nm, 40 mW
	488 nm, 50 mW
	561 nm, 100 mW
Flow cell dimensions	200 × 250 μm

Table 9.3: Optical parameters of the MACSQuant Instrument

Fluidics

Minimum uptake volume	1 μL (25μL recommended for volumetric applications)
Sample flow rate	25, 50, or 100 μL/min plus automated flow rate to maintain 500, 1000, or 2000 events/second
Sheath fluid consumption	2-10 mL/min
Maximal event rate	up to 15,000 events/second
Absolute cell count performance	accuracy +/- 7%; reproducibility (CV) <5%
Sample carryover	0.01%

Table 9.4: Fluidic parameters of the MACSQuant Instrument

Instrument variants

Model	Channel	Laser [nm]	Filter [nm]
MACSQuant Analyzer 10	FSC	488	488/10
	SSC	488	488/10
	V1	405	450/50
	V2	405	525/50
	B1	488	525/50
	B2	488	585/40
	B3	488	655-730
	B4	488	750 LP
	R1	640	655-730
R2	640	750 LP	
MACSQuant Analyzer 16	FSC	488	488/10
	SSC	405	405/10
	V1	405	450/50
	V2	405	525/50
	V3	405	579/34
	V4	405	615/20
	V5	405	667/30
	B1	488	525/50
	B2	488	579/34
	B3	488	615/20
	B4	488	667/30
	B5	488	725/40
	B6	488	785/62
	R1	640	667/30
R2	640	725/40	
R3	640	785/62	
MACSQuant VYB	FSC	561	561/4
	SSC	561	561/4
	V1	405	452/45
	V2	405	525/50
	B1	488	525/50
	B2	488	593-650
	Y1	561	586/15
	Y2	561	615/20
	Y3	561	661/20
Y4	561	740 LP	

Table 9.5: Instrument variants

9.2 Excess volume of the MACSQuant Analyzer 10, VYB, and Analyzer 16

The MACSQuant Analyzer 10, VYB, and Analyzer 16 pick up some additional volume (excess volume) depending on the selected **Mode** in the **Experiment** tab. Make sure to provide enough sample volume as the excess volume is additionally picked up to the programmed uptake volume.

Excess volume	
Mode	Additionally picked up volume (excess volume)
Screen	19.9 μ L
Fast	9.9 μ L
Standard	9.9 μ L
Extended	9.9 μ L
Enrich	9.9 μ L
EnrichS	9.9 μ L
EnrichS2	9.9 μ L

Table 9.6: Excess volume of the MACSQuant Analyzer 10, VYB, and Analyzer 16

9.3 FCC compliance statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications of the instrument, unless expressly approved by Miltenyi Biotec, may void your authority to operate the instrument pursuant to FCC 47 CFR.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

9.4 EC/EU Declaration of Conformity

This declaration applies to the MACSQuant Analyzer 16, MACSQuant Analyzer 10, and MACSQuant VYB. The legislative requirements are equivalent for all three MACSQuant Instruments.

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

This declaration relates exclusively to the machinery in the state in which it was placed on the market, and excludes components which are added and/or operations carried out subsequently by the final user.

The declaration of conformity refers to the machinery identified as follows:

Description: Laboratory equipment
Model: MACSQuant Analyzer 10, MACSQuant VYB and MACSQuant Analyzer 16

The machinery complies with all essential requirements of the following directives:

2006/42/EC Machinery
2011/65/EU Restriction of the use of certain hazardous substances in electrical & electronic equipment
2014/30/EU Electromagnetic compatibility

The machinery is in conformity with the following harmonized standards:

EN 60825-1:2014
EN 61010-1:2010
EN 61010-2-081:2015
EN 61326-1:2013

Person authorized to compile the relevant technical documentation:

Dr. Bernd Schröder
Global Head Regulatory Affairs
Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

Additional international conformity (MACSQuant Analyzer 10, MACSQuant VYB, and MACSQuant Analyzer 16)

The MACSQuant Instruments are in conformity with:

- IEC 61010-1
- IEC 61010-2-081
- UL 61010-1
- UL 61010-2-081
- CAN/CSA-C22.2 No. 61010-1

- CAN/CSA-C22.2 No. 61010-2-081
- IEC 60825-1
- CAN/CSA-E60825-1

9.5 UK Declaration of Conformity

This declaration applies to the MACSQuant Analyzer 16, MACSQuant Analyzer 10, and MACSQuant VYB. The legislative requirements are equivalent for all three MACSQuant Instruments.

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Miltenyi Biotec B.V. & Co. KG
 Friedrich-Ebert-Straße 68
 51429 Bergisch Gladbach
 Germany

This declaration relates exclusively to the machinery in the state in which it was placed on the market, and excludes components which are added and/or operations carried out subsequently by the final user.

The declaration of conformity refers to the machinery identified as follows:

Description: Laboratory equipment
 Model: MACSQuant Analyzer 10, MACSQuant VYB, and MACSQuant Analyzer 16

The machinery complies with all essential requirements of the following legislations:

Supply of machinery (Safety) regulations 2008
 The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012
 Electromagnetic Compatibility Regulations 2016

The machinery is in conformity with the following UK designated standards:

EN 60825-1:2014
 EN 61010-1:2010
 EN 61010-2-081:2015
 EN 61326-1:2013

Person authorized to compile the relevant technical documentation:

Dr. Bernd Schröder
 Global Head Regulatory Affairs
 Miltenyi Biotec B.V. & Co. KG
 Friedrich-Ebert-Straße 68
 51429 Bergisch Gladbach
 Germany

Additional international conformity (MACSQuant Analyzer 10, MACSQuant VYB, and MACSQuant Analyzer 16)

The MACSQuant Instruments are in conformity with:

- IEC 61010-1
- IEC 61010-2-081
- UL 61010-1
- UL 61010-2-081
- CAN/CSA-C22.2 No. 61010-1
- CAN/CSA-C22.2 No. 61010-2-081
- IEC 60825-1
- CAN/CSA-E60825-1

9.6 Technical data and specifications of the MACS MiniSampler Plus

The MACS MiniSampler Plus is labeled as a protection class III device and must be plugged into the connector of the Instrument labeled with **External CAN**.

For other safety considerations, refer to the product label, or visit www.miltenyibiotec.com.

Design and specifications are subject to change without notice.

MACS MiniSampler Plus	
Size without lid	182 mm × 148 mm × 47 mm
Size with lid	280 mm × 153 mm × 172 mm
Weight	1.5 kg
Input voltage	24 VDC
Power	19.2 W
Sub D9 interface with shielding	Pins 1, 4 and 8: NC Pin 2: CAN-L Pins 3 and 6: GND Pins 5 and 9: +24 VDC/2A Pin 7: CAN-H
Working temperature	+ 20 °C to + 25 °C
Working environment	Indoor use only
Humidity	0% to 85% relative humidity, non-condensing
Altitude	Max. 2000 m

Table 9.7: Technical data of the MACS MiniSampler Plus

9.7 EC/EU Declaration of Conformity

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

This declaration relates exclusively to the machinery in the state in which it was placed on the market, and excludes components which are added and/or operations carried out subsequently by the final user.

The declaration of conformity refers to the machinery identified as follows:

Description: Laboratory equipment
Model: MACS MiniSampler Plus

The machinery complies with all essential requirements of the following directives:

2006/42/EC Machinery
2011/65/EU Restriction of the use of certain hazardous substances in electrical & electronic equipment
2014/30/EU Electromagnetic compatibility

The machinery is in conformity with the following harmonized standards:

EN 61010-1:2010
EN 61010-2-081:2015
EN 61326-1:2013

Person authorized to compile the relevant technical documentation:

Dr. Bernd Schröder
Global Head Regulatory Affairs
Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

9.8 UK Declaration of Conformity

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

This declaration relates exclusively to the machinery in the state in which it was placed on the market, and excludes components which are added and/or operations carried out subsequently by the final user.

The machinery complies with all essential requirements of the following legislations:

Description: Laboratory equipment
Model: MACS MiniSampler Plus

The machinery complies with all essential requirements of the following legislations:

Supply of machinery (Safety) regulations 2008
The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012
Electromagnetic Compatibility Regulations 2016

The machinery is in conformity with the following UK designated standards:

EN 61010-1:2010
EN 61010-2-081:2015
EN 61326-1:2013

Person authorized to compile the relevant technical documentation:

Dr. Bernd Schröder
Global Head Regulatory Affairs
Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

10

Technical support

For technical support, contact your local Miltenyi Biotec representative or Miltenyi Biotec Technical Support at Miltenyi Biotec headquarters:

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone: +49 2204 8306 3803
technicalsupport@miltenyi.com

Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

11

Legal notes

11.1 Limited warranty

Except as stated in a specific warranty statement which may accompany this product or as otherwise agreed in writing by an authorized representative of Miltenyi Biotec, Miltenyi Biotec's warranty to you, the original purchaser and end user ("you" or "your"), with respect to the product accompanied by this limited warranty shall be subject to the following provisions and the general terms and conditions of sale of the company within the Miltenyi Biotec group which supplied the product in effect at the date of purchase. Those terms and conditions of sale may vary by country and region. Nothing in this document should be construed as constituting an additional warranty.

Miltenyi Biotec warrants that this product will operate or perform substantially in conformance with Miltenyi Biotec's published specifications and be free from material defects in material and workmanship, when subjected to normal, proper, and intended usage by properly trained personnel, for the period of time set forth in the product documentation or package inserts accompanying the product (the "Warranty Period").

Miltenyi Biotec agrees, during the Warranty Period, to repair or replace, at Miltenyi Biotec's option, the defective product so as to cause the same to operate in substantial conformance with said published specifications; provided that you shall (a) promptly notify Miltenyi Biotec in writing upon the discovery of any nonconformity or defect, which notice shall include the product model and serial number (if applicable) and details of the warranty claim; and (b) return the nonconforming or defective product to Miltenyi Biotec, freight prepaid, only after receipt of a Return Material Authorization ("RMA") from Miltenyi Biotec, which may include biohazard decontamination procedures and other product-specific handling instructions, if applicable.

Miltenyi Biotec shall have no obligation to make repairs, replacements, or corrections to the product or any component thereof required, in whole or in part, as the result of (i) normal wear and tear, (ii) improper handling, installation, operation, storage, service, maintenance, or repair, (iii) failure to follow the instructions, cautions, warnings, and notes set forth in the product documentation provided with the product or provided by Miltenyi Biotec from time to time, (iv) abnormal use, misuse, neglect, abuse, mishandling, misapplication, modification, or alteration of the product, (v) use of the product in a manner for which it was not designed, (vi) causes external to the product such as, but not limited to, power failure or electrical power surges, (vii) use of the product in combination with equipment, accessories, consumables, or software not supplied or approved by Miltenyi Biotec, or (viii) accident, disaster, or acts of God. ANY INSTALLATION, MAINTENANCE, REPAIR, SERVICE, OR ALTERATION TO OR OF, OR OTHER TAMPERING WITH, THE PRODUCT PERFORMED BY ANY PERSON OR ENTITY OTHER THAN MILTENYI BIOTEC AUTHORIZED PERSONNEL WITHOUT MILTENYI BIOTEC'S PRIOR WRITTEN APPROVAL, OR ANY USE OF REPLACEMENT PARTS NOT SUPPLIED BY MILTENYI BIOTEC, SHALL IMMEDIATELY VOID AND CANCEL ALL WARRANTIES WITH RESPECT TO THE AFFECTED PRODUCT.

Miltenyi Biotec's warranty does not cover products sold AS IS or WITH ALL FAULTS, or which had its serial number defaced, altered, or removed, or any consumables or parts identified as being supplied by a third party.

Miltenyi Biotec must be informed promptly if a claim is made under this limited warranty. If a material or manufacturing defect occurs within the Warranty Period, Miltenyi Biotec will take the appropriate steps, at Miltenyi Biotec's option, to make repairs, replacements, or corrections to the product or any component thereof as may be required to restore the full usability of your product. If Miltenyi Biotec determines that a product for which you have requested warranty services is not covered by the warranty hereunder, you shall pay or reimburse Miltenyi Biotec for all costs of investigating and responding to such request at Miltenyi Biotec's then prevailing time and materials rates. If Miltenyi Biotec provides repair services or replacement parts that are not covered by this warranty, you shall pay Miltenyi Biotec therefor at Miltenyi Biotec's then prevailing time and materials rates.

THE OBLIGATIONS CREATED BY THIS WARRANTY STATEMENT TO REPAIR OR REPLACE A DEFECTIVE PRODUCT ARE EXCLUSIVE AND SHALL BE THE SOLE REMEDY OF BUYER IN THE EVENT OF A DEFECTIVE PRODUCT. EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY STATEMENT, MILTENYI BIOTEC HEREBY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, ORAL OR WRITTEN, WITH RESPECT TO THE PRODUCT, INCLUDING WITHOUT LIMITATION ALL IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. MILTENYI BIOTEC DOES NOT WARRANT THAT THE PRODUCT IS ERROR-FREE OR WILL ACCOMPLISH ANY PARTICULAR RESULT.

UNDER NO CIRCUMSTANCES WILL MILTENYI BIOTEC BE LIABLE FOR ANY LOSS OF USE, INTERRUPTION OF BUSINESS, OR ANY INDIRECT, SPECIAL, INCIDENTAL, PUNITIVE, OR CONSEQUENTIAL DAMAGES OF ANY KIND (INCLUDING LOST PROFITS) REGARDLESS OF THE FORM OF ACTION WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE), STRICT PRODUCT LIABILITY, OR OTHERWISE, EVEN IF MILTENYI BIOTEC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. IN NO EVENT SHALL THE TOTAL LIABILITY OF MILTENYI BIOTEC HEREUNDER EXCEED THE GREATER OF 100.00 EUROS OR THE AMOUNT YOU ACTUALLY PAID FOR THE PRODUCT GIVING RISE TO SUCH LIABILITY, REGARDLESS OF THE CAUSE OF ACTION, IN CONTRACT, TORT, STRICT LIABILITY, OR OTHERWISE.

NOT ALL JURISDICTIONS ALLOW SUCH LIMITATIONS OF DAMAGES SO THE FOREGOING LIMITATIONS MAY NOT APPLY TO YOU. This warranty statement gives you specific legal rights and you may have other rights, which may vary from jurisdiction to jurisdiction.

11.2 Trademarks

autoMACS, MACS, MACSima (MACSQuant/MACSima Storage Solution), MACSQuant, MACSQuantify, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only.

12

Index

A

Accessories **71**

Acquisition mode **83**

Air in system **109**

 Sheath particle filter **109**

B

Blockage

 Remove **102**

Bottle illumination **95**

C

Calibration

 Chill 5 Rack **87**

 Chill 96 Rack **88**

 Custom racks **89**

 MACSPlex Filter Plate **88**

 Needle arm **85**

 Rack detection **90**

 Racks and plates **87**

 Single tube rack **87**

 Single tube rack detection **90**

 Universal Reagent Rack **89**

Cleaning

 Hot bleach **102**

 Instrument **101**

 Programs **101**

 Uptake needle **104**

 Washing station **104**

 Weekly cleaning routine **102**

Cleaning programs **101**

Cleaning recommendations **102**

Column leakage **109**

Complementary components **71**

Consumables **70**

D

Decontamination **99**

Disinfection

 Fluidics system **99**

 Surface **99**

E

Error messages **107**

Exchange

 Buffer and solution bottles **99**

 Hydrophobic air filter **100**

 MACSQuant Columns **101**

 Waste bottle **100**

H

Hardware monitor **93**

Lasers and detectors **94**

Hot bleach **102**

Hydrophobic air filter

Exchange **100**

Install **73**

Location **73**

I

Install

Fluid bottles **71**

Hydrophobic air filter **73**

MACS MiniSampler Plus **76**

MACSQuant Column **74**

Needle guard **76**

Single tube rack **73**

Universal Reagent Rack **79**

Webcam **81**

Instrument monitoring

Lasers and detectors **94**

LED colors **95**

Intended use **65**

K

Key features **65**

L

Laser excitation **113**

Lasers and detectors **94**

Leakage **108**

Column **109**

Syringe **109**

Long-term storage **104**

Loose connections **109**

M

MACS MiniSampler Plus

Install **76**

MACSQuant Column

Exchange **101**

Install **74**

Maintenance

Cleaning **101**

Decontamination **99**

Exchange buffer bottles **99**

N

Needle guard

Install **76**

No events **102**

O

Optics **113**

P

Priming the instrument **83**

R

Remote support **97**

Remove

Blockage **102**

Replace

Check valve **108**

S

Service contracts **71**

Sheath particle filter

Release air **109**

Shutdown

Automatic **84**

Manual **84**

Single tube rack

Install **73**

Software

Start **83**

Startup **83**

Surface disinfection **99**

Switch off **83**

Switch on **83**

Syringe pump leakage **109**

T

Troubleshooting **106-107**

U

Uptake needle

Cleaning **104**

V

View

Hardware monitor **93**

W

Waste bottle

Exchange **100**

Webcam **81**

Weekly cleaning **102**

Miltenyi Biotec

Germany/Austria
Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macsde@miltenyi.com

USA/Canada
Miltenyi Biotec, Inc.
1201 Clopper Road
Gaithersburg, MD 20878
USA
Phone 800 FOR MACS
Phone +1 866 811 4466
Fax +1 877 591 1060
macsus@miltenyi.com

Australia
Miltenyi Biotec Australia Pty. Ltd.
Unit 11, 2 Eden Park Drive
Macquarie Park NSW 2113
Australia
Phone +61 2 8877 7400
Fax +61 2 9889 5044
macsau@miltenyi.com

Benelux
Miltenyi Biotec B.V.
Dellaertweg 9C
2316 WZ Leiden
The Netherlands
macsnl@miltenyi.com
**Customer service for:
The Netherlands**
Phone 0800 4020120
Fax 0800 4020100

Belgium
Phone 0800 94016
Fax 0800 99626

Luxembourg
Phone 800 24971
Fax 800 24984

China
Miltenyi Biotec Technology &
Trading (Shanghai) Co., Ltd.
Room A401, 4/F
No. 1077, Zhangheng Road
Pudong New Area
201203 Shanghai, P.R. China
Phone +86 21 6084 0210
Fax +86 21 6235 0953
macscn@miltenyi.com.cn

France
Miltenyi Biotec SAS
10 rue Mercoeur
75011 Paris, France
Phone +33 1 56 98 16 16
macsfr@miltenyi.com

Hong Kong
Miltenyi Biotec Hong Kong Ltd.
Unit 301, Lakeside 1
No. 8 Science Park West Avenue
Hong Kong Science Park
Pak Shek Kok, New Territories
Hong Kong
Phone +852 3751 6698
Fax +852 3619 5772
macshk@miltenyi.com.hk

India
Miltenyi Biotec India Pvt. Ltd.
Vatika Business Centre,
Floor No. 6
Divya Sree Omega
Kondapur, Serilingampally
K.V. Rangareddy
Telangana 500084
India
Phone +91 040 45175910
macsin@miltenyi.com

Italy
Miltenyi Biotec S.r.l.
Via Paolo Nanni Costa, 30
40133 Bologna, Italy
Phone +39 051 6 460 411
Fax +39 051 6 460 499
macsit@miltenyi.com

Japan
Miltenyi Biotec K.K.
NEX-Eitai Building 5F
16-10 Fuyuki, Koto-ku,
Tokyo 135-0041, Japan
Phone +81 3 5646 8910
Fax +81 3 5646 8911
macsjp@miltenyi.com

Nordics and Baltics
Miltenyi Biotec Norden AB
Medicon Village
Scheeletorget 1
223 81 Lund, Sweden
macsse@miltenyi.com
Customer service for:

Sweden
Phone 0200 111 800
Fax +46 280 72 99

Denmark
Phone 80 20 30 10
Fax +46 46 280 72 99

**Norway, Finland, Iceland,
and Baltic countries**
Phone +46 46 280 72 80
Fax +46 46 280 72 99

Singapore
Miltenyi Biotec Asia Pacific Pte Ltd.
438B Alexandra Road, Block B
Alexandra Technopark
#06-01
Singapore 119968
Phone +65 6238 8183
Fax +65 6238 0302
macssg@miltenyi.com

South Korea
Miltenyi Biotec Korea Co., Ltd.
Donggeuk 7F,
562 Nonhyeon-ro
Gangnam-gu
Seoul 06136, South Korea
Phone +82 2 555 1988
Fax +82 2 555 8890
macskr@miltenyi.com

Spain
Miltenyi Biotec S.L.
C/Virgilio 2, Edificio II, Planta -1
28223 Pozuelo de Alarcón, Madrid
Spain
Phone +34 91 512 12 90
Fax +34 91 512 12 91
macses@miltenyi.com

Switzerland
Miltenyi Biotec Swiss AG
Soodstrasse 52
8134 Adliswil
Switzerland
Phone +41 32 623 08 47
Fax +49 2204 85197
macsch@miltenyi.com

United Kingdom
Miltenyi Biotec Ltd.
Almac House, Church Lane
Bisley, Surrey GU24 9DR, UK
Phone +44 1483 799 800
Fax +44 1483 799 811
macsuk@miltenyi.com

www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. autoMACS, MACS, MACSima (MACSQuant/MACSima Storage Solution), MACSQuant, MACSQuantify, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide.

Copyright © 2025 Miltenyi Biotec and/or its affiliates. All rights reserved.