

Contents

1. Description
 - 1.1 Background information
 - 1.2 Applications
 - 1.3 Reagent and instrument requirements
2. Protocols
 - 2.1 Immunofluorescent staining of nucleated cells, e.g., PBMCs
 - 2.2 Immunofluorescent staining and lysis of whole blood or leukapheresis products (lyse/no wash)
 - 2.3 Immunofluorescent staining and lysis of whole blood or leukapheresis products (lyse/wash)
 - 2.4 Flow cytometric data acquisition with the MACSQuant® Analyzer 10 using an Express Mode
3. Examples of immunofluorescent staining with the StainExpress Immune Cell Composition Cocktail

1. Description

This product is for research use only.

Components 24 tubes **StainExpress Immune Cell Composition Cocktail, human:**

Dried cocktail of fluorochrome-conjugated recombinant engineered REAfinity™ Antibodies (isotype: recombinant human IgG1) containing:

CD45-VioBlue® (clone: REA747),
 CD4-VioGreen™ (clone: REA623),
 CD3-FITC (clone: REA613),
 CD56-PE (clone: REA196),
 CD16-PE (clone: REA423),
 CD19-PE-Vio® 770 (clone: REA675),
 CD14-APC (clone: REA599),
 CD8-APC-Vio 770 (clone: REA734),
 7-AAD Staining Solution.

2 StainExpress Immune Cell Composition Compensation Sets, human:

Each set contains seven dried single antibody tubes of fluorochrome-conjugated recombinant engineered REAfinity Antibodies (isotype: recombinant human IgG1) for compensation controls:

CD45-VioBlue (clone: REA747),
 CD4-VioGreen (clone: REA623),
 CD3-FITC (clone: REA613),
 CD16-PE (clone: REA423),
 CD19-PE-Vio 770 (clone: REA675),

CD14-APC (clone: REA599),
 CD8-APC-Vio 770 (clone: REA734).

- Capacity** 24 tests, one test for up to 10⁶ total cells.
- Product format** Antibodies are supplied in a dry format containing stabilizer.
- Storage** Store at dry conditions in a closed pouch. Store protected from light at 19–25 °C. The expiration date is indicated on the pouch label.

1.1 Background information

The StainExpress Immune Cell Composition Cocktail simplifies flow cytometric evaluation of cell fractions for immunofluorescent staining of whole blood, leukapheresis products, peripheral mononuclear cells (PBMCs), or other single cell suspensions from human tissue.

The cocktail has been designed for the reliable identification of human monocytes, neutrophils, eosinophils, and T, B, and NK lymphocyte populations as well as CD4⁺, CD8⁺, and NKT cell subsets in human blood.

For flow cytometric analysis use a flow cytometer equipped with a red (638 nm), a blue (488 nm), and a violet (405 nm) laser, for example, the MACSQuant Analyzer 16.

1.2 Applications

- Evaluation of leukocyte subsets in whole blood, leukapheresis products, PBMCs, or other single cell suspensions from human tissue.
- Determination of cell counts and population subsets during CAR T cell manufacturing procedures using the MACSQuant Analyzer.

1.3 Reagent and instrument requirements

- **Buffer:** Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS® BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (2–8 °C).

▲ **Note:** EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). BSA can be replaced by other proteins such as human serum albumin, human serum, or fetal bovine serum (FBS). Buffers or media containing Ca²⁺ or Mg²⁺ are not recommended for use.

- Flow cytometer, e.g., MACSQuant Analyzer 16 (# 130-109-803) or MACSQuant Analyzer 10 (# 130-096-343)
- ▲ **Note:** The MACSQuant VYB cannot be used.
- (Optional) MACS Comp Bead Kit, anti-REA (# 130-104-693), for compensation of the fluorescence spillover from fluorochrome-conjugated antibodies.
- (Optional) 10× Red Blood Cell Lysis Solution (# 130-094-183)
- (Optional) MACS MiniSampler Plus (# 130-105-745)

- (Optional) Chill 5 Rack (# 130-092-951)
- (Optional) CAR T Cell Express Mode Package (# 160-002-376)

2. Protocols

2.1 Immunofluorescent staining of nucleated cells, e.g., PBMCs

▲ Volumes given below are for up to 10^6 nucleated cells. When working with fewer than 10^6 cells, use the same volumes as indicated. When working with higher cell numbers, use multiple tubes accordingly (e.g. for 2×10^6 nucleated cells, use two tubes).

1. Determine cell number.
2. Adjust cell concentration to up to 10^6 nucleated cells per 100 μL using buffer.

▲ **Note:** If necessary, centrifuge cell suspension at $300 \times g$ for 5 minutes, aspirate supernatant completely, and resuspend up to 10^6 nucleated cells per 100 μL of buffer.
3. Add 100 μL of cell suspension to one tube of StainExpress Immune Cell Composition Cocktail.
4. Incubate for 10–20 seconds, then vortex at high speed for 5 seconds and incubate for 10 minutes in the dark at 19–25 °C.

▲ **Note:** Working at lower temperatures requires increased incubation times.
5. (Optional) Wash cells by adding 1 mL of buffer and centrifuge at $300 \times g$ for 5 minutes. Aspirate supernatant completely.
6. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry.

▲ **Note:** Store samples at 2–8 °C protected from light until analysis. Acquire samples within 1 hour after staining.
7. Proceed to flow cytometric analysis.

2.2 Immunofluorescent staining and lysis of whole blood or leukapheresis products (lyse/no wash)

1. Dilute 10 \times Red Blood Cell Lysis Solution 1:10 with double-distilled water (ddH₂O), for example, dilute 1 mL of 10 \times Red Blood Cell Lysis Solution with 9 mL of ddH₂O.

▲ **Note:** Do not dilute with deionized water. Store prepared 1 \times Red Blood Cell Lysis Solution at room temperature. Discard unused solution at the end of the day.
2. Add 100 μL of whole blood or leukapheresis product to one StainExpress Immune Cell Composition Cocktail tube.
3. Incubate for 10–20 seconds, then vortex at high speed for 5 seconds and incubate for 10 minutes in the dark at 19–25 °C.

▲ **Note:** Working at lower temperatures requires increased incubation times.
4. Add 1900 μL of 1 \times Red Blood Cell Lysis Solution and immediately vortex thoroughly for 3 seconds. Incubate for 15 minutes in the dark at 19–25 °C.

▲ **Note:** Store samples at 2–8 °C protected from light until analysis. Acquire samples within 1 hour after staining.
5. Proceed to flow cytometric analysis.

2.3 Immunofluorescent staining and lysis of whole blood or leukapheresis products (lyse/wash)

1. Dilute 10 \times Red Blood Cell Lysis Solution 1:10 with double-distilled water (ddH₂O), for example, dilute 1 mL of 10 \times Red Blood Cell Lysis Solution with 9 mL of ddH₂O.

▲ **Note:** Do not dilute with deionized water. Store prepared 1 \times Red Blood Cell Lysis Solution at room temperature. Discard unused solution at the end of the day.
2. Add 100 μL of whole blood or leukapheresis product to one StainExpress Immune Cell Composition Cocktail tube.
3. Incubate for 10–20 seconds, then vortex at high speed for 5 seconds and incubate for 10 minutes in the dark at 19–25 °C.

▲ **Note:** Working at lower temperatures requires increased incubation times.
4. Add 1900 μL of 1 \times Red Blood Cell Lysis Solution and immediately vortex thoroughly for 3 seconds. Incubate for 15 minutes in the dark at 19–25 °C.
5. Centrifuge at $300 \times g$ for 5 minutes. Aspirate supernatant completely.
6. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry.

▲ **Note:** Store samples at 2–8 °C protected from light until analysis. Acquire samples within 1 hour after staining.
7. Proceed to flow cytometric analysis.

2.4 Flow cytometric data acquisition with the MACSQuant Analyzer 16 using an Express Mode

▲ Please refer to the MACSQuant Instrument user manual and software guide for detailed information on using the MACSQuant Analyzer.

▲ Please refer to the data sheet of the MACS Comp Bead Kit, anti-REA, when using beads for compensation.

1. Prepare and prime the MACSQuant Analyzer. Make sure the calibration and instrument settings of the instrument have been optimized for acquisition of the StainExpress Immune Cell Composition Cocktail.
2. For optimal compensation, prepare single stainings of suitable beads or cells with all single antibody tubes from one pouch of the matching StainExpress Compensation Set. If using beads, add one full drop of MACS Comp Beads – anti-REA and one full drop of MACS Comp Beads – blank directly to each compensation tube.

▲ **Note:** For optimal compensation always use the provided StainExpress Compensation Set matching the StainExpress Cocktails used. To ensure they match, check for identical order number-related lot numbers indicated on the pouches next to the expiration date.

▲ **Note:** One full drop of beads is approximately 50 μL .

▲ **Note:** For compensation on cells add 100 μL cell suspension instead of beads directly to each compensation tube.
3. Incubate for 10–20 seconds, then vortex at high speed for 5 seconds and incubate for 10 minutes in the dark at 19–25 °C.


▲ **Note:** Working at lower temperatures requires increased incubation times.
4. Dilute each sample by adding 600 μL of MACSQuant Running Buffer. Mix well.

5. If using the MACS Comp Bead Kit, anti-REA, follow the protocol for compensation set up of the MACSQuant Analyzer.
▲ **Note:** For automated compensation choose the “CompensationMultiColor” Express Mode.
6. Choose appropriate voltage settings for forward scatter (FSC) and side scatter (SSC).
7. Define an appropriate threshold, based on FSC versus SSC, for the exclusion of debris from the data acquisition.
8. Select the corresponding Express Mode.
9. Start flow cytometric data acquisition.


3. Examples of immunofluorescent staining with the StainExpress Immune Cell Composition Cocktail

Whole blood from a healthy donor was stained with the StainExpress Immune Cell Composition Cocktail. Staining was carried out for 10 minutes at room temperature (19–25 °C). Subsequently, red blood cells were lysed by incubation with 1× Red Blood Cell Lysis Solution at room temperature for 15 minutes. Cells were analyzed by flow cytometry using the MACSQuant Analyzer 10.

To exclude debris, a gate was set on FSC versus SSC encompassing all cells (A). To exclude residual erythrocytes and to identify leukocytes, CD45 was used to gate on CD45⁺ leukocytes (B), and dead cells were excluded by 7-AAD (C). CD3⁺ cells were identified (D) and discriminated into T cells and NKT cells based on expression of CD56 (E). T cells were further divided into CD4⁺ and CD8⁺ T cells (F). Among CD3⁻ cells, monocytes were defined by CD14 expression and B cells by CD19 expression (G). The remaining CD14⁻/CD19⁻ cells were further divided into SSC^{high}/CD16⁻ eosinophils, SSC^{high}/CD16⁺ neutrophils, and SSC^{low}/CD56⁺/CD16⁺ cells (H).



Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.



Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, MACSQuant, the Miltenyi Biotec logo, REAfinity, StainExpress, Vio, VioBlue, and VioGreen are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2022 Miltenyi Biotec and/or its affiliates. All rights reserved.