

## Contents

1. Description
  - 1.1 Principle of the MACS dPCR LENTiCheck Kit (VSV-G), human
  - 1.2 Background information
2. Protocol

## 1. Description

This product is for research use only.

<b>Components</b>	<p><b>80 µL MACS dPCR LENTiCheck (VSV-G) Primers</b> (20 µM each)</p> <p><b>80 µL MACS dPCR LENTiCheck (VSV-G) Probe</b> (5.5 µM): Probe is labeled with FAM dye.</p> <p><b>100 µL MACS dPCR LENTiCheck (VSV-G) Control</b> (10<sup>4</sup> copies/µL)</p>
<b>Capacity</b>	<p>1.6 nmol of each primer and 0.44 nmol of the probe</p> <p>▲ <b>Note:</b> Numbers of reactions depend on used protocol. Please refer to the protocol section for further information.</p>
<b>Product format</b>	<p>MACS dPCR LENTiCheck (VSV-G) Primers and Probe are supplied in nuclease-free water.</p> <p>MACS dPCR LENTiCheck (VSV-G) Control is supplied in a buffered solution of 1 mM TRIS-Cl and 0.1 mM EDTA, at pH 8.0, using nuclease-free water and poly-A-RNA carrier.</p>
<b>Storage</b>	<p>Store protected from light at -20 °C. The expiration date is indicated on the vial label.</p>

### 1.1 Principle of the MACS dPCR LENTiCheck Kit (VSV-G), human

The dPCR replication-competent lentivirus (RCL) assay is intended to detect the lentiviral envelope gene sequence (vesicular stomatitis virus G glycoprotein (VSV-G) gene) of potential RCLs that might arise, e.g., during unintended recombination of plasmids used in vector production or due to recombination of the viral vector with viral sequences in cells after transduction.

The MACS dPCR LENTiCheck Kit (VSV-G), human provides specific primers to amplify an amplicon of the lentiviral envelope VSV-G gene sequence. The envelope gene (VSV-G) is thought to be the most suitable target for detecting RCLs as the incorporation of an envelope gene sequence is required to generate a replication-competent virus.

A VSV-G-specific probe comprising a fluorescent dye (FAM) and quencher moiety is additionally used. During the PCR amplification, the probe binds to the VSV-G gene and gets hydrolyzed by the 5'-3'-exonuclease activity of the PCR polymerase. Thereby, the fluorescent dye is separated from the quencher leading to the emission of a fluorescence signal. The number of positive events

determined by digital PCR is used to calculate the detected copies of the respective gene in the sample.

The MACS dPCR LENTiCheck (VSV-G) Control is used to determine the linearity and sensitivity of the VSV-G gene detection. Additional positive and negative controls are used to ensure proper detection of the VSV-G gene. No signals should be gained for the viral vector or the transduced cells under investigation.

The primers included in the MACS dPCR LENTiCheck Kit (VSV-G), human are only compatible with the lentiviral vectors provided by Lentigen and Miltenyi Biotec.

### 1.2 Background information

During the production of lentiviral vectors as well as in cells transduced with lentiviral vectors, there is the risk of an interplay between the viral vector and viral properties gained from other potentially present viruses or virus genes. Such an interplay might give rise to RCLs, also named replication-competent retroviruses (RCRs). The MACS dPCR LENTiCheck Kit (VSV-G), human has been developed to detect the gene encoding for the pseudotyping envelope protein VSV-G. The VSV-G gene is used to express the envelope protein during the production of lentiviral vectors to enable the packaging of the vector. However, VSV-G should not be present in the viral vector. Therefore, the viral vector as well as the transduced cells are intended to be free of VSV-G gene. The presence of the VSV-G gene can be detected using the MACS dPCR LENTiCheck Kit (VSV-G), human. The absence of the VSV-G gene in *ex vivo* transduced cells according to the dPCR assay can be used for lot release testing, in particular, when time constraints are present.

## 2. Protocol

▲ The primers included in the MACS dPCR LENTiCheck Kit (VSV-G), human are only compatible with the lentiviral vectors provided by Lentigen and Miltenyi Biotec.

▲ Detailed protocols for the MACS LENTiCheck assay using digital PCR systems, e.g. BIO-RAD® QX200™/QX600™/QX ONE Droplet digital PCR™ systems or QIAGEN® QIAcuity® One, Four, or Eight Systems, can be found at [www.miltenyibiotec.com/applicationprotocols](http://www.miltenyibiotec.com/applicationprotocols).

	Number of reactions and volumes
BIO-RAD QX200/QX600/QX ONE Droplet digital PCR systems	80 reactions (1.75 mL reaction volume)
QIAGEN QIAcuity Nanoplate 8.5k	91 reactions (1.1 mL reaction volume)

**Table 1:** Examples of reaction volumes. Contact Miltenyi Biotec Global Technical Support for non-listed digital PCR (dPCR) instruments .

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for local Miltenyi Biotec Technical Support contact information.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

MACS and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2025 Miltenyi Biotec and/or its affiliates. All rights reserved.