

Contents

1. Description
 - 1.1 Background information
 - 1.2 Reagent and instrument requirements
2. Cell preparation

1. Description

This product is for research use only.

Components	100 mL MACSQuant® Tyto® Running Buffer or 6× 100 mL MACSQuant® Tyto® Running Buffer
Capacity	100 mL for up to 10 sorts or 6× 100 mL for 60 sorts with up to 10 mL per sort.
Product format	The filtered (0.22 µm) MACSQuant® Tyto® Running Buffer (pH 7.4) contains low concentration of salts and stabilizer.
Storage	Store the MACSQuant Tyto Running Buffer protected from light at 2–8 °C. The expiration date is indicated on the bottle label.

1.1 Background information

The MACSQuant Tyto Running Buffer has been developed for the optimal processing of cells in combination with the MACSQuant Tyto Cartridges and the MACSQuant Tyto cell sorter.

1.2 Reagent and instrument requirements

- MACSQuant Tyto (# 130-103-931)
- MACSQuant Tyto Cartridges (# 130-106-088)
- Pre-Separation Filters (20 µm) (# 130-101-812)
- Washing buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, and 0.5% bovine serum albumin (BSA). Pass the buffer through a 0.22 µm filter for sterile filtration. Store buffer at 2–8 °C.

2. Cell preparation

▲ Refer to the MACSQuant Tyto user manual for instructions on how to use the instrument.

1. Stain the cells with the fluorochrome-conjugated antibodies of interest according to the respective antibody protocol.
▲ **Note:** If the staining buffer contains EDTA please wash the cells once with washing buffer to remove EDTA.
2. Wash cells by adding washing buffer (10–40 mL) and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
▲ **Note:** The washing buffer volume is depending on the cell number. Please refer to the respective antibody protocol.
3. Resuspend the cell pellet in MACSQuant Tyto Running Buffer (0.5–10 mL) depending on cell number.
4. Pass cells through a 20 µm nylon mesh (Pre-Separation Filters (20 µm), # 130-101-812) to remove cell clumps which may clog the cartridge. Moisten filter with buffer before use.
5. Add the sample into a MACSQuant Tyto Cartridge. Please refer to the data sheet of the MACSQuant Tyto Cartridges for detailed information.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

MACSQuant, the Miltenyi Biotec logo, and Tyto are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.