English | en

1 Intended purpose

The 7-AAD Staining Solution is intended for identification and exclusion of interfering non-viable/degenerating cells in flow cytometry by semi-quantitative, non-automated fluorescent detection of those cells in human whole blood. The reagent serves as an aid to diagnosis to evaluate cell viability in mono- or multiparametric analyses of human leukocytes in flow cytometry in the context of providing information on physiological or pathological state and is to be used by laboratory professionals.

2 Technical data and background information

7-AAD (7-amino-actinomycin D) is a fluorescent dye that intercalates into double-stranded DNA (cytosine and guanine rich regions). It is excluded from viable cells but can penetrate cell membranes of non-viable or degenerating cells.

Technical data

Excitation laser	Blue (488 nm)
Fluorescence	Em _{Max} ^a 647 nm
Product formulation	The ready-to-use 7-AAD Staining Solution is supplied in phosphate-buffered saline, at a concentration of 262.5 µg/mL (0.03%) ^b and containing 1.04% methanol and 0.04% sodium azide
Quantity supplied per vial	100 tests in 0.2 mL (2 µL/test)

^a Emission maximum

^b Values in % refer to %w/w

Table 1: Technical data 7-AAD Staining Solution

3 Important information

Interpretation of results is the full responsibility of the user and analysis results obtained with the reagent should never be used as the sole basis for classifying disease states.

4 Warnings and precautions

WARNING

Risk of infection

All biological specimens and all materials that come into contact with blood and blood products must be treated as infectious material. Regulations for the handling and disposal of potentially infectious material must be observed.¹

Risk of explosion

Sodium azide may react with lead and copper plumbing to form highly explosive buildups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide buildups in plumbing.

Risk of damage to environment

Reagent contains components that are (aqua-)toxic in pure form. National and local regulations for disposal must be observed.

Risk of poisoning

Reagent contains 7-AAD, methanol and sodium azide. In pure form, 7-AAD is a potentially carcinogenic and reproductive toxin; methanol and sodium azide are toxic. Although this mixture is not classified as dangerous at product concentration, safety guidelines must be observed (notably: use of gloves, gown and protective glasses). Never

pipet by mouth and avoid any contact with the skin, mucosa, and eyes. Keep away from heat. Do not breathe in vapors and aerosols.

Hazardous components

1 – <5% methanol

Relevant H and EUH statements

H225 Highly flammable liquid and vapor.

H301 Toxic if swallowed.

H311 Toxic in contact with skin.

H331 Toxic if inhaled.

H370 Causes damage to organs.

EUH210 Safety data sheet available on request.

Table 2: Relevant H and EUH statements

To request the Safety Data Sheet (SDS), contact Miltényi Biotec Technical Support (see chapter 10 "Notice to user").

CAUTION

Risk of delay of results.

Any deviation from instruction for use and protocol can cause impairment of performance. Always observe the instructions for use and protocol.

Consult instructions for use. Strictly adhere to the protocol and only use the reagent as described.

Store and handle at +2 to +8 °C. Do not freeze.

Minimize exposure to light during handling and storage.

For in-use stability at +2 to +8 °C storage temperature refer to the use-by date printed on the vial label.

Do not use after the use-by date.

Do not use the reagent if the vial is damaged or if signs of leakage are observed. Avoid ingress of bacteria into the reagent vial.

- Callihan, D. R. et al. (2014) Protection of Laboratory Workers from Occupationally Acquired Infections.; Clinical and Laboratory Standards Institute (CLSI) 4th ed. CLSI M29-A4.

5 Evidence of deterioration

Any change in the physical appearance of the reagent may indicate deterioration and the reagent should not be used. In case of unexpected results, repeat the measurement or contact Miltényi Biotec Technical Support (see chapter 10 "Notice to user").

6 Materials required but not provided

For material required but not provided, the manufacturer's recommendations and safety regulations must be observed.

- Disposable capped polystyrene tubes, 12×75 mm
- PBS/EDTA/BSA (PEB) buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting e.g., MACS® BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (+2 to +8 °C).

- Red blood cell (RBC) lysing solution: e.g., Red Blood Cell Lysis Solution (10×), Miltenyi Biotec (# 130-094-183), or equivalent
- Double-distilled water (ddH₂O)
- Micropipettes with disposable tips, capable of dispensing volumes of 2 µL up to 2 mL, respectively
- Low speed centrifuge: minimum speed 300×g, with 12×75 mm tube carriers
- Vortex mixer
- Flow cytometer with appropriate laser and filter setting, e.g., MACSQuant® Dx (# 170-082-000). For excitation and emission data for 7-AAD Staining Solution, see chapter 2 "Technical data and background information".
- For multiparametric use: additional immunofluorescent reagents

7 Protocol

For all handling, consideration of good laboratory practice (GLP) regulations is recommended.

7.1 Principle of method

The incubation of human whole blood samples with the provided reagent leads to fluorescent staining of non-viable/degenerating cells. Analysis of the sample is performed in a flow cytometer at a single-cell level. The analysis is based on the detection of characteristic light emission patterns emitted by the fluorescently labeled cell upon excitation with laser light. The collected data can be processed and analyzed using flow cytometry software.

7.2 Sample requirements and limitations

Storage conditions and limitations for whole blood samples should be considered prior to collection and analysis. For collection of patient samples European and national legislation must be followed.

- Venous blood must be taken using sterile tubes containing EDTA as anticoagulant.
- For monoparametric use: whole blood should be stained within 72 hours.
- For multiparametric use: whole blood should be stained within 24 hours.
- Cell count of white blood cells should not exceed 10,000 cells/µL

Note: If necessary, dilute cell sample with buffer.

7.3 Limitations

For multiparametric analysis, it is strongly recommended to check inter-compatibility of the reagent with other immunofluorescent reagents. In particular, be aware of incompatibilities with reagents that have a similar emission spectrum (e.g., PerCP or PE-Cy[®]5 reagents). See chapter 2 "Technical data and background information" for the specifications of the flow cytometer.

Specimens with impaired integrity (hemolytic, icteric, lipemic sample or sample containing cell aggregates) might lead to changes in performance. Hence, such samples shall be rejected.

7.4 Protocol for staining

Note: Do not use reagents containing permeating or fixating agents during the procedure in order to prevent false positive or misleading results.

1. If necessary, prepare a working solution for lysis of red blood cells according to manufacturer's recommendation.
2. Add 2 µL of 7-AAD Staining Solution to 100 µL of cell sample. Optional for multiparametric use: Add the recommended amount of immunofluorescent reagent(s).
3. Mix well and incubate for 15 minutes in the dark at room temperature (+20 to +25 °C).
Note: Higher temperatures and/or longer incubation times may lead to non-specific cell labeling.

4. Perform lysis of red blood cells by following the recommendations of the lysis reagent used.
If using Red Blood Cell Lysis Solution (10×), Miltenyi Biotec (# 130-094-183), add 2 mL of the working solution (1×), vortex immediately and incubate for 10 minutes at room temperature in the dark. Check for successful lysis (clear, red solution) by optical inspection of the tube before continuation. Extend lysis time to up to 20 minutes if necessary.
5. Centrifuge at 300×g for 10 minutes. Remove supernatant.
6. Wash cells by adding 1–2 mL of buffer, centrifuge at 300×g for 10 minutes. Remove supernatant.
7. Resuspend cell pellet in a suitable amount of buffer and proceed to flow cytometric analysis. Store samples at +2 to +8 °C and minimize exposure to light until analysis.
Note: Time until analysis should not exceed 2 hours.

Quality control:

For diagnostic flow cytometry, it is recommended to regularly run control samples from healthy adult specimens or commercially available whole blood controls.

8 Performance characteristics

7-AAD Staining Solution does not have a medical decision point and will not have clinical performance characteristics such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio, or expected values in normal and affected populations.

The values obtained from healthy individuals may vary from laboratory to laboratory. Hence, for multiparametric use, each laboratory should establish its own normal reference ranges for the populations under study.

Analytical performance parameters were tested according to CLSI-H62² recommendations, where applicable. EDTA-whole blood samples from presumed healthy donors (N) were collected and stored up to 96 hours for repeated testing. Samples were acquired within 2 hours post staining on a MACSQuant Dx flow cytometer.

8.1 Analytical specificity

7-AAD is a well-documented viability marker which binds specifically to cytosine and guanine bases of the DNA double strand of non-viable cells.^{3,4,5} Combined with its spectral properties, it is well-suited for flow cytometry analysis.^{6,7} Furthermore, trueness was tested to confirm the specificity of the reagent.

8.2 Accuracy

Precision

Frequency (%pos) of 7-AAD+ subpopulation was tested for repeatability using ≥10 EDTA-whole blood samples over ≥3 days of testing and for reproducibility using ≥5 EDTA-whole blood samples over ≥3 days of testing. Each sample was run in a 4-run cross design (two operators, two instruments, same sample tested in the morning and afternoon). Statistical analysis was conducted based on variance components analysis using JMP[®] version 15.2.

	N	Mean %pos	SD ^a	SD between operator	SD between instrument	% CV ^b
Repeatability	14	41.6	3.09	2.35	0.13	8.59
Reproducibility	11	0.47	0.16	0.08	0.05	33.2

a Standard deviation
b Coefficient of variation

Table 3: Repeatability and reproducibility of 7-AAD Staining Solution

Trueness

Frequency (%pos) of 7-AAD+ subpopulation detected by 7-AAD Staining Solution was tested by comparing the results with a reference reagent on ≥10 whole blood samples.

N	Mean %pos	Mean delta	Mean bias
15	32.0	-1.40	-2.03

Table 4: Trueness of 7-AAD Staining Solution

- Litwin, V. *et al.* (2021) Validation of Assays performed by flow cytometry. Clinical and Laboratory Standards Institute (CLSI). 1st ed. CLSI H62.
- Graves, D. E. and Wadkins, R. M. (1989) 7-Azidoactinomycin D: a novel probe for examining actinomycin D-DNA interactions. *J Biol Chem.* 5;264(13):7262-6.
- Zelenin, A. V. *et al.* (1984) 7-Amino-actinomycin D as a specific fluorophore for DNA content analysis by laser flow cytometry. *Cytometry* 5(4):348-54.
- Gill, J. E. *et al.* (1975) 7-Amino-actinomycin D as a cytochemical probe. I. Spectral properties. *J. Histochem. Cytochem.*, 23, 793-799.
- Schmid, I. *et al.* (1992) Dead cell discrimination with 7-aminoactinomycin D in combination with dual color immunofluorescence in single laser flow cytometry. *Cytometry* 13 (2): 204–208.
- Cossarizza, A. *et al.* (2017) Guidelines for the use of flow cytometry and cell sorting in immunological studies. *Eur. J. Immunol.* 47, 1584–1797.

9 Glossary of symbols

-  In-use stability
-  Contains following number of tests
-  Phone
-  E-Mail
-  Website

10 Notice to user

Any serious incident that has occurred in relation to this product should be reported to Miltenyi Biotec B.V. & Co. KG – using the contact information provided – and the competent authority of the member state in which the user of this product is established.

Miltenyi Biotec B.V. & Co. KG

Friedrich-Ebert-Str. 68
51429 Bergisch Gladbach
Germany

Miltenyi Biotec Technical Support

+49 2204 8306-3803
technicalsupport@miltenyi.com
www.miltenyibiotec.com

11 Document revision history

Revision No.	Date of issue	Modifications since last revision
01	2023-12	N/A, new issue
02	2024-02	Correction of headline format

Table 5: Document revision history

12 Legal notes

autoMACS, MACS, MACSQuant, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. Cy is a registered trademark of GE Healthcare UK Limited. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only. Copyright © 2024 Miltenyi Biotec and/or its affiliates. All rights reserved.

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2 Litwin, V. *et al.* (2021) Validation of Assays performed by flow cytometry. Clinical and Laboratory Standards Institute (CLSI). Erstaufgabe CLSI H62.

3 Graves, D. E. and Wadkins, R. M. (1989) 7-Azidoactinomycin D: a novel probe for examining actinomycin D-DNA interactions. J Biol Chem. 5;264(13):7262-6.

4 Zelenin, A. V. *et al.* (1984) 7-Amino-actinomycin D as a specific fluorophore for DNA content analysis by laser flow cytometry. Cytometry 5(4):348-54.

5 Gill, J. E. *et al.* (1975) 7-Amino-actinomycin D as a cytochemical probe. I. Spectral properties. J. Histochem. Cytochem., 23, 793-799.

6 Schmid, I. *et al.* (1992) Dead cell discrimination with 7-aminoactinomycin D in combination with dual color immunofluorescence in single laser flow cytometry. Cytometry 13 (2): 204–208.

7 Cossarizza, A. *et al.* (2017) Guidelines for the use of flow cytometry and cell sorting in immunological studies. Eur. J. Immunol. 47, 1584–1797.

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