

Contents

1. Description
 - 1.1 Principle of the MACS® Separation
 - 1.2 Background information
 - 1.3 Applications
 - 1.4 Reagent and instrument requirements
2. Protocol
 - 2.1 Sample preparation
 - 2.2 Magnetic labeling of non-PDCs
 - 2.3 Magnetic separation: Depletion of non-PDCs
 - 2.4 Magnetic labeling of PDCs
 - 2.5 Magnetic separation: Positive selection of PDCs
3. Example of a separation using the Diamond Plasmacytoid Dendritic Cell Isolation Kit II
4. References

1. Description


Components	2 mL Non-PDC Biotin-Antibody Cocktail II, human: Cocktail of biotin-conjugated monoclonal antibodies against antigens that are not expressed by plasmacytoid dendritic cells. 2 mL Non-PDC MicroBead Cocktail II, human: Cocktail of MicroBead-conjugated monoclonal antibodies against antigens that are not expressed by plasmacytoid dendritic cells and MicroBeads conjugated to a monoclonal anti-biotin antibody (isotype: mouse IgG1). 2 mL CD304 (BDCA-4/Neuropilin-1) Diamond MicroBeads, human: MicroBeads conjugated to monoclonal anti-CD304 (BDCA-4/Neuropilin-1) antibodies (isotype: mouse IgG1).
Capacity	For 2×10 ⁹ total cells, up to 20 separations.
Product format	All components are supplied in buffer containing stabilizer and 0.05% sodium azide.
Storage	Store protected from light at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label.

1.1 Principle of the MACS® Separation

The Diamond Plasmacytoid Dendritic Cell Isolation Kit II has been developed for the isolation of extremely pure plasmacytoid dendritic cells (PDCs) from peripheral blood mononuclear cells

(PBMCs). The expected purity is almost 100%.

The isolation of PDCs is performed in a two-step procedure. First, the non-PDCs are indirectly magnetically labeled with a cocktail of biotin-conjugated antibodies against lineage-specific antigens and Anti-Biotin MicroBeads. Additionally, non-PDCs are directly magnetically labeled with a cocktail of MicroBead-conjugated antibodies against antigens that are not expressed on PDCs. The labeled cells are subsequently depleted by separation over a MACS® Column, which is placed in the magnetic field of a MACS Separator. In the second step, the pre-enriched PDCs are labeled with PDC-specific CD304 (BDCA-4/Neuropilin-1) Diamond MicroBeads and isolated by positive selection fraction by separation over a MACS Column, which is placed in the magnetic field of a MACS Separator. After removing the column from the magnetic field, the magnetically retained PDCs can be eluted as the positively selected cell fraction. To increase the purity, the positively selected cell fraction containing the PDCs must be separated over a second column.



1.2 Background information

PDCs are one of three subsets of dendritic cells originally identified in human peripheral blood. They are also known as plasmacytoid T cells, plasmacytoid monocytes, lymphoid dendritic cells, IFNα/β-producing cells (IPCs), or type 2 pre-dendritic cells (pDC2).^{1,2} Upon viral infection, they produce high amounts of type I interferons, which block viral replication and stimulate innate and adaptive immune responses. In culture, they mature into potent antigen-presenting cells, after exposure to IL-3 alone or in combination with an appropriate stimulus. In healthy donors, PDCs represent about 0.4% of total PBMCs. Apart from blood, immature PDCs have been found in human lymphoid tissue and in inflammatory sites, e.g., skin of systemic lupus erythematosus (SLE)³ or psoriasis vulgaris patients⁴.

In blood and bone marrow, PDCs are identified as being CD303 (BDCA-2)⁺, CD304 (BDCA-4/Neuropilin-1)⁺, CD123⁺, CD11c⁻. Unlike CD123, CD303 (BDCA-2) and CD304 (BDCA-4/Neuropilin-1) are exclusively expressed on PDCs.^{5–8} Further, they

are CD4⁺, CD45RA⁺, CD141 (BDCA-3)^{dim}, CD1c (BDCA-1)⁺, CD2⁻, lack expression of lineage markers (CD3, CD14, CD16, CD19, CD20, CD56), and express neither myeloid markers, e.g., CD13 and CD33, nor Fc receptors such as CD16, CD64, or FcεRI.⁹

Binding of antibodies to CD304 (BDCA-4/Neuropilin-1) does not have a substantial effect on IFN type I production in PDCs after stimulation.^{1,9,10}

▲ The BDCA-2 and BDCA-4 (Neuropilin-1) antigens have been assigned to CD303 and CD304, respectively, at the 8th International Workshop in Human Leukocyte Differentiation Antigens in Adelaide, Australia, in 2004.

1.3 Applications

- Isolation of PDCs from human PBMCs for phenotypic, functional or molecular analyses, e.g., studies on expression of Toll-like receptors^{8,11-13}, chemokine receptors^{7,12,14}, or new antigens, and on dendritic cell activation⁸, migration⁷, cytokine production^{1,8,12}, and T cell polarization^{1,8,15}.

1.4 Reagent and instrument requirements

- Buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (2–8 °C). Degas buffer before use, as air bubbles could block the column.

▲ **Note:** EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). BSA can be replaced by other proteins such as human serum albumin, human serum, or fetal bovine serum (FBS). Buffers or media containing Ca²⁺ or Mg²⁺ are not recommended for use.

- MACS Columns and MACS Separators: Depletion of non-PDCs can be performed on an LD Column. The subsequent positive selection of PDCs can be performed on two MS Columns. Positive selection and depletion can also be performed by using the autoMACS Pro or the autoMACS Separator.

Column	Max. number of labeled cells	Max. number of total cells	Separator
Positive selection			
MS	10 ⁷	2×10 ⁸	MiniMACS, OctoMACS, VarioMACS, SuperMACS II
Depletion			
LD	10 ⁸	5×10 ⁸	MidiMACS, QuadroMACS, VarioMACS, SuperMACS II
Positive selection or depletion			
autoMACS	2×10 ⁸	4×10 ⁹	autoMACS Pro, autoMACS

▲ **Note:** Column adapters are required to insert certain columns into the VarioMACS™ or SuperMACS™ II Separators. For details refer to the respective MACS Separator data sheet.

- (Optional) Fluorochrome-conjugated antibodies for flow cytometric analysis, e.g., CD303 (BDCA-2)-FITC (# 130-090-510), CD303 (BDCA-2)-PE (# 130-090-511), CD303 (BDCA-2)-APC (# 130-090-905), CD123-FITC (# 130-090-897), CD123-PE (# 130-090-899), or CD123-APC (# 130-090-901). For more information about antibodies refer to www.miltenyibiotec.com/antibodies.
- (Optional) Propidium Iodide Solution (# 130-093-233) or 7-AAD for flow cytometric exclusion of dead cells.
- (Optional) Dead Cell Removal Kit (# 130-090-101) for the depletion of dead cells.

- (Optional) Pre-Separation Filters, 30 μm (# 130-041-407) to remove cell clumps.

2. Protocol

2.1 Sample preparation

When working with anticoagulated peripheral blood or buffy coat, peripheral blood mononuclear cells (PBMCs) should be isolated by density gradient centrifugation, for example, using Ficoll-Paque™.

▲ **Note:** To remove platelets after density gradient separation, resuspend cell pellet in buffer and centrifuge at 200×g for 10–15 minutes at 20 °C. Carefully aspirate supernatant. Repeat washing step.

When working with tissues or lysed blood, prepare a single-cell suspension using standard methods.

For details see the protocols section at www.miltenyibiotec.com/protocols.

▲ Dead cells may bind non-specifically to MACS MicroBeads. To remove dead cells, we recommend using density gradient centrifugation or the Dead Cell Removal Kit (# 130-090-101).



2.2 Magnetic labeling of non-PDCs

▲ Work fast, keep cells cold, and use pre-cooled solutions. This will prevent capping of antibodies on the cell surface and non-specific cell labeling.

▲ Volumes for magnetic labeling given below are for up to 10⁸ total cells. When working with fewer than 10⁸ cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10⁸ total cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ For optimal performance it is important to obtain a single-cell suspension before magnetic labeling. Pass cells through 30 μm nylon mesh (Pre-Separation Filters, 30 μm, # 130-041-407) to remove cell clumps which may clog the column. Moisten filter with buffer before use.

▲ The recommended incubation temperature is 2–8 °C. Working on ice may require increased incubation times. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling.

- Determine cell number.
- Centrifuge cell suspension at 300×g for 10 minutes. Aspirate supernatant completely.
- Resuspend cell pellet in 400 μL of buffer per 10⁸ total cells.
- Add 100 μL of Non-PDC Biotin-Antibody Cocktail II per 10⁸ total cells.
- Mix well and incubate for 10 minutes in the refrigerator (2–8 °C).
- Wash cells by adding 5–10 mL of buffer per 10⁸ cells and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
- Resuspend cell pellet in 400 μL of buffer per 10⁸ total cells.
- Add 100 μL of Non-PDC MicroBead Cocktail II per 10⁸ total cells.
- Mix well and incubate for 15 minutes at 2–8 °C.

10. Wash cells by adding 5–10 mL of buffer per 10^8 cells and centrifuge at $300\times g$ for 10 minutes. Aspirate supernatant completely.
11. Resuspend cells in a final volume of 500 μL per 10^8 total cells.
12. Proceed to magnetic separation (2.3).



2.3 Magnetic separation: Depletion of non-PDCs

- ▲ Choose an appropriate MACS Column and MACS Separator according to the number of total cells and the number of PDCs. For details refer to table in section 1.4.
- ▲ Always wait until the column reservoir is empty before proceeding to the next step.

Depletion with LD Columns

1. Place LD Column in the magnetic field of a suitable MACS Separator. For details refer to LD Column data sheet.
2. Prepare column by rinsing with 2 mL of buffer.
3. Apply cell suspension onto the column.
4. Collect unlabeled cells that pass through and wash column with 2×1 mL of buffer. Collect total effluent; this is the unlabeled pre-enriched PDC fraction. Perform washing steps by adding buffer two times. Only add new buffer when the column reservoir is empty.
5. Proceed to 2.4 for the labeling of PDCs.

Depletion with the autoMACS® Pro Separator or the autoMACS® Separator

- ▲ Refer to the respective user manual for instructions on how to use the autoMACS® Pro Separator or the autoMACS Separator.
- ▲ Buffers used for operating the autoMACS Pro Separator or the autoMACS Separator should have a temperature of $\geq 10^\circ\text{C}$.
- ▲ Program choice depends on the isolation strategy, the strength of magnetic labeling, and the frequency of magnetically labeled cells. For details refer to the section describing the cell separation programs in the respective user manual.

Magnetic separation with the autoMACS® Pro Separator

1. Prepare and prime the instrument.
2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube in row A of the tube rack and the fraction collection tubes in rows B and C.
3. For a standard separation choose the following program:
Depletion: Depletes
Collect negative fraction in row B of the tube rack.
4. Proceed to 2.4 for the labeling of PDCs.

Magnetic separation with the autoMACS® Separator

1. Prepare and prime the instrument.
2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube at the uptake port and the fraction collection tubes at port neg1 and port pos1.

3. For a standard separation choose the following program:

Depletion: Depletes

Collect negative fraction from outlet port neg1.

4. Proceed to 2.4 for the labeling of PDCs.



2.4 Magnetic labeling of PDCs

▲ Volumes for magnetic labeling given below are for an initial starting cell number of up to 10^8 total cells. For higher initial cell numbers, scale up all volumes accordingly.

1. Centrifuge cell suspension at $300\times g$ for 10 minutes. Aspirate supernatant completely.
2. Resuspend cells directly in 100 μL CD304 (BDCA-4/Neuropilin-1) Diamond MicroBeads.
3. Mix well and incubate for 15 minutes in the dark in the refrigerator ($2-8^\circ\text{C}$).
4. Wash cells by adding 1–2 mL of buffer and centrifuge at $300\times g$ for 10 minutes. Aspirate supernatant completely.
5. Resuspend up to 10^8 cells in 500 μL of buffer.
▲ **Note:** For higher cell numbers, scale up buffer volume accordingly.
6. Proceed to magnetic separation (2.5).



2.5 Magnetic separation: Positive selection of PDCs

Positive selection with MS Columns

▲ To achieve highest purities, perform two consecutive column runs.

1. Place MS Column in the magnetic field of a suitable MACS Separator. For details refer to MS Column data sheet.
2. Prepare column by rinsing with 2 mL of buffer.
3. Apply cell suspension onto the column. Collect flow-through containing unlabeled cells.
4. Wash column with 2×1 mL of buffer. Collect unlabeled cells that pass through and combine with the effluent from step 3.
▲ **Note:** Perform washing steps by adding buffer aliquots only when the column reservoir is empty.
5. Remove column from the separator and place it on a suitable collection tube.
▲ **Note:** To perform a second column run, you may elute the cells directly from the first onto the second, equilibrated column instead of a collection tube.
6. Pipette 1 mL of buffer onto the column. Immediately flush out the magnetically labeled cells by firmly pushing the plunger into the column.
7. To increase purity of PDCs, the eluted fraction must be enriched over a second MS Column. Repeat the magnetic separation procedure as described in steps 1 to 6 by using a new column.

Positive selection with the autoMACS® Pro Separator or the autoMACS® Separator

Magnetic separation with the autoMACS® Pro Separator

1. Prepare and prime the instrument.

2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube in row A of the tube rack and the fraction collection tubes in rows B and C.
3. For a standard separation choose the following program:

Positive selection: Possel


Collect positive fraction in row C of the tube rack. This is the enriched PDC fraction.

Magnetic separation with the autoMACS® Separator

1. Prepare and prime the instrument.
2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube at the uptake port and the fraction collection tubes at port neg1 and port pos1.
3. For a standard separation choose the following program:



Positive selection: Possel

Collect positive fraction from outlet port pos1. This is the enriched PDC fraction.



3. Example of a separation using the Diamond Plasmacytoid Dendritic Cell Isolation Kit II

PDCs were isolated from human PBMCs by using the Diamond Plasmacytoid Dendritic Cell Isolation Kit II, an LD and two MS Columns, a MidiMACS™ and a MiniMACS™ Separator. The cells were fluorescently stained with CD123-PE and CD303 (BDCA-2)-APC and analyzed by flow cytometry using the MACSQuant® Analyzer. Cell debris and dead cells were excluded from the analysis based on scatter signals and propidium iodide fluorescence.



Isolated PDCs

4. References

1. Dzionek, A. *et al.* (2002) Plasmacytoid dendritic cells: from specific surface markers to specific cellular functions. *Hum. Immunol.* 63: 1133–1148.
2. McKenna, K. *et al.* (2005) Plasmacytoid dendritic cells: linking innate and adaptive immunity. *J. Virol.* 79: 17–27.
3. Rönnblom, L. and Alm, G.V. (2001) A pivotal role for the natural interferon α -producing cells (plasmacytoid dendritic cells) in the pathogenesis of lupus. *J. Exp. Med.* 194: F59–F63.
4. Gilliet, M. *et al.* (2004) Psoriasis triggered by toll-like receptor 7 agonist imiquimod in the presence of dermal plasmacytoid dendritic cell precursors. *Arch. Dermatol.* 140: 1490–1495.
5. Dzionek, A. *et al.* (2000) BDCA-2, BDCA-4 and BDCA-4: three markers for distinct subsets of dendritic cells in human peripheral blood. *J. Immunol.* 165: 6037–6046.
6. Grabbe, S. *et al.* (2000) Dendritic cells: multi-lineal and multi-functional. *Immunol. Today* 21: 431–433.
7. Penna, G. *et al.* (2001) Cutting edge: selective usage of chemokine receptors by plasmacytoid dendritic cells. *J. Immunol.* 167: 1862–1866.
8. Krug, A. *et al.* (2001) Toll-like receptor reveals CpG DNA as a unique microbial stimulus for plasmacytoid dendritic cells which synergizes with CD40 ligand to induce high amounts of IL-12. *Eur. J. Immunol.* 31: 3026–3037.
9. Dzionek, A. *et al.* (2001) BDCA-2, a novel plasmacytoid dendritic cell-specific type II C-type lectin, mediates antigen capture and is a potent inhibitor of interferon α/β induction. *J. Exp. Med.* 194: 1823–1834.
10. Colonna, M. *et al.* (2004) Plasmacytoid dendritic cells in immunity. *Nature Immunol.* 5: 1219–1226.
11. Hornung, V. *et al.* (2002) Quantitative expression of Toll-like receptor 1–10 mRNA in cellular subsets of human peripheral blood mononuclear cells and sensitivity to CpG oligodeoxynucleotides. *J. Immunol.* 168: 4531–4537.
12. Jarrossay, D. *et al.* (2001) Specialization and complementarity in microbial molecule recognition by human myeloid and plasmacytoid dendritic cells. *Eur. J. Immunol.* 31: 3388–3393.
13. Matsumoto, M. *et al.* (2003) Subcellular localization of Toll-like receptor 3 in human dendritic cells. *J. Immunol.* 171: 3154–3162.
14. Penna, G. *et al.* (2002) Cutting edge: differential chemokine production by myeloid and plasmacytoid dendritic cells. *J. Immunol.* 169: 6673–6676.
15. Kaser, A. *et al.* (2004) Interleukin-18 attracts plasmacytoid dendritic cells (DC2s) and promotes Th1 induction by DC2s through IL-18 receptor expression. *Blood* 103: 648–655.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, MACSQuant, MidiMACS, the Miltenyi Biotec logo, MiniMACS, OctoMACS, QuadroMACS, SuperMACS, and VarioMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Ficoll-Paque is a trademark of GE Healthcare companies.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.