

FeraSpin™ Series: FeraSpin™ XL

MRI contrast agent for
pre-clinical imaging

1 vial (5×100 µL injections)

130-095-172

5 vials (25×100 µL injections)

130-095-173

Contents

1. Description
 - 1.1 Background information
 - 1.2 Applications
 - 1.3 Physico-chemical properties
 - 1.4 Requirements
2. Protocol
 - 2.1 Preparation
 - 2.2 Injection
 - 2.3 Imaging
3. References
4. Related products

1. Description

| | |
|----------------|---|
| Components | 850 µL FeraSpin™ XL MRI contrast agent (superparamagnetic iron oxide nanoparticles) or 5×850 µL FeraSpin™ XL MRI contrast agent (superparamagnetic iron oxide nanoparticles). |
| Capacity | 5×100 µL injections or 25×100 µL injections. |
| Product format | FeraSpin XL is supplied as a sterile isotonic solution with an iron concentration of 10 mM. |
| Appearance | Clear, amber liquid. |
| Storage | Store at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label. |

For laboratory and animal research use only. Warning: Not for human or animal therapeutic or diagnostic use. Make sure to comply with all laws and regulations governing research on animals.

1.1 Background information

FeraSpin XL is a nanoparticulate superparamagnetic iron oxide contrast agent specifically formulated for pre-clinical magnetic resonance imaging (MRI).

FeraSpin XL belongs to the FeraSpin Series, which encompasses size-selected, narrowly distributed nanoparticulate agents with well-defined particle sizes derived from FeraSpin R (for details refer to the product portfolio at www.miltenyibiotec.com/viscover). The FeraSpin Series includes FeraSpin XS, FeraSpin S, FeraSpin M, FeraSpin L, FeraSpin XL, and FeraSpin XXL. Their identical composition allows for exclusive size-dependent studies and selection of the most suitable contrast agent for the intended application.

FeraSpin XS, S, M, L, XL, and XXL are agents of high relaxivity. They enhance the contrast in T2- and T2*-weighted MRI due to a shortening of the spin-spin relaxation time (T2) and increase the signal intensity in T1-weighted MRI due to a shortening of the spin-lattice relaxation time (T1). The T2-effect increases with increasing particle size whereas the T1-effect increases with decreasing particle size. On accumulation in cells, the T1-effect diminishes and the T2-effect increases.

Upon intravenous injection, all agents of this series circulate in the bloodstream and are taken up by macrophages. They accumulate in the liver and spleen and are degraded within a few days with their iron being transferred into the physiological iron stores. The macrophage uptake varies in dependence of particle size. Increased uptake by the Kupffer cells (macrophages of the liver) with increasing particle size leads to a rapid accumulation in the liver and spleen and a short blood circulation time. With decreasing particle size, the uptake by the Kupffer cells is reduced leading to a prolonged circulation time and increased uptake by other macrophages. The following schematic diagram shows the characteristics of the FeraSpin contrast agents in dependence of their particle size.




Figure 1: Characteristics of the FeraSpin Series in dependence of particle size.

1.2 Applications

The FeraSpin Series is an innovative research toolbox that offers solutions for a wide range of imaging applications.

It provides the possibility for selection of the MR contrast properties, the pharmacokinetic properties as well as the iron oxide load and, thus, allows a tailoring of the contrast agent to customer needs.

FeraSpin XS, S, M, L, XL, and XXL are indicated for use in MRI of small animals, for example, mice. They can be used in various applications, such as *ex vivo* cell labeling or *in vivo* macrophage labeling for inflammation imaging. Principally, all agents of this series can be used to facilitate the visualization of the liver and spleen as well as visualization of the vasculature.

▲ **Note:** For liver and spleen imaging the use of FeraSpin R is recommended. In applications where a long blood circulation time and a strong T1-effect are favorable, for example in MR angiography, the use of FeraSpin XS is most suited. For more details refer to the product portfolio at www.miltenyibiotec.com/viscover.

▲ **Note:** The contrast agents of this series are provided with equal composition for reasons of comparability. For custom-tailored concentrations please contact the customer support.

1.3 Physico-chemical properties

| FeraSpin | Mean particle size [#] | Relaxivity, r2/r1 ^{##} |
|----------|---------------------------------|---------------------------------|
| XS | 10–20 nm | 3–5 |
| S | 20–30 nm | 5–9 |
| M | 30–40 nm | 9–22 |
| L | 40–50 nm | 22–32 |
| XL | 50–60 nm | 32–39 |
| XXL | 60–70 nm | 39–46 |

[#] hydrodynamic diameter, ^{##} in water, 37 °C, 1.41 T




Figure 2: Schematic diagram of small, medium, and large-sized nanoparticles of the FeraSpin Series.

1.4 Requirements

- Sterile syringes and needles (27G–30G)

▲ **Note:** To allow sufficient volume for 5×100 µL injections per vial, the syringe/needle dead volume should be kept below 70 µL. Tip: Use insulin or tuberculin syringes.

- 70% ethanol

2. Protocol

2.1 Preparation

▲ Read the entire protocol before starting.

▲ **Tip:** For optimum device settings perform initial studies in a suitable imaging phantom.

▲ The contrast agent is ready for injection as provided.

▲ Standard animal-handling procedures and local regulations must be followed.

▲ The dosing varies in dependence of the intended application as well as the selected FeraSpin contrast agent and, thus, has to be adapted accordingly. For a mouse weighing 20–30 g an injection volume of 100 µL corresponds to a dose of 40 µmol Fe/kg body weight (for a 25 g mouse).

2.2 Injection

▲ FeraSpin XL contains no preservatives. Avoid microbial contamination and discard any unused material after 24 hours.

1. Vortex the vial to ensure thorough mixing.
2. Disinfect the septum with 70% ethanol. Let septum dry.
3. Warm the mouse tail to dilate the veins and enhance their visibility.
4. Inject FeraSpin XL via the lateral tail vein of the mouse.

2.3 Imaging

▲ Imaging can be performed in a multitude of devices at all commonly used field strengths including high-field MRI.

▲ FeraSpin XL can be detected by T1- as well as T2- and T2*-weighted sequences.

▲ The time interval between injection and imaging depends on the application. For applications involving imaging of the vasculature begin imaging immediately after injection. For other applications imaging over an extended time period after injection, for example 24 hours, is recommended.

Find examples of FeraSpin-enhanced MR images at www.miltenyibiotec.com/viscover.

3. References

1. Metz, S. *et al.* (2004) Capacity of human monocytes to phagocytose approved iron oxide MR contrast agents *in vitro*. *Eur. Radiol.* 14: 1851–1858.
2. Golovko, D. M. *et al.* (2009) Accelerated stem cell labeling with ferucarbotran and protamine. *Eur. Radiol.* 20: 640–648.
3. Allkemper, T. *et al.* (2002) Contrast-enhanced blood-pool MR angiography with optimized iron oxides: effect of size and dose on vascular contrast enhancement in rabbits. *Radiology* 223: 432–438.
4. Bremer, C. *et al.* (1999) RES-specific imaging of the liver and spleen with iron oxide particles designed for blood pool MR-angiography. *J. Magn. Reson. Imaging* 10: 461–467.
5. Briley-Saebo, K. C. *et al.* (2006) Clearance of iron oxide particles in rat liver. effect of hydrated particle size and coating material on liver metabolism. *Invest. Radiol.* 41: 560–570.
6. Tsuda, N. *et al.* (2005) Potential of superparamagnetic iron oxide in the differential diagnosis of metastasis and inflammation in bone marrow: experimental study. *Invest. Radiol.* 40: 676–681.
7. Persigehl, P. *et al.* (2007) Antiangiogenic tumor treatment: early noninvasive monitoring with USPIO-enhanced MR imaging in mice. *Radiology* 244: 449–456.
8. Frericks, B. B. *et al.* (2009) Magnetic resonance imaging of experimental inflammatory bowel disease: quantitative and qualitative analyses with histopathologic correlation in a rat model using the ultrasmall iron oxide SHU 555 C. *Invest. Radiol.* 44: 23–30.
9. Simon, G. H. *et al.* (2006) Ultrasmall superparamagnetic iron oxide-enhanced magnetic resonance imaging of antigen-induced arthritis: a comparative study between SHU555 C, Ferumoxtran-10, and Ferumoxytol. *Invest Radiol* 41: 45–51.

4. Related products

| | |
|---------------|------------------------------|
| FeraSpin™ R | # 130-095-138, # 130-095-139 |
| FeraSpin™ XS | # 130-095-140, # 130-095-141 |
| FeraSpin™ S | # 130-095-166, # 130-095-167 |
| FeraSpin™ M | # 130-095-168, # 130-095-169 |
| FeraSpin™ L | # 130-095-170, # 130-095-171 |
| FeraSpin™ XXL | # 130-095-174, # 130-095-175 |
| GadoSpin™ M | # 130-095-134, # 130-095-135 |
| GadoSpin™ P | # 130-095-136, # 130-095-137 |
| GadoSpin™ F | # 130-095-162, # 130-095-163 |
| GadoSpin™ D | # 130-095-164, # 130-095-165 |

A comprehensive product portfolio for the imaging modalities MRI, CT, US, OI, SPECT, and PET is available at www.miltenyibiotec.com/viscover.

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

MACS is a registered trademark of Miltenyi Biotec GmbH.

FeraSpin and GadoSpin are trademarks of nanoPET GmbH.

Manufacturer: nanoPET Pharma GmbH, Berlin, Germany, exclusively for Miltenyi Biotec GmbH.

Copyright © 2010 Miltenyi Biotec GmbH. All rights reserved.