

Contents

1. Description
 - 1.1 Background information
 - 1.2 Applications
 - 1.3 Reagent and instrument requirements
2. Protocols
 - 2.1 Immunofluorescent staining of nucleated cells, e.g., PBMCs, using the StainExpress CD19 CAR T Monitoring Cocktail, optional in combination with additional antibodies
 - 2.2 Immunofluorescent staining and lysis of whole blood or leukapheresis products using the StainExpress CD19 CAR T Monitoring Cocktail, optional in combination with additional antibodies
 - 2.3 Flow cytometric data acquisition with the MACSQuant® Analyzer 16
3. Examples of immunofluorescent staining with the StainExpress CD19 CAR T Monitoring Cocktail

1. Description

This product is for research use only.

Components 12 tubes StainExpress CD19 CAR T Monitoring Cocktail, human:

Dried cocktail of fluorochrome-conjugated recombinant engineered REAfinity® Antibodies (isotype: recombinant human IgG1) containing: CD45-VioGreen™ (clone: REA747), CD4-APC-Vio® 770 (clone: REA623), CD3-FITC (clone: REA613), CD19 CAR FMC63-PE (clone: REA1297), CD8-VioBlue® (clone: REA734), 7-AAD Staining Solution.

2 StainExpress CD19 CAR T Monitoring Compensation Sets, human:

Each set consists of five tubes containing dried single antibody fluorochrome-conjugated recombinant engineered REAfinity Antibodies (isotype: recombinant human IgG1) for compensation controls: CD45-VioGreen (clone: REA747), CD4-APC-Vio 770 (clone: REA623), CD3-FITC (clone: REA613), CD19 CAR FMC63-PE (clone: REA1297), CD8-VioBlue (clone: REA734).

2 StainExpress CD19 CAR T Monitoring FMO Sets, human:

Each set contains six dried cocktails of fluorochrome-conjugated recombinant engineered REAfinity Antibodies (isotype: recombinant human IgG1) containing:

CD45-VioGreen (clone: REA747),
CD4-APC-Vio 770 (clone: REA623),
CD3-FITC (clone: REA613),
CD8-VioBlue (clone: REA734),
7-AAD Staining Solution.

Capacity 12 tests, one test for up to 10⁶ total cells.

Product format Antibodies are supplied in a dry format containing stabilizer.

Storage Store at dry conditions in a closed pouch. Store protected from light at +19 to +25°C. The expiration date is indicated on the pouch label.

1.1 Background information

The StainExpress CD19 CAR T Monitoring Cocktail, human, offers a pre-formulated antibody backbone panel designed to simplify and standardize the characterization of CD19 CAR T cells using flow cytometry. This comprehensive panel includes the CD19 CAR FMC63 Idiotypic Antibody, REAfinity™, specifically developed to detect genetically modified T cells engineered to express a chimeric antigen receptor (CAR) containing a single-chain variable fragment (scFv) derived from the mouse anti-human CD19 antibody clone FMC63.

The panel can be expanded with additional liquid markers in available channels, enabling customization for phenotyping, such as assessing the differentiation state of CD19 CAR T cells and conducting other analyses. The cocktail is optimized to evaluate the expansion and persistence of CD19 CAR T cells in human peripheral blood.

1.2 Applications

- Identification and analysis of CD19 CAR⁺ T cells by flow cytometry.
- Evaluation of the expansion and persistence of CD19 CAR⁺ T cells in patient monitoring research.

1.3 Reagent and instrument requirements

- Buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS® BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (+2 to +8 °C).

▲ **Note:** EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). BSA can be replaced by other proteins such as human serum albumin, human serum, or fetal bovine serum (FBS). Buffers or media containing Ca^{2+} or Mg^{2+} are not recommended for use.

- Flow cytometer, e.g., MACSQuant Analyzer 16 (# 130-109-803) or MACSQuant Analyzer 10 (# 130-096-343).

▲ **Note:** The MACSQuant VYB cannot be used.

- MACSQuant Running Buffer (# 130-092-747)
- (Optional) MACS Comp Bead Kit, anti-REA (# 130-104-693), for compensation of the fluorescence spillover from fluorochrome-conjugated antibodies.
- (Optional) MACS MiniSampler Plus (# 130-105-745).
- (Optional) Chill 5 Rack (# 130-092-951).
- (Optional) Red Blood Cell Lysis Solution (10×) (# 130-094-183)
- (Optional) StainExpress Immune Cell Composition Cocktail, human (# 130-127-637)
- (Optional) Markers for customization to expand the StainExpress CD19 CAR T Monitoring Cocktail by incorporating liquid antibodies into available channels. Recommended antibodies include:

Antibody	Clone
CD14	REA599
CD19	REA675
CD45RO	REA611
CD197 (CCR7)	REA546

2. Protocols

2.1 Immunofluorescent staining of nucleated cells, e.g., PBMCs, using the StainExpress CD19 CAR T Monitoring Cocktail, optional in combination with additional antibodies

▲ Volumes given below are for up to 10^6 nucleated cells. When working with fewer than 10^6 cells, use the same volumes as indicated. When working with higher cell numbers, use multiple tubes accordingly (e.g. for 2×10^6 nucleated cells, use two tubes).

1. Determine cell number.
2. Adjust cell concentration to up to 10^6 nucleated cells per 100 μL using buffer.
▲ **Note:** If necessary, centrifuge cell suspension at $300 \times g$ for 5 minutes, aspirate supernatant completely, and resuspend up to 10^6 nucleated cells per 100 μL of buffer.
3. Add 100 μL of cell suspension to one tube of StainExpress CD19 CAR T Monitoring Cocktail.
4. (Optional) Add suitable amounts of other liquid antibodies for available channels (refer to chapter 1.3, table 1) according to manufacturer's recommendation.
5. Incubate for 10–20 seconds, then vortex at high speed for 5 seconds and incubate for 10 minutes in the dark at room temperature (+19 to +25 °C).

▲ **Note:** Working at lower temperatures requires increased incubation times.

6. Wash cells by adding 1 mL of buffer and centrifuge at $300 \times g$ for 5 minutes. Aspirate supernatant completely.
7. (Optional for fixation) Add 250 μL of buffer and 250 μL of Inside Fix to the cells and incubate for 20 minutes in the dark at room temperature (+19 to +25 °C).
8. (Optional for fixation) Add up to 2 mL of buffer.
9. (Optional for fixation) Centrifuge cells at $300 \times g$ for 5 minutes at room temperature (+19 to +25 °C). Aspirate supernatant completely.
10. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry.
▲ **Note:** Store samples at +2 to +8 °C protected from light until analysis. Acquire samples within 1 hour after staining.
11. (Optional) As a negative control, fluorescence-minus-one (FMO) controls can be prepared by adding a sample into a tube from the matching StainExpress FMO Set. These tubes do not contain the CD19 CAR FMC63 Idiotypic Antibody, PE, REAfinity. These controls should be treated exactly the same as the fully stained samples to ensure accurate gating of CD19 CAR⁺ T cell populations
12. Proceed to flow cytometric analysis.

2.2 Immunofluorescent staining and lysis of whole blood or leukapheresis products using the StainExpress CD19 CAR T Monitoring Cocktail, optional in combination with additional antibodies

▲ (Optional) To determine cell counts, prepare a replicate sample using a no wash protocol. For optimal convenience and comprehensive cellular composition information, run the StainExpress Immune Cell Composition Cocktail on a MACSQuant Analyzer with precise volumetric counting. Calculate cell count/mL by multiplying the number of cells/mL obtained with the StainExpress Cocktail by the dilution factor (e.g. $\times 20$). To determine the number of PE-positive cells, such as CAR⁺ T cells, use the CD19 CAR T Monitoring Cocktail. For example, if 50% of the CD3⁺ T cells are PE-positive and the measure is 1×10^6 CD3⁺ T cells/mL, then 5×10^5 CAR⁺ T cells/mL are in the sample.

1. Dilute 10× Red Blood Cell Lysis Solution 1:10 with double-distilled water (ddH_2O), for example, dilute 1 mL of 10× Red Blood Cell Lysis Solution with 9 mL of ddH_2O .
▲ **Note:** Do not dilute with deionized water. Store prepared 1× Red Blood Cell Lysis Solution at room temperature. Discard unused solution at the end of the day.
2. Add 100 μL of whole blood or leukapheresis product to one StainExpress CD19 CAR T Monitoring Cocktail tube.
3. (Optional) Add suitable amounts of other liquid antibodies for available channels (refer to chapter 1.3, table 1) according to manufacturer's recommendation.
4. Incubate for 10–20 seconds, then vortex at high speed for 5 seconds and incubate for 10 minutes in the dark at room temperature (+19 to +25 °C).
5. Add 1900 μL of 1× Red Blood Cell Lysis Solution and immediately vortex thoroughly for 3 seconds. Incubate for 15 minutes in the dark at room temperature (+19 to +25 °C).
6. Centrifuge at $300 \times g$ for 5 minutes. Aspirate supernatant completely.

7. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry.
▲ **Note:** Store samples at +2 to +8 °C protected from light until analysis. Acquire samples within 1 hour after staining.
8. (Optional for fixation) Add 250 µL of buffer and 250 µL of Inside Fix to the cells and incubate for 20 minutes in the dark at room temperature (+19 to +25 °C).
9. (Optional for fixation) Add up to 2 mL of buffer.
10. (Optional for fixation) Centrifuge cells at 300×g for 5 minutes at room temperature (+19 to +25 °C). Aspirate supernatant completely.
11. (Optional) As a negative control, fluorescence-minus-one (FMO) controls can be prepared by adding a sample into a tube from the matching StainExpress FMO Set. These tubes do not contain the CD19 CAR FMC63 Idiotypic Antibody, PE, REAfinity. These controls should be treated exactly the same as the fully stained samples to ensure accurate gating of CD19 CAR⁺ T cell populations.
12. Proceed to flow cytometric analysis.

2.3 Flow cytometric data acquisition with the MACSQuant Analyzer 16

▲ Please refer to the MACSQuant Instrument user manual and software guide for detailed information on using the MACSQuant Analyzer.

▲ Please refer to the data sheet of the MACS Comp Bead Kit, anti-REA, when using beads for compensation.

1. Prepare and prime the MACSQuant Analyzer. Make sure the calibration and instrument settings of the instrument have been optimized for acquisition of the StainExpress CD19 CAR T Monitoring Cocktail.
2. For optimal compensation, prepare single stainings of suitable beads or cells with all single antibody tubes from one pouch of the matching StainExpress Compensation Set. If using beads, add one full drop of MACS Comp Beads – anti-REA and one full drop of MACS Comp Beads – blank directly to each compensation tube.
(Optional) If additional liquid markers have been added in available channels, prepare an additional tube for optimal compensation of the used antibody conjugate. For example, if using CD19 Antibody, anti-human, APC, REAfinity, prepare an additional tube with 98 µL buffer and add 2 µL of CD19 Antibody, anti-human, APC, REAfinity.
(Optional) Prepare an additional tube for optimal compensation of 7-AAD. If using beads, add 500 µL of MACSQuant Running Buffer and one drop of MACS Comp Beads - blank. Use this sample to set the compensation values for 7-AAD, ensuring precise correction for any signal spillover into the 7-AAD channel.

▲ **Note:** For optimal compensation always use the provided StainExpress Compensation Set matching the StainExpress Cocktails used. To ensure they match, check for identical order number-related lot numbers indicated on the pouches next to the expiration date.

▲ **Note:** One full drop of beads is approximately 50 µL.

▲ **Note:** For compensation on cells add 100 µL cell suspension instead of beads directly to each compensation tube.


3. Incubate for 10–20 seconds, then vortex at high speed for 5 seconds and incubate for 10 minutes in the dark at room temperature (+19 to +25 °C).
▲ **Note:** Working at lower temperatures requires increased incubation times.
4. Dilute each sample by adding 1 mL of MACSQuant Running Buffer. Mix well.
5. If using the MACS Comp Bead Kit, anti-REA, follow the protocol for compensation set up of the MACSQuant Analyzer.
▲ **Note:** For automated compensation choose the “CompensationMultiColor” Express Mode.
6. Choose appropriate voltage settings for forward scatter (FSC) and side scatter (SSC).
7. Define an appropriate threshold, based on FSC versus SSC, for the exclusion of debris from the data acquisition.
8. Start flow cytometric data acquisition.


3. Examples of immunofluorescent staining with the StainExpress CD19 CAR T Monitoring Cocktail

CAR T cells, generated within 12 days using the CliniMACS Prodigy® T Cell Transduction process, were spiked into healthy donor blood samples and were stained with the StainExpress CD19 CAR T Monitoring Cocktail.

Staining was carried out at +19 to +25 °C for 10 minutes. Cells were analyzed by flow cytometry using the MACSQuant Analyzer 16.

To exclude debris, a gate was set on FSC versus SSC encompassing lymphocytes and monocytes (A). Single cells were identified and doublets excluded by gating on the FSC-A and FSC-H plot, selecting the population with uniform scatter characteristics (B). To exclude residual erythrocytes and to identify leukocytes, CD45 was used to gate on CD45⁺ leukocytes (C). Dead cells were excluded by gating on 7-AAD⁻ cells (D). CD3⁺ cells were identified (E) and discriminated in CAR⁺ and CAR⁻ cells (F). CAR⁺ cells were further divided into CD4⁺ and CD8⁺ T cells (G), as were CAR⁻ cells (H).





Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, CliniMACS Prodigy, MACS, MACSQuant, the Miltenyi Biotec logo, REAfinity, StainExpress, Vio, VioBlue, and VioGreen are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2024 Miltenyi Biotec and/or its affiliates. All rights reserved.