Miltényi Biotec

CryoMACS® Freezing Bags

Instructions for use

CryoMACS® Freezing Bags

Instructions for use

Issued: 2026-06

35210/07

For U.S. CAUTION:

Federal (USA) law restricts this device to sale by or on the order of a physician.

Miltenyi Biotec B.V. & Co. KG

Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

Miltenyi Biotec Technical Support

☎ +49 2204 8306-3803

✉ technicalsupport@miltenyi.com

🏠 www.miltenyibiotec.com

Table of contents

en	Instructions for use	4
ar	إرشادات الاستخدام	7
bg	Инструкции за употреба	8
cs	Návod k použití	10
da	Brugsanvisning	12
de	Gebrauchsanweisung	14
el	Οδηγίες χρήσης	16
es	Instrucciones de uso	18
fi	Käyttöohje	20
fr	Notice d'utilisation	22
he	הוראות שימוש	25
hr	Upute za upotrebu	26
hu	Használati útmutató	28
it	Istruzioni per l'uso	30
lt	Naudojimo instrukcija	32
lv	Lietošanas pamācība	34
nl	Gebruiksaanwijzing	36
no	Bruksanvisning	38
pl	Instrukcja używania	40
pt	Instruções de utilização	42
ro	Instrucțiuni de utilizare	44
ru	Инструкции по использованию	46
sk	Návod na použitie	48
sl	Navodila za uporabo	50
sr	Uputstvo za upotrebu	52
sv	Bruksanvisning	54
tr	Kullanma talimatları	56

1 Contents

Each CryoMACS® Freezing Bag consists of one freezing bag with access ports as the primary containment for hematopoietic progenitor cells and one overwrap bag as a secondary containment.

Part number (P/N)	Catalog number (REF)	Product	Nominal volume	Recommended fill volume
74400	200-074-400	CryoMACS Freezing Bag 50	50 mL	10 to 30 mL
74401	200-074-401	CryoMACS Freezing Bag 250	250 mL	30 to 70 mL
74402	200-074-402	CryoMACS Freezing Bag 500	500 mL	55 to 100 mL
74403	200-074-403	CryoMACS Freezing Bag 750	750 mL	80 to 190 mL
74404	200-074-404	CryoMACS Freezing Bag 1000	1 000 mL	125 to 270 mL

The fill volume recommendations are based on freezing in a horizontal position with the use of a protective metal cassette with an internal thickness of 4 mm. If any other storage form is used, the optimal fill volume and freezing conditions must be validated by the user.

2 Intended purpose

The CryoMACS Freezing Bags are intended for a single cycle of freezing, storage (down to -196°C [-321°F]), and subsequent thawing (at $+37^{\circ}\text{C}$ [$+99^{\circ}\text{F}$]) of hematopoietic progenitor cells.

3 Precautions

 All materials that have been in contact with blood or blood components, including hematopoietic progenitor cells, must be treated and disposed as biohazardous materials. Therefore, standard hospital or institutional requirements and national legislations must be observed.

Before application to humans, the suitability of the stored hematopoietic progenitor cells must be demonstrated regarding the indication, quality, and quantity. For the use of the stored hematopoietic progenitor cells in humans the national legislation and regulations must be followed.

Any use of disinfectants on the freezing bag and the overwrap bag of the CryoMACS Freezing Bag is exclusively within the responsibility of the user and risks must be assessed within the user's own risk management.

Full performance of this product requires correct operation by the user after filling of the product according to the instructions for use. Pay special attention to use of the overwrap bag and metal cassettes as well as drying of the product and removal of residual air.

The user of the CryoMACS Freezing Bag is responsible for the selection of adequate cryopreservation protocols for the target components. Cryo-preservation using these bags should be performed by trained operators only, such as healthcare professionals in a clinical environment. Thus, CryoMACS Freezing Bags are not intended to be used by lay-persons. Handle frozen bags with care.

4 Warnings

Do not use after the use-by date indicated on the product.

Do not use if package is damaged. Use product only if package is undamaged and sealed and if port and tubing closures are in place. Use only intact bags without any signs of leakage. Use only the provided labeling pouch for identification of contents. Use only permanent sealing for tubing closure. Administer DMSO only at a concentration up to 20% (v/v) to the bag. DMSO can cause degradation of plastic materials.

Do not re-use. Re-use causes contamination of the interior fluid path. The product cannot be sterilized after use. The fluid path can be compromised by the thermal stress of re-use and the irreversible breach of the bag closures.

Do not connect the product directly to the patient.

5 Storage

Store the product at $+5^{\circ}\text{C}$ to $+35^{\circ}\text{C}$ ($+41^{\circ}\text{F}$ to $+95^{\circ}\text{F}$) until usage.

Keep away from sunlight.

6 Further information

Sterilized using irradiation.

The depicted symbol indicates the single sterile barrier of the product with protective packaging outside.

Any serious incident that has occurred in relation to this product should be reported to Miltenyi Biotec B.V. & Co. KG – using the contact information provided – and the competent authority of the member state in which the user of this product is established.

7 Handling instructions

Follow these handling instructions for the product and institutional protocols for filling, freezing (including freezing parameters), storage, and thawing. Follow aseptic working procedures for fluid transfer.

Labeling: Place an identification tag into the label pocket of the freezing bag and spot seal intermittently along the edge. Do not seal the pocket completely.

Fluid transfer (in): Use only the filling assembly for transfer of fluids into the bag. Do not use the spike ports. Use only the injection port for transfer of fluids containing DMSO. Attach transfer devices to either of the luer connectors, the injection port or by sterile docking to the PVC tubing (inner/outer diameter 2.6/4.1 mm) distal to the injection port. Observe the instructions for use supplied by the manufacturer of the sterile docking device. Carefully remove any residual air from the freezing bag after filling.

⚠ WARNING

Risk of loss of sterility and integrity. Inadequate sealing procedures may result in loss of sterility and integrity.

- Pay careful attention to the compatibility of sealing methods with ethyl vinyl acetate (EVA) material.
- Adjust dielectric or radio frequency (RF) sealing procedures accordingly.
- After filling, close the bag by sealing the EVA tubing.
- Test the seal integrity and remove the tubing distal to the seal.

Freezing: Ensure during this procedure that all surfaces of the bags and the cassette are completely dry before freezing, because ice particles on the surface may damage the bags. Place the freezing bag into the overwrap bag and vacuum seal. Place the doubled bag into a protective metal cassette. Avoid bending, folding, or kinks of the bag when closing the cassette. Seal off the tubing proximal to the freezing bag before freezing. Use of a controlled rate freezer is recommended.

Storage: Storage should be performed in the vapor phase of liquid nitrogen. Temperature fluctuations during storage should be avoided.

⚠ WARNING

Risk of loss of function and contents. If the freezing bag remains in the metal cassette for an extended period after being taken out of the nitrogen tank, it may freeze onto the metal cassette causing loss of function and contents.

- Use the overwrap bag as secondary containment during frozen storage.
- Immediately remove the freezing bag from the metal cassette after retrieval from the nitrogen tank.

Thawing: If the bag is stored in the liquid phase of nitrogen, place the bag into the vapor phase of nitrogen for a minimum of 4 hours before removal from the nitrogen tank and thawing. Remove the bag with its overwrap from the cassette and inspect for leaks or ruptures before thawing in a water bath at +37 °C (+99 °F).

⚠ WARNING

Risk of rupture and loss of content. Air or nitrogen in the freezing bag or overwrap bag can rapidly expand, causing rupture and loss of content.

- Remove air after filling and vacuum-seal the overwrap bag.

Fluid transfer (out): Completely dry the bag surface, remove the overwrap bag, twist off the top part of one of the spike ports and carefully access the bag contents by puncturing the port seal with the spike of a sterile transfusion assembly.

8 Description of components

The components of a CryoMACS Freezing Bag and the filling assembly are described in the drawing shown on the fold-out page.

CryoMACS Freezing Bag components	
No.	
1	Overwrap bag
2	Label pouch
3	Spike ports
4	Sealing area
5	EVA tubing
Filling assembly components	
No.	
6	Injection port
7	PVC tubing
8	Roller clamp
9	Female Luer lock connectors
10	Male Luer lock connector

9 Glossary of symbols

An explanation of the symbols used for the labeling of the products can be found in the glossary in the column for the respective language 'English | en' on page 58.

10 Legal notes

CryoMACS and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. Copyright © 2026 Miltenyi Biotec and/or its affiliates. All rights reserved.

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CryoMACS und das Miltenyi Biotec-Logo sind eingetragene Warenzeichen oder Warenzeichen der Miltenyi Biotec B.V. & Co. KG und/oder ihrer verbundenen Unternehmen in verschiedenen Ländern weltweit. Copyright © 2026 Miltenyi Biotec und/oder verbundene Unternehmen. Alle Rechte vorbehalten.
CryoMACS y el logotipo de Miltenyi Biotec son marcas comerciales registradas o marcas comerciales de Miltenyi Biotec B.V. & Co. KG y/o de sus filiales en distintos países del mundo. Copyright © 2026 Miltenyi Biotec y/o sus filiales. Todos los derechos reservados.
CryoMACS ja Miltenyi Biotec-logo ovat rekisteröityjä tuotemerkkejä tai Miltenyi Biotec B.V. & Co. KG:n ja/tai sen osakkuusyhtiöiden tuotemerkkejä eri maissa ympäri maailmaa. Copyright © 2026 Miltenyi Biotec ja/tai sen osakkuusyhtiöt. Kaikki oikeudet pidätetään.
CryoMACS et le logo Miltenyi Biotec sont des marques ou marques déposées de Miltenyi Biotec B.V. & Co. KG et/ou de ses filiales dans différents pays du monde. Copyright © 2026 Miltenyi Biotec et/ou ses filiales. Tous droits réservés.
CryoMACS e il logo Miltenyi Biotec sono marchi registrati o marchi depositati di Miltenyi Biotec B.V. & Co. KG e/o sue affiliate in vari paesi in tutto il mondo. Copyright © 2026 Miltenyi Biotec e/o sue affiliate. Tutti i diritti riservati.
CryoMACS en het Miltenyi Biotec-logo zijn geregistreerde handelsmerken of handelsmerken van Miltenyi Biotec B.V. & Co. KG en/of diens filialen in verschillende landen wereldwijd. Copyright © 2026 Miltenyi Biotec en/of diens filialen. Alle rechten voorbehouden.