

## Contents

1. Description
2. Background information
3. Applications
  - 3.1 General applications
  - 3.2 Specific applications
4. Instructions for use
  - 4.1 Recommended concentrations
  - 4.2 Reconstitution protocol
5. References

## 1. Description

<b>Components</b>	<b>130-100-107:</b> 200 µg lyophilized ODN 1 mL 1× TE Buffer or <b>130-100-278:</b> 1 mg lyophilized ODN 1 mL 1× TE Buffer.
<b>Description</b>	B-class control oligodeoxyribonucleotide (ODN).
<b>Product format</b>	Lyophilized product without carrier protein or preservatives.
<b>Sequence</b>	dT*dG*dC*dT*dG*dC*dT*dT*dT*dG*dT* dG*dC*dT*dT*dT*dG*dT*dG*dC*dT*dT* * Phosphorothioate backbone
<b>Endotoxin level</b>	Low endotoxin (<1 EU/mg) as determined by kinetic Limulus Amebocyte Lysate (LAL) assay.
<b>Storage</b>	Store lyophilized product at -20 °C. Upon reconstitution, aliquots should be stored at -20 °C and are stable for 6 months. Avoid repeated freeze-thaw cycles. The expiration date is indicated on the vial label.

## 2. Background information

TLR9 is a prominent member of the toll-like-receptor (TLR) family recognizing pathogen-associated molecular patterns. TLR9 recognizes specifically unmethylated CpG motifs in bacterial DNA leading to activation of immune cells<sup>1,2</sup>.

These effects can be mimicked by short synthetical ODNs containing unmethylated CpG motifs<sup>3</sup>.

Control compounds for CpG ODNs share the sequence with their stimulatory counterparts but contain GpC dinucleotides instead of the CpG motifs and thus control for the CpG specificity of the effects.

## 3. Applications

- ODN 2006 Control (ODN 2137) can be used as a sequence control for B-class ODN 2006 (# 130-100-105, # 130-100-106).

## 4. Instructions for use

### 4.1 Recommended concentrations

Recommended concentrations for use are

for human and murine immune cells: 0.05–2 µM

for recombinant cell lines: 0.05–10 µM

### 4.2 Reconstitution protocol

1. Spin down pellet.
2. a) For 200 µg lyophilized ODN:  
To obtain a 200 µM solution resuspend pellet in 130 µL of 1× TE Buffer.  
▲ **Note:** Alternatively, PBS or water can be used for reconstitution.  
b) For 1 mg lyophilized ODN:  
To obtain a 200 µM solution resuspend pellet in 650 µL of 1× TE Buffer.  
▲ **Note:** Alternatively, PBS or water can be used for reconstitution.
3. Vortex and incubate overnight at 4 °C.
4. Store aliquots at -20 °C.

## 5. References

1. Hemmi, H. *et al.* (2000) A Toll-like receptor recognizes bacterial DNA. *Nature* 408: 740–745.
2. Krieg, A. M. *et al.* (1995) CpG motifs in bacterial DNA trigger direct B-cell activation. *Nature* 374: 546–549.
3. Bauer, S. *et al.* (2001) Human TLR9 confers responsiveness to bacterial DNA via species-specific CpG motif recognition. *Proc. Natl. Acad. Sci. U.S.A.* 98: 9237–9242.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty").

Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

The Miltenyi Biotec logo is a registered trademark or trademark of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.