



Miltenyi Biotec

Miltenyi Biotec B.V. & Co. KG

Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

Phone: +49 2204 8306-0
E-Mail: macs@miltenyibiotec.de
Web: www.miltenyibiotec.com

Responsible Department: Technical Support
Phone: +49 2204 8306-3803
E-Mail: technicalsupport@miltenyi.com

Printdate: 24.11.2025

Product ID: Lung Dissociation Kit, mouse
Product code: 130-095-927

Statement

A safety datasheet is provided for the following component(s) :

120-075-999: Enzyme D 80-160 mg

If the product contains other components, the following is valid:

It is hereby stated that the respective product poses no health risk or physical and environmental hazard. It is not classified as hazardous based on our evaluations following the standards of the European Regulation (EC) No. 1272/2008 on classification, labelling and packaging of substances and mixtures. A safety data sheet is not required according to Regulation (EC) No. 1907/2006 (REACH).

Despite the classification as non-hazardous, we strongly recommend reviewing the data sheets or instructions for use as well as applying appropriate safety measures like using personal protective equipment and considering disposal legislation and local regulations when handling chemicals and reagents.

Please contact technicalsupport@miltenyi.com for any further questions or additional safety information.

Safety Data Sheet

according to Regulation (EC) No 1907/2006

Enzyme D 80-160 mg

Revision date: 25.01.2024

Product code: 120-075-999

Page 1 of 8

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Enzyme D 80-160 mg

Further trade names

Material-No. 110-008-993

Substance name: protease

CAS No: 9001-92-7

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Scientific research and development

1.3. Details of the supplier of the safety data sheet

Company name: Miltenyi Biotec B.V. & Co. KG
Street: Friedrich-Ebert-Str. 68
Place: D-51429 Bergisch Gladbach
Telephone: +49 2204 8306-0
E-mail (Contact person): technicalsupport@miltenyi.com
Internet: www.miltenyibiotec.com
Responsible Department: Technical Support
Tel.: +49 2204 8306 3803

1.4. Emergency telephone number:

GBK GmbH +49 (0)6132-84463 (24h)

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Regulation (EC) No 1272/2008

Skin Irrit. 2; H315
Eye Irrit. 2; H319
Resp. Sens. 1; H334
STOT SE 3; H335


Full text of hazard statements: see SECTION 16.

2.2. Label elements

Regulation (EC) No 1272/2008

Signal word: Danger

Pictograms:



Hazard statements

H315 Causes skin irritation.
H319 Causes serious eye irritation.
H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335 May cause respiratory irritation.

Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

Safety Data Sheet

according to Regulation (EC) No 1907/2006

Enzyme D 80-160 mg

Revision date: 25.01.2024

Product code: 120-075-999

Page 2 of 8

P302+P352 IF ON SKIN: Wash with plenty of water.
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P342+P311 If experiencing respiratory symptoms: Call a POISON CENTER/doctor.

2.3. Other hazards

none

SECTION 3: Composition/information on ingredients

3.1. Substances

Hazardous components

CAS No	Chemical name	Quantity		
	EC No	Index No	REACH No	
	Classification (Regulation (EC) No 1272/2008)			
9001-92-7	protease			100 %
	Skin Irrit. 2, Eye Irrit. 2, Resp. Sens. 1, STOT SE 3; H315 H319 H334 H335			

Full text of H and EUH statements: see section 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

After inhalation

IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a doctor if you feel unwell.

After contact with skin

Water and soap. In case of skin irritation, seek medical treatment. Change contaminated clothing.

After contact with eyes

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink plenty of water. Seek medical advice immediately. Call a doctor if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

May cause allergy or asthma symptoms or breathing difficulties if inhaled.

May cause sensitization by skin contact.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings. Water spray jet, alcohol resistant foam, Dry extinguishing powder, Carbon dioxide (CO₂)

Unsuitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixture

Non-flammable. In case of fire may be liberated: Carbon monoxide, Carbon dioxide (CO₂).

5.3. Advice for firefighters

Safety Data Sheet

according to Regulation (EC) No 1907/2006

Enzyme D 80-160 mg

Revision date: 25.01.2024

Product code: 120-075-999

Page 3 of 8

Wear a self-contained breathing apparatus and chemical protective clothing.

Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water. Move undamaged containers from immediate hazard area if it can be done safely.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General advice

Wear personal protection equipment. Do not breathe gas/vapour/aerosol. Provide adequate ventilation.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Other information

Take up mechanically, placing in appropriate containers for disposal. Wash with plenty of water.

6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Avoid contact with skin and eyes. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear personal protection equipment. Wash hands thoroughly after handling.

Advice on general occupational hygiene

Wash hands before breaks and after work. Wear personal protection equipment.

Further information on handling

Avoid dust formation. Avoid the formation of aerosol. Provide adequate ventilation.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

The following must be prevented: UV-radiation/sunlight.

Hints on joint storage

No data available

7.3. Specific end use(s)

Scientific research and development

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Additional advice on limit values

Does not contain substances above concentration limits fixing an occupational exposure limit.

8.2. Exposure controls



Safety Data Sheet

according to Regulation (EC) No 1907/2006

Enzyme D 80-160 mg

Revision date: 25.01.2024

Product code: 120-075-999

Page 4 of 8

Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust. Use extractor hood (laboratory).

If local exhaust ventilation is not possible or not sufficient, the entire working area should be ventilated by technical means.

Individual protection measures, such as personal protective equipment

Eye/face protection

Suitable eye protection: Eye glasses with side protection EN 166

Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Tested protective gloves are to be worn: Single-use gloves. German Industry Norms (DIN) / European Norms (EN): EN ISO 374

Suitable material: NBR (Nitrile rubber). ($\geq 0,11$ mm, > 480 min.)

Recommended protective gloves brand: Supreno, Manufacturer. Microfex Cooperation

Or comparable articles from other companies.

Skin protection

Personal protection equipment

Wash hands before breaks and after work.

Respiratory protection

In case of inadequate ventilation wear respiratory protection. Respiratory protection necessary at: aerosol or mist generation. (P3)

Environmental exposure controls

Do not allow to enter into surface water or drains.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state:	solid	
Colour:	No information available.	
Odour:	No information available.	
Melting point/freezing point:		not determined
Boiling point or initial boiling point and boiling range:		not determined
Flammability:		not determined not applicable
Lower explosion limits:		not determined
Upper explosion limits:		not determined
Flash point:		not determined
Decomposition temperature:		not determined
pH-Value:	No information available.	
Viscosity / kinematic:		not determined
Water solubility:		not determined
Solubility in other solvents		not determined
Partition coefficient n-octanol/water:		not determined
Vapour pressure:		not determined
Density:		not determined
Relative vapour density:		not determined
Particle characteristics:		Nanoform: none

Safety Data Sheet

according to Regulation (EC) No 1907/2006

Enzyme D 80-160 mg

Revision date: 25.01.2024

Product code: 120-075-999

Page 5 of 8

9.2. Other information

Information with regard to physical hazard classes

Self-ignition temperature

Solid: not determined

Gas: not applicable

Other safety characteristics

Evaporation rate: not determined

Viscosity / dynamic: not determined

Further Information

No information available.

SECTION 10: Stability and reactivity

10.1. Reactivity

No information available.

10.2. Chemical stability

The product is chemically stable under recommended conditions of storage, use and temperature.

10.3. Possibility of hazardous reactions

No information available.

10.4. Conditions to avoid

May cause decomposition by long-term light influence.

10.5. Incompatible materials

Oxidising agent, strong

10.6. Hazardous decomposition products

No known hazardous decomposition products.

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Toxicokinetics, metabolism and distribution

Toxicological data are not available.

Acute toxicity

Based on available data, the classification criteria are not met.

Irritation and corrosivity

Causes skin irritation.

Causes serious eye irritation.

Sensitising effects

May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

May cause respiratory irritation.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Specific effects in experiment on an animal

No information available.

Safety Data Sheet

according to Regulation (EC) No 1907/2006

Enzyme D 80-160 mg

Revision date: 25.01.2024

Product code: 120-075-999

Page 6 of 8

Practical experience

No information available.

11.2. Information on other hazards

Endocrine disrupting properties

none

SECTION 12: Ecological information

12.1. Toxicity

Toxicological data are not available.

12.2. Persistence and degradability

No information available.

12.3. Bioaccumulative potential

No information available.

12.4. Mobility in soil

No information available.

12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of REACH, annex XIII.

12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

12.7. Other adverse effects

No information available.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations

Dispose of waste according to applicable legislation.

Contaminated packaging

Dispose of waste according to applicable legislation.

SECTION 14: Transport information

Land transport (ADR/RID)

14.1. UN number or ID number:

-

14.2. UN proper shipping name:

Not a hazardous material with respect to these transportation regulations.

14.3. Transport hazard class(es):

-

14.4. Packing group:

-

Inland waterways transport (ADN)

14.1. UN number or ID number:

-

14.2. UN proper shipping name:

Not a hazardous material with respect to these transportation regulations.

14.3. Transport hazard class(es):

-

14.4. Packing group:

-

Marine transport (IMDG)

14.1. UN number or ID number:

-

14.2. UN proper shipping name:

Not a hazardous material with respect to these transportation regulations.

14.3. Transport hazard class(es):

-

14.4. Packing group:

-

Air transport (ICAO-TI/IATA-DGR)

14.1. UN number or ID number:

-

14.2. UN proper shipping name:

Not a hazardous material with respect to these transportation regulations.

Safety Data Sheet

according to Regulation (EC) No 1907/2006

Enzyme D 80-160 mg

Revision date: 25.01.2024

Product code: 120-075-999

Page 7 of 8

14.3. Transport hazard class(es): -

14.4. Packing group: -

14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

14.6. Special precautions for user

No information available.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU regulatory information

Information according to 2012/18/EU (SEVESO III): Not subject to 2012/18/EU (SEVESO III)

Additional information

Safety Data Sheet according to Regulation (EC) No. 1907/2006 (REACH)

National regulatory information

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers.

Water hazard class (D): 1 - slightly hazardous to water

15.2. Chemical safety assessment

For this substance a chemical safety assessment has not been carried out.

SECTION 16: Other information

Changes

SECTION 1: Identification of the substance/mixture and of the company/undertaking

Safety Data Sheet

according to Regulation (EC) No 1907/2006

Enzyme D 80-160 mg

Revision date: 25.01.2024

Product code: 120-075-999

Page 8 of 8

Abbreviations and acronyms

CLP: Classification, labelling and Packaging
REACH: Registration, Evaluation and Authorization of Chemicals
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals
UN: United Nations
CAS: Chemical Abstracts Service
DNEL: Derived No Effect Level
DMEL: Derived Minimal Effect Level
PNEC: Predicted No Effect Concentration
ATE: Acute toxicity estimate
LC50: Lethal concentration, 50%
LD50: Lethal dose, 50%
LL50: Lethal loading, 50%
EL50: Effect loading, 50%
EC50: Effective Concentration 50%
ErC50: Effective Concentration 50%, growth rate
NOEC: No Observed Effect Concentration
BCF: Bio-concentration factor
PBT: persistent, bioaccumulative, toxic
vPvB: very persistent, very bioaccumulative
ADR: Accord européen sur le transport des marchandises dangereuses par Route
(European Agreement concerning the International Carriage of Dangerous Goods by Road)
RID: Regulations concerning the international carriage of dangerous goods by rail
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)
IMDG: International Maritime Code for Dangerous Goods
EmS: Emergency Schedules
MFAG: Medical First Aid Guide
IATA: International Air Transport Association
ICAO: International Civil Aviation Organization
MARPOL: International Convention for the Prevention of Marine Pollution from Ships
IBC: Intermediate Bulk Container
SVHC: Substance of Very High Concern
For abbreviations and acronyms, see: ECHA Guidance on information requirements and chemical safety assessment, chapter R.20 (Table of terms and abbreviations).
Skin Irrit: Skin irritation
Eye Irrit: Eye irritation
Resp. Sens: Respiratory sensitisation
STOT SE: Specific target organ toxicity - single exposure

Relevant H and EUH statements (number and full text)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.

Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights.