



Miltenyi Biotec

# CliniMACS® Workbench

Roll scalability into your cell manufacturing concept

The CliniMACS Workbench frees the CliniMACS Prodigy® Platform from a benchtop to be positioned optimally within a cleanroom. Designed for use in a GMP environment, the CliniMACS Workbench facilitates efficient manufacturing scale-out with a small footprint, an integrated uninterruptible power supply, and storage space for accessories and consumables.

- Deploy multiple CliniMACS Prodigy Instruments in one cleanroom for easy manufacturing scale out
- Count on safe and reliable operation thanks to smart design with uninterruptible power supply
- Keep all modules, connections, and supplies on one durable, mobile unit

► [miltenyibiotec.com/CliniMACSWorkbench](https://miltenyibiotec.com/CliniMACSWorkbench)



## Compact and smart

The streamlined design concentrates everything you need to run the CliniMACS Prodigy Platform into a compact and sturdy unit. The central supply simplifies power, gas, alarm connection, and data exchange via ceiling sockets\* to keep the cleanroom floor clear of cables and supply lines. Plus, an uninterruptible power supply (UPS) safeguards workflows from costly interruptions and downtime. The CliniMACS Workbench is the perfect complement to maximize the impact of the CliniMACS Prodigy Platform.



**Figure 1:** The CliniMACS Workbench features a central supply to feed the CliniMACS Prodigy from ceiling sockets and a UPS.

## Fit for the ballroom concept

On the CliniMACS Workbench, the CliniMACS Prodigy Platform is easily integrated into your cell manufacturing cleanroom. With small footprint and sturdy wheelbase, the CliniMACS Workbench supports a ballroom manufacturing concept. Automated, closed system cell processing reduces your cleanroom requirements and facilitates parallel manufacturing of clinical cell products on multiple instruments.



**Figure 2:** Use your cleanroom space efficiently: movable manufacturing units with a small footprint.

## One platform, two workbench options

The CliniMACS Workbench is the perfect fit for the CliniMACS Prodigy and the CliniMACS Formulation Unit, allowing for easy cell processing and automated formulation of your cell products. To add non-viral cell transfection to your workflows with the CliniMACS Electroporator, the CliniMACS Workbench Large provides space for a seamless integration and instrument connection. Expand your cell manufacturing possibilities while minimizing the instrument footprint.



**Figure 3:** Two options: The CliniMACS Workbench and the CliniMACS Workbench Large.

VISIT



Create modern spaces for cell manufacturing.  
To learn more about the CliniMACS Workbench visit:

► [miltenyibiotec.com/clinimacsworkbench](https://www.miltenyibiotec.com/clinimacsworkbench)

\*Stäubli CombiTac Ceiling Socket, wired (150-001-914) must be purchased separately and installed in the facility by the user.

Miltenyi Biotec B.V. & Co. KG | Phone +49 2204 8306-0 | Fax +49 2204 85197 | [macsde@miltenyi.com](mailto:macsde@miltenyi.com) | [www.miltenyibiotec.com](https://www.miltenyibiotec.com)

Miltenyi Biotec provides products and services worldwide. Visit [www.miltenyibiotec.com/local](https://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. For further regulatory notes, refer to [miltenyibiotec.com/regulatory-notes](https://www.miltenyibiotec.com/regulatory-notes). CliniMACS, CliniMACS Prodigy, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. Copyright © 2023 Miltenyi Biotec and/or its affiliates. All rights reserved.