

Contents

1. Description
 - 1.1 Background information
 - 1.2 Applications
2. Protocol
 - 2.1 Preparation of complete medium

1. Description

This product is for research use only.

Components	500 mL StemMACS PSC-Brew Basal Medium XF, human 10 mL StemMACS PSC-Brew Supplement XF (50×), human
Specifications	Osmolality: 300–340 mOsmol/kg pH: 7.2–7.6
Quality control	Functionality assay: Maintenance of hPSCs, morphology, and surface phenotype (TRA-1-60 ⁺ , SSEA-4 ⁺) over five continuous passages. Low endotoxin level by Limulus Amoebocyte Lysate (LAL) assay. Tested negative for mycoplasma.
Storage	Store StemMACS PSC-Brew Basal Medium XF, human protected from light at 2–8 °C. Do not freeze. Upon arrival, store StemMACS PSC-Brew Supplement XF (50×), human at –20 °C. Aliquots of the supplemented complete medium can be stored at –20 °C for up to 2 months. Avoid repeated freeze-thaw cycles. Once thawed, aliquots should be kept at +2 to +8 °C and be used within 2 weeks. The expiration dates are indicated on the vial labels.

1.1 Background information

StemMACS PSC-Brew XF, human is a new generation cell culture medium for highly efficient maintenance and expansion of human pluripotent stem cells (hPSC). It has a xeno-free and serum-free formulation that ensures a reduced variability and standardized performance. The carefully balanced composition and the presence of stable FGF-2 ensure constant exposure levels and provide flexibility with feeding schedules.

StemMACS PSC-Brew XF is compatible with both defined and undefined cell attachment matrices, e.g., MACSmatrix Laminin 511 or Corning® Matrigel®. This allows the maintenance of a highly consistent pluripotent phenotype even in low density cultures. hPSCs grown in StemMACS PSC-Brew XF show a typical morphology, expression of pluripotency-associated markers and retain the ability to differentiate into cell types that derive from the three embryonic germ layers. Moreover, the maintenance of a stable karyotype in long-term cultivation is supported.

StemMACS PSC-Brew XF allows rapid culture re-initiation of hPSC cultures after cryopreservation and enables cluster passaging without the addition of any inhibitor.

StemMACS PSC-Brew XF is fully compatible with StemMACS PSC-Support XF and, when combined, can be used in challenging applications that would normally stress the cells, such as gene editing approaches or after reprogramming.

1.2 Applications

- Highly efficient maintenance and expansion of hPSCs under xeno- and feeder-free conditions.

2. Protocol

▲ Further protocols are available at <https://www.miltenyibiotec.com/applications/pluripotent-stem-cells.html>.

2.1 Preparation of complete medium

▲ Before StemMACS PSC-Brew XF can be used in cell culture, the two kit components need to be mixed according to the following protocol.

▲ Recommended feeding volumes for every day feeding are 2 mL/well of a 6-well plate. It is recommended to split the cells at a confluency of 80–90%.

1. Thaw StemMACS PSC-Brew Supplement XF (50×) at +2 to +8 °C prior to use.
2. Add 10 mL StemMACS PSC-Brew Supplement XF (50×) to 500 mL StemMACS PSC-Brew Basal Medium XF. Mix well. The medium is ready-to-use now. Use the complete medium within 2 weeks when stored at +2 to +8 °C.
3. For longer storage, prepare aliquots and store at –20 °C for up to 2 months. Thaw aliquots of complete medium overnight at +2 to +8 °C. Once thawed, keep aliquots at +2 to +8 °C and use within 2 weeks.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

The Miltenyi Biotec logo and StemMACS are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2025 Miltenyi Biotec and/or its affiliates. All rights reserved.